

## Bacto® Neisseria Meningitidis Antisera

Neisseria Meningitidis Antiserum Group A · Neisseria Meningitidis Antiserum Group B · Neisseria Meningitidis Antiserum Group C  
Neisseria Meningitidis Antiserum Group D · Neisseria Meningitidis Antiserum Group W135 · Neisseria Meningitidis Antiserum Group X · Neisseria Meningitidis Antiserum Group Y · Neisseria Meningitidis Antiserum Group Z · Neisseria Meningitidis Antiserum Group Z' · Neisseria Meningitidis Antiserum Poly (Groups A, B, C, D)  
Neisseria Meningitidis Antiserum Poly 2 (Groups X, Y, Z)

### Intended Use

Bacto Neisseria Meningitidis Antisera are used in the slide agglutination test for serotyping *Neisseria meningitidis*.

### Summary and Explanation

*Neisseria meningitidis* is found in the oropharynx and nasopharynx of humans. Because the organism survives poorly in the environment, humans are the primary reservoir. In asymptomatic persons, the carrier state lasts for various periods, usually several weeks. The microorganism is transmitted from person-to-person by direct contact with respiratory secretions or airborne droplets.<sup>1,2</sup>

In some colonized persons, the organism spreads from the nasopharynx through the bloodstream to produce meningococemia, meningitis or both. Meningococemia is characterized by a petechial or purpuric skin rash. In fulminating infections (Waterhouse-Friderichsen syndrome), widespread coagulation and fulminant sepsis occur, resulting in shock and, usually, death.<sup>3</sup> Persons with inherited complement deficiencies are at greater risk for acquiring systemic meningococcal infections and may experience repeated episodes.<sup>2</sup>

Typical human specimens for isolating the organism are cerebrospinal fluid (CSF), blood, skin lesions (in cases where petechiae occur) and nasopharyngeal swabs. *N. meningitidis* occurs in the cervix and vagina of females and can cause serious pelvic disease. Other sources for the organism are the anal canal and, in males, the urethra.

*N. meningitidis* is divided into serological groups based on the presence of either capsular or outer membrane protein antigens. Among the currently recognized groups are A, B, C, D, 29E, H, I, K, L, X, Y, Z, Z' and W135. Groups A, B, C, Y, and W135 are most frequently implicated in systemic disease.<sup>4</sup> Classically, group A and C strains cause epidemic meningococcal disease.<sup>3</sup> Group B strains have been associated with sporadic infections. Other serogroups are sporadically isolated from carriers and patients with disease.

*N. meningitidis* are gram-negative cocci, usually occurring in pairs called diplococci. They are strict aerobes and produce the enzyme, cytochrome oxidase. The growth of *N. meningitidis* is enhanced by a CO<sub>2</sub>-enriched atmosphere.

The Quellung reaction (capsular swelling) has been performed for serotyping *N. meningitidis*. However, capsules have not been demonstrated in strains of serogroup B organisms. In addition, the Quellung reaction is very nonspecific. *N. meningitidis* should be defined serologically by the slide agglutination test rather than by the Quellung reaction.

### Principles of the Procedure

Identification of *N. meningitidis* includes isolation of the microorganism, biochemical identification and serological confirmation.

### User Quality Control

#### Identity Specifications

Neisseria Meningitidis Antiserum Poly  
Neisseria Meningitidis Antiserum Poly 2  
Neisseria Meningitidis Antiserum A  
Neisseria Meningitidis Antiserum B  
Neisseria Meningitidis Antiserum C  
Neisseria Meningitidis Antiserum D  
Neisseria Meningitidis Antiserum X  
Neisseria Meningitidis Antiserum Y  
Neisseria Meningitidis Antiserum Z  
Neisseria Meningitidis Antiserum Z'  
Neisseria Meningitidis Antiserum W135

Lyophilized Appearance: Light gold to amber button to powdered cake.

Rehydrated Appearance: Light gold to amber liquid.

#### Performance Response

Rehydrate Neisseria Meningitidis Antisera per label directions. Test as described (see Test Procedure). Known positive and negative control cultures must give appropriate reactions.

Serological confirmation involves the reaction in which the microorganism (antigen) reacts with its corresponding antibody. This *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction is rapid, has at least a 3+ reaction, does not dissociate (high avidity), and binds (high affinity).

Because a microorganism (antigen) may agglutinate with an antibody produced in response to another species, heterologous reactions are possible. These are characterized as weak in strength or slow in formation. Such unexpected and, perhaps, unpredictable reactions may lead to some confusion in serological identification. Therefore, a positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism.

Homologous reactions are rapid and strong. Heterologous reactions are slow and weak.

## Reagents

Neisseria Meningitidis Antisera are lyophilized, polyclonal rabbit antisera containing approximately 0.02% Thimerosal as a preservative. Neisseria Meningitidis Antisera Poly and Group D are absorbed for detection of Group D; Neisseria Meningitidis Antisera Poly 2, Z', W135, A, B, C, X, Y and Z are not absorbed for detection of Group D, which is rarely isolated.

Neisseria Meningitidis Antisera detect the following antigenic groups:

ANTISERUM	ANTIGENIC GROUP(S) DETECTED
Poly	A, B, C, D
Poly 2	X, Y, Z
W135	W135
A	A
B	B
C	C
D	D
X	X
Y	Y
Z	Z, Z'
Z'	Z'

When rehydrated and used as described, each 1 ml vial of Neisseria Meningitidis Antiserum contains sufficient reagent for 20 slide tests.

## Precautions

1. For In Vitro Diagnostic Use.
2. The Packaging of This Product Contains Dry Natural Rubber.
3. Follow proper established laboratory procedure in handling and disposing of infectious materials.

## Storage

Store lyophilized and rehydrated Neisseria Meningitidis Antisera at 2-8°C.

## Expiration Date

The expiration date applies to the product in its intact container when stored as directed. Do not use a product if it fails to meet specifications for identity and performance.

## Procedure

### Materials Provided

Neisseria Meningitidis Antiserum Poly  
Neisseria Meningitidis Antiserum Poly 2

Neisseria Meningitidis Antiserum A  
Neisseria Meningitidis Antiserum B  
Neisseria Meningitidis Antiserum C  
Neisseria Meningitidis Antiserum D  
Neisseria Meningitidis Antiserum X  
Neisseria Meningitidis Antiserum Y  
Neisseria Meningitidis Antiserum Z  
Neisseria Meningitidis Antiserum Z'  
Neisseria Meningitidis Antiserum W135

### Materials Required But Not Provided

Agglutination slides  
Applicator sticks  
Sterile distilled or deionized water  
Sterile 0.85% NaCl solution

### Reagent Preparation

**Neisseria Meningitidis Antisera:** To rehydrate, add 1 ml sterile distilled or deionized water and rotate gently to completely dissolve the contents.

Equilibrate all materials to room temperature prior to performing the tests. Ensure that all glassware and pipettes are clean and free of detergent residues.

### Specimen Collection and Preparation

*N. meningitidis* can be recovered on blood agar or chocolate agar. Determine that the test organism has the following characteristics of *N. meningitidis*:

<i>Morphology:</i>	grayish-white, opaque, smooth, butyrous, non-pigmented
<i>Gram stain:</i>	gram-negative diplococcus
<i>Oxidase:</i>	positive
<i>Catalase:</i>	positive
<i>ONPG reaction:</i>	negative
<i>Nitrate reduction:</i>	negative
<i>Glucose:</i>	positive
<i>Maltose:</i>	positive
<i>Sucrose:</i>	negative
<i>Fructose:</i>	negative
<i>Lactose:</i>	negative

Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as *N. meningitidis*. For more detailed information on the biochemical identification of *N. meningitidis*, consult appropriate references.<sup>3,5</sup> After these criteria are met, serological identification can proceed.

### Testing the Isolate for Autoagglutination

1. From the test culture on chocolate agar, transfer a loopful of growth to a drop of sterile 0.85% NaCl solution on a clean slide and emulsify the organism.
2. Rotate the slide for one minute and then observe for agglutination. If agglutination (autoagglutination) occurs, the culture is rough and cannot be tested. Subculture to chocolate agar, incubate, and test the organism again as described in steps 1 and 2. If no agglutination occurs, proceed with testing the organism.

### Choosing Antisera to Test

1. Test the organism first with Neisseria Meningitidis Antisera Poly, Poly 2 and Group W135.
2. Depending on the reaction, continue testing as follows.
 

<b>If agglutination occurs with</b>	<b>Test with</b>
Neisseria Meningitidis Antiserum: Poly	Neisseria Meningitidis Antiserum: Groups A, B, C, D
Poly 2	Groups X, Y, Z, Z' (See NOTE.)
Group W135	No further testing is required.

NOTE: *N. meningitidis* Group Z' organisms may agglutinate monospecific Neisseria Meningitidis Antiserum Group Z. However, *N. meningitidis* Group Z organisms will not agglutinate Neisseria Meningitidis Antiserum Group Z'. The expected agglutination reactions of Neisseria Meningitidis Antiserum Groups Z' and Z with test organisms are:

Test Organism	Neisseria Meningitidis Antiserum	
	Group Z'	Group Z
<i>N. meningitidis</i> Group Z'	3+	+
<i>N. meningitidis</i> Group Z	-	3+

### Test Procedure

1. **Neisseria Meningitidis Antiserum:** Dispense 1 drop of the antiserum to be tested on an agglutination slide.
2. **Test isolate:** Transfer a loopful of growth to the drop of antiserum and mix thoroughly.
3. Rotate the slide for one minute and read for agglutination.
4. Repeat this procedure for known positive and negative cultures.

### Results

1. Read and record results as follows.
 

4+	100% agglutination; background is clear to slightly hazy.
3+	75% agglutination; background is slightly cloudy.
2+	50% agglutination; background is moderately cloudy.
1+	25% agglutination; background is cloudy.
-	No agglutination.
2. **Positive control:** Should produce 3+ or greater agglutination.  
**Negative control:** Should produce no agglutination.  
**Test isolate:** A positive test result is defined as agglutination of 3+ or greater within one minute.

### Limitations of the Procedure

1. Correct interpretation of serological reactions depends on culture purity, as well as morphological characteristics and biochemical reactions that are consistent with identification of the microorganism as *N. meningitidis*.
2. Serological methods alone cannot identify the isolate as *N. meningitidis*. Organisms unrelated to *Neisseria*, yet capable of causing meningitis, and other species of *Neisseria* can cross-react with meningococcal antisera. Cultural isolation must precede serological examination.
3. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent a smooth suspension of the microorganism or may cause evaporation or precipitation of the test mixture. False-positive reactions may occur.

4. Rough culture isolates occur and will agglutinate spontaneously, causing agglutination of the negative control (autoagglutination). Smooth colonies must be selected and tested in serological procedures.
5. *N. meningitidis* Group A and *N. meningitidis* Group C may cross-react due to the presence of common capsular polysaccharides.
6. Group Z' meningococci may agglutinate group Z antiserum. Group Z meningococci will not agglutinate group Z' antiserum.
7. Neisseria Meningitidis Antisera have been tested using undiluted cultures taken from agar media. These antisera have not been tested using antigen suspensions in NaCl solution or other diluents. If the user employs a variation of the recommended procedure, each lot of antiserum must be tested with known control cultures to verify that expected reactions are obtained under the modified procedure.
8. Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.
9. A rehydrated Neisseria Meningitidis Antiserum that is cloudy or develops a precipitate during use should be discarded.

### References

1. **Given, K. F., B. W. Thomas, and A. G. Johnston.** 1977. Isolation of *Neisseria meningitidis* from the urethra, cervix, and anal canal: further observations. *Br. J. Vener. Dis.* **53**:109-112.
2. **Janda, W. M., M. Bohnhoff, J. A. Morello, and S. A. Lerner.** 1980. Prevalence and site-pathogen studies of *Neisseria meningitidis* and *N. gonorrhoeae* in homosexual men. *JAMA* **244**:2060-2064.
3. **Knapp, J. S., and R. J. Rice.** 1995. *Neisseria* and *Branhamella*, p. 324-340. In P. R. Murray, E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.), *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
4. **Zollinger, W. D., B. L. Brandt, and E. C. Tramont.** 1986. Immune response to *Neisseria meningitidis*, p. 346-352. In N. R. Rose, H. Friedman, and J. L. Fahey (ed.), *Manual of clinical laboratory immunology*, 3rd ed. American Society for Microbiology, Washington, D.C.
5. **Pezzo, M.** 1994. Aerobic bacteriology, p. 1.0.1-1.20.47. In H. D. Isenberg (ed.), *Clinical microbiology procedures handbook*, vol. 1. American Society for Microbiology, Washington, D.C.

### Packaging

Neisseria Meningitidis Antiserum Poly	1 ml	2232-50
Neisseria Meningitidis Antiserum Poly 2	1 ml	2910-50
Neisseria Meningitidis Antiserum Group A	1 ml	2228-50
Neisseria Meningitidis Antiserum Group B	1 ml	2229-50
Neisseria Meningitidis Antiserum Group C	1 ml	2230-50
Neisseria Meningitidis Antiserum Group D	1 ml	2231-50
Neisseria Meningitidis Antiserum Group X	1 ml	2880-50
Neisseria Meningitidis Antiserum Group Y	1 ml	2881-50
Neisseria Meningitidis Antiserum Group Z	1 ml	2891-50
Neisseria Meningitidis Antiserum Group Z'	1 ml	2252-50
Neisseria Meningitidis Antiserum Group W135	1 ml	2253-50

# Bacto® Proteus Antigens and Antisera · The Weil-Felix Test

## Proteus OX2 Antigen (Slide) · Proteus OX2 Antigen (Tube)

## Proteus OX19 Antigen (Slide) · Proteus OX19 Antigen (Tube)

## Proteus OXK Antigen (Slide) · Proteus OXK Antigen (Tube)

## Proteus OX2 Antiserum · Proteus OX19 Antiserum

## Proteus OXK Antiserum · Febrile Negative Control

### Intended Use

Bacto Proteus OX2, OX19 and OXK Antigens (Slide) and (Tube) are used for detecting antibodies by the slide and tube agglutination tests.

Bacto Proteus OX2, OX19 and OXK Antisera are used in the quality control testing of Proteus OX2, OX19 and OXK Antigens in slide and tube agglutination tests.

### User Quality Control

#### Identity Specifications

**Proteus OX2 Antigen (Slide)**

**Proteus OX19 Antigen (Slide)**

**Proteus OXK Antigen (Slide)**

Appearance: Turquoise-blue-violet suspension.

**Proteus OX2 Antigen (Tube)**

**Proteus OX19 Antigen (Tube)**

**Proteus OXK Antigen (Tube)**

Appearance: Light gray to white suspension.

**Proteus OX2 Antiserum**

**Proteus OX19 Antiserum**

**Proteus OXK Antiserum**

Lyophilized Appearance: Light gold to amber, button to powdered cake.

Rehydrated Appearance: Light gold to amber, clear liquid.

#### Febrile Negative Control

Lyophilized Appearance: Colorless to light gold, button to powdered cake.

Rehydrated Appearance: Colorless to light gold, clear liquid.

#### Performance Response

Rehydrate Proteus OX Antisera and Febrile Negative Control per label directions. Perform the slide or tube agglutination test using Proteus OX Antigen (Slide) or (Tube). Both positive and negative controls are diluted in the same proportion as a patient serum and processed in the same manner following procedures for the rapid slide test or the macroscopic tube test (see Test Procedure).

An antigen is considered satisfactory if it does not agglutinate with the negative control, and if it reacts 2+ or greater at a titer of 1:160 or more with the positive control.

### Summary and Explanation

Rickettsiae cause a variety of human diseases that share symptoms such as chills, fever, malaise and myalgia. These symptoms occur suddenly, usually within 3 to 14 days after exposure. Patients frequently have a rash and may have mild pulmonary symptoms. The rickettsiae are obligate intracellular bacteria and multiply within arthropods (lice, ticks, fleas, etc.) which may serve as the vectors of infection.

The spotted fevers are caused by species of *Rickettsia*, with Rocky Mountain spotted fever caused by *Rickettsia rickettsii* being well known. Epidemic and murine typhus are caused by *R. prowazekii* and *R. typhi*, respectively. Scrub typhus is caused by *Orientalis tsutsugamushi*. For a complete discussion of the rickettsiae, consult an appropriate reference.<sup>1-5</sup>

Because rickettsial diseases develop as a febrile illness, patient diagnosis has frequently involved measurements of antibody response. The Weil-Felix test became popular in the 1920's after it was observed that certain strains of *Proteus* would agglutinate early-convalescent-phase sera from patients with suspected rickettsial disease.<sup>6</sup> *Proteus* antigens (OX2, OX19 and OXK) will cross-react in predictable patterns, although the reactions are not highly sensitive or specific. *Rickettsiae* are pathogenic microorganisms that, upon invasion, produce a fever in their host. *Proteus* antigens are often called "Febrile Antigens" because they are used to detect the response to a rickettsial infection.

The human immune response to a particular microorganism results in measurable antibody production that in some cases can help in completing the patient's clinical diagnosis. In blood samples, the antibody titer during the initial phase of the infection (acute) is compared to the antibody titer 7 to 14 days later (convalescent). Antibody titers that are high initially in the acute phase or an acute or convalescent pair of samples that shows an increase in antibody titer are helpful in the diagnosis of disease.

Diagnosis of the cause of febrile disease cannot be based solely on analysis of serum samples for antibody response. Many factors may affect measurable antibody levels. For example, the patient's immune response can be affected by age, immune status, general state of health and previous immunizations.

Certain organisms may share cross-reacting antigens leading to the production of heterologous antibodies. These heterologous antibodies may react with one or more antigens in an antibody test procedure resulting in low-level antibody titers that may not, as a single result, suggest disease. The Weil-Felix test is not specific for rickettsial diseases.

## Principles of the Procedure

Agglutination tests involving the use of *Proteus* antigens determine the presence of antibodies that react with the test antigen. The serological procedure involves serially diluting the patient serum and then adding a standard volume of an antigen. The end point of the test is the last dilution of the serum that shows a specific amount of agglutination. The end point converted to a dilution of the serum is called the patient's antibody "titer."

## Reagents

### Antigens

- Proteus Antigens** are ready to use, nonmotile strains of the organisms listed below. Proteus Antigen (Slide) contains 20% glycerin. Each vial of Proteus Antigen (Slide) contains sufficient reagent for 33 slide tests. Each vial of Proteus Antigen (Tube) contains sufficient reagent for 6 tube tests.

Proteus OX2 Antigen (Slide) and (Tube) - *Proteus vulgaris* OX2  
 Proteus OX19 Antigens (Slide) and (Tube) - *Proteus vulgaris* OX19  
 Proteus OXK Antigen (Slide) and (Tube) - *Proteus mirabilis* OXK

- Concentration of Antigen:** Antigen density may vary because it is adjusted for optimum performance when standardized with hyperimmune sera obtained from laboratory animals.

Variation in color intensity is normal and will not affect test performance.

- Proteus antigens contain the following preservative(s):  
**Proteus OX2, OX19 and OXK Antigens (Slide):** 0.5% formaldehyde, and approximately 0.002% crystal violet and 0.005% brilliant green.  
**Proteus OX2, OX19 and OXK Antigens (Tube):** 0.25% formaldehyde.

### Antisera

- Proteus Antisera** are lyophilized, polyclonal rabbit antisera containing approximately 0.04% Thimerosal as a preservative. Each vial contains sufficient reagent for 19 slide tests or 30 tube tests.
- Febrile Negative Control is a standard protein solution containing 0.02% Thimerosal as a preservative. Each vial of Febrile Negative Control contains sufficient reagent for 32 slide tests.

## Precautions

- For In Vitro Diagnostic Use.
- Observe universal blood and body fluid precautions in the handling and disposing of specimens.<sup>7,8</sup>
- Proteus OX2 Antigen (Slide)**  
**Proteus OX2 Antigen (Tube)**  
**Proteus OX19 Antigen (Slide)**  
**Proteus OX19 Antigen (Tube)**  
**Proteus OXK Antigen (Slide)**  
**Proteus OXK Antigen (Tube)**

POSSIBLE RISK OF IRREVERSIBLE EFFECTS. (US) Avoid contact with skin and eyes. Do not breathe mist. Wear suitable protective clothing. Keep container tightly closed. TARGET ORGAN(S): Eyes, Kidneys, Lungs, Skin.

**FIRST AID:** In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. If inhaled, remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Seek medical advice. If swallowed seek medical advice immediately and show this container or label.

- Proteus OX2 Antiserum**  
**Proteus OX19 Antiserum**  
**Proteus OXK Antiserum**

The Packaging of This Product Contains Dry Natural Rubber.

- Follow proper established laboratory procedure in handling and disposing of infectious materials.
- Proteus Antigens are not intended for use in the immunization of humans or animals.

## Storage

Store Proteus OX2, OX19 and OXK Antigens (Slide) and (Tube) at 2-8°C.

Store lyophilized and rehydrated Proteus OX2, OX19 and OXK Antisera at 2-8°C.

Store lyophilized and rehydrated Febrile Negative Control at 2-8°C.

## Expiration Date

The expiration date applies to the product in its intact container when stored as directed. Do not use a product if it fails to meet specifications for identity and performance.

## Procedure

### Materials Provided

Proteus OX2 Antigen (Slide)  
 Proteus OX2 Antigen (Tube)  
 Proteus OX19 Antigen (Slide)  
 Proteus OX19 Antigen (Tube)  
 Proteus OXK Antigen (Slide)  
 Proteus OXK Antigen (Tube)  
 Proteus OX2 Antiserum  
 Proteus OX19 Antiserum  
 Proteus OXK Antiserum  
 Febrile Negative Control

### Materials Required But Not Provided

#### Slide test

Agglutination slides, 5 squares, 1" each  
 Applicator sticks  
 Sterile 0.85% NaCl solution  
 Sterile distilled or deionized water  
 Serological pipettes, 0.2 ml

#### Tube Test

Culture tubes 12 x 75 mm and rack  
 Waterbath, 35-37°C  
 Serological pipettes, 1 ml and 5 ml  
 Sterile 0.85% NaCl solution  
 Sterile distilled or deionized water

### Reagent Preparation

Proteus OX2, OX19 and OXK Antigens (Slide) and (Tube) are ready to use.

Equilibrate all materials to room temperature before performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergent.

**Proteus OX2, OX19 and OXK Antisera:** To rehydrate, add 3 ml sterile 0.85% NaCl solution and rotate gently to completely dissolve the contents. The rehydrated antiserum is considered a 1:2 working dilution.

**Febrile Negative Control:** To rehydrate, add 5 ml sterile distilled or deionized water and rotate gently to completely dissolve the contents.

### Specimen Collection and Preparation

Collect a blood specimen by aseptic venipuncture. Serum is required for the test. Store serum specimens at room temperature for no longer than 4 hours; for prolonged storage, keep at 2-8°C for up to 5 days or maintain below -20°C.

Serum specimens must be clear, free of hemolysis and show no visible evidence of bacterial contamination (turbidity, hemolysis or particulate matter). Consult appropriate references for more information on collection of specimens.<sup>9,10</sup> Serum specimens must not be heated. Heat may inactivate or destroy certain antibodies.

### Test Procedure

#### Slide Test

Use the slide test only as a screening test. Confirm positive results with the tube test.

1. **Test serum:** Using a 0.2 ml serological pipette, dispense 0.08, 0.04, 0.02, 0.01 and 0.005 ml of serum into a row of squares on the agglutination slide.
2. **Positive control:** Using a 0.2 ml serological pipette, dispense 0.08, 0.04, 0.02, 0.01 and 0.005 ml of Proteus Antiserum into a row of squares on the agglutination slide.
3. **Negative control:** Using a 0.2 ml serological pipette, dispense 0.08, 0.04, 0.02, 0.01 and 0.005 ml of Febrile Negative Control into a row of squares on the agglutination slide.
4. **Proteus Antigen:** Shake the vial of antigen well to ensure a smooth, uniform suspension. Add one drop (approximately 35 µl) of antigen to each drop of diluted test serum, positive control and negative control.
5. Mix each row of test sera and control sera, using a separate applicator stick for each row. Start with the most dilute mixture (0.005 ml) and work to the most concentrated (0.08 ml).
6. Rotate the slide for 1 minute and read for agglutination.
7. The final dilutions in squares 1-5 correspond with tube dilutions of 1:20, 1:40, 1:80, 1:160, 1:320, respectively.

#### Results

1. Read and record results as follows.
  - 4+ 100% agglutination; background is clear to slightly hazy.
  - 3+ 75% agglutination; background is slightly cloudy.
  - 2+ 50% agglutination; background is moderately cloudy.
  - 1+ 25% agglutination; background is cloudy.
  - No agglutination.

2. **Positive control:** Should show 2+ or greater agglutination at 1:160.
3. **Negative control:** Should show no agglutination.
4. If results for either the positive or negative controls are not as described, the test is invalid and results cannot be read.
5. **Test specimens:** The serum titer is that dilution which shows 2+ or greater agglutination.
6. The slide test is a screening test, only; results must be confirmed with the tube test.

#### Tube Test

1. Prepare a row of 8 culture tubes (12 x 75 ml) for each test serum, including a row for the appropriate Proteus Antiserum.
2. **Sterile 0.85% NaCl solution:** Dispense 0.9 ml in the first tube of each row and 0.5 ml in the remaining tubes.
3. **Test serum:** Using a 1 ml serological pipette, add 0.1 ml of the serum in the first tube in the row and mix thoroughly. Transfer 0.5 ml from tube 1 to tube 2 and mix thoroughly. In like manner, continue transferring 0.5 ml through tube 7, discarding 0.5 ml from tube 7 after mixing. Tube 8 is the antigen control tube and contains only sterile 0.85% NaCl solution.
4. **Positive control:** Using a 1 ml serological pipette, add 0.1 ml of the appropriate Proteus Antiserum to the first tube in the row and mix thoroughly. Transfer 0.5 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.5 ml through tube 7, discarding 0.5 ml from tube 7 after mixing. Tube 8 is the antigen control tube and contains only sterile 0.85% NaCl solution.
5. **Proteus Antigen:** Shake the vial of antigen to ensure a smooth, uniform suspension. Add 0.5 ml of the antigen to each of the 8 tubes in each row and shake the rack to mix the suspensions. Final dilutions in tubes 1-7 are 1:20, 1:40, 1:80, 1:160, 1:320, 1:640 and 1:1280, respectively.
6. Incubate in a waterbath at 35-37°C for 2 hours; then refrigerate at 2-8°C for 22 ± 2 hours.
7. Remove from incubation. Avoid excessive shaking before reading the reactions either when the tubes are incubated or when removing them from the incubation.
8. Read and record the results.

#### Results

1. Read and record results as follows.
  - 4+ 100% agglutination; background is clear to slightly hazy.
  - 3+ 75% agglutination; background is slightly cloudy.
  - 2+ 50% agglutination; background is moderately cloudy.
  - 1+ 25% agglutination; background is cloudy.
  - No agglutination.
2. **Positive control:** Should show a 2+ or greater agglutination at 1:160.
3. **Antigen control (tube 8 of each row):** Should show no agglutination.
4. If results of the positive control or antigen control are not as described, the test is invalid and results cannot be read.
5. For each test serum, the serum titer is that dilution which shows 2+ or greater agglutination.

**Interpretation<sup>1</sup>**

Compare results:

DISEASE	AGENT	PROTEUS OX2	PROTEUS OX19	PROTEUS OXK
Epidemic Typhus*	<i>R. prowazekii</i>	+	+	-
Murine Typhus*	<i>R. typhi</i>	+	+	-
Scrub Typhus	<i>O. tsutsugamushi</i>	-	-	+
Rocky Mountain Spotted Fever**	<i>R. rickettsii</i>	+	+	-
Other Spotted Fevers**	<i>Rickettsia</i> sp.	+	+	-

\*In cases of epidemic and murine typhus, the strength of the antibody agglutination with Proteus OX19 is usually stronger (4+) than the agglutination with Proteus OX2 (2+).

\*\*In cases of spotted fevers, antibodies may agglutinate either or both strains of Proteus OX19 or OX2, and the strength of agglutination may vary from 1+ to 4+.

For a single serum specimen, a titer of 1:160 is suggestive of infection.

A pair of serum specimens (acute and convalescent) showing a two-dilution difference in the titers is a significant increase in antibody level and is suggestive of infection. A one dilution difference is within the limits of laboratory error.

**Limitations of the Procedure**

- The slide test is for screening only and results should be confirmed by performing the tube test. The slide test dilutions are made to detect a prozone reaction and do not represent true quantitation of the antibody. A serum specimen with a prozone reaction shows no agglutination because of excessively high antibody concentrations. To avoid this occurrence, all 5 serum dilutions (slide test) should be run.
- The detection of antibodies in serum specimens may complete the clinical picture of a patient having infection. However, the isolation of the causative agent from patient specimens may be required. A definitive diagnosis must be made by a physician based on patient history, physical examination and data from all laboratory tests.
- Cross-reacting heterologous antibodies are responsible for many low titer reactions. Infections with other organisms, vaccinations and past history of disease may result in low level of antibody titers. Antimicrobial therapy may suppress antibody production.  
The Weil-Felix test is not specific for rickettsial diseases. *Rickettsia* species cause cross-reacting antibodies, and infections with *Proteus* species can also cause cross-reacting antibodies.
- While a single serum specimen showing a titer of 1:160 suggests infection, it is not diagnostic.
- To test for a significant rise in antibody titer, at least two specimens are necessary, an acute specimen obtained at time of initial symptoms and a convalescent specimen obtained 7 to 14 days later. A two-dilution difference in titer is a significant increase in antibody level and is suggestive of infection.
- The Weil-Felix test does not differentiate between epidemic and murine typhus.
- Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.
- Exposure of the antigen reagents to temperatures below 2°C can result in autoagglutination. Antigens must be smooth, uniform suspensions. Examine antigen vials for agglutination before use. Suspensions with agglutination are not usable and should be discarded.

- Rehydrated Proteus OX2, OX19 and OXK Antisera that is cloudy or has a precipitate during use should be discarded.

**References**

- Eisemann, C. S., and J. V. Osterman.** 1986. *Rickettsiae*, p. 593-599. In N. R. Rose, H. Friedman, and J. L. Fahey, (ed.), Manual of clinical laboratory immunology, 3rd ed. American Society for Microbiology, Washington, D.C.
- McDade, J. E.** 1991. *Rickettsiae*, p. 1036-1044. In A. Balows (ed.), Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
- Olson, J. G., and J. E. McDade.** 1995. *Rickettsia* and *Coxiella*, p. 678-685. In P. R. Murray, E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
- Miller, L. E., H. R. Ludke, J. E. Peacock, and R. H. Tomar.** 1991. Manual of laboratory immunology, 2nd ed. Lea & Febiger.
- Turgeon, M. L.** 1990. Immunology and serology in laboratory medicine. The C. V. Mosby Company, St. Louis, MO.
- Weil, E., and A. Felix.** 1916. Zur serologischen Diagnosis des Fleckfiebers. Wien. Klin. Wochenschr. **29**:33-35.
- Centers for Disease Control.** 1988. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. Morbidity and Mortality Weekly Reports **37**:377-382, 387-388.
- Occupational Safety and Health Administration, U.S. Department of Labor.** 1991. 29 CFR, part 1910. Occupational exposure to bloodborne pathogens; final rule. Federal Register **56**:64175-64182.
- Pezlo, M.** 1992. Aerobic bacteriology, p. 1.0.1-1.20.47. In H. D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
- Miller, J. M., and H. T. Holmes.** 1995. Specimen collection, transport and storage. In P. R. Murray, E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.

**Packaging**

Proteus OX2 Antigen (Slide)	5 ml	2243-56
Proteus OX2 Antigen (Tube)	25 ml	2248-65
Proteus OX19 Antigen (Slide)	5 ml	2234-56
Proteus OX19 Antigen (Tube)	25 ml	2247-65
Proteus OXK Antigen (Slide)	5 ml	2244-56
Proteus OXK Antigen (Tube)	25 ml	2249-65
Proteus OX2 Antiserum	3 ml	2245-47
Proteus OX19 Antiserum	3 ml	2235-47
Proteus OXK Antiserum	3 ml	2246-47
Febrile Negative Control	5 ml	3239-56

## Bacto® QC Antigens Salmonella

QC Antigen Salmonella O Group A · QC Antigen Salmonella O Group B · QC Antigen Salmonella O Group C<sub>1</sub> · QC Antigen Salmonella O Group C<sub>2</sub> · QC Antigen Salmonella O Group D · QC Antigen Salmonella O Group E<sub>1</sub> · QC Antigen Salmonella O Group E<sub>2</sub> · QC Antigen Salmonella O Group E<sub>4</sub> · QC Antigen Salmonella O Group F · QC Antigen Salmonella O Group G<sub>1</sub> · QC Antigen Salmonella O Group H · QC Antigen Salmonella O Group I · QC Antigen Salmonella Vi  
 Febrile Negative Control

### Intended Use

Bacto QC Antigens Salmonella are used in the quality control testing of Salmonella Antisera by the slide agglutination test.

### User Quality Control

#### Identity Specifications

**QC Antigens Salmonella O Groups A, B, C<sub>1</sub>, C<sub>2</sub>, D, E<sub>1</sub>, E<sub>2</sub>, E<sub>4</sub>, F, G<sub>1</sub>, H, I and Salmonella Vi**

Appearance: Liquid, light gray to white, may settle on standing.

#### Febrile Negative Control

Lyophilized appearance: Colorless to light gold, button to powdered cake.

Rehydrated appearance: Colorless to light gold, clear liquid.

#### Cultural Response

Rehydrate the Salmonella antiserum per label directions. Perform the slide agglutination test using an appropriate QC Antigen Salmonella as the homologous (positive) or heterologous (negative) control. The homologous control should produce 3+ or greater agglutination. The heterologous control should not produce agglutination. Infrequently, a +/- reaction will occur.

The following chart lists the identifying antigen(s) of various Salmonella Antisera and the recommended homologous QC Antigen(s) Salmonella (positive control). To demonstrate a heterologous (negative control) reaction, use a QC Antigen Salmonella that contains antigens unrelated to those in the homologous control.

*continued on following page*

### Summary and Explanation

*Salmonella* species cause a variety of human diseases called salmonellosis. The range of disease is from mild self-limiting gastroenteritis to a more severe form, possibly with bacteremia to typhoid fever, which can be severe and life-threatening. Severe disease and bacteremia are associated primarily with *S. choleraesuis*, *S. paratyphi* A and *S. typhi*, while most of the other 2300 or more strains are associated with gastroenteritis. The severity of the diarrheal disease depends on the virulence of the strain and the condition of the human host.

Salmonellae are found in nature and occur in the intestinal tract of many animals, both wild and domestic. The microorganism can spread to man through environmental contact or from eating contaminated meat or vegetable food products.

The genus *Salmonella* is in the family *Enterobacteriaceae*. Salmonellae are facultatively anaerobic, gram-negative bacilli that typically are oxidase negative, lactose negative, H<sub>2</sub>S positive and produce gas.

Serotypes of *Salmonella* are defined based on the antigenic structure of both the somatic or cell wall (O) antigens and the flagellar (H) antigens. The antigenic formula provides the O antigen(s) first, followed by the H antigen(s). In characterizing serotypes of *Salmonella*, the somatic O heat-stable antigens are identified first and are numbered 1-67 using Arabic numerals. The numbers are not completely continuous because certain strains were reclassified to other genera and the antigenic Arabic numbers were deleted from the scheme.

Serogroups, which represent the organization of the *Salmonella* strains based on the antigen(s) shared in common, are designated by the letters A-Z. After exhausting the alphabet, the serogroups were numbered beginning with the numeral 51 (the serogroup Z organism having antigen number 50). While one somatic antigen identifies each serogroup, certain other antigens may be shared among several serogroups.

The use of Salmonella antisera in the serological identification of *Salmonella* requires the use of quality control test suspensions to verify that the antisera are performing as expected. Most laboratories are required to test antisera with positive and negative controls prior to use.<sup>1,2</sup> QC Antigens Salmonella are designed as homologous controls for testing the efficacy of the *Salmonella* grouping antisera employed in routine laboratory procedures.

## Principles of the Procedure

Serological procedures that confirm the identification of an organism are usually agglutination reactions. Agglutination reactions may be either homologous or heterologous. Homologous reactions occur between a microorganism (antigen) and the corresponding antibody. These reactions occur rapidly and are strong. Heterologous reactions occur when a microorganism (antigen) reacts with an antibody

produced in response to some other species or serotype. These reactions occur slowly and are weak.

Heterologous reactions may be unexpected and unpredictable and may lead to confusion in serological identification. Therefore, only strongly positive homologous agglutination reactions should be regarded as significant.

## Reagents

QC ANTIGEN SALMONELLA	ORGANISM USED FOR ANTIGEN PREPARATION	HOMOLOGOUS IDENTIFYING ANTIGEN(S)
O Group A	<i>S. paratyphi</i> A var. Durazzo, Factors 2, 12	2
O Group B	<i>S. typhimurium</i> Factors <u>1</u> , 4, [5], 12	4, 5
O Group C <sub>1</sub>	<i>S. choleraesuis</i> factors 6, 7	7
O Group C <sub>2</sub>	<i>S. newport</i> factors 6, 8	8
O Group D	<i>S. gallinarum</i> factors <u>1</u> , 9, 12	9
O Group E <sub>1</sub>	<i>S. anatum</i> factors 3, 10	10
O Group E <sub>2</sub>	<i>S. newington</i> factors 3, <u>15</u>	15
O Group E <sub>4</sub>	<i>S. senftenberg</i> factors 1, 3, 19	19
O Group F	<i>S. rubislaw</i> factor 11	11
O Group G <sub>1</sub>	<i>S. poona</i> factors [1], 13, 22, [36], [37]	22
O Group H	<i>S. carrau</i> factors 6, 14, 24	14
O Group I	<i>S. hvittingfoss</i> factor 16	16
Vi	<i>Citrobacter ballerup</i> O29	Vi

Note: Brackets [ ] indicate that the antigen may be absent.

Underlining indicates that the O antigen has been lysogenized in that strain.

These antigen suspensions are ready to use. QC Antigens Salmonella O are preserved with 0.5% phenol USP; QC Antigen Salmonella Vi contains 0.01% Thimerosal. When used as described, each vial of QC Antigen Salmonella has sufficient reagent for 20 slide tests.

Febrile Negative Control is a lyophilized standard protein solution, containing approximately 0.04% Thimerosal as a preservative. When used as described, each vial of Febrile Negative Control has sufficient reagent for 100 slide tests.

## Precautions

- For In Vitro Diagnostic use.
- QC Antigens Salmonella**  
The Packaging of This Product Contains Dry Natural Rubber.
- Follow proper established laboratory procedure in handling and disposing of infectious materials.
- QC Antigens Salmonella are not to be used for immunization of humans or animals.

## Storage

Store QC Antigens Salmonella at 2-8°C.

Store lyophilized and rehydrated Febrile Negative Control at 2-8°C.

## Expiration Date

The expiration date applies to the product in its intact container when stored as directed. Do not use a product if it fails to meet specifications for identity and performance.

### User Quality Control cont.

SALMONELLA ANTISERUM	QC ANTIGEN SALMONELLA HOMOLOGOUS CONTROL(S)
Poly A-I & Vi	A, B, D, E <sub>1</sub> , E <sub>2</sub> , E <sub>4</sub> , F, G <sub>1</sub> , H, I, Vi
Poly A	A, B, D, E <sub>1</sub> , E <sub>2</sub> , E <sub>4</sub>
Poly B	C <sub>1</sub> , C <sub>2</sub> , F, G <sub>1</sub> , H
Poly C	I
Group A Factors 1, 2, 12	A
Group B Factors 1, 4, 5, 12	B
Group B Factors 1, 4, 12, 27	B
Group C <sub>1</sub> Factors 6, 7	C <sub>1</sub>
Group C <sub>2</sub> Factors 6, 8	C <sub>2</sub>
Group D <sub>1</sub> Factors 1, 9, 12	D
Group E Factors 1, 3, 10, 15, 19, 34	E <sub>1</sub> , E <sub>2</sub> , E <sub>4</sub>
Group E <sub>1</sub> Factors 3, 10	E <sub>1</sub>
Group E <sub>2</sub> Factors 3, 15	E <sub>2</sub>
Group E <sub>4</sub> Factors 1, 3, 19	E <sub>4</sub>
Group F Factor 11	F
Group G Factors 13, 22, 23, (36), (37)	G <sub>1</sub>
Group G <sub>1</sub> Factors 13, 22, (36), (37)	G <sub>1</sub>
Group H Factors 1, 6, 14, 24, 25	H
Group I Factor 16	I
Vi	Vi
Factor 2	A
Factor 4	B
Factors 4, 5	B
Factor 5	B
Factor 7	C <sub>1</sub>
Factor 8	C <sub>2</sub>
Factor 9	D
Factor 10	E <sub>1</sub>
Factor 15	E <sub>2</sub>
Factor 19	E <sub>4</sub>
Factor 22	G <sub>1</sub>
Factor 14	H

Note: Parentheses ( ) indicate that the antigen is poorly developed or agglutinates weakly. For a complete and current explanation of the classification of *Salmonella*, consult appropriate references.<sup>3,4,5</sup>

## Procedure

### Materials Provided

QC Antigens Salmonella  
Febrile Negative Control

### Materials Required But Not Provided

Agglutination slides  
Applicator sticks  
Sterile distilled or deionized water

### Reagent Preparation

QC Antigens Salmonella are ready to use.

Equilibrate all materials to room temperature before performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergents.

**Febrile Negative Control:** To rehydrate, add 5 ml sterile distilled or deionized water and rotate gently to completely dissolve the contents.

Before using QC Antigens Salmonella, examine the Salmonella antisera (Poly, Group or Factor) chosen for use. The antisera must meet all product specifications.

### Test Procedure

1. **Positive control:** Dispense 1 drop (35 Fl) of the Salmonella Antiserum to be tested on an agglutination slide. Add 1 drop of the QC Antigen Salmonella chosen as the positive control and mix thoroughly.
2. **Negative control:** Dispense 1 drop of Febrile Negative Control on the agglutination slide. Add 1 drop of the QC Antigen Salmonella chosen as the positive control and mix thoroughly.
3. Rotate the slide for 1 minute and read for agglutination. Results must be read within 1 minute.

### Results

1. Read and record results as follows.
  - 4+ 100% agglutination; background is clear to slightly hazy.
  - 3+ 75% agglutination; background is slightly cloudy.
  - 2+ 50% agglutination; background is moderately cloudy.
  - 1+ 25% agglutination; background is cloudy.
  - No agglutination.
2. **Positive control:** Should show 3+ or greater agglutination.
3. **Negative control:** Should show no agglutination. Rarely, a +/- reaction is possible.

### Limitations of the Procedure

1. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False-positive reactions may occur.
2. QC Antigens Salmonella will react with their corresponding homologous Salmonella polyvalent, group or factor antiserum. Some single factor antisera may give weaker reactions than polyvalent or grouping antisera due to specificity of the single factor antiserum for the identifying antigen(s).

3. The density of QC Antigens Salmonella is adjusted so that they give negative reactions with heterologous Salmonella Group Antisera. However, an antiserum may have a particular avidity for a given factor and agglutinate with another QC Antigen Salmonella having the common factor. For example, Salmonella O Group E<sub>1</sub> Antiserum Factors 3, 10 will agglutinate QC Antigen Salmonella O Group E<sub>1</sub> (3, 10) but may also agglutinate QC Antigen Salmonella O Group E<sub>2</sub> (3, 15) if the antiserum is of high titer for factor 3.
4. Exposure to temperatures below 2°C can result in autoagglutination. Antigens must be smooth, uniform suspensions. Examine the antigen vial for agglutination before use. Suspensions with agglutination are not usable and should be discarded.

### References

1. Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D. C.
2. Isenberg, H. D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D. C.
3. Ewing, W. H. (ed.). 1986. Edwards and Ewing's identification of Enterobacteriaceae, 4th ed. Elsevier Science Publishing Co., Inc., New York, NY.
4. McWhorter-Murlin, A. C., and F. W. Hickman-Brenner. 1994. Identification and serotyping of *Salmonella* and an update of the Kauffmann-White Scheme. U. S. Dept. Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, Atlanta, GA.
5. Popoff, M. Y., and L. LeMinor. 1997. Antigenic formulas of the *Salmonella* serovars. WHO Collaborating Centre for Reference and Research on *Salmonella*. Institut Pasteur, Paris, France.

### Packaging

QC Antigen Salmonella O Group A	1 ml	2130-50
QC Antigen Salmonella O Group B	1 ml	2131-50
QC Antigen Salmonella O Group C <sub>1</sub>	1 ml	2132-50
QC Antigen Salmonella O Group C <sub>2</sub>	1 ml	2133-50
QC Antigen Salmonella O Group D	1 ml	2134-50
QC Antigen Salmonella O Group E <sub>1</sub>	1 ml	2135-50
QC Antigen Salmonella O Group E <sub>2</sub>	1 ml	2136-50
QC Antigen Salmonella O Group E <sub>4</sub>	1 ml	2137-50
QC Antigen Salmonella O Group F	1 ml	2138-50
QC Antigen Salmonella O Group G <sub>1</sub>	1 ml	2139-50
QC Antigen Salmonella O Group H	1 ml	2140-50
QC Antigen Salmonella O Group I	1 ml	2141-50
QC Antigen Salmonella Vi	1 ml	2142-50
Febrile Negative Control	5 ml	3239-56

## Bacto® QC Antigens Shigella

### QC Antigen Shigella Group A · QC Antigen Shigella Group A<sub>1</sub> QC Antigen Shigella Group B · QC Antigen Shigella Group C QC Antigen Shigella Group C<sub>1</sub> · QC Antigen Shigella Group C<sub>2</sub> QC Antigen Shigella Group D · QC Antigen Alkalescens-Dispar Group 1

#### Intended Use

Bacto QC Antigens Shigella and Bacto QC Antigen Alkalescens-Dispar Group 1 are used in the quality control testing of Shigella Antisera Poly and Alkalescens-Dispar Antiserum Poly by the slide agglutination test.

#### User Quality Control

##### Identity Specifications

###### QC Antigens Shigella

Appearance: Light gray to white suspension.

###### QC Antigen Alkalescens-Dispar Group 1

Appearance: Light gray to white suspension.

###### Febrile Negative Control

Lyophilized appearance: Colorless to light gold, button to powdered cake.

Rehydrated appearance: Colorless to light gold, clear liquid.

##### Performance Response

Rehydrate Shigella Antiserum Poly and Alkalescens-Dispar Antiserum Poly per label directions. Perform the slide agglutination test using an appropriate QC Antigens Shigella or Alkalescens-Dispar Group 1.

The following chart lists the QC Antigens Shigella or QC Antigen Alkalescens-Dispar Group 1 recommended as the homologous (positive) control antigen. The homologous control antigen has certain identifying antigen(s) in common with the antiserum.

ANTISERUM	QC ANTIGEN HOMOLOGOUS CONTROL
Shigella Antiserum Poly Group A	Shigella Group A
Shigella Antiserum Poly Group A <sub>1</sub>	Shigella Group A <sub>1</sub>
Shigella Antiserum Poly Group B	Shigella Group B
Shigella Antiserum Poly Group C	Shigella Group C
Shigella Antiserum Poly Group C <sub>1</sub>	Shigella Group C <sub>1</sub>
Shigella Antiserum Poly Group C <sub>2</sub>	Shigella Group C <sub>2</sub>
Shigella Antiserum Poly Group D	Shigella Group D
Alkalescens-Dispar Antiserum Poly	Alkalescens-Dispar Group 1

#### Summary and Explanation

*Shigella* species cause the diarrheal disease known as shigellosis (classic bacillary dysentery) in humans. The range of illness is from mild diarrhea to severe dysentery characterized by abdominal cramps and frequent passage of bloody, mucoid stools. While the disease is usually self-limiting, it can be life threatening to the young, the elderly and malnourished persons. *Shigella* species are carried primarily in humans and are not generally distributed in nature. While transmission is usually direct person-to-person or through contaminated water supplies, food borne outbreaks do occur.

The genus *Shigella* belongs to the family Enterobacteriaceae. *Shigella* are facultatively anaerobic, gram-negative bacilli that typically are oxidase negative, lactose negative, H<sub>2</sub>S negative and do not produce gas. *Shigella* and *Escherichia* are genetically related. Certain strains of *E. coli* may resemble *Shigella* biochemically because both can be lactose negative, nonmotile and non-gas-producing. These anaerogenic, nonmotile types have historically been called the Alkalescens-Dispar group and are presently classified as *E. coli*.<sup>1-6</sup>

Serological testing with polyvalent and group specific antisera should be used to confirm the identification of isolates that are morphologically and biochemically identified as *Shigella* species. *Shigella* are nonmotile, so serological identification is based on somatic ("O") antigens. However, some strains have envelope antigens that prevent agglutination in somatic antisera. Heating the suspension at 100°C for 15-60 minutes destroys these interfering antigens. The four named species or serotypes of *Shigella* are *S. dysenteriae* (10 serovars), *S. flexneri* (six serovars), *S. boydii* (15 serovars) and *S. sonnei*. For a complete and current explanation of the classification of *Shigella*, consult appropriate references.<sup>1</sup>

Shigella Antisera Poly and Alkalescens-Dispar Antiserum Poly are used in the serological identification of *Shigella* species and the Alkalescens-Dispar (A-D) Group. QC Antigens Shigella and Alkalescens-Dispar Group 1 are designed as positive controls for testing the efficacy of the *Shigella* grouping antisera used in laboratory procedures.

QC Antigens Shigella and QC Antigen Alkalescens-Dispar Group 1 may also be used as negative controls by using a heterologous antigen (possessing no common antigen) with a given test serum. However, cross reactivity may occur. Consult appropriate references for further details on cross reactivity.<sup>1</sup>

## Principles of the Procedure

Serological procedures that confirm the identification of an organism are usually agglutination reactions. Agglutination reactions may be either homologous or heterologous. Homologous reactions occur between a microorganism (antigen) and the corresponding antibody. These reactions occur rapidly and are strong. Heterologous reactions occur when a microorganism (antigen) reacts with an antibody produced in response to another species or serotype. These reactions occur slowly and are weak.

Heterologous reactions may be unexpected and unpredictable and may lead to confusion in serological identification. Therefore, only strongly positive homologous agglutination reactions should be regarded as significant.

## Reagents

QC ANTIGEN	ORGANISM IDENTITY
Shigella Group A	<i>Shigella dysenteriae</i> type 1
Shigella Group A <sub>1</sub>	<i>Shigella dysenteriae</i> type 8
Shigella Group B	<i>Shigella flexneri</i> type 6
Shigella Group C	<i>Shigella boydii</i> type 3
Shigella Group C <sub>1</sub>	<i>Shigella boydii</i> type 8
Shigella Group C <sub>2</sub>	<i>Shigella boydii</i> type 12
Shigella Group D	<i>Shigella sonnei</i>
Alkalescens-Dispar Group 1	<i>E. coli</i> A-D Group 1

QC Antigens Shigella and QC Antigen Alkalescens-Dispar Group 1 contain killed whole organisms preserved in 0.5% formaldehyde. They are ready to use.

When used as described, each vial of antigen contains sufficient reagent for 20 slide tests.

## Precautions

- For In Vitro Diagnostic Use.
- QC Antigen Shigella Group A**  
**QC Antigen Shigella Group A<sub>1</sub>**  
**QC Antigen Shigella Group B**  
**QC Antigen Shigella Group C**  
**QC Antigen Shigella Group C<sub>1</sub>**  
**QC Antigen Shigella Group C<sub>2</sub>**  
**QC Antigen Shigella Group D**  
**QC Antigen Alkalescens-Dispar Group 1**

POSSIBLE RISK OF IRREVERSIBLE EFFECTS. (US) Avoid contact with skin and eyes. Do not breathe mist. Wear suitable protective clothing. Keep container tightly closed. TARGET ORGAN(S): Eyes, Kidneys, Lungs, Skin.

FIRST AID: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. If inhaled, remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Seek medical advice. If swallowed seek medical advice immediately and show this container or label.

The Packaging of This Product Contains Dry Natural Rubber.

- Follow proper established laboratory procedure in handling and disposing of infectious materials.

- QC Antigens Shigella and QC Antigen Alkalescens-Dispar Group 1 are not to be used for immunization of humans or animals.

## Storage

Store QC Antigens Shigella and QC Antigen Alkalescens-Dispar Group 1 at 2-8°C. Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

## Expiration Date

The expiration date applies to the product in its intact container when stored as directed. Do not use a product if it fails to meet specifications for identity and performance.

## Procedure

### Materials Provided

QC Antigens Shigella  
QC Antigen Alkalescens-Dispar Group 1

### Materials Required But Not Provided

Febrile Negative Control  
Agglutination slides  
Applicator sticks  
Sterile 0.85% NaCl solution

### Reagent Preparation

QC Antigens Shigella and QC Antigen Alkalescens-Dispar Group 1 are ready to use.

Equilibrate all materials to room temperature prior to performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergents.

### Test Procedure

- Positive control:** Dispense 1 drop (35 µl) of the Shigella Antiserum or Alkalescens-Dispar Antiserum Poly to be tested on an agglutination slide. Add 1 drop of the appropriate QC Antigen Shigella or QC Antigen Alkalescens-Dispar Group 1 chosen as the positive control and mix thoroughly.
- Negative control:** Dispense 1 drop of sterile 0.85% NaCl solution or Febrile Negative Control on the agglutination slide. Add 1 drop of the appropriate QC Antigen Shigella or QC Antigen Alkalescens-Dispar Group 1 and mix thoroughly.
- Rotate the slide for 1 minute and read for agglutination.

### Results

- Read and record results as follows:
  - 4+ 100% agglutination; background is clear to slightly hazy.
  - 3+ 75% agglutination; background is slightly cloudy.
  - 2+ 50% agglutination; background is moderately cloudy.
  - 1+ 25% agglutination; background is cloudy.
  - No agglutination.
- Positive control:** Should show 3+ or greater agglutination.
- Negative control:** Should show no agglutination. Rough reactions can occur. If so, repeat the test using Febrile Negative Control.

## Limitations of the Procedure

1. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent a smooth suspension of the microorganism or may cause evaporation or precipitation of the test mixture. False-positive reactions can occur.
2. Exposure to temperatures below 2°C can cause autoagglutination. Antigens must be smooth uniform suspensions. Examine antigen vials for agglutination before use. Suspensions with agglutination are not usable and should be discarded.
3. Allow the QC Antigens, the antisera and all equipment used to be at room temperature at the time of testing. The test reagents, if cold, may cause false-negative reactions.
4. Shake the antigen well before use to suspend the organisms.

## References

1. **Ewing, W.H. (ed.)**. 1986. Edwards and Ewing's identification of *Enterobacteriaceae*, 4th ed. Elsevier Science Publishing Co., Inc., New York, NY.
2. **Gray, L. D.** 1995. *Escherichia, Salmonella, Shigella and Yersinia*, p. 450-456. In P. R. Murray, E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
3. **Baron, E. J., L. R. Peterson, and S. M. Finegold.** 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis, MO.

4. **Pezzlo, M. (ed.)**. 1994. Aerobic Bacteriology, p. 1.0.1-1.20.47. In H. D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D. C.
5. **Andrews, W. H., G. A. June, and P. S. Sherrod.** 1995. Shigella, p. 6.01-6.06. In FDA Bacteriological Analytical Manual, 8th ed. AOAC International, Gaithersburg, MD.
6. **Smith, J. L.** 1992. *Shigella*, p. 423-431. In C. Vanderzant and D. F. Splittstoesser (eds.), Compendium of methods for the microbiological examination of foods, 3rd ed. American Public Health Association, Washington, D.C.

## Packaging

QC Antigen Shigella Group A	1 ml	2100-50
QC Antigen Shigella Group A <sub>1</sub>	1 ml	2101-50
QC Antigen Shigella Group B	1 ml	2102-50
QC Antigen Shigella Group C	1 ml	2103-50
QC Antigen Shigella Group C <sub>1</sub>	1 ml	2104-50
QC Antigen Shigella Group C <sub>2</sub>	1 ml	2105-50
QC Antigen Shigella Group D	1 ml	2106-50
QC Antigen Alkaescens-Dispar Group 1	1 ml	2116-50
Febrile Negative Control	5 ml	3239-56

# Bacto® Salmonella Antisera

## Salmonella O Antisera · Salmonella H Antisera · Salmonella H Antisera Spicer-Edwards

### Intended Use

Bacto Salmonella O Antisera are used in agglutination tests for the identification of *Salmonella* by somatic (O) antigens.

Bacto Salmonella H Antisera are used in tube agglutination tests for the identification of *Salmonella* by flagellar (H) antigens.

Bacto Salmonella H Antisera Spicer-Edwards are used in tube agglutination tests for screening and identifying the most commonly encountered salmonellae by flagellar (H) antigens.

### User Quality Control

#### Identity Specifications

##### Salmonella O Antisera

Lyophilized Appearance: Light gold to amber, button to powdered cake.

Rehydrated appearance: Light gold to amber, clear liquid.

##### Salmonella H Antisera

Lyophilized Appearance: Light gold to amber, button to powdered cake.

Rehydrated appearance: Light gold to amber, clear liquid.

*continued on following page*

### Summary and Explanation

*Salmonella* species cause a variety of human diseases called salmonellosis. The range of disease is from mild self-limiting gastroenteritis to more severe forms, possibly with bacteremia or typhoid fever, which can be life-threatening. Severe disease and bacteremia are associated primarily with three serovars of *S. enterica* subsp. *enterica* (Choleraesuis, Paratyphi A and Typhi) while most of the other 2,300 or more strains are associated with gastroenteritis. The severity of the diarrheal disease depends upon the virulence of the strain and the condition of the human host.

*Salmonella* is found in nature and occurs in the intestinal tract of many animals, both wild and domestic. The microorganism can spread to man from contact with the environment or from eating meat or vegetable food products.

**User Quality Control cont.****Performance Response**

Rehydrate Salmonella O, Salmonella Vi, and Salmonella H Antisera per label directions. Perform the slide agglutination test using appropriate Salmonella O and Vi Antisera and QC Antigens Salmonella O Groups A -I and Vi.

The chart below includes the QC Antigens Salmonella O recommended as homologous (positive) control antigens. (The homologous control antigen has certain identifying antigen(s) in common with the antiserum.) For a negative (heterologous) antigen control, use a QC Antigen Salmonella containing antigens unrelated to those in the homologous control.

**Homologous control:** Should show 3+ or greater agglutination.

**Negative control:** Should show no agglutination. Rarely, a +/- reaction is possible.

For Salmonella H Antisera, maintain stock cultures of known serological identification, and prepare antigen positive and negative controls by using known serotypes and following the procedure described above in Tube Test Preparation.

SALMONELLA O ANTISERUM	QC ANTIGEN SALMONELLA O HOMOLOGOUS CONTROL
Poly A-I & Vi	Groups A,B,C <sub>1</sub> ,C <sub>2</sub> ,D,E <sub>1</sub> ,E <sub>2</sub> ,E <sub>4</sub> ,F,G <sub>1</sub> ,H,I,Vi
Poly A	Groups A,B,D,E <sub>1</sub> ,E <sub>2</sub> ,E <sub>4</sub>
Poly B	Groups C <sub>1</sub> ,C <sub>2</sub> ,F,G <sub>1</sub> ,H
Poly C	Group I
Group A Factors 1,2,12	Group A
Group B Factors 1,4,5,12	Group B
Group B Factors 1,4,12,27	Group B
Group C1 Factors 6,7	Group C <sub>1</sub>
Group C2 Factors 6,8	Group C <sub>2</sub>
Group D1 Factors 1,9,12	Group D
Group E Factors 1,3,10,15,19,34	Groups E <sub>1</sub> ,E <sub>2</sub> ,E <sub>4</sub>
Group E1 Factors 3,10	Group E <sub>1</sub>
Group E2 Factors 3,15	Group E <sub>2</sub>
Group E4 Factors 1,3,19	Group E <sub>4</sub>
Group F Factor 11	Group F
Group G Factors 13,22,23,(36),(37)	Group G <sub>1</sub>
Group G1 Factors 13,22,(36),(37)	Group G <sub>1</sub>
Group H Factors 1,6,14,24,25	Group H
Group I Factor 16	Group I
Vi	Group Vi
Factor 2	Group A
Factor 4	Group B
Factors 4,5	Group B
Factor 5	Group B
Factor 7	Group C <sub>1</sub>
Factor 8	Group C <sub>2</sub>
Factor 9	Group D
Factor 10	Group E <sub>1</sub>
Factor 15	Group E <sub>2</sub>
Factor 19	Group E <sub>4</sub>
Factor 22	Group G <sub>1</sub>
Factor 14	Group H

Note: Parentheses ( ) enclosing the designation for an antigen indicate that the antigen is poorly developed or agglutinates weakly agglutinates. For a complete and current explanation of the classification of *Salmonella*, consult appropriate references.<sup>1,2,3,8,11</sup>

**Table 1.** Differentiation of the genus *Salmonella* from other genera.<sup>1</sup>

Test	<i>Salmonella</i>	<i>Citrobacter</i>			<i>Edwardsiella</i>
		<i>amolonaticus</i>	<i>diversus</i>	<i>freundii</i>	
Indole production	-	+	+	-	d
Citrate, Simmons	+	[+]	+	+	-
H <sub>2</sub> S production	+	-	-	[+]	d
Urea hydrolysis	-	[+]	d	d	-
Lysine decarboxylase	+	-	-	-	+
Ornithine decarboxylase	+	+	+	[-]	[+]
D-Adonitol, acid production	-	-	+	-	-
L-Arabinose, acid production	+	+	+	+	-
L-Rhamnose, acid production	+	+	+	+	-
D-Sorbitol, acid production	+	+	+	+	-
D-Xylose, acid production	+	+	+	+	-
Acetate utilization	+	[+]	[+]	[+]	-

+ 90-100% positive  
 [+] 76-89% positive  
 d 26-75% positive  
 [-] 11-25% positive  
 - 0-10% positive

All *Salmonella* serovars belong to two species: *S. bongori*, which contains 18 serovars, and *S. enterica*, which contains the remaining 2,300 or more serovars divided among six subspecies.<sup>2,3</sup>

The six subspecies of *S. enterica* are:

- S. enterica* subsp. *enterica* (I or 1)
- S. enterica* subsp. *salamae* (II or 2)
- S. enterica* subsp. *arizonae* (IIIa or 3a)
- S. enterica* subsp. *diarizonae* (IIIb or 3b)
- S. enterica* subsp. *houtenae* (IV or 4)
- S. enterica* subsp. *indica* (VI or 6)

(The legitimate species name for *S. enterica* is *S. choleraesuis*. However, this name may be confused with the serotype named "choleraesuis." At the International Congress for Microbiology in 1986, the International Subcommittee for *Enterobacteriaceae* agreed to adopt the species name, *S. enterica*.<sup>4</sup> LeMinor and Popoff<sup>5</sup> published a request to the Judicial Commission to use *S. enterica* as a species name. The Commission ruled that *S. choleraesuis* is the legitimate name.<sup>6,7</sup> *S. enterica* is used in many countries and is favorably accepted as the species name.<sup>8,9</sup> The Centers for Disease Control has adopted this designation until the problem of naming this species is resolved.<sup>2</sup>)

Nomenclature and classification of these bacteria are constantly changing.<sup>1</sup> *Salmonella* and the former *Arizona* should be considered a single genus, *Salmonella*.<sup>9</sup> It is recommended that laboratories report the names of *Salmonella* serovars for the subspecies *enterica*. The serovar names are no longer italicized and the first letter is capitalized. For example, the strain that used to be identified as *Salmonella typhimurium* is now known as *Salmonella* Typhimurium.

Serovars of other subspecies of *S. enterica* (except some in the subspecies *salamae* and *houtenae*) and those of *S. bongori* are not named and are designated by their antigenic formula. For the most recent information on nomenclature, consult appropriate references.<sup>1,7,10-13</sup>

Results are for 48-hours incubation. Tests were performed at 35-37°C.

Serotypes of *Salmonella* are defined based on the antigenic structure of both somatic or cell wall (O) antigens and flagellar (H) antigens. The antigenic formula lists the O antigen(s) first, followed by the H antigen(s). The major antigens are separated by colons and the components of the antigens separated by commas. For example, the antigenic formula for *Salmonella* Typhimurium is *Salmonella* 1,4,5,12:i:1,2. This means that the strain has O antigen factors 1,4,5 and 12, the flagella phase 1 antigen i, and flagella phase 2 antigens 1 and 2.

**Table 2a.** Differentiation of *Salmonella* species, subspecies, and some serovars.<sup>1,2</sup>

Test	<i>Salmonella enterica</i>					
	<i>Salmonella bongori</i>	subsp. <i>arizonae</i>	subsp. <i>enterica</i>	subsp. <i>diaizonae</i>	subsp. <i>houtenae</i>	subsp. <i>indica</i>
Citrate, Simmons	+	+	+	+	+	[+]
H <sub>2</sub> S production	+	+	+	+	+	+
Lysine decarboxylase	+	+	+	+	+	+
Ornithine decarboxylase	+	+	+	+	+	+
Motility	+	+	+	+	+	+
KCN, growth	+	-	-	-	+	-
Malonate utilization	-	+	-	+	-	-
D-Glucose, gas	[+]	+	+	+	+	+
L-arabinose, acid	+	+	+	+	+	+
Dulcitol, acid	+	-	+	-	-	d
Lactose, acid	-	[-]	-	[+]	-	[-]
Maltose, acid	+	+	+	+	+	+
Melibiose, acid	[+]	+	+	+	+	[+]
L-Rhamnose, acid	+	+	+	+	+	+
D-Sorbitol	+	+	+	+	+	-
Trehalose, acid	+	+	+	+	+	+
D-Xylose, acid	+	+	+	+	+	+
Mucate, acid	+	+	+	d	-	+
Tartrate, Jordans	-	-	+	[-]	d	+
ONPG	+	+	-	+	-	D

+ 90-100% positive  
 [+] 76-89% positive  
 d 26-75% positive

Complete identification of *Salmonella* requires cultural isolation, biochemical characterization and serotyping. However well-defined the serology of *Salmonella*, the use of serological procedures does not supersede cultural isolation and biochemical characterization. Any serological results obtained before biochemical identification must be considered as presumptive identification only. Consult to appropriate references for complete identification of *Salmonella*.<sup>1,2,3,8,11-14</sup>

**Characterizing the Serotypes of *Salmonella***

**Salmonella O Antigens:** The somatic (O) heat-stable antigens are identified first. The O antigens are numbered 1-67 using Arabic numerals. The numbers are not completely continuous because certain strains were reclassified to other genera and the antigenic Arabic numbers were deleted from the schema.

**Table 2b.** Differentiation of *Salmonella* species, subspecies and some serovars.<sup>1,2</sup>

Test	subsp. <i>salamae</i>	<i>Salmonella enterica</i> subsp. <i>enterica</i>				
		serovar Choleraesuis	serovar Gallinarum	Serovar Paratyphi A	serovar Pullorum	Serovar Typhi
Citrate, Simmons	+	[-]	-	-	-	-
H <sub>2</sub> S production	+	d	+	-	+	+
Lysine decarboxylase	+	+	+	-	+	+
Ornithine decarboxylase	+	+	-	+	+	-
Motility	+	+	-	+	-	+
KCN, growth	-	-	-	-	-	-
Malonate utilization	+	-	-	-	-	-
D-Glucose, gas	+	+	-	+	+	-
L-arabinose, acid	+	-	[+]	+	+	-
Dulcitol, acid	+	-	+	+	-	-
Lactose, acid	-	-	-	-	-	-
Maltose, acid	+	+	+	+	-	+
Melibiose, acid	-	d	-	+	-	+
L-Rhamnose, acid	+	+	-	+	+	-
D-Sorbitol	+	[+]	-	+	[-]	+
Trehalose, acid	+	-	d	+	[+]	+
D-Xylose, acid	+	+	d	-	[+]	[+]
Mucate, acid	+	-	d	-	-	-
Tartrate, Jordans	d	[+]	+	-	-	+
ONPG	[-]	-	-	-	-	-

[-] 11-25% positive  
 0 10% positive

Serogroups represent the organization of *Salmonella* strains based on the antigen(s) shared in common and are designated by the letters A-Z. After exhausting the alphabet, the serogroups were numbered beginning with No. 51 (the serogroup Z organism having antigen No. 50). While one somatic antigen identifies each serogroup, certain other antigens may be shared among several serogroups.

Most organisms contain antigens in common that will cause cross-reactions in an unabsorbed or “partially absorbed” antiserum. One somatic antigen identifies a serogroup and is shared in common by all members of a given serogroup. For example, serogroup A is represented by three members, *Salmonella* Paratyphi A (somatic antigens 1,2,12), *Salmonella* Kiel (somatic antigens 1,2,12), and *Salmonella* Nitra (somatic antigens 2,12). All three members of this serogroup contain antigens 2 and 12 in common. Serogroup B is represented by many organisms consisting having different combinations of somatic antigens 1,4,5,12 and 27. Serogroup D organisms contain somatic antigens 1,9,12, etc.

In the above example, all three serogroups A, B, and D contain antigens 1 and 12. An antiserum prepared from a 1,2,12 culture, if not absorbed, will react with cultures of serogroups B and D in varying degrees depending on the concentration of the commonly shared 1 and 12 factors. This must be taken into consideration when choosing an antiserum to be used in the examination of the salmonellae.

Several different antisera are available. Some represent group antigens. Others are single factor sera, which should be used when testing for an identifiable antigen in a given serogroup. Such a single factor serum is not called a “group” serum, though it contains the group identifiable agglutinin. (It has been recommended by the CDC that the term “group” be applied only to those sera possessing all the major agglutinins found in that group.)

In unabsorbed antisera, cross-reactions occur if strains sharing some “like” antigens are tested, even when they are in separate serogroups based on the major group antigen(s) they possess. Unabsorbed antisera are available as group antisera containing all factors in that group.

In absorbed antisera, cross reactions are less likely and are weaker. Absorbed antisera are available as factor specific antisera.

**Flagellar *Salmonella* H Antigens:** The flagellar (H) antigens are heat labile and are usually associated with motility. Cultures are ordinarily flagellated and actively motile, although flagellated cultures can be nonmotile. H antigen characterization is done after the serogroup of the strain is determined. The H antigens of *Salmonella* are designated by letters of the alphabet, a-z, followed by z, z<sub>1</sub>, z<sub>2</sub>, etc., and by Arabic numerals. H antigens exist in 2 phases, phase 1 and phase 2. Phase 1 antigens are expressed in letters a-z, etc., and the phase 2 antigens are most often expressed in Arabic numerals. Older cultures may express both phases of a diphasic serotype, but recent clinical isolates more often express only one phase. Phase reversal may be necessary to isolate both phases of a diphasic culture. Consult an appropriate reference for more detailed information.<sup>8</sup>

A pure H antiserum cannot be prepared without some somatic content. However, since H antigens are highly antigenic, the serum derived from motile cultures may be used at a dilution that reduces somatic agglutination below the detection level.

***Salmonella* Vi Antigen:** The Vi Antigen is a heat-labile envelope antigen that may surround a cell wall and mask somatic antigen activity. Microorganisms having the Vi Antigen will not agglutinate in O antisera.

### Using *Salmonella* Antisera

***Salmonella* O Antisera:** The recommended serological Identification scheme begins with *Salmonella* O Antisera Poly A through Poly G, which contain the following:

SALMONELLA POLY GROUP ANTISERA	SOMATIC GROUPS PRESENT
Bacto <i>Salmonella</i> O Antiserum Poly A	A,B,D,E <sub>1</sub> ,(E <sub>2</sub> ,E <sub>3</sub> ),*E <sub>4</sub> ,L
Bacto <i>Salmonella</i> O Antiserum Poly B	C <sub>1</sub> ,C <sub>2</sub> ,F,G,H
Bacto <i>Salmonella</i> O Antiserum Poly C	I,J,K,M,N,O
Bacto <i>Salmonella</i> O Antiserum Poly D	P,Q,R,S,T,U
Bacto <i>Salmonella</i> O Antiserum Poly E	V,W,X,Y,Z
Bacto <i>Salmonella</i> O Antiserum Poly F	51–55
Bacto <i>Salmonella</i> O Antiserum Poly G	56–61

\*Strains of groups E<sub>2</sub> and E<sub>3</sub> are lysogenized by phage 15, then by phage 34. These strains are now classified into group E<sub>1</sub>.<sup>3</sup>

If agglutination occurs, use individual *Salmonella* O Group Factor Antisera to determine the specific serogroup to which the isolate belongs. For efficiency, test first with individual *Salmonella* O Group Factor Antisera.

If agglutination does not occur with Poly A or B, test the isolate with *Salmonella* O Antiserum Vi. If positive, heat and retest. If agglutination does not occur with *Salmonella* O Antiserum Vi, the isolate is not likely to be *Salmonella*. Results should be examined. If questions exist, the isolate should be sent to a reference laboratory.

If agglutination does not occur with Poly C, D, E, F, and G, the isolate is not likely to be *Salmonella*.

**Table 3.** Schema for using *Salmonella* O Antisera Poly Groups A, B, C, D, E, F and G.

Test with	<i>Salmonella</i> O antisera Poly Groups A, B, C, D, E, F and G			
Test Result	+	– with Poly A or B		– with Poly C, D, E, F and G
Test with	Individual <i>Salmonella</i> O Antisera	Vi Antiserum		
Test Result	+ with one <i>Salmonella</i> O Antiserum (required)	+	–	
Test Conclusion or Next Action	Determine the <i>Salmonella</i> H Antigen	Heat and retest with individual <i>Salmonella</i> O Antisera	Test isolate is not a <i>Salmonella</i>	Test isolate is not a <i>Salmonella</i>

***Salmonella* O Antiserum Poly A-I and Vi:** This antiserum detects factors 1-16, 19, 22-25, 34 and Vi. This combination of factors represents the most frequently isolated Groups A-I and the Vi antigens and is used to screen possible *Salmonella* isolates.

A positive reaction indicates that further serological testing is needed to identify the isolate using *Salmonella* O Group Factor Antisera. The most common serogroups are B, D and C<sub>1</sub>. For efficiency, first use the *Salmonella* O Group Factor Antisera for these serogroups.

If the isolate is positive with *Salmonella* O Antiserum Poly A-I and Vi but negative with Poly A-Poly G, test the isolate with *Salmonella* Vi Antiserum. If positive with *Salmonella* Vi Antiserum, heat and retest using individual *Salmonella* O Antisera. If negative with *Salmonella*

Vi Antiserum, the isolate is not likely to be Salmonella. Results should be examined. If questions exist, the isolate should be sent to a reference laboratory.

A negative reaction indicates the isolate is not in serogroups A-I. If the biochemical reactions are consistent with *Salmonella*, a serogroup other than A-I is possible. Further testing with antisera for other serogroup antigens is necessary.

**Table 4.** Schema for using Salmonella O Antiserum Poly A-I & Vi.

Test with	Salmonella O Antisera Poly A-I and Vi			
Test Result	+		-	
Test with	Individual Salmonella O Antisera			
Test Result	+	-		
Test with	Salmonella Vi Antiserum			
Test Result		+	-	
Test Conclusion or Next Action	Determine the Salmonella H Antigen	Heat and retest with individual Salmonella O Antisera	Test isolate is not a <i>Salmonella</i>	May be a <i>Salmonella</i> detectable by use of Salmonella O Antisera Poly C, D, E, F or G

**Salmonella O Group Factor Antisera and Single Factor Antisera:**

Use selected Salmonella O Group Factor Antisera. Cross reactions may occur between serogroups that share O antigens. Consider this partial list of Salmonella O Group Factor Antisera as an example:

- Salmonella O Antiserum Group A Factors 1, 2, 12
- Salmonella O Antiserum Group B Factors 1, 4, 5, 12
- Salmonella O Antiserum Group B Factors 1, 4, 12, 27

Factors 1 and 12 occur in combination with other antigens and may cause cross-reactions. The strength of the reactions will help in interpretation. Rapidly forming 3+ or greater agglutination indicates a homologous reactions.

Use selected Salmonella O Factor Antisera. Absorbed antisera specific for an identifiable antigen in a given serogroup is used to identify the isolate further. In the example above, Salmonella O Factor Antisera could be used:

- Salmonella O Antiserum Factor 2
- Salmonella O Antiserum Factor 4
- Salmonella O Antiserum Factors 4, 5
- Salmonella O antiserum Factor 5

**Polyvalent Salmonella H Antisera:** Further identification of a *Salmonella* isolate includes the characterization of the flagellar antigens. Agglutination with the following Polyvalent H Antisera can be done:

SALMONELLA POLY GROUP ANTISERA	FLAGELLAR ANTIGENS PRESENT
Salmonella H Antiserum Poly a-z	Groups EN,G,L,Z <sub>4</sub> , 1 complexes and a-k,r-z,Z <sub>6</sub> ,Z <sub>10</sub> ,Z <sub>29</sub>
Salmonella H Antiserum Poly A	Groups a,b,c,d,i,Z <sub>10</sub> ,Z <sub>29</sub>
Salmonella H Antiserum Poly B	Groups eh,en,enz,Z <sub>15</sub> , G complex
Salmonella H Antiserum Poly C	Groups k,l,r,y,z,Z <sub>4</sub>
Salmonella H Antiserum Poly D	Groups z <sub>35</sub> ,Z <sub>36</sub> ,Z <sub>37</sub> ,Z <sub>38</sub> ,Z <sub>39</sub> ,Z <sub>41</sub> ,Z <sub>42</sub>
Salmonella H Antiserum Poly E	1 complex Z <sub>6</sub>

Absorbed H antisera specific for single antigens or a complex of antigens can be used to identify the isolate further.

**Unabsorbed and Absorbed Salmonella H Antisera:** Complete identification of a *Salmonella* isolate involves analysis of phase 1 and phase 2 antigens using H antisera. For the complex pattern of analysis and procedures, consult appropriate references.<sup>8</sup>

**Salmonella H Antisera Spicer-Edwards:** Salmonella H Antisera Spicer-Edwards is used for screening and identifying the most commonly encountered *Salmonella* using a combination of polyvalent and single complex antisera.

**Table 5.** Identification of *Salmonella* H using Salmonella H Antisera Spicer-Edwards.

H Antigen(s)	Salmonella H Antisera Spicer-Edwards			
	1	2	3	4
a	+	+	+	-
b	+	+	-	+
c	+	+	-	-
d	+	-	+	+
e,h	+	+	+	-
G Complex*	+	-	-	+
i	+	-	-	-
k	-	+	+	+
r	-	+	-	+
y	-	+	-	-
z	-	-	+	+
Z <sub>4</sub> Complex**	-	-	+	-
Z <sub>10</sub>	-	-	-	+
Z <sub>29</sub>	-	+	+	-

**Table 6.** Identification of *Salmonella* H using Salmonella H Antisera.

H Antigen(s)	Salmonella H Antisera
e,n,x, e,n,Z <sub>15</sub>	EN Complex
I,v I,w I,Z <sub>13</sub> I,Z <sub>28</sub>	L Complex
1,2 1,5 1,6 1,7	1 Complex

\* The G complex component of Salmonella H Antisera Spicer-Edwards 1 and 4 reacts with antigens f,g; f,g,s; f,g,t; g,m; g,m,q; g,m,s; g,m,t; g,m,t; g,p; g,p,s; g,p,u; g,q; g,s,t; g,t; m,p,t,u; and m,t.

\*\* The Z<sub>4</sub> Complex component reacts with z<sub>4</sub>,z<sub>23</sub>; z<sub>4</sub>,z<sub>24</sub>; and z<sub>4</sub>, z<sub>32</sub>.

Note that no antigen is positive with all four Salmonella H Antisera Spicer-Edwards Any antigen that reacts with all four sera should be checked for smoothness.

### Extent of Serological Identification Necessary

Complete serological characterization of *Salmonella* is not required for successful detection of the microorganism when it occurs as a pathogen. The use of adequate isolation procedures and differential biochemical tests is of primary importance. Possible *Salmonella* isolates can be presumptively identified with a minimum of serological identification. Isolates can be sent to laboratories that perform the level of testing necessary to completely identify the microorganism.

For a further discussion of the serological identification of *Salmonella*, consult appropriate references.<sup>1,2,3,8,11</sup>

### Principles of the Procedure

Identification of *Salmonella* species includes both biochemical and serological identification. Serological confirmation involves the procedure in which the microorganism (antigen), reacts with its corresponding antibody. This *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction is rapid, does not dissociate (high avidity) and bonds strongly (high affinity).

Because a microorganism (antigen) may agglutinate with an antibody produced in response to another species, heterologous reactions are possible. These are characterized as weak in strength or slow in formation. Such unexpected and perhaps unpredictable reactions may lead to some confusion in serological identification. Therefore, a positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism.

Agglutination of the somatic antigen in the slide test appears as a firm granular clumping. Homologous reactions are rapid and strong (3+). Heterologous reactions are slow and weak.

Agglutination of the flagellar antigens in the tube test appears as a loose flocculation that can easily be resuspended.

### Reagents

Salmonella O, H, and Vi Antisera are lyophilized, polyclonal rabbit antisera containing approximately 0.04% Thimerosal as a preservative.

**Salmonella O Poly Antisera** are polyvalent antisera. Each antiserum is specific for certain serogroup antigens. When properly rehydrated and used as recommended, each vial of Salmonella O or Vi Antisera contains sufficient reagent for 60 tests. Salmonella O Antisera Poly A-I and Vi is prepared with representative strains of these serogroups and is not absorbed. It may cross-react with other antisera because of shared common O antigens.

Salmonella O Group Antisera are specific for the major factors present in the serogroup. Salmonella O Factor Antisera are specific for the factors of the individual serogroups. When using Salmonella O Group Antisera, cross-reactions are possible because serogroups may share non-major group antigens. Salmonella O Factor Antisera are absorbed as necessary to render each antiserum as specific as practical without reducing the homologous reactions to an unsatisfactory level.

**Salmonella H Poly Antisera** are polyvalent antisera. Each antiserum is specific for certain flagellar antigens. Each vial of Salmonella H Antiserum contains sufficient reagent to perform between 150-1500 tests, depending on the antiserum used. Salmonella H Antisera are either absorbed or unabsorbed specifically for either phase 1 or phase 2 antigens. Salmonella H Antisera Spicer-Edwards are pooled, polyvalent antisera and additional adjunctive antisera to identify the more commonly occurring H antigens.

### Precautions

1. For In Vitro Diagnostic Use.
2. The Packaging of This Product Contains Dry Natural Rubber.
3. Follow proper established laboratory procedure in handling and disposing of infectious materials.

### Storage

Store lyophilized and rehydrated Salmonella O, H and Vi antisera at 2-8°C.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products. Discard any antiserum that becomes cloudy during storage.

### Expiration Date

The expiration date applies to the product in its intact container when stored as directed. Do not use a product if it fails to meet specifications for identity and performance.

### Procedure

#### Materials Provided

Salmonella O Antisera  
Salmonella H Antisera  
Salmonella Vi Antiserum  
(See Packaging.)

#### Materials Required But Not Provided

##### Slide Test

0.85% NaCl solution, sterile  
Agglutination slides with 1 inch squares  
Applicator sticks  
Boiling waterbath  
Centrifuge  
QC Antigens Salmonella

##### Tube Test

0.85% NaCl solution, sterile  
Culture tubes, 12 x 75 mm, and rack  
Waterbath, 50 ± 2°C  
Serological pipettes, 1 ml  
Formaldehyde

#### Reagent Preparation

Equilibrate all materials to room temperature prior to performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergents.

**Salmonella O, H and Vi Antisera:** To rehydrate, add 3 ml of sterile 0.85% NaCl solution and rotate gently to completely dissolve the contents. Rehydrated antisera are considered a 1:2 dilution. Subsequent Salmonella O Antisera dilutions are based on this as a starting dilution. The H antisera are further diluted for use.

#### Specimen Collection and Preparation

**Clinical specimens:** *Salmonella* can be recovered from selective differential media such as Hektoen Enteric Agar or XLD agar. For specific recommendations, consult appropriate references.<sup>11,12</sup> Determine that a pure culture of the microorganism has been obtained and

that biochemical test reactions are consistent with the identification of the organism as a *Salmonella* species. After these criteria are met, serological identification can be performed.

**Food samples:** *Salmonella* can be recovered when samples are processed to recover injured microorganisms and prevent overgrowth of competing microorganisms. Consult appropriate references for recommended procedures for isolation of *Salmonella* from foods.<sup>13,14</sup> Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as a *Salmonella* species. After these criteria have been met, serological identification can be performed.

### Slide Test Procedure

#### Salmonella O and Vi Antisera

Use this procedure to test the isolate with each selected antiserum.

1. **Salmonella Antiserum:** Add Dispense 1 drop (3  $\mu$ l) of each antiserum to be tested on an agglutination slide.
2. **Negative control:** Dispense 1 drop of 0.85% sterile NaCl solution on an agglutination slide.
3. **Test isolate:** From a solid agar medium, transfer a portion of a loopful of an isolated colony to each of the two reaction areas above and mix thoroughly.
4. **Positive control:** Dispense 1 drop of each Salmonella O Antiserum to be tested on an agglutination slide. Add 1 drop of an appropriate QC Antigen Salmonella.
5. Rotate the slides for 1 minute and read for agglutination. Results must be read within 1 minute.

### Slide Test Results

1. Read and record results as follows:
  - 4+ 100% agglutination; background is clear to slightly hazy.
  - 3+ 75% agglutination; background is slightly cloudy.
  - 2+ 50% agglutination; background is moderately cloudy.
  - 1+ 25% agglutination; background is cloudy.
  - No agglutination.
2. **Positive control:** Should show a 3+ or greater agglutination.
3. **Negative control:** Should show no agglutination. If agglutination occurs, the culture is rough and cannot be tested. Subculture to a non-inhibitory medium, incubate and test the organism again.
4. **Test isolates:** 3+ or greater agglutination is a positive result.
5. A partial (less than 3+) or delayed agglutination reaction should be considered negative.
6. If a positive reaction occurs with Salmonella Vi Antiserum, follow this procedure:
  - Prepare a dense suspension of the isolate from an agar medium in 3-5 ml of sterile 0.85% NaCl solution.
  - Heat in a boiling waterbath for 30-60 minutes and cool. The suspension should not precipitate after heating. If this occurs, select another colony for testing.
  - Centrifuge at 1,000 rpm for 10-15 minutes.
  - Aspirate and discard the supernatant.
  - Resuspend the sediment in 0.5 ml sterile 0.85% NaCl solution.
  - Retest a drop of the sediment with Salmonella O Group Antisera, as outlined.

7. If the heated culture continues to react with Salmonella Vi Antiserum and not with the Salmonella O Antisera, the isolate may not be *Salmonella*. Test the isolate further to determine if it is correctly identified.
8. If an H antigen identification is required, proceed to the next section.
9. When a negative reaction is obtained with Salmonella O Antiserum Poly A-I and Vi in the above procedure, the organism is presumptively negative for *Salmonella* that belong to serogroups A-I. Biochemical tests should be performed to confirm this negative result. If biochemical tests prove the organism to be a *Salmonella*, a serogroup beyond serogroup A-I is probably involved.
 

If the organism reacts with Poly A-I and Vi but does not react with the specific somatic antisera, it should be checked with Salmonella Vi Antiserum by the above procedure.

### Tube Test Preparation

1. **0.6% formalized saline:** Prepare by adding 6 ml formaldehyde per 1,000 ml of sterile 0.85% NaCl solution.
2. **Test organism:** It is often necessary to increase the motility of the test organism. To accomplish this, make several consecutive transfers in Motility GI Medium.
  - Inoculate the tube slightly below the surface of the medium using the stab method.
  - Incubate at 35-37°C for 18-20 hours.
  - Transfer only those organisms that have migrated to the bottom of the tube.
  - When the organism successfully travels 50-60 mm through the medium in 18-20 hours, it is ready for use.
  - An infusion broth such as Veal Infusion Broth is recommended for cultivating the motile *Salmonella* prior to testing. It should be inoculated and incubated at 35°C for 24 hours. Brain Heart Infusion Broth may be used with incubation at 35°C for 4-6 hours. If Tryptic Soy Broth is used, incubate at 35°C for 24 hours.
  - Prepare the test organism suspension by using equal volumes of broth culture and 0.6% formalized saline. The final density of this test suspension should be that of a McFarland Barium Sulfate Standard No. 3.
3. **Positive control:** Commercially prepared QC Salmonella H antigens are not available. The user must maintain stock cultures of known serological identification for use in quality control. Prepare the antigen by using known serotypes and following the procedure described above. (See Test organism, above.)
4. **Salmonella H Antisera:** Rehydrated antisera is considered a 1:2 working dilution. Prepare dilutions as follows and use on the day prepared. Discard any unused portion.
  - Most **Salmonella H Antisera:** A 1:1,000 final dilution is used. Prepare by adding 0.1 ml rehydrated antiserum to 25 ml of 0.85% NaCl solution.
  - **Salmonella H Antisera x, z<sub>13</sub>, z<sub>15</sub> and z<sub>28</sub>:** A 1:500 final dilution is used. Prepare by adding 0.1 ml rehydrated antiserum to 12.5 ml of 0.85% NaCl solution.
  - **Salmonella H Antiserum Poly a-z:** A 1:100 final dilution is used. Prepare by adding 0.1 ml rehydrated antiserum to 2.5 ml of 0.85% NaCl solution.

### Tube Test Procedure

#### Salmonella H Antiserum:

1. Prepare a 12 x 75 ml culture tube for each organism to be tested.
2. **Diluted antiserum:** Dispense Add 0.5 ml of diluted antiserum in each tube.
3. **Test isolate:** Add 0.5 ml to the appropriate tube.
4. **Positive control:** Add 0.5 ml of antigen positive control to a tube containing 0.5 ml of antiserum.
5. **Negative control:** Add 0.5 ml of 0.85% NaCl solution to a tube containing 0.5 ml of test isolate.
6. Incubate all tubes in a waterbath at  $50 \pm 2^\circ\text{C}$  for 1 hour.
7. Read for flocculation (agglutination).
8. Repeat the Tube Test using a phase-reversed test organism. (See the procedure for phase reversal below.)

#### Phase Reversal:

1. Prepare Motility GI Medium phase reversal medium according to directions.
2. Prepare the antiserum opposite to the phase desired. For example, incubating *Salmonella* Typhimurium phase 1[i] in GI Motility Medium containing i antiserum allows the growth and spread of *S. Typhimurium* phase 2 [1,2].
3. Add 1 ml of a 1:10 dilution of antiserum to 25 ml of sterile GI Motility Medium and mix well. Pour into a sterile Petri dish and allow to solidify.
4. Inoculate by punching the edge of the solidified medium.
5. Incubate at  $35\text{--}37^\circ\text{C}$  for 24 hours.
6. Transfer growth from the spreading edge opposite the inoculation site to a liquid medium for testing according to steps under Tube Test Procedure – Salmonella H Antisera.
7. If motility is not acceptable, pass through GI Motility Medium again.

### Tube Test Procedure

#### Salmonella H Antiserum Spicer-Edwards

1. Prepare the test organism and the 1:2 antiserum dilution as described above in Tube Test Preparation.
2. **Final 1:1,000 dilution of antiserum:** Prepare by adding 0.1 ml of rehydrated antiserum (1:2 working dilution) to 25 ml of 0.85% NaCl solution.
3. Prepare 4 culture tubes (12 x 75 ml) for each test organism.
4. **Salmonella H Antisera Spicer-Edwards 1-4:** Add 0.5 ml of the diluted antiserum to the culture tubes.
5. **Test organism:** Add 0.5 ml to each tube.
6. Incubate tubes in waterbath at  $50 \pm 2^\circ\text{C}$  for 1 hour.
7. Remove from the waterbath. Avoid excessive shaking when the tubes are in the waterbath or when removing them from the waterbath prior to reading the reactions.
8. Read for flocculation (agglutination).

### Tube Test Results

Compare results with the flocculation (agglutination) patterns for the Spicer-Edwards schemae. (Table 5).

### Limitations of the Procedure

1. Complete O and H antigen characterization of a *Salmonella* isolate is required for final identification. Due to the complexity of the laboratory procedures, identification with polyvalent antisera may be sufficient for most laboratories.
2. Possible *Salmonella* isolates having inconsistencies in biochemical reactions and O and H antigen tests should be referred to a reference laboratory for further testing.
3. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent making a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False-positive reactions may occur.
4. Rough culture isolates do occur and will agglutinate spontaneously, causing agglutination of the negative control reaction (autoagglutination). Smooth colonies must be selected and tested in serological procedures.
5. In the slide agglutination procedure for O antigen testing, it is recommended that several colonies be tested and that unabsorbed polyvalent antisera be used followed by absorbed single factor antisera. For example, colonies of a 1,2,12 culture on an agar plate will have varying degrees of each antigen. A 1,2,12 antiserum absorbed of 1 and 12 antibodies will be highly specific but will show weak or no agglutination with colonies that have less of antigen 2 and more of antigens 1 and 12. Using unabsorbed Salmonella O Antiserum Group Factors 1,2,12 to test several suspicious colonies on a plate followed by testing with absorbed Salmonella O Antiserum Factor 2 gives the needed balance of sensitivity and specificity.
6. The slide agglutination antisera (Salmonella O Antisera) have been prepared for use in identifying cultures already defined biochemically. Such cultures are taken from an agar medium using a bacterial loop. A portion of the isolated colony is emulsified in a drop of antiserum. It is recommended that more than one colony be tested. The sera have been tested and absorbed using this method. They have not been tested employing an antigen suspension in NaCl solution or alcohol-treated cultures. If variations in the recommended procedures are to be used, the investigator is advised to test each lot of antiserum with known positive control cultures to ensure its proper homologous and heterologous reactions under their test conditions.
7. Agglutination reactions of 3+ or greater are interpreted as positive reactions. Cross-reactions resulting in a 1+ or 2+ agglutination are likely since there are somatic antigens shared among different groups as non-major group antigens.
8. The tube agglutination technique is recommended for H antigen testing because cross-reactions with somatic antigens may occur at the dilutions used in the slide technique.
9. No attempt has been made to absorb or test for O antibodies in H antisera.
10. In the tube test, make certain that the proper dilution is prepared for a given antiserum. Various dilutions are used for various sera. The information is given under Tube Test Preparation. Also, it is important in this test to use the recommended time and temperature of incubation. Make certain that the waterbath is in a location free of mechanical vibration.

11. There may exist common antigens between various "O" serogroups of Salmonella. As an example, Salmonella O Antiserum Poly A contains, among others, agglutinins for factor 1, since cultures possessing factor 1 were used in immunization. It may be expected that this polyvalent antiserum will react with cultures other than those contained in "O" serogroups A, B, D, E, and L due to the common 1 antigen (those organisms in Group G<sub>1</sub>, G<sub>2</sub>, H, R, T, etc., which contain factor 1).
12. Salmonella O Antiserum Poly A-I & Vi has been prepared with representative members of those somatic groups and has not been absorbed. It is obvious that this serum may and will react with higher O groups of Salmonella.

## References

1. **Holt, J. G., N. R. Krieg, P. H. Sneath, J. T. Staley, and S. T. Williams.** 1994. Bergey's manual of determinative bacteriology, 9th ed. Williams & Wilkins, Baltimore, MD.
2. **McWhorter-Murlin, A. C., and F. W. Hickman-Brenner.** 1994. Identification and serotyping of *Salmonella* and an update of the Kauffmann-White Scheme. Centers for Disease Control and Prevention, Atlanta, GA.
3. **Popoff, M. Y., and L. LeMinor.** 1997. Antigenic Formulas of the *Salmonella* Serovars. WHO Collaborating Centre for Reference and Research on *Salmonella*. Institut Pasteur, Paris, France.
4. **Penner, J. L.** 1988. International committee on systematic bacteriology taxonomic subcommittee on *Enterobacteriaceae*. Int. J. Syst. Bacteriol. **38**:223-224.
5. **LeMinor, L., and M. Y. Popoff.** 1987. Request for an opinion. Designation of *Salmonella enterica* sp. nov., nom. rev., as the type and only species of the genus *Salmonella*. Int. J. Syst. Bacteriol. **37**:465-468.
6. **Wayne, L. G.** 1991. Judicial Commission of the International Committee on Systematic Bacteriology. Int. J. Syst. Bacteriol. **41**:185-187.
7. **Wayne, L. G.** 1994. Actions of the Judicial Commission of the International Committee on Systematic Bacteriology on requests for opinions published between January 1985 and July 1993. Int. J. Syst. Bacteriol. **44**:177.
8. **Ewing, W.H.** 1986. Edwards and Ewing's Identification of *Enterobacteriaceae*, 4th ed. Elsevier Science Publishing Co., Inc., New York, NY.
9. **Old, D. C.** 1992. Nomenclature of *Salmonella*. J. Med. Microbiol. **37**:361-363.
10. **Farmer III, J. J., III, A. C. McWhorter, D. J. Brenner, and G. D. Morris.** 1984. The *Salmonella-Arizona* group of *Enterobacteriaceae*: nomenclature, classification and reporting. Clin. Microbiol. Newsl. **6**:63-66.
11. **Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.).** 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
12. **Isenberg, H. D. (ed.)** 1992. Clinical microbiology procedures handbook, vol. 2. American Society for Microbiology, Washington, D.C.
13. **Andrews, W. H., G. A. June, P. Sherrod, T. S. Hammack, and R. M. Amaguana.** 1995. Food and drug administration bacteriological analytical manual, 8th ed. AOAC International, Gaithersburg, MD.
14. **Russell, S. F., J. D'Aoust, W. H. Andrews, and J. S. Bailey.** 1992. *Salmonella*. In C. Vanderzant, C., and Splittstoesser, D.F. (eds.), Compendium of methods for the microbiological examination of foods, 3rd ed. American Public Health Association, Washington, D.C.

## Packaging

Salmonella H Antiserum a	3 ml	2820-47
Salmonella H Antiserum b	3 ml	2821-47
Salmonella H Antiserum c	3 ml	2822-47
Salmonella H Antiserum d	3 ml	2823-47
Salmonella H Antiserum eh	3 ml	2273-47
Salmonella H Antiserum f	3 ml	2544-47
Salmonella H Antiserum h	3 ml	2545-47
Salmonella H Antiserum I	3 ml	2824-47
Salmonella H Antiserum k	3 ml	2274-47
Salmonella H Antiserum m	3 ml	2546-47
Salmonella H Antiserum p	3 ml	2548-47
Salmonella H Antiserum r	3 ml	2275-47
Salmonella H Antiserum s	3 ml	2550-47
Salmonella H Antiserum t	3 ml	2551-47
Salmonella H Antiserum w	3 ml	2554-47
Salmonella H Antiserum x	3 ml	2555-47
Salmonella H Antiserum y	3 ml	2276-47
Salmonella H Antiserum z	3 ml	2277-47
Salmonella H Antiserum Z <sub>6</sub>	3 ml	2473-47
Salmonella H Antiserum Z <sub>10</sub>	3 ml	2279-47
Salmonella H Antiserum Z <sub>13</sub>	3 ml	2556-47
Salmonella H Antiserum Z <sub>15</sub>	3 ml	2557-47
Salmonella H Antiserum Z <sub>23</sub>	3 ml	2558-47
Salmonella H Antiserum Z <sub>28</sub>	3 ml	2561-47
Salmonella H Antiserum Z <sub>29</sub>	3 ml	2280-47
Salmonella H Antiserum Z <sub>32</sub>	3 ml	2562-47
Salmonella H Antiserum EN Complex	3 ml	2270-47
Salmonella H Antiserum G Complex	3 ml	2269-47
Salmonella H Antiserum L Complex	3 ml	2271-47
Salmonella H Antiserum Z <sub>4</sub> Complex	3 ml	2278-47
Salmonella H Antiserum Poly a-z	3 ml	2406-47
Salmonella H Antiserum Poly A (a,b,c,d,i,Z <sub>10</sub> ,Z <sub>29</sub> )	3 ml	2539-47
Salmonella H Antiserum Poly B (eh,en,enz <sub>15</sub> , and G Complex)	3 ml	2540-47
Salmonella H Antiserum Poly C (k,l,r,y,Z <sub>1</sub> ,Z <sub>4</sub> )	3 ml	2541-47
Salmonella H Antiserum Poly D (Z <sub>35</sub> ,Z <sub>36</sub> ,Z <sub>37</sub> ,Z <sub>38</sub> ,Z <sub>39</sub> ,Z <sub>41</sub> ,Z <sub>42</sub> )	3 ml	2542-47
Salmonella H Antiserum Poly E (I Complex, z <sub>6</sub> )	3 ml	2543-47
Salmonella H Antiserum Single Factor 2	3 ml	2474-47

Salmonella H Antiserum Single Factor 5	3 ml	2475-47	Salmonella O Antiserum Group E Factors 1,3,10,15,19,34	3 ml	2819-47
Salmonella H Antiserum Single Factor 6	3 ml	2476-47	Salmonella O Antiserum Group E <sub>1</sub> Factors 3,10	3 ml	2952-47
Salmonella H Antiserum Single Factor 7	3 ml	2477-47	Salmonella O Antiserum Group E <sub>2</sub> Factors 3,15	3 ml	2954-47
Salmonella H Antiserum Spicer-Edwards 1	3 ml	2265-47	Salmonella O Antiserum Group E <sub>3</sub> Factors (3),(15),34	3 ml	3018-47
Salmonella H Antiserum Spicer-Edwards 2	3 ml	2266-47	Salmonella O Antiserum Group E <sub>4</sub> Factors 1,3,19	3 ml	3019-47
Salmonella H Antiserum Spicer-Edwards 3	3 ml	2267-47	Salmonella O Antiserum Group F Factor 11	3 ml	2260-47
Salmonella H Antiserum Spicer-Edwards 4	3 ml	2268-47	Salmonella O Antiserum Group G Factors 13,22,23,(36),(37)	3 ml	3029-47
Salmonella H Antiserum 1 Complex	3 ml	2272-47	Salmonella O Antiserum Group G <sub>1</sub> Factors 13,22,(36)	3 ml	2261-47
Salmonella O Antiserum Factor 2	3 ml	2814-47	Salmonella O Antiserum Group G <sub>2</sub> Factors 1,13,23,(37)	3 ml	3020-47
Salmonella O Antiserum Factor 4	3 ml	2659-47	Salmonella O Antiserum Group H Factors 1,6,14,24,25,47	3 ml	2262-47
Salmonella O Antiserum Factors 4,5	3 ml	2815-47	Salmonella O Antiserum Group I Factor 16	3 ml	2263-47
Salmonella O Antiserum Factor 5	3 ml	2660-47	Salmonella O Antiserum Group J Factor 17	3 ml	2517-47
Salmonella O Antiserum Factor 7	3 ml	2816-47	Salmonella O Antiserum Group K Factor 18	3 ml	2518-47
Salmonella O Antiserum Factor 8	3 ml	2817-47	Salmonella O Antiserum Group L Factor 21	3 ml	2519-47
Salmonella O Antiserum Factor 9	3 ml	2818-47	Salmonella O Antiserum Group M Factor 28	3 ml	2520-47
Salmonella O Antiserum Factor 10	3 ml	2257-47	Salmonella O Antiserum Group N Factor 30	3 ml	2521-47
Salmonella O Antiserum Factor 12	3 ml	2779-47	Salmonella O Antiserum Group O Factor 35	3 ml	2522-47
Salmonella O Antiserum Factor 14	3 ml	2661-47	Salmonella O Antiserum Poly A-I & Vi	3 ml	2264-47
Salmonella O Antiserum Factor 15	3 ml	2258-47	Salmonella O Antiserum Poly A (A,B,D,E <sub>1</sub> ,E <sub>2</sub> ,E <sub>3</sub> ,E <sub>4</sub> , & L)	3 ml	2534-47
Salmonella O Antiserum Factor 19	3 ml	2259-47	Salmonella O Antiserum Poly B (C <sub>1</sub> ,C <sub>2</sub> ,F,G, & H)	3 ml	2535-47
Salmonella O Antiserum Factor 20	3 ml	2662-47	Salmonella O Antiserum Poly C (I,J,K,M,N,& O)	3 ml	2536-47
Salmonella O Antiserum Factor 22	3 ml	2663-47	Salmonella O Antiserum Poly D (P,Q,R,S,T,& U)	3 ml	2537-47
Salmonella O Antiserum Factor 23	3 ml	2664-47	Salmonella O Antiserum Poly E (V,W,X,Y, & Z)	3 ml	2538-47
Salmonella O Antiserum Factor 25	3 ml	2666-47	Salmonella O Antiserum Poly F (Groups 51-55)	3 ml	2645-47
Salmonella O Antiserum Factor 27	3 ml	2667-47	Salmonella O Antiserum Poly G (Groups 56-61)	3 ml	2646-47
Salmonella O Antiserum Factor 34	3 ml	2512-47	Salmonella O Antiserum Vi	3 ml	2827-47
Salmonella O Antiserum Group A Factors 1,2,12	3 ml	2947-47			
Salmonella O Antiserum Group B Factors 1,4,5,12	3 ml	2948-47			
Salmonella O Antiserum Group B Factors 1,4,12,27	3 ml	2973-47			
Salmonella O Antiserum Group C <sub>1</sub> Factors 6,7	3 ml	2949-47			
Salmonella O Antiserum Group C <sub>2</sub> Factors 6,8	3 ml	2950-47			
Salmonella O Antiserum Group C <sub>3</sub> Factors (8), 20	3 ml	3016-47			
Salmonella O Antiserum Group D <sub>1</sub> Factors 1,9,12	3 ml	2951-47			
Salmonella O Antiserum Group D <sub>2</sub> Factors (9),46	3 ml	3017-47			

# Salmonella, Antigenic Scheme

## Update of the Kauffmann-White Scheme

The Centers for Disease Control<sup>1</sup> has modified the Kauffmann-White<sup>2</sup> Antigenic Scheme originally proposed by Ewing.<sup>3</sup> The updated scheme is used with Difco *Salmonella* Antisera as an aid in the serological identification of *Salmonella*.

All of the *Salmonella* serovars belong to two species: *S. bongori* containing 18 serovars and *S. enterica* containing the remaining 2300+ serovars divided among 6 subspecies.<sup>1,4</sup>

The six subspecies of *S. enterica* are:

- S. enterica* subsp. *enterica* (I or 1)
- S. enterica* subsp. *salamae* (II or 2)
- S. enterica* subsp. *arizonae* (IIIa or 3a)
- S. enterica* subsp. *diarizonae* (IIIb or 3b)
- S. enterica* subsp. *houtenae* (IV or 4)
- S. enterica* subsp. *indica* (VI or 6)

The legitimate species name for the above subspecies is *S. choleraesuis*. However, this name may be confused with the serotype named “*choleraesuis*.” At the International Congress for Microbiology in 1986, the International Subcommittee for *Enterobacteriaceae* agreed to adopt the species name “*S. enterica*.”<sup>5</sup> LeMinor and Popoff<sup>6</sup> published a request for the use of *S. enterica* as a species name to the Judicial Commission. The Judicial Commission ruled that *S. choleraesuis* is the legitimate name.<sup>7,8</sup> *S. enterica* is used in many countries and is favorably accepted as the species name.<sup>3,9</sup> The Centers for Disease Control has adopted this designation until the problem of naming this species is resolved.<sup>1</sup>

Nomenclature and classification of these bacteria are ever changing.<sup>13</sup> *Salmonella* and the former *Arizona* should be considered a single genus - *Salmonella*.<sup>9</sup> It is recommended that laboratories report names of *Salmonella* serovars for the subspecies *enterica*. These serovar names are no longer italicized.<sup>1,4</sup> The first letter of the serovar name begins with a capital letter. For example, the strain that used to be identified as “*Salmonella typhimurium*” is now known as “*Salmonella* Typhimurium.”

Serovars of other subspecies of *S. enterica* (except some in subspecies *salamae* and *houtenae*) and those of *S. bongori* are not named, and are designated by their antigenic formula. For the most recent information on nomenclature, consult appropriate references.<sup>1,4,8,10-14</sup>

Serotypes of *Salmonella* are defined based on the antigenic structure of both somatic or cell wall (O) antigens and flagellar (H) antigens. The O antigen groups were first designated by letters of the alphabet. Additional antigens were later delineated. Since each letter of the alphabet was already used to describe an O antigen group, numbers 51 to 67 were used to describe the next O antigen groups.

O ANTIGEN GROUP	O ANTIGEN FACTORS PRESENT
A	1,2,12
B	4,12; 1,4,5,12; or 1,4,12,27
C <sub>1</sub>	6,7,[Vi] or 6,7,14
C <sub>2</sub>	6,8
C <sub>3</sub>	8; or 8,20
D <sub>1</sub>	1,9,12
D <sub>2</sub>	9,46
D <sub>3</sub>	1,9,12,46,27
E <sub>1</sub>	3,10
E <sub>2</sub>	3,15
E <sub>3</sub>	3,15,34
E <sub>4</sub>	1,3,19
F	11
G	13,22 or 13,23
H	6,14; 6,14,24; or 1,6,14,25
I	16
J	17
K	18
L	21
M	28
N	30
O	35
P	38
Q	30
R	40
S	41
T	42
U	43
V	44
W	45
X	47
Y	48
Z	50
51	51
52	52
53	53
54	54
55	55
56	56
57	57
58	58
59	59
60	60
61	61
62	62
63	63
65	65
66	66
67	67

When writing an antigenic formula, list the O antigen(s) first followed by the H antigen(s). Separate the major antigens by colons and the components of the antigens by commas. For example, the antigenic formula for *Salmonella* Typhimurium is *Salmonella* 1,4,5,12:i:1,2. This means that the strain has O antigen factors 1,4,5 and 12, the flagella phase 1 antigen i, and flagellar phase 2 antigens 1 and 2.

Complete identification of *Salmonella* requires cultural isolation, biochemical characterization and serotyping. Any serological results obtained before biochemical identification must be considered as presumptive identification only. Refer to appropriate references for complete identification of *Salmonella*.<sup>1,3,4,10,12-15</sup>

The Kauffmann-White Scheme is presented in two forms. Appendix A contains a list of *Salmonella* Serotypes by O Group. Appendix B contains an alphabetical list of *Salmonella* Serotypes.

## References

1. **McWhorter-Murlin, A. C. and F. W. Hickman-Brenner.** 1994. Identification and Serotyping of *Salmonella* and an update of the Kauffmann-White Scheme. Centers for Disease Control and Prevention, Atlanta, GA.
2. **Kauffmann, F.** 1969. *Enterobacteriaceae*, 2nd ed. Munksgaard, Copenhagen.
3. **Ewing, W. H.** 1986. Edwards and Ewing's identification of *Enterobacteriaceae*, 4th ed. Elsevier Science Publishing Co. Inc., New York, NY.
4. **Popoff, M. Y. and L. LeMinor.** 1997. Antigenic Formulas of the *Salmonella* Serovars. WHO Collaborating Centre for Reference and Research on *Salmonella*. Institut Pasteur, Paris, France.
5. **Penner, J. L.** 1988. International committee on systematic bacteriology taxonomic subcommittee on *Enterobacteriaceae*. Int. J. Syst. Bacteriol. **38**:223-224.
6. **LeMinor, L. and M. Y. Popoff.** 1987. Request for an opinion. Designation of *Salmonella enterica* sp. nov., nom. rev., as the type and only species of the genus *Salmonella*. Int. J. Syst. Bacteriol. **37**:465-468.
7. **Wayne, L. G.** 1991. Judicial Commission of the International Committee on Systematic Bacteriology. Int. J. Syst. Bacteriol. **41**:185-187.
8. **Wayne, L. G.** 1994. Actions of the Judicial Commission of the International Committee on Systematic Bacteriology on requests for opinions published between January 1985 and July 1993. Int. J. Syst. Bacteriol. **44**:177.
9. **Old, D. C.** 1992. Nomenclature of *Salmonella*. J. Med. Microbiol. **37**:361-363.
10. **Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.).** 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
11. **Farmer, J. J., III, A. C. McWhorter, D. J. Brenner, and G. D. Morris.** 1984. The *Salmonella*-*Arizona* group of *Enterobacteriaceae*: nomenclature, classification and reporting. Clin. Microbiol. Newsl. **6**:63-66.
12. **Isenberg, H. D. (ed.)** 1992. Clinical microbiology procedures handbook, vol. 2. American Society for Microbiology, Washington, D.C.
13. **Holt, J. G., N. R. Krieg, P. H. Sneath, J. T. Staley, S. T. Williams.** 1994. Bergey's manual of determinative bacteriology, 9th ed. Williams & Wilkins, Baltimore, MD.
14. **Andrews, W. H., G. A. June, P. Sherrod, T. S. Hammack, and R. M. Amaguana.** 1995. Food and drug administration bacteriological analytical manual, 8th edition. AOAC International, Gaithersburg, MD.
15. **Russell, S. F., J. D'Aoust, W. H. Andrews, and J. S. Bailey.** 1992. *Salmonella*. In Vanderzant, C. and D. F. Splittstoesser, (eds.), Compendium of methods for the microbiological examination of foods, 3rd ed. American Public Health Association, Washington, D.C.
16. **Popoff, M. Y., and L. LeMinor.** 1992. Antigenic formulas of the *Salmonella* serovars, 6<sup>th</sup> revision. WHO Collaborating Centre for Reference and Research on *Salmonella*. Pasteur Institute, Paris, France.
17. **Rohde, R.** 1979. Serological integration of all known *Arizona* species into the Kauffmann-White scheme. Zentralbl. Bakteriol. Hyg, I. Abt. Orig. A. **243**:148-176.
18. **LeMinor, L.** 1984. Genus III. *Salmonella*. In N. R. Krieg, (ed.), Bergey's manual of systematic bacteriology, vol. 1. The Williams and Wilkins Co., Baltimore, MD.

## Appendix A

### Kauffmann-White Scheme

#### List of *Salmonella* Serotypes by O Group (Updated 1994)

Appendix A contains a list of *Salmonella* Serotypes by O Group. The serotypes are sorted by O group first, then by Phase I and Phase 2 of the H antigens. The z antigens do not appear in the correct numerical order - z<sub>4</sub> will be listed after z<sub>10</sub> because the 1 of 10 is read first, and will appear after z<sub>29</sub>.

#### Key:

**IP** Institut Pasteur. See reference #16 in Kauffmann-White Scheme monograph published by WHO Collaborating Centre for Reference and Research on *Salmonella*.

**Ar.** "Arizona" antigenic formula

**Rohde** Refer to reference #17 in Kauffmann-White Scheme monograph by R. Rohde. He incorporated all known *Arizona* serotypes into the Kauffmann-White Scheme.

**Bergey** Refer to reference #18 in Kauffmann-White Scheme monograph by L. LeMinor for "Bergey's Manual of Systematic Bacteriology."

#### Underlined

**Numbers** Numbers that are underlined in a serotype represent somatic factors determined by phage conversion. They are present if the culture is lysogenized by the corresponding converting phage.

[ ] O or H factors may be present or absent without relation to phage conversion.

( ) O or H factor is weakly agglutinable.

**R phases** Abnormal specificities of H antigens that were described by Kauffmann. They are uncommon.

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	A	Paratyphi A	<u>1,2,12</u>	a	[1,5]	
I	A	Nitra	2,12	g,m	-	
I	A	Kiel	<u>1,2,12</u>	g,p	-	
I	A	Koessen	2,12	l,v	1,5	
II	B		4,12	-	1,6	
I	B	Abortusequi	4,12	-	e,n,x	
I	B	Kisangani	<u>1,4,[5],12</u>	a	1,2	
I	B	Fulica	4,[5],12	a	1,5	
I	B	Hessarek	4,12, <u>27</u>	a	1,5	
I	B	Arechavaleta	4,[5],12	a	[1,7]	
I	B	Bispejerg	<u>1,4,[5],12</u>	a	e,n,x	
II	B		<u>1,4,[5],12,<u>27</u></u>	a	e,n,x	
II	B	Makoma	<u>1,4,[5],12,<u>27</u></u>	a	[e,n,x]	
I	B	Tinda	<u>1,4,5,<u>27</u></u>	a	e,n,z <sub>15</sub>	
I	B	Huettwilen	<u>1,4,12</u>	a	l,w	
II	B		<u>1,4,12,<u>27</u></u>	a	z <sub>39</sub>	
I	B	Nakuru	<u>1,4,12,27</u>	a	z <sub>6</sub>	
I	B	Schleissheim	4,12, <u>27</u>	b	-	
I	B	Java	<u>1,4,5,12</u>	b	[1,2], (tartrate +)	IP calls Java, Paratyphi B var. Java. Java is often monophasic in the U.S. May possess H phase Rz <sub>33</sub> .
I	B	Limete	<u>1,4,12,<u>27</u></u>	b	1,5	
II	B		4,12	b	1,5	
I	B	Canada	4,12, <u>27</u>	b	1,6	
I	B	Uppsala	4,12, <u>27</u>	b	1,7	
I	B	Abony	1,4,5,12	b	e,n,x	IP combined Sladun (1,4,12,27:b:e,n,x) with Abony to form Abony <u>1,4,[5],12,27:b:e,n,x</u> (gelatin neg.).
I	B	Abortusbovis	<u>1,4,12,27</u>	b	e,n,x	Gelatin pos., mucate pos. <sup>1-4 days</sup> . Abortusbovis was combined with Abony (1,4,5,12:b:e,n,x), gelatin neg. The name Abortusbovis was dropped.
I	B	Sladun	<u>1,4,12,27</u>	b	e,n,x	IP combined Sladun with Abony (1,3,4,12:b:e,n,x) to form Abony <u>1,4,[5],12,27:b:e,n,x</u> . Sladun is now called Abony var. O 27+. The name Sladun has been dropped.
II	B	Sofia	<u>1,4,12,27</u>	b	[e,n,x]	
I	B	Wagenia	<u>1,4,12,27</u>	b	e,n,z <sub>15</sub>	
I	B	Wien	<u>1,4,12,27</u>	b	l,w	
I	B	Abortuscanis	4,5,12	b	Rz <sub>5</sub>	Abortuscanis was combined with Paratyphi B ( <u>1,4,[5],12:b:1,2</u> ) in 1938 and the name Abortuscanis was dropped.
I	B	Tripoli	<u>1,4,12,27</u>	b	z <sub>6</sub>	
I	B	Paratyphi B	1,4,[5],12	[b]	[1,2]	Paratyphi B is tartrate neg.; Paratyphi B var. Java (CDC calls this S. ser. Java) is often monophasic ( <u>1,4,5,12:b:-</u> ) and is tartrate pos. Paratyphi B and Java may possess H phase Rz <sub>33</sub> .
I	B	Legon	<u>1,4,12,27</u>	c	1,5	
I	B	Abortusovis	4,12	c	1,6	
I	B	Altendorf	4,12	c	1,7	IP combined Womba (4,12,27:c:1,7) with Altendorf to form Altendorf <u>4,12,27:c:1,7</u> .
I	B	Womba	4,12, <u>27</u>	c	1,7	IP combined Womba with Altendorf (4,12:c:1,7) to form Altendorf <u>4,12,27:c:1,7</u> . Womba is now called Altendorf var. 027+. The name Womba has been dropped.

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	B	Bissau	4,12	c	e,n,x	
I	B	Jericho	<u>1,4,12,27</u>	c	e,n,z <sub>15</sub>	
I	B	Hallfold	<u>1,4,12,27</u>	c	l,w	
I	B	Bury	4,12,27	c	z <sub>6</sub>	
I	B	Cairo	<u>1,4,12,27</u>	d	1,2	IP combined Cairo with Stanley (4,5,12:d:1,2) to form Stanley <u>1,4,[5],12,27</u> :d:1,2. The name Cairo has been dropped.
I	B	Stanley	4,5,12	d	1,2	IP combined Cairo ( <u>1,4,12,27</u> :d:1,2) with Stanley to form Stanley <u>1,4,[5],12,27</u> :d:1,2
I	B	Eppendorf	<u>1,4,12,27</u>	d	1,5	
I	B	Brezany	<u>1,4,12,27</u>	d	1,6	
I	B	Schwarzengrund	<u>1,4,12,27</u>	d	1,7	
I	B	Sarajane	<u>1,4,[5],12,27</u>	d	e,n,x	
II	B	Kluetjenfelde	4,12	d	e,n,x	
I	B	Duisburg	<u>1,4,12,27</u>	d	e,n,z <sub>15</sub>	
I	B	Mons	<u>1,4,12,27</u>	d	1,w	
I	B	Ayinde	<u>1,4,12,27</u>	d	z <sub>6</sub>	
I	B	Salinatis	4,12	d,e,h	d,e,n,z <sub>15</sub>	IP states that Salinatis was combined with Duisburg ( <u>1,4,12,27</u> :d:e,n,z <sub>15</sub> ). This is incorrect; IP should have stated that it was combined with Sandiego (4,[5],12:e,h:e,n,z <sub>15</sub> ), because Salinatis loses the d and becomes Sandiego.
I	B	Saintpaul	<u>1,4,[5],12</u>	e,h	1,2	
I	B	Reading	<u>1,4,[5],12</u>	e,h	[1,5]	
I	B	Eko	4,12	e,h	1,6	
I	B	Kaapstad	4,12	e,h	1,7	
I	B	Chester	1,4,[5],12	e,h	e,n,x	
I	B	Sandiego	4,[5],12	[e,h]	e,n,z <sub>15</sub>	
II	B		4,12	e,n,x	1,2,7	
II	B	Makumira	<u>1,4,12,27</u>	e,n,x	1,[5],7	
II	B		4,12	(f),g	–	Not in IP book
I	B	Derby	<u>1,4,[5],12</u>	f,g	[1,2]	
I	B	Agona	<u>1,4,[5],12</u>	f,g,s	[1,2]	
II	B		4,[5],12	f,g,t	z <sub>6</sub> ,z <sub>42</sub>	
I	B	Essen	4,12	g,m	–	
I	B	Hato	4,[5],12	g,m,s	–	
I	B	California	4,12	g,m,t	–	
II	B		4,12	g,m,t	–	IP calls this monophasic var. of Bechuana.
II	B		4,12	g,m,t	z <sub>39</sub>	
II	B	Caledon	<u>1,4,12,27</u>	g,[m],[s],t	e,n,x	
II	B	Bechuana	<u>1,4,12,27</u>	g,[m],t	[1,5]	
I	B	Joenkoeping	4,5,12	g,s,t	–	IP combined Joenkoeping with Kingston ( <u>1,4,12,27</u> :g,s,t:-) to form Kingston <u>1,4,[5],12,27</u> :g,s,t:-. The name Joenkoeping has been dropped.
I	B	Kingston	<u>1,4,12,27</u>	g,s,t	[1,2]	IP combined Kingston with Joenkoeping (4,5,12:g,s,t:-) to form Kingston <u>1,4,[5],12,27</u> :g,s,t:[1,2]. Kingston may possess H phase Rz <sub>27</sub> or Rz <sub>43</sub> .
I	B	Budapest	<u>1,4,12,27</u>	g,t	–	
I	B	Travis	4,[5],12	g,z <sub>51</sub>	1,7	
I	B	Tennyson	4,5,12	g,z <sub>51</sub>	e,n,z <sub>15</sub>	
II	B		4,12	g,z <sub>62</sub>	–	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	B	Typhimurium var. Copenhagen	<u>1</u> ,4,12	i	1,2	
I	B	Typhimurium	<u>1</u> ,4,5,12	i	1,2,[7]	
I	B	Lagos	<u>1</u> ,4,[5],12	i	1,5	
I	B	Agama	4,12	i	1,6	
I	B	Farsta	4,12	i	e,n,x	
I	B	Tsevie	4,12	i	e,n,z <sub>15</sub>	
I	B	Gloucester	<u>1</u> ,4,12, <u>27</u>	i	1,w	
II	B		4,12, <u>27</u>	i	z <sub>35</sub>	
I	B	Tumodi	<u>1</u> ,4,12	i	z <sub>6</sub>	
I	B	Massenia	<u>1</u> ,4,12, <u>27</u>	k	1,5	
I	B	Neumuenster	<u>1</u> ,4,12, <u>27</u>	k	1,6	
II	B		<u>1</u> ,4,12, <u>27</u>	k	1,6	
I	B	Ljubljana	4,12, <u>27</u>	k	e,n,x	
I	B	Texas	4,[5],12	k	e,n,z <sub>15</sub>	
I	B	Fyris	4,[5],12	l,v	1,2	
I	B	Azteca	4,[5],12, <u>27</u>	l,v	1,5	
I	B	Bredeney	<u>1</u> ,4,12, <u>27</u>	l,v	1,7	Bredeney may possess H phase Rl,z <sub>40</sub> instead of l,v.
I	B	Kimuenza	<u>1</u> ,4,12, <u>27</u>	l,v	e,n,x	
II	B		<u>1</u> ,4,12, <u>27</u>	l,v	e,n,x	
I	B	Brandenburg	<u>1</u> ,4,[5],12, <u>27</u>	l,v	e,n,z <sub>15</sub>	
II	B		<u>1</u> ,4,12, <u>27</u>	l,v	z <sub>39</sub>	
I	B	Clackamas	4,12	l,v,[z <sub>13</sub> ]	1,6	IP does not get z <sub>13</sub>
I	B	Mono	4,12	l,w	1,5	
I	B	Togo	4,12	l,w	1,6	
II	B	Kilwa	4,12	l,w	e,n,x	
I	B	Ayton	<u>1</u> ,4,12, <u>27</u>	l,w	z <sub>6</sub>	
I	B	Haduna	4,12	l,z <sub>13</sub> ,[z <sub>28</sub> ]	1,6	
I	B	Kubacha	<u>1</u> ,4,12, <u>27</u>	l,z <sub>13</sub> ,z <sub>28</sub>	1,7	
I	B	Kano	<u>1</u> ,4,12, <u>27</u>	l,z <sub>13</sub> ,z <sub>28</sub>	e,n,x	
I	B	Vom	<u>1</u> ,4,12, <u>27</u>	l,[z <sub>13</sub> ],[z <sub>28</sub> ]	e,n,z <sub>15</sub>	
I	B	Tyresoe	4,12	l,[z <sub>13</sub> ],z <sub>28</sub>	1,5	
I	B	Kunduchi	<u>1</u> ,4,[5],12, <u>27</u>	l,[z <sub>13</sub> ],[z <sub>28</sub> ]	[1,2]	
II	B		4,12	l,z <sub>28</sub>	–	
I	B	Reinickendorf	4,12	l,z <sub>28</sub>	e,n,x	
I	B	Banana	<u>1</u> ,4,[5],12	m,t	[1,5]	Banana may possess H phase Rz <sub>45</sub> .
I	B	Madras	4,[5],12	m,t	e,n,z <sub>15</sub>	
I	B	Heidelberg	<u>1</u> ,4,[5],12	r	1,2	
I	B	Bradford	4,12, <u>27</u>	r	1,5	
I	B	Winneba	4,12	r	1,6	
I	B	Remo	<u>1</u> ,4,12, <u>27</u>	r	1,7	
I	B	Bochum	4,[5],12	r	1,w	
I	B	Southampton	<u>1</u> ,4,12, <u>27</u>	r	z <sub>6</sub>	
I	B	Africana	4,12	r,i	1,w	
I	B	Drogana	<u>1</u> ,4,12, <u>27</u>	r,(i)	e,n,z <sub>15</sub>	IP calls Drogana r,i.
I	B	Coeln	4,[5],12	y	1,2	
I	B	Trachau	4,12, <u>27</u>	y	1,5	
I	B	Finaghy	4,12	y	1,6	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	B	Teddington	<u>1,4,12,27</u>	y	1,7	
I	B	Ball	<u>1,4,12,27</u>	y	e,n,x	IP combined Ball with Ruki (4,5,12:y:e,n,x) and Dalat (4,5,27:y:e,n,x) to form Ball <u>1,4,[5],12,27:y:e,n,x</u> .
I	B	Ruki	4,5,12	y	e,n,x	IP combined Ruki with Ball ( <u>1,4,12,27:y:e,n,x</u> ) and Dalat (4,5,27:y:e,n,x) to form Ball <u>1,4,[5],12,27:y:e,n,x</u> . The name Ruki has been dropped.
I	B	Dalat	4,5,27	y	e,n,x	Dalat was combined with Ball. The name Dalat has been dropped.
I	B	Jos	<u>1,4,12,27</u>	y	e,n,z <sub>15</sub>	
I	B	Kamoru	4,12, <u>27</u>	y	z <sub>6</sub>	
I	B	Shubra	4,[5],12	z	1,2	
I	B	Kiambu	4,12	z	1,5	
II	B		<u>1,4,12,27</u>	z	1,5	
I	B	Loubomo	4,12	z	1,6	
I	B	Indiana	<u>1,4,12</u>	z	1,7	
II	B		4,12	z	1,7	
I	B	Neftenbach	4,12	z	e,n,x	
II	B	Nordenham	<u>1,4,12,27</u>	z	e,n,x	
I	B	Koenigstuhl	<u>1,4,[5],12</u>	z	e,n,z <sub>15</sub>	
I	B	Preston	<u>1,4,12</u>	z	1,w	
II	B		4,12	z	z <sub>39</sub>	
I	B	Entebbe	<u>1,4,12,27</u>	z	z <sub>6</sub>	
I	B	Haifa	<u>1,4,[5],12</u>	z <sub>10</sub>	1,2	
I	B	Ituri	<u>1,4,12</u>	z <sub>10</sub>	1,5	
I	B	Tudu	4,12	z <sub>10</sub>	1,6	
I	B	Albert	4,12	z <sub>10</sub>	e,n,x	
I	B	Tokoin	4,12	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	B	Mura	<u>1,4,12</u>	z <sub>10</sub>	1,w	
I	B	Vellore	<u>1,4,12,27</u>	z <sub>10</sub>	z <sub>35</sub>	
I	B	Fortune	<u>1,4,12,27</u>	z <sub>10</sub>	z <sub>6</sub>	
I	B	Brancaster	<u>1,4,12,27</u>	z <sub>29</sub>	–	
II	B	Helsinki	<u>1,4,12</u>	z <sub>29</sub>	[e,n,x]	
I	B	Pasing	4,12	z <sub>35</sub>	1,5	
I	B	Tafo	<u>1,4,12,27</u>	z <sub>35</sub>	1,7	
I	B	Yaounde	<u>1,4,12,27</u>	z <sub>35</sub>	e,n,z <sub>15</sub>	
I	B	Sloterdijk	<u>1,4,12,27</u>	z <sub>35</sub>	z <sub>6</sub>	
I	B	Tejas	4,12	z <sub>36</sub>	–	
I	B	Wilhelmsburg	<u>1,4,[5],12,27</u>	z <sub>38</sub>	[e,n,z <sub>15</sub> ]	
II	B	Durbanville	<u>1,4,12,27</u>	[z <sub>39</sub> ]	1,[5],7	
I	B	Jaja	4,12, <u>27</u>	z <sub>4</sub> ,z <sub>23</sub>	–	IP combined Jaja with Stanleyville ( <u>1,4,[5],12:z<sub>4</sub>,z<sub>23</sub>:1,5</u> ) to form Stanleyville <u>1,4,[5],12,27:z<sub>4</sub>,z<sub>23</sub>:1,5</u> . Jaja is now called Stanleyville var. O 27+. The name Jaja has been dropped.
I	B	Stanleyville	<u>1,4,[5],12</u>	z <sub>4</sub> ,z <sub>23</sub>	[1,2]	IP combined Jaja (4,12, <u>27:z<sub>4</sub>,z<sub>23</sub>:–</u> ) with Stanleyville to form Stanleyville <u>1,4,[5],12,27:z<sub>4</sub>,z<sub>23</sub>:1,2</u> .
I	B	Vuadens	4,12, <u>27</u>	z <sub>4</sub> ,z <sub>23</sub>	z <sub>6</sub>	
I	B	Kalamu	4,[5],12	z <sub>4</sub> ,z <sub>24</sub>	[1,5]	
I	B	Thayngen	<u>1,4,12,27</u>	z <sub>41</sub>	1,(2),5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	B	Maska	<u>1,4,12,27</u>	Z <sub>41</sub>	e,n,Z <sub>15</sub>	
II	C <sub>1</sub>		6,7	–	1,6	
I	C <sub>1</sub>	Sanjuan	6,7	a	1,5	
II	C <sub>1</sub>		6,7, <u>14</u>	a	1,5	
I	C <sub>1</sub>	Umhlali	6,7	a	1,6	
I	C <sub>1</sub>	Austin	6,7	a	1,7	
I	C <sub>1</sub>	Oslo	6,7, <u>14</u>	a	e,n,x	
I	C <sub>1</sub>	Denver	6,7	a	e,n,Z <sub>15</sub>	
I	C <sub>1</sub>	Coleypark	6,7, <u>14</u>	a	1,w	
II	C <sub>1</sub>	Calvinia	6,7	a	Z <sub>42</sub>	
I	C <sub>1</sub>	Damman	6,7	a	Z <sub>6</sub>	
II	C <sub>1</sub>		6,7	a	Z <sub>6</sub>	
I	C <sub>1</sub>	Nissii	6,7, <u>14</u>	b	–	Nissii was combined with Nienstedten (6,7, <u>14</u> :b:l,w) as a monophasic variant of Nienstedten. Nienstedten is now called a variant of Ohio by IP. The name Nissii has been dropped.
I	C <sub>1</sub>	Brazzaville	6,7	b	1,2	
I	C <sub>1</sub>	Edinburg	6,7	b	1,5	
I	C <sub>1</sub>	Adime	6,7	b	1,6	
I	C <sub>1</sub>	Koumra	6,7	b	1,7	
I	C <sub>1</sub>	Lockleaze	6,7, <u>14</u>	b	e,n,x	
II	C <sub>1</sub>	Bloemfontein	6,7	b	[e,n,x]:Z <sub>42</sub>	
I	C <sub>1</sub>	Georgia	6,7	b	e,n,Z <sub>15</sub>	
I	C <sub>1</sub>	Ohio	6,7	b	1,w	IP combined Nienstedten (6,7, <u>14</u> :b:[1,w]) with Ohio to form Ohio 6,7, <u>14</u> :b:[1,w]. Ohio may possess H phase R <sub>Z59</sub> .
I	C <sub>1</sub>	Nienstedten	6,7, <u>14</u>	b	[1,w]	Nienstedten was combined with Nissii (6,7, <u>14</u> :b:-) and called Nienstedten; then IP combined Nienstedten with Ohio (6,7:b:l,w) to form Ohio (6,7, <u>14</u> :b:[1,w]). Nienstedten is now called Ohio var. O 14+ by IP.
I	C <sub>1</sub>	Kotte	6,7	b	Z <sub>35</sub>	
II	C <sub>1</sub>		6,7	b	Z <sub>39</sub>	
I	C <sub>1</sub>	Leopoldville	6,7, <u>14</u>	b	Z <sub>6</sub>	
I	C <sub>1</sub>	Hissar	6,7, <u>14</u>	c	1,2	
I	C <sub>1</sub>	Choleraesuis	6,7	c	1,5	H <sub>2</sub> S negative
I	C <sub>1</sub>	Decatur	6,7	c	1,5	IP has dropped Decatur and calls it dulcitol positive, mucate positive variant of Choleraesuis.
I	C <sub>1</sub>	Paratyphi C	6,7,[Vi]	c	1,5	
I	C <sub>1</sub>	Typhisuis	6,7	c	1,5	Typhisuis is a bioserotype found in pigs. It is like Choleraesuis except tartrate negative.
I	C <sub>1</sub>	Choleraesuis var. Kunzendorf		6,7	c	[1,5] H <sub>2</sub> S positive
I	C <sub>1</sub>	Birkenhead	6,7	c	1,6	
I	C <sub>1</sub>	Schwabach	6,7	c	1,7	
I	C <sub>1</sub>	Namibia	6,7	c	e,n,x	
I	C <sub>1</sub>	Kaduna	6,7, <u>14</u>	c	e,n,Z <sub>15</sub>	
I	C <sub>1</sub>	Kisii	6,7	d	1,2	
I	C <sub>1</sub>	Kambole	6,7	d	1,[2],7	
I	C <sub>1</sub>	Isangi	6,7, <u>14</u>	d	1,5	
I	C <sub>1</sub>	Mission	6,7	d	1,5	Mission was combined with Isangi 6,7, <u>14</u> :d:1,5. The name Mission has been dropped.
I	C <sub>1</sub>	Kivu	6,7	d	1,6	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	C <sub>1</sub>	Amersfoort	6,7	d	e,n,x	IP combined Omderman (6,7,14:d:e,n,x) with Amersfoort to form Amersfoort 6,7,14:d:e,n,x.
I	C <sub>1</sub>	Omderman	6,7,14	d	e,n,x	IP combined Omderman with Amersfoort (6,7:d:e,n,x) to form Amersfoort 6,7,14:d:e,n,x. Omderman is now called Amersfoort var. O 14+ by IP.
I	C <sub>1</sub>	Gombe	6,7,14	d	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Eimsbuettel	6,7,14	d	l,w	IP combined Eimsbuettel with Livingstone (6,7:d:l,w) to form Livingstone 6,7,14:d:l,w. Eimsbuettel is now called Livingstone var. O 14+ by IP.
I	C <sub>1</sub>	Livingstone	6,7	d	l,w	IP combined Eimsbuettel (6,7,14:d:l,w) with Livingstone to form Livingstone 6,7,14:d:l,w.
I	C <sub>1</sub>	Wil	6,7	d	l,z <sub>13</sub> ,z <sub>28</sub>	
II	C <sub>1</sub>		6,7	d	z <sub>42</sub>	
I	C <sub>1</sub>	Nieuwerk	6,7,14	d	z <sub>6</sub>	
I	C <sub>1</sub>	Larochelle	6,7	e,h	1,2	
I	C <sub>1</sub>	Lomita	6,7	e,h	1,5	
I	C <sub>1</sub>	Norwich	6,7	e,h	1,6	
I	C <sub>1</sub>	Nola	6,7	e,h	1,7	
I	C <sub>1</sub>	Braenderup	6,7,14	e,h	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Ardwick	6,7,14	f,g	-	IP combined Ardwick with Rissen (6,7:f,g:-) to form Rissen 6,7,14:f,g:-. Ardwick is now called Rissen var. O 14+ by IP.
I	C <sub>1</sub>	Rissen	6,7	f,g	-	IP combined Ardwick (6,7,14:f,g:-) with Rissen for form Rissen 6,7,14:f,g:-.
I	C <sub>1</sub>	Eingedi	6,7	f,g,t	1,2,7	
I	C <sub>1</sub>	Afula	6,7	f,g,t	e,n,x	
II	C <sub>1</sub>		6,7	(g),m,[s],t	[1,5]	
I	C <sub>1</sub>	Montevideo	6,7,14	g,m,[p],s	[1,2,7]	
II	C <sub>1</sub>		6,7	g,m,[s],t	e,n,x	
I	C <sub>1</sub>	Othmarschen	6,7,14	g,m,[t]	-	
II	C <sub>1</sub>		6,7	g,[m],s,t	[z <sub>42</sub> ]	
I	C <sub>1</sub>	Menston	6,7	g,s,[t]	[1,6]	
I	C <sub>1</sub>	Riggil	6,7	g,t	-	
II	C <sub>1</sub>		6,7	g,t	e,n,x:z <sub>42</sub>	
IV	C <sub>1</sub>		6,7	g,z <sub>51</sub>	-	
I	C <sub>1</sub>	Alamo	6,7	g,z <sub>51</sub>	1,5	
I	C <sub>1</sub>	Augustenborg	6,7,14	i	1,2	
I	C <sub>1</sub>	Oritamerin	6,7	i	1,5	
I	C <sub>1</sub>	Garoli	6,7	i	1,6	
I	C <sub>1</sub>	Lika	6,7	i	1,7	
I	C <sub>1</sub>	Athinai	6,7	i	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Norton	6,7	i	l,w	
I	C <sub>1</sub>	Stuttgart	6,7,14	i	z <sub>6</sub>	
I	C <sub>1</sub>	Galiema	6,7,14	k	1,2	
I	C <sub>1</sub>	Daytona	6,7	k	1,6	
I	C <sub>1</sub>	Baiboukoum	6,7	k	1,7	
I	C <sub>1</sub>	Singapore	6,7	k	e,n,x	
I	C <sub>1</sub>	Escanaba	6,7	k	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Cardiff	6,7	k	R1,10	IP combined Cardiff that contains H phase R1,10 with Thompson (6,7,14:k:1,5).

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	C <sub>1</sub>		6,7	k	[z <sub>6</sub> ]	
IIIa	C <sub>1</sub>		6,7	(k)	z:[z <sub>55</sub> ]	(Ar. 27:22:31:[37])
I	C <sub>1</sub>	Thompson	6,7,14	[k]	[1,5]	IP combined Cardiff that contains H phase R1,10 (6,7:k:R1,10) with Thompson.
I	C <sub>1</sub>	Concord	6,7	l,v	1,2	
I	C <sub>1</sub>	Irumu	6,7	l,v	1,5	
I	C <sub>1</sub>	Mkamba	6,7	l,v	1,6	
I	C <sub>1</sub>	Kortrijk	6,7	l,v	1,7	
I	C <sub>1</sub>	Bonn	6,7	l,v	e,n,x	
I	C <sub>1</sub>	Potsdam	6,7,14	l,v	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Coromandel	6,7	l,v	z <sub>35</sub>	
IIIb	C <sub>1</sub>		6,7	l,v	z <sub>53</sub>	(Ar. 27:23:25)
I	C <sub>1</sub>	Gdansk	6,7	l,v	z <sub>6</sub>	IP combined Gelsenkirchen (6,7,14:l,v:z <sub>6</sub> ) with Gdansk to form Gdansk 6,7,14:l,v:z <sub>6</sub> .
I	C <sub>1</sub>	Gelsenkirchen	6,7,14	l,v	z <sub>6</sub>	IP combined Gelsenkirchen with Gdansk (6,7:l,v:z <sub>6</sub> ) to form Gdansk 6,7,14:l,v:z <sub>6</sub> . Gelsenkirchen is now called Gdansk var. O 14+ by IP.
I	C <sub>1</sub>	Gabon	6,7	l,w	1,2	
I	C <sub>1</sub>	Colorado	6,7	l,w	1,5	
II	C <sub>1</sub>		6,7	l,w	1,5,7	
II	C <sub>1</sub>		6,7	l,w	z <sub>42</sub>	
I	C <sub>1</sub>	Nessziona	6,7	l,z13	1,5	
I	C <sub>1</sub>	Kenya	6,7	l,z13	e,n,x	
I	C <sub>1</sub>	Strathcona	6,7	l,z13,z28	1,7	
I	C <sub>1</sub>	Makiso	6,7	l,z13,z28	z <sub>6</sub>	
I	C <sub>1</sub>	Neukoelin	6,7	l,z13,[z28]	e,n,z <sub>15</sub>	
II	C <sub>1</sub>	Heilbron	6,7	l,z28	1,5:[z <sub>42</sub> ]	
II	C <sub>1</sub>		6,7	l,z28	e,n,x	
II	C <sub>1</sub>		6,7	l,z28	z <sub>6</sub>	
I	C <sub>1</sub>	Haelsingborg	6,7	m,p,t,[u]	–	
I	C <sub>1</sub>	Oranienburg	6,7	m,t	–	IP combined Theilallee (6,7,14:m,t:-) with Oranienburg to form Oranienburg 6,7,14:m,t:-. Oranienburg may possess H phase Rz <sub>57</sub> .
I	C <sub>1</sub>	Theilallee	6,7,14	m,t	–	IP combined Theilallee with Oranienburg (6,7:m,t:-) to form Oranienburg 6,7,14:m,t:-. Theilallee is now called Oranienburg var. O 14+ by IP.
II	C <sub>1</sub>		6,7	m,t	–	
I	C <sub>1</sub>	Winston	6,7	m,t	1,6	
I	C <sub>1</sub>	Oakey	6,7	m,t	z <sub>64</sub>	
I	C <sub>1</sub>	Virchow	6,7	r	1,2	
I	C <sub>1</sub>	Infantis	6,7,14	r	1,5	Infantis may possess H phase Rz <sub>49</sub> .
I	C <sub>1</sub>	Nigeria	6,7	r	1,6	
I	C <sub>1</sub>	Colindale	6,7	r	1,7	
I	C <sub>1</sub>	Papuana	6,7	r	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Grampian	6,7	r	1,w	
I	C <sub>1</sub>	Richmond	6,7	y	1,2	
I	C <sub>1</sub>	Bareilly	6,7,14	y	1,5	
I	C <sub>1</sub>	Oyonnax	6,7	y	1,6	
I	C <sub>1</sub>	Gatow	6,7	y	1,7	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	C <sub>1</sub>	Hartford	6,7	y	e,n,x	Hartford may possess H phase R <sub>Z50</sub> or R <sub>Z67</sub> .
I	C <sub>1</sub>	Mikawasima	6,7,14	y	e,n,z <sub>15</sub>	Mikawasima may possess H phase R <sub>Z47</sub> or R <sub>Z50</sub> .
I	C <sub>1</sub>	Chile	6,7	z	1,2	
I	C <sub>1</sub>	Poitiers	6,7	z	1,5	
II	C <sub>1</sub>	Tosamanga	6,7	z	1,5	
I	C <sub>1</sub>	Oakland	6,7	z	1,6,[7]	
I	C <sub>1</sub>	Cayar	6,7	z	e,n,x	
I	C <sub>1</sub>	Businga	6,7	z	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Bruck	6,7	z	1,w	
II	C <sub>1</sub>		6,7	z	Z <sub>39</sub>	
II	C <sub>1</sub>	Oysterbeds	6,7	z	Z <sub>42</sub>	
II	C <sub>1</sub>		6,7	z	Z <sub>6</sub>	
I	C <sub>1</sub>	Menden	6,7	Z <sub>10</sub>	1,2	
I	C <sub>1</sub>	Inganda	6,7	Z <sub>10</sub>	1,5	
I	C <sub>1</sub>	Eschweiler	6,7	Z <sub>10</sub>	1,6	
I	C <sub>1</sub>	Ngili	6,7	Z <sub>10</sub>	1,7	
I	C <sub>1</sub>	Djugu	6,7	Z <sub>10</sub>	e,n,x	
I	C <sub>1</sub>	Mbandaka	6,7,14	Z <sub>10</sub>	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Jerusalem	6,7,14	Z <sub>10</sub>	1,w	
I	C <sub>1</sub>	Omuna	6,7	Z <sub>10</sub>	Z <sub>35</sub>	
II	C <sub>1</sub>		6,7	Z <sub>10</sub>	Z <sub>35</sub>	
I	C <sub>1</sub>	Redba	6,7	Z <sub>10</sub>	Z <sub>6</sub>	
II	C <sub>1</sub>		6,7	Z <sub>29</sub>	-	
I	C <sub>1</sub>	Tennessee	6,7,14	Z <sub>29</sub>	[1,2,7]	
I	C <sub>1</sub>	Tienba	6,7	Z <sub>35</sub>	1,6	
I	C <sub>1</sub>	Palime	6,7	Z <sub>35</sub>	e,n,z <sub>15</sub>	
IV	C <sub>1</sub>	Argentina	6,7	Z <sub>36</sub>	-	
I	C <sub>1</sub>	Tampico	6,7	Z <sub>36</sub>	e,n,z <sub>15</sub>	
II	C <sub>1</sub>	Bacongo	6,7	Z <sub>36</sub>	Z <sub>42</sub>	
I	C <sub>1</sub>	Bornum	6,7,14	Z <sub>38</sub>	-	IP combined Bornum with Lille (6,7:z <sub>38</sub> :-) to form Lille 6,7,14:z <sub>38</sub> :-. Bornum is now called Lille var. O 14+ by IP.
I	C <sub>1</sub>	Lille	6,7	Z <sub>38</sub>	-	IP combined Lille with Bornum (6,7,14:z <sub>38</sub> :-) to form Lille 6,7,14:z <sub>38</sub> :-.
I	C <sub>1</sub>	Rumford	6,7	Z <sub>38</sub>	1,2	
II	C <sub>1</sub>	Gilbert	6,7	Z <sub>39</sub>	1,5,7	
IV	C <sub>1</sub>	Roterberg	6,7	Z <sub>4</sub> ,Z <sub>23</sub>	-	
I	C <sub>1</sub>	Obogu	6,7	Z <sub>4</sub> ,Z <sub>23</sub>	1,5	
I	C <sub>1</sub>	Planckendael	6,7	Z <sub>4</sub> ,Z <sub>23</sub>	1,6	
I	C <sub>1</sub>	Goma	6,7	Z <sub>4</sub> ,Z <sub>23</sub>	Z <sub>6</sub>	
I	C <sub>1</sub>	Aequatoria	6,7	Z <sub>4</sub> ,Z <sub>23</sub>	[e,n,z <sub>15</sub> ]	
I	C <sub>1</sub>	Somone	6,7	Z <sub>4</sub> ,Z <sub>24</sub>	-	
IV	C <sub>1</sub>	Kralendyk	6,7	Z <sub>4</sub> ,Z <sub>24</sub>	-	
II	C <sub>1</sub>		6,7	Z <sub>4</sub> ,Z <sub>24</sub>	Z <sub>42</sub>	
VI	C <sub>1</sub>		6,7	z41	1,7	
I	C <sub>1</sub>	Hillsborough	6,7	z41	1,w	
I	C <sub>1</sub>	Tamilnadu	6,7	Z <sub>41</sub>	Z <sub>35</sub>	
II	C <sub>1</sub>	Sullivan	6,7	Z <sub>42</sub>	1,7	
II	C <sub>1</sub>		6,7	Z <sub>42</sub>	e,n,x:1,6	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	C <sub>1</sub>	Bulovka	6,7	Z <sub>44</sub>	–	
II	C <sub>1</sub>	Cape	6,7	Z <sub>6</sub>	1,7	
I	C <sub>2</sub>	Newport var. Puerto Rico	6,8	–	1,2	
II	C <sub>2</sub>		6,8	–	1,5,7	
I	C <sub>2</sub>	Valdosta	6,8	a	1,2	
I	C <sub>2</sub>	Doncaster	6,8	a	1,5	
I	C <sub>2</sub>	Curacao	6,8	a	1,6	
I	C <sub>2</sub>	Nordufer	6,8	a	1,7	
I	C <sub>2</sub>	Narashino	6,8	a	e,n,x	
II	C <sub>2</sub>		6,8	a	e,n,x	
I	C <sub>2</sub>	Leith	6,8	a	e,n,Z <sub>15</sub>	
II	C <sub>2</sub>		6,8	a	Z <sub>39</sub>	
II	C <sub>2</sub>	Tulear	6,8	a	Z <sub>52</sub>	
I	C <sub>2</sub>	Be	8,20	a	[Z <sub>6</sub> ]	
I	C <sub>2</sub>	Djelfa	8	b	1,2	
I	C <sub>2</sub>	Skansen	6,8	b	1,2	
I	C <sub>2</sub>	Korbol	8,20	b	1,5	
I	C <sub>2</sub>	Nagoya	6,8	b	1,5	
II	C <sub>2</sub>		6,8	b	1,5	
I	C <sub>2</sub>	Stourbridge	6,8	b	1,6	
I	C <sub>2</sub>	Eboko	6,8	b	1,7	
I	C <sub>2</sub>	Sanga	8	b	1,7	
I	C <sub>2</sub>	Gatuni	6,8	b	e,n,x	
I	C <sub>2</sub>	Konstanz	8	b	e,n,x	
I	C <sub>2</sub>	Presov	6,8	b	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Shipley	8,20	b	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Bukuru	6,8	b	1,w	
I	C <sub>2</sub>	Banalia	6,8	b	Z <sub>6</sub>	
I	C <sub>2</sub>	Tounouma	8,20	b	Z <sub>6</sub>	
I	C <sub>2</sub>	Wingrove	6,8	c	1,2	
I	C <sub>2</sub>	Utah	6,8	c	1,5	
I	C <sub>2</sub>	Bronx	6,8	c	1,6	
I	C <sub>2</sub>	Belfast	6,8	c	1,7	
I	C <sub>2</sub>	Belem	6,8	c	e,n,x	
I	C <sub>2</sub>	Santiago	8,20	c	e,n,x	
I	C <sub>2</sub>	Quiniela	6,8	c	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Alexanderpolder	8	c	1,w	
I	C <sub>2</sub>	Tado	8,20	c	Z <sub>6</sub>	
I	C <sub>2</sub>	Mexicana	6,8	d	1,2	Mexicana was combined with Muenchen. The name Mexicana has been dropped.
I	C <sub>2</sub>	Muenchen	6,8	d	1,2	
I	C <sub>2</sub>	Virginia	8	d	[1,2]	
I	C <sub>2</sub>	Manhattan	6,8	d	1,5	
I	C <sub>2</sub>	Yovokome	8,20	d	1,5	
I	C <sub>2</sub>	Dunkwa	6,8	d	1,7	
I	C <sub>2</sub>	Portanigra	8,20	d	1,7	
I	C <sub>2</sub>	Sterrenbos	6,8	d	e,n,x	
I	C <sub>2</sub>	Herston	6,8	d	e,n,Z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	C <sub>2</sub>	Labadi	8, <u>20</u>	d	Z <sub>6</sub>	
II	C <sub>2</sub>		6,8	d	Z <sub>6</sub> :Z <sub>42</sub>	
I	C <sub>2</sub>	Bardo	8	e,h	1,2	
I	C <sub>2</sub>	Newport	6,8, <u>20</u>	e,h	1,2	Newport may possess H phase RZ <sub>50</sub> or RZ <sub>58</sub> or RZ <sub>78</sub> or R1,12
I	C <sub>2</sub>	Ferruch	8	e,h	1,5	
I	C <sub>2</sub>	Kottbus	6,8	e,h	1,5	
I	C <sub>2</sub>	Cremieu	6,8	e,h	1,[6]	
I	C <sub>2</sub>	Atakpame	8, <u>20</u>	e,h	1,7	
I	C <sub>2</sub>	Tshiongwe	6,8	e,h	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Rehovot	8, <u>20</u>	e,h	Z <sub>6</sub>	
I	C <sub>2</sub>	Sadow	6,8	f,g	e,n,Z <sub>15</sub>	
II	C <sub>2</sub>		6,8	f,g,m,t	[e,n,x]	
I	C <sub>2</sub>	Emek	8, <u>20</u>	g,m,s	-	
I	C <sub>2</sub>	Chincol	6,8	g,m,[s]	[e,n,x]	
I	C <sub>2</sub>	Reubeuss	8, <u>20</u>	g,m,t	-	
II	C <sub>2</sub>		6,8	g,m,t	1,7	
I	C <sub>2</sub>	Alminko	8, <u>20</u>	g,s,t	-	
I	C <sub>2</sub>	Nanergou	6,8	g,s,t	-	
I	C <sub>2</sub>	Lindenburg	6,8	i	1,2	
I	C <sub>2</sub>	Bargny	8, <u>20</u>	i	1,5	
I	C <sub>2</sub>	Takoradi	6,8	i	1,5	
I	C <sub>2</sub>	Warnow	6,8	i	1,6	
I	C <sub>2</sub>	Malmoe	6,8	i	1,7	
I	C <sub>2</sub>	Bonariensis	6,8	i	e,n,x	
I	C <sub>2</sub>	Aba	6,8	i	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Cyprus	6,8	i	1,w	
I	C <sub>2</sub>	Magherafelt	8, <u>20</u>	i	1,w	
I	C <sub>2</sub>	Kentucky	8, <u>20</u>	i	Z <sub>6</sub>	
I	C <sub>2</sub>	Kallo	6,8	k	1,2	
I	C <sub>2</sub>	Blockley	6,8	k	1,5	Blockley may possess H phase RZ <sub>58</sub>
I	C <sub>2</sub>	Haardt	8	k	1,5	
I	C <sub>2</sub>	Schwerin	6,8	k	e,n,x	
I	C <sub>2</sub>	Charlottenburg	6,8	k	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Litchfield	6,8	l,v	1,2	
I	C <sub>2</sub>	Pakistan	8	l,v	1,2	
I	C <sub>2</sub>	Loanda	6,8	l,v	1,5	
I	C <sub>2</sub>	Manchester	6,8	l,v	1,7	
I	C <sub>2</sub>	Holcomb	6,8	l,v	e,n,x	
II	C <sub>2</sub>		6,8	l,v	e,n,x	
I	C <sub>2</sub>	Edmonton	6,8	l,v	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Amherstiana	8	l,(v)	1,6	
I	C <sub>2</sub>	Fayed	6,8	l,w	1,2	
II	C <sub>2</sub>		6,8	l,w	Z <sub>6</sub> :Z <sub>42</sub>	
I	C <sub>2</sub>	Hiduddify	6,8	l,Z <sub>13</sub> ,Z <sub>28</sub>	1,5	
II	C <sub>2</sub>		6,8	l,Z <sub>28</sub>	e,n,x	
I	C <sub>2</sub>	Breukelen	6,8	l,Z <sub>13</sub> ,[Z <sub>28</sub> ]	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Bassa	6,8	m,t	-	
I	C <sub>2</sub>	Yokoe	8, <u>20</u>	m,t	-	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	C <sub>2</sub>	Baragwanath	6,8	m,t	1,5	
II	C <sub>2</sub>	Germiston	6,8	m,t	e,n,x	
I	C <sub>2</sub>	Bsilla	6,8	r	1,2	
I	C <sub>2</sub>	Hindmarsh	8,20	r	1,5	
I	C <sub>2</sub>	Akanji	6,8	r	1,7	
I	C <sub>2</sub>	Noya	8	r	1,7	
I	C <sub>2</sub>	Goldcoast	6,8	r	l,w	
I	C <sub>2</sub>	Pikine	8,20	r	z <sub>6</sub>	Pikine was combined with Altona (8,20:r,[i]:z <sub>6</sub> ). The name Pikine has been dropped.
I	C <sub>2</sub>	Cocody	8,20	r,i	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Hidalgo	6,8	r,i	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Bovismorbificans	6,8	r,[i]	1,5	
I	C <sub>2</sub>	Brikama	8,20	r,[i]	l,w	
I	C <sub>2</sub>	Altona	8,20	r,[i]	z <sub>6</sub>	Pikine (8,20:r:z <sub>6</sub> ) was combined with Altona and called Altona.
I	C <sub>2</sub>	Giza	8,20	y	1,2	
I	C <sub>2</sub>	Brunei	8,20	y	1,5	
I	C <sub>2</sub>	Tananarive	6,8	y	1,5	
I	C <sub>2</sub>	Bulgaria	6,8	y	1,6	
II	C <sub>2</sub>		6,8	y	1,6:z <sub>42</sub>	
I	C <sub>2</sub>	Alagbon	8	y	1,7	
I	C <sub>2</sub>	Inchpark	6,8	y	1,7	
I	C <sub>2</sub>	Daarle	6,8	y	e,n,x	
I	C <sub>2</sub>	Sunnycove	8	y	e,n,x	
I	C <sub>2</sub>	Praha	6,8	y	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Benue	6,8	y	l,w	
I	C <sub>2</sub>	Sindelfingen	8,20	y	l,w	
I	C <sub>2</sub>	Kralingen	8,20	y	z <sub>6</sub>	
I	C <sub>2</sub>	Mowanjum	6,8	z	1,5	
II	C <sub>2</sub>		6,8	z	1,5	
I	C <sub>2</sub>	Kalumburu	6,8	z	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Phaliron	8	z	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Kuru	6,8	z	l,w	
I	C <sub>2</sub>	Daula	8,20	z	z <sub>6</sub>	
I	C <sub>2</sub>	Bazenheid	8,20	z <sub>10</sub>	1,2	
I	C <sub>2</sub>	Zerifin	6,8	z <sub>10</sub>	1,2	
I	C <sub>2</sub>	Mapo	6,8	z <sub>10</sub>	1,5	
I	C <sub>2</sub>	Paris	8,20	z <sub>10</sub>	1,5	
I	C <sub>2</sub>	Cleveland	6,8	z <sub>10</sub>	1,7	
I	C <sub>2</sub>	Hadar	6,8	z <sub>10</sub>	e,n,x	
I	C <sub>2</sub>	Istanbul	8	z <sub>10</sub>	e,n,x	
I	C <sub>2</sub>	Chomedey	8,20	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Glostrup	6,8	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Remiremont	8,20	z <sub>10</sub>	l,w	
I	C <sub>2</sub>	Molade	8,20	z <sub>10</sub>	z <sub>6</sub>	
I	C <sub>2</sub>	Wippra	6,8	z <sub>10</sub>	z <sub>6</sub>	
II	C <sub>2</sub>		6,8	z <sub>29</sub>	1,5	
II	C <sub>2</sub>		6,8	z <sub>29</sub>	e,n,x	
II	C <sub>2</sub>		8	z <sub>29</sub>	e,n,x:z <sub>42</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	C <sub>2</sub>	Tamale	8,20	Z <sub>29</sub>	[e,n,z <sub>15</sub> ]	
I	C <sub>2</sub>	Uno	6,8	Z <sub>29</sub>	[e,n,z <sub>15</sub> ]	
I	C <sub>2</sub>	Kolda	8,20	Z <sub>35</sub>	1,2	
I	C <sub>2</sub>	Yarm	6,8	Z <sub>35</sub>	1,2	
I	C <sub>2</sub>	Angers	8,20	Z <sub>35</sub>	Z <sub>6</sub>	
I	C <sub>2</sub>	Apeyeme	8,20	Z <sub>38</sub>	–	
I	C <sub>2</sub>	Bellevue	8	Z <sub>4</sub> ,Z <sub>23</sub>	1,7	
I	C <sub>2</sub>	Lezennes	6,8	Z <sub>4</sub> ,Z <sub>23</sub>	1,7	
I	C <sub>2</sub>	Breda	6,8	Z <sub>4</sub> ,Z <sub>23</sub>	e,n,x	
I	C <sub>2</sub>	Chailey	6,8	Z <sub>4</sub> ,Z <sub>23</sub>	[e,n,z <sub>15</sub> ]	
I	C <sub>2</sub>	Dabou	8,20	Z <sub>4</sub> ,Z <sub>23</sub>	l,w	
I	C <sub>2</sub>	Corvallis	8,20	Z <sub>4</sub> ,Z <sub>23</sub>	[z <sub>6</sub> ]	
I	C <sub>2</sub>	Albany	8,20	Z <sub>4</sub> ,Z <sub>24</sub>	–	Albany may possess H phase RZ <sub>45</sub> .
I	C <sub>2</sub>	Duesseldorf	6,8	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	C <sub>2</sub>	Tallahassee	6,8	Z <sub>4</sub> ,Z <sub>32</sub>	–	
I	C <sub>2</sub>	Diogoye	8,20	Z <sub>41</sub>	Z <sub>6</sub>	
I	C <sub>2</sub>	Aesch	6,8	Z <sub>60</sub>	1,2	
I	D <sub>1</sub>	Gallinarum	1,9,12	–	–	Gallinarum must be identified biochemically.
I	D <sub>1</sub>	Pullorum	1,9,12	–	–	IP combined Pullorum with Gallinarum (1,9,12:-:-). They must be identified biochemically.
I	D <sub>1</sub>	Miami	1,9,12	a	1,5	Miami must be differentiated from Sendai with biochemical tests. Miami is pos. for H <sub>2</sub> S, citrate, and tartrate; Sendai is neg.
I	D <sub>1</sub>	Sendai	1,9,12	a	1,5	Sendai must be differentiated from Miami with biochemical tests. Sendai is neg. for H <sub>2</sub> S, citrate, and tartrate; Miami is pos.
II	D <sub>1</sub>		9,12	a	1,5	
I	D <sub>1</sub>	Os	9,12	a	1,6	
I	D <sub>1</sub>	Saarbruecken	1,9,12	a	1,7	
I	D <sub>1</sub>	Lomalinda	1,9,12	a	e,n,x	
II	D <sub>1</sub>		1,9,12	a	e,n,x	
I	D <sub>1</sub>	Durban	9,12	a	e,n,z <sub>15</sub>	
II	D <sub>1</sub>		9,12	a	Z <sub>39</sub>	
II	D <sub>1</sub>		1,9,12	a	Z <sub>42</sub>	
I	D <sub>1</sub>	Onarimon	1,9,12	b	1,2	
I	D <sub>1</sub>	Frintrop	1,9,12	b	1,5	
II	D <sub>1</sub>	Mjimwema	1,9,12	b	e,n,x	
II	D <sub>1</sub>	Suederelbe	1,9,12	b	Z <sub>39</sub>	
II	D <sub>1</sub>	Blankenese	1,9,12	b	Z <sub>6</sub>	
I	D <sub>1</sub>	Goeteborg	9,12	c	1,5	
I	D <sub>1</sub>	Ipeko	9,12	c	1,6	
I	D <sub>1</sub>	Elokate	9,12	c	1,7	
I	D <sub>1</sub>	Alabama	9,12	c	e,n,z <sub>15</sub>	
I	D <sub>1</sub>	Ridge	9,12	c	Z <sub>6</sub>	
I	D <sub>1</sub>	Typhi	9,12,[Vi]	d	–	Typhi may possess H phase Rj or RZ <sub>66</sub> .
I	D <sub>1</sub>	Ndolo	1,9,12	d	1,5	
I	D <sub>1</sub>	Tarshyne	9,12	d	1,6	
I	D <sub>1</sub>	Eschberg	9,12	d	1,7	
II	D <sub>1</sub>	Rhodesiense	9,12	d	e,n,x	
I	D <sub>1</sub>	Bangui	9,12	d	e,n,z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	D <sub>1</sub>	Jaffna	<u>1</u> ,9,12	d	Z <sub>35</sub>	
II	D <sub>1</sub>		9,12	d	Z <sub>39</sub>	
I	D <sub>1</sub>	Zega	9,12	d	Z <sub>6</sub>	
I	D <sub>1</sub>	Bournemouth	9,12	e,h	1,2	
I	D <sub>1</sub>	Eastbourne	<u>1</u> ,9,12	e,h	1,5	
I	D <sub>1</sub>	Westafrica	9,12	e,h	1,7	
I	D <sub>1</sub>	Israel	9,12	e,h	e,n,Z <sub>15</sub>	
II	D <sub>1</sub>	Lindrick	9,12	e,n,x	1,[5],7	
II	D <sub>1</sub>		9,12	e,n,x	1,6	
I	D <sub>1</sub>	Enteritidis	<u>1</u> ,9,12	[f],g,m,[p],[t]	[1,7]	
I	D <sub>1</sub>	Berta	<u>1</u> ,9,12	[f],g,t	–	
I	D <sub>1</sub>	Blegdam	9,12	g,m,q	–	
II	D <sub>1</sub>	Muizenberg	9,12	g,m,s,t	1,5	IP combined Muizenberg with Hamburg ( <u>1</u> ,9,12:g,t:-) and Manica ( <u>1</u> ,9,12:g,m,s,t:Z <sub>42</sub> ) to form S. II <u>1</u> ,9,12:g,[m],[s],t:[1,5,7]:[Z <sub>42</sub> ].
II	D <sub>1</sub>	Kuilsrivier	<u>1</u> ,9,12	g,m,s,t	e,n,x	
II	D <sub>1</sub>	Manica	<u>1</u> ,9,12	g,m,s,t	Z <sub>42</sub>	IP combined Manica with Hamburg ( <u>1</u> ,9,12:g,t:-) and Muizenberg (9,12:g,m,s,t:1,5) to form S. II <u>1</u> ,9,12:g,[m],[s],t:[1,5,7]:[Z <sub>42</sub> ].
II	D <sub>1</sub>		<u>1</u> ,9,12	g,m,[s],t	[1,5,7]:[Z <sub>42</sub> ]	
I	D <sub>1</sub>	Dublin	<u>1</u> ,9,12,[Vi]	g,p	–	
I	D <sub>1</sub>	Naestved	<u>1</u> ,9,12	g,p,s	–	
I	D <sub>1</sub>	Rostock	<u>1</u> ,9,12	g,p,u	–	
I	D <sub>1</sub>	Moscow	9,12	g,q	–	
II	D <sub>1</sub>	Neasden	9,12	g,s,t	e,n,x	
II	D <sub>1</sub>	Hamburg	<u>1</u> ,9,12	g,t	–	IP combined Hamburg with Manica ( <u>1</u> ,9,12:g,m,s,t:Z <sub>42</sub> ) and Muizenberg (9,12:g,m,s,t:1,5) to form S. II <u>1</u> ,9,12:g,m,[s],t:[1,5,7]:[Z <sub>42</sub> ].
I	D <sub>1</sub>	Newmexico	9,12	g,Z <sub>51</sub>	1,5	
II	D <sub>1</sub>		<u>1</u> ,9,12	g,Z <sub>62</sub>	–	
I	D <sub>1</sub>	Antarctica	9,12	g,Z <sub>63</sub>	–	
I	D <sub>1</sub>	Seremban	9,12	i	1,5	
I	D <sub>1</sub>	Claibornei	<u>1</u> ,9,12	k	1,5	
I	D <sub>1</sub>	Goverdhan	9,12	k	1,6	
I	D <sub>1</sub>	Mendoza	9,12	l,v	1,2	
I	D <sub>1</sub>	Panama	<u>1</u> ,9,12	l,v	1,5	Panama may possess H phase R1,11
I	D <sub>1</sub>	Kapemba	9,12	l,v	1,7	
I	D <sub>1</sub>	Zaiman	9,12	l,v	e,n,x	
II	D <sub>1</sub>		9,12	l,v	e,n,x	
I	D <sub>1</sub>	Goettingen	9,12	l,v	e,n,Z <sub>15</sub>	
I	D <sub>1</sub>	Italiana	9,12	l,v	R1,11	IP combined Italiana that contains H phase R1,11 with Panama ( <u>1</u> ,9,12:l,v:1,5). The name Italiana has been dropped.
II	D <sub>1</sub>		9,12	l,v	Z <sub>39</sub>	
I	D <sub>1</sub>	Victoria	<u>1</u> ,9,12	l,w	1,5	
II	D <sub>1</sub>	Daressalaam	<u>1</u> ,9,12	l,w	e,n,x	
I	D <sub>1</sub>	Itami	9,12	l,Z <sub>13</sub>	1,5	
I	D <sub>1</sub>	Miyazaki	9,12	l,Z <sub>13</sub>	1,7	
I	D <sub>1</sub>	Napoli	<u>1</u> ,9,12	l,Z <sub>13</sub>	e,n,x	
I	D <sub>1</sub>	Javiana	<u>1</u> ,9,12	l,Z <sub>28</sub>	1,5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	D <sub>1</sub>		9,12	1,z <sub>28</sub>	1,5:[z <sub>42</sub> ]	
I	D <sub>1</sub>	Kotu	9,12	1,z <sub>28</sub>	1,6	
II	D <sub>1</sub>		9,12	1,z <sub>28</sub>	e,n,x	
II	D <sub>1</sub>		9,12	m,t	–	
I	D <sub>1</sub>	Pensacola	<u>1</u> ,9,12	m,t	[1,2]	
II	D <sub>1</sub>		<u>1</u> ,9,12	m,t	1,5	
II	D <sub>1</sub>		<u>1</u> ,9,12	m,t	z <sub>39</sub>	
I	D <sub>1</sub>	Jamaica	9,12	r	1,5	
I	D <sub>1</sub>	Camberwell	9,12	r	1,7	
I	D <sub>1</sub>	Campinense	9,12	r	e,n,z <sub>15</sub>	
I	D <sub>1</sub>	Lome	9,12	r	z <sub>6</sub>	
I	D <sub>1</sub>	Powell	9,12	y	1,7	
I	D <sub>1</sub>	Lawndale	<u>1</u> ,9,12	z	1,5	
I	D <sub>1</sub>	Kimpese	9,12	z	1,6	
II	D <sub>1</sub>	Stellenbosch	<u>1</u> ,9,12	z	1,7	
II	D <sub>1</sub>	Hueningen	9,12	z	z <sub>39</sub>	
II	D <sub>1</sub>	Angola	<u>1</u> ,9,12	z	z <sub>6</sub>	
I	D <sub>1</sub>	Portland	9,12	z <sub>10</sub>	1,5	
I	D <sub>1</sub>	Ruanda	9,12	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	D <sub>1</sub>	Treguier	<u>1</u> ,9,12	z <sub>10</sub>	z <sub>6</sub>	
II	D <sub>1</sub>	Canastel	9,12	z <sub>29</sub>	[1,5]	
II	D <sub>1</sub>		<u>1</u> ,9,12	z <sub>29</sub>	e,n,x	
I	D <sub>1</sub>	Penarth	9,12	z <sub>35</sub>	z <sub>6</sub>	
I	D <sub>1</sub>	Elomrane	<u>1</u> ,9,12	z <sub>38</sub>	–	
II	D <sub>1</sub>	Wynberg	<u>1</u> ,9,12	z <sub>39</sub>	1,7	
I	D <sub>1</sub>	Wangata	<u>1</u> ,9,12	z <sub>4</sub> ,z <sub>23</sub>	[1,7]	
I	D <sub>1</sub>	Natal	9,12	z <sub>4</sub> ,z <sub>24</sub>	–	
I	D <sub>1</sub>	Ottawa	<u>1</u> ,9,12	z <sub>41</sub>	1,5	
II	D <sub>1</sub>		<u>1</u> ,9,12	z <sub>42</sub>	1,[5],7	
I	D <sub>1</sub>	Franken	<u>1</u> ,9,12	z <sub>60</sub>	z <sub>67</sub>	
I	D <sub>2</sub>	Baildon	9,46	a	e,n,x	
I	D <sub>2</sub>	Doba	9,46	a	e,n,z <sub>15</sub>	
I	D <sub>2</sub>	Cheltenham	9,46	b	1,5	
I	D <sub>2</sub>	Zadar	9,46	b	1,6	
I	D <sub>2</sub>	Worb	9,46	b	e,n,x	
II	D <sub>2</sub>	Lundby	9,46	b	e,n,x	
I	D <sub>2</sub>	Bamboye	9,46	b	1,w	
I	D <sub>2</sub>	Kolar	9,46	b	z <sub>35</sub>	
I	D <sub>2</sub>	Linguere	9,46	b	z <sub>6</sub>	
I	D <sub>2</sub>	Itutaba	9,46	c	z <sub>6</sub>	
I	D <sub>2</sub>	Ontario	9,46	d	1,5	
I	D <sub>2</sub>	Quentin	9,46	d	1,6	
I	D <sub>2</sub>	Strasbourg	9,46	d	1,7	
I	D <sub>2</sub>	Olten	9,46	d	e,n,z <sub>15</sub>	
I	D <sub>2</sub>	Plymouth	9,46	d	z <sub>6</sub>	
I	D <sub>2</sub>	Bergedorf	9,46	e,h	1,2	
I	D <sub>2</sub>	Waedenswil	9,46	e,h	1,5	
I	D <sub>2</sub>	Guerin	9,46	e,h	z <sub>6</sub>	
II	D <sub>2</sub>		9,46	e,n,x	1,5,7	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	D <sub>2</sub>	Wernigerode	9,46	f,g	–	
I	D <sub>2</sub>	Hillingdon	9,46	g,m	–	
I	D <sub>2</sub>	Macclesfield	9,46	g,m,s,t	1,(2),7	
II	D <sub>2</sub>	Duivenhoks	9,46	g,[m],[s],t	[e,n,x]	
I	D <sub>2</sub>	Gateshead	9,46	g,s,t	–	
II	D <sub>2</sub>		9,46	g,Z <sub>62</sub>	–	
I	D <sub>2</sub>	Mathura	9,46	i	e,n,Z <sub>15</sub>	
I	D <sub>2</sub>	Potto	9,46	i	Z <sub>6</sub>	
I	D <sub>2</sub>	Marylebone	9,46	k	1,2	
I	D <sub>2</sub>	Cochin	9,46	k	1,5	
I	D <sub>2</sub>	Ceyco	9,46	k	Z <sub>35</sub>	
I	D <sub>2</sub>	India	9,46	l,v	1,5	
I	D <sub>2</sub>	Geraldton	9,46	l,v	1,6	
I	D <sub>2</sub>	Toronto	9,46	l,v	e,n,x;[Z <sub>44</sub> ]	
I	D <sub>2</sub>	Ackwepe	9,46	l,w	–	
I	D <sub>2</sub>	Sangalkam	9,46	m,t	–	
II	D <sub>2</sub>		9,46	m,t	e,n,x	
I	D <sub>2</sub>	Deckstein	9,46	r	1,7	
I	D <sub>2</sub>	Shoreditch	9,46	r	e,n,Z <sub>15</sub>	
I	D <sub>2</sub>	Sokode	9,46	r	Z <sub>6</sub>	
I	D <sub>2</sub>	Benin	9,46	y	1,7	
I	D <sub>2</sub>	Irchel	9,46	y	e,n,x	
I	D <sub>2</sub>	Nantes	9,46	y	l,w	
I	D <sub>2</sub>	Mayday	9,46	y	Z <sub>6</sub>	
II	D <sub>2</sub>		9,46	z	1,5	
II	D <sub>2</sub>	Haarlem	9,46	z	e,n,x	
I	D <sub>2</sub>	Bambylor	9,46	z	e,n,Z <sub>15</sub>	
I	D <sub>2</sub>	Lishabi	9,46	Z <sub>10</sub>	1,7	
I	D <sub>2</sub>	Inglis	9,46	Z <sub>10</sub>	e,n,x	
I	D <sub>2</sub>	Mahina	9,46	Z <sub>10</sub>	e,n,Z <sub>15</sub>	
II	D <sub>2</sub>		9,46	Z <sub>10</sub>	Z <sub>39</sub>	
I	D <sub>2</sub>	Louisiana	9,46	Z <sub>10</sub>	Z <sub>6</sub>	
II	D <sub>2</sub>		9,46	Z <sub>10</sub>	Z <sub>6</sub>	
I	D <sub>2</sub>	Ouakam	9,46	Z <sub>29</sub>	–	
I	D <sub>2</sub>	Hillegersberg	9,46	Z <sub>35</sub>	1,5	
I	D <sub>2</sub>	Basingstoke	9,46	Z <sub>35</sub>	e,n,Z <sub>15</sub>	
I	D <sub>2</sub>	Trimdon	9,46	Z <sub>35</sub>	Z <sub>6</sub>	
I	D <sub>2</sub>	Fresno	9,46	Z <sub>38</sub>	–	
II	D <sub>2</sub>		9,46	Z <sub>39</sub>	1,7	
I	D <sub>2</sub>	Ekotedo	9,46	Z <sub>4</sub> ,Z <sub>23</sub>	–	
I	D <sub>2</sub>	Ngaparou	9,46	Z <sub>4</sub> ,Z <sub>24</sub>	–	
II	D <sub>2</sub>	Maarssen	9,46	Z <sub>4</sub> ,Z <sub>24</sub>	Z <sub>39</sub> ;Z <sub>42</sub>	
I	D <sub>2</sub>	Wuppertal	9,46	Z <sub>41</sub>	–	
II	D <sub>3</sub>	Zuerich	1,9,12,46,27	c	Z <sub>39</sub>	
II	D <sub>3</sub>		9,12,46,27	g,t	e,n,x	
II	D <sub>3</sub>		1,9,12,46,27	l,Z <sub>13</sub> ,Z <sub>28</sub>	Z <sub>39</sub>	
II	D <sub>3</sub>		1,9,12,46,27	y	Z <sub>39</sub>	
II	D <sub>3</sub>		1,9,12,46,27	Z <sub>10</sub>	1,5	
II	D <sub>3</sub>		1,9,12,46,27	Z <sub>10</sub>	e,n,x	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	D <sub>3</sub>		1,9,12,46,27	Z <sub>10</sub>	Z <sub>39</sub>	
II	D <sub>3</sub>		1,9,12,46,27	Z <sub>4</sub> ,Z <sub>24</sub>	[1,5]	
I	E <sub>1</sub>	Aminatu	3,10	a	1,2	
I	E <sub>1</sub>	Goelzau	3,10	a	1,5	IP combined Clichy (3,15:a:1,5) with Goelzau to form Goelzau 3,10,[15]:a:1,5.
I	E <sub>1</sub>	Oxford	3,10	a	1,7	IP combined Khartoum (3,15,34:a:1,7) with Oxford to form Oxford 3,10,[15],[15,34]:a:1,7.
I	E <sub>1</sub>	Masembe	3,10	a	e,n,x	Masembe may possess H phase Rz <sub>5</sub> .
II	E <sub>1</sub>	Matroosfontein	3,10	a	e,n,x	
I	E <sub>1</sub>	Galil	3,10	a	e,n,z <sub>15</sub>	
II	E <sub>1</sub>		3,10	a	Z <sub>39</sub>	
I	E <sub>1</sub>	Kalina	3,10	b	1,2	
I	E <sub>1</sub>	Butantan	3,10	b	1,5	IP combined Rosenthal (3,15:b:1,5) and unnamed 3,15,34:b:1,5 with Butantan to form Butantan 3,10,[15],[15,34]:b:1,5.
I	E <sub>1</sub>	Allerton	3,10	b	1,6	
I	E <sub>1</sub>	Huvudsta	3,10	b	1,7	
I	E <sub>1</sub>	Benfica	3,10	b	e,n,x	
II	E <sub>1</sub>		3,10	b	e,n,x	
I	E <sub>1</sub>	Yaba	3,10,[15]	b	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Epicrates	3,10	b	1,w	
I	E <sub>1</sub>	Westminster	3,10,[15]	b	Z <sub>35</sub>	CDC has no 3,10:b:z <sub>35</sub> .
II	E <sub>1</sub>		3,10	b	Z <sub>39</sub>	
I	E <sub>1</sub>	Wilmington	3,10	b	Z <sub>6</sub>	
I	E <sub>1</sub>	Asylanta	3,10	c	1,2	
I	E <sub>1</sub>	Gbadago	3,10,[15]	c	1,5	
I	E <sub>1</sub>	Ikayi	3,10,[15]	c	1,6	Ikayi Var. O 15+ was described after E <sub>1</sub> and E <sub>2</sub> were combined.
I	E <sub>1</sub>	Pramiso	3,10	c	1,7	
I	E <sub>1</sub>	Agege	3,10	c	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Anderlecht	3,10	c	1,w	
I	E <sub>1</sub>	Okefoko	3,10	c	Z <sub>6</sub>	
I	E <sub>1</sub>	Stormont	3,10	d	1,2	
I	E <sub>1</sub>	Shangani	3,10	d	1,5	IP combined Pankow (3,15:d:1,5) with Shangani to form Shangani 3,10,[15]:d:1,5.
I	E <sub>1</sub>	Lekke	3,10	d	1,6	
I	E <sub>1</sub>	Onireke	3,10	d	1,7	
I	E <sub>1</sub>	Souza	3,10	d	e,n,x	IP combined Eschersheim (3,15:d:e,n,x) with Souza to form Souza 3,10,[15]:d:e,n,x.
II	E <sub>1</sub>		3,10	d	e,n,x	
I	E <sub>1</sub>	Madjorio	3,10	d	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Birmingham	3,10	d	1,w	
I	E <sub>1</sub>	Maron	3,10	d	Z <sub>35</sub>	
I	E <sub>1</sub>	Weybridge	3,10	d	Z <sub>6</sub>	
I	E <sub>1</sub>	Vejle	3,10	e,h	1,2	IP combined Goerlitz (3,15:e,h:1,2) with Vejle to form Vejle 3,10,[15]:e,h:1,2.
I	E <sub>1</sub>	Muenster	3,10	e,h	1,5	IP combined Muenster with Newhaw (3,15:e,h:1,5) and Arkansas (3,15,34:e,h:1,5) to form Muenster 3,10,[15],[15,34]:e,h:1,5.
I	E <sub>1</sub>	Anatum	3,10	e,h	1,6	IP combined Newington (3,15:e,h:1,6) and Minneapolis (3,15,34:e,h:1,6) with Anatum to form Anatum 3,10,[15],[15,34]:e,h:1,6.

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>1</sub>	Nyborg	3,10	e,h	1,7	IP combined Selandia (3,15:e,h:1,5) with Nyborg to form Nyborg 3,10,[15]:e,h:1,7.
I	E <sub>1</sub>	Newlands	3,10,[15,34]	e,h	e,n,x	
I	E <sub>1</sub>	Lamberhurst	3,10	e,h	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Meleagridis	3,10	e,h	l,w	IP combined Meleagridis with Cambridge (3,15:e,h:l,w) and Wildwood (3,15,34:e,h:l,w) to form Meleagridis 3,10,[15],[15,34]:e,h:l,w.
I	E <sub>1</sub>	Sekondi	3,10	e,h	z <sub>6</sub>	
II	E <sub>1</sub>	Chudleigh	3,10	e,n,x	1,7	
I	E <sub>1</sub>	Alfort	3,10	f,g	e,n,x	
I	E <sub>1</sub>	Regent	3,10	f,g,[s]	[1,6]	
I	E <sub>1</sub>	Suberu	3,10	g,m	–	
I	E <sub>1</sub>	Amsterdam	3,10	g,m,s	–	IP combined Drypool (3,15:g,m,s) and Drypool var. O 34+ with Amsterdam to form Amsterdam 3,10,[15],[15,34]:g,m,s:-.
II	E <sub>1</sub>	Parow	3,10,[15]	g,m,s,t	–	
I	E <sub>1</sub>	Westhampton	3,10	g,s,t	–	IP combined Halmstad (3,15:g,s,t:-) and Canoga (3,15,34:g,s,t:-) with Westhampton to form Westhampton 3,10,[15],[15,34]:g,s,t:-. Westhampton may possess H phase Rz <sub>37</sub> or Rz <sub>43</sub> or Rz <sub>45</sub> .
II	E <sub>1</sub>	Islington	3,10	g,t	–	
I	E <sub>1</sub>	Bloomsbury	3,10	g,t	1,5	
I	E <sub>1</sub>	Cuckmere	3,10	i	1,2	
I	E <sub>1</sub>	Amounderness	3,10	i	1,5	
I	E <sub>1</sub>	Truro	3,10	i	1,7	
I	E <sub>1</sub>	Bessi	3,10	i	e,n,x	
I	E <sub>1</sub>	Falkensee	3,10	i	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Hoboken	3,10	i	l,w	
I	E <sub>1</sub>	Yeerongpilly	3,10	i	z <sub>6</sub>	
I	E <sub>1</sub>	Wimborne	3,10	k	1,2	
I	E <sub>1</sub>	Zanzibar	3,10,[15]	k	1,5	
I	E <sub>1</sub>	Serrekunda	3,10	k	1,7	
I	E <sub>1</sub>	Yundum	3,10	k	e,n,x	
I	E <sub>1</sub>	Marienthal	3,10	k	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Newrochelle	3,10	k	l,w	
I	E <sub>1</sub>	Nchanga	3,10	l,v	1,2	IP combined Nchanga with Nancy (3,15:l,v:1,2) to form Nchanga 3,10,[15]:l,v:1,2.
I	E <sub>1</sub>	Sinstorf	3,10	l,v	1,5	
I	E <sub>1</sub>	London	3,10	l,v	1,6	IP combined London with Portsmouth (3,15:l,v:1,6) to form London 3,10,[15]:l,v:1,6.
I	E <sub>1</sub>	Give	3,10	l,v	1,7	IP combined Newbrunswick (3,15:l,v:1,7) and Menhaden (3,15,34:l,v:1,7) with Give to form Give 3,10,[15],[15,34]:[d],l,v:[d],1,7. Give may possess H phase d; Rl,z <sub>40</sub> ; or Rz <sub>77</sub> .
II	E <sub>1</sub>		3,10	l,v	e,n,x	
I	E <sub>1</sub>	Ruzizi	3,10	l,v	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Sinchew	3,10	l,v	z <sub>35</sub>	
II	E <sub>1</sub>	Fuhlsbuettel	3,10	l,v	z <sub>6</sub>	
I	E <sub>1</sub>	Assinie	3,10	l,w	z <sub>6</sub>	Assinie may possess H phase Rz <sub>45</sub> .
I	E <sub>1</sub>	Freiburg	3,10	l,z <sub>13</sub>	1,2	
I	E <sub>1</sub>	Uganda	3,10,15	l,z <sub>13</sub>	1,5	
I	E <sub>1</sub>	Fallowfield	3,10	l,z <sub>13</sub> ,z <sub>28</sub>	e,n,z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>1</sub>	Hoghton	3,10	1,z <sub>13</sub> ,z <sub>28</sub>	z <sub>6</sub>	
II	E <sub>1</sub>		3,10	1,z <sub>28</sub>	1,5	
I	E <sub>1</sub>	Joal	3,10	1,z <sub>28</sub>	1,7	
I	E <sub>1</sub>	Lamin	3,10	1,z <sub>28</sub>	e,n,x	
II	E <sub>1</sub>	Westpark	3,10	1,z <sub>28</sub>	e,n,x	
II	E <sub>1</sub>		3,10	1,z <sub>28</sub>	z <sub>39</sub>	
II	E <sub>1</sub>		3,10	m,t	1,5	
I	E <sub>1</sub>	Southbank	3,10,15,34	m,t	[1,6]	
II	E <sub>1</sub>	Stikland	3,10	m,t	e,n,x	
I	E <sub>1</sub>	Ughelli	3,10	r	1,5	
I	E <sub>1</sub>	Elisabethville	3,10	r	1,7	
I	E <sub>1</sub>	Simi	3,10	r	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Weltevreden	3,10	r	z <sub>6</sub>	IP combined Lanka 3,15:r:z <sub>6</sub> with Weltevreden to form Weltevreden 3,10,[15]:r:z <sub>6</sub> .
I	E <sub>1</sub>	Seegefeld	3,10	r,[i]	1,2	
I	E <sub>1</sub>	Dumfries	3,10	r,[i]	1,6	
I	E <sub>1</sub>	Rutgers	3,10	R1,z <sub>40</sub>	1,7	Rutgers has been dropped from the scheme and the H phase R1,z <sub>40</sub> is now considered an R phase of Give.
I	E <sub>1</sub>	Amager	3,10	y	1,2	IP combined Tuebingen (3,15:y:1,2) with Amager to form Amager 3,10,[15]:y:1,2. Amager may possess H phase Rz <sub>45</sub> .
I	E <sub>1</sub>	Orion	3,10	y	1,5	IP combined Binza (3,15:y:1,5) and Thomasville (3,15,34:y:1,5) with Orion to form Orion 3,10,[15],[15,34]:y:1,5.
I	E <sub>1</sub>	Mokola	3,10	y	1,7	
I	E <sub>1</sub>	Ohlstedt	3,10	y	e,n,x	
I	E <sub>1</sub>	Bolton	3,10	y	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Langensalza	3,10	y	1,w	
I	E <sub>1</sub>	Stockholm	3,10	y	z <sub>6</sub>	IP combined Tournai (3,15:y:z <sub>6</sub> ) with Stockholm to form Stockholm 3,10,[15]:y:z <sub>6</sub> .
I	E <sub>1</sub>	Fufu	3,10	z	1,5	
II	E <sub>1</sub>	Alexander	3,10	z	1,5	
I	E <sub>1</sub>	Harleystreet	3,10	z	1,6	
I	E <sub>1</sub>	Huddinge	3,10	z	1,7	
II	E <sub>1</sub>	Finchley	3,10	z	e,n,x	
I	E <sub>1</sub>	Clerkenwell	3,10	z	1,w	
II	E <sub>1</sub>	Tafelbaai	3,10	z	z <sub>39</sub>	
I	E <sub>1</sub>	Landwasser	3,10	z	z <sub>6</sub>	
I	E <sub>1</sub>	Okerara	3,10	z <sub>10</sub>	1,2	
I	E <sub>1</sub>	Lexington	3,10	z <sub>10</sub>	1,5	IP combined Lexington with Manila (3,15:z <sub>10</sub> :1,5) and Illinois (3,15,34:z <sub>10</sub> :1,5) to form Lexington 3,10,[15],[15,34]:z <sub>10</sub> :1,5. Lexington may possess H phase Rz <sub>49</sub> .
I	E <sub>1</sub>	Harrisonburg	3,10,[15],[15,34]	z <sub>10</sub>	1,6	
I	E <sub>1</sub>	Coquilhatville	3,10	z <sub>10</sub>	1,7	
I	E <sub>1</sub>	Kristianstad	3,10	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Biafra	3,10	z <sub>10</sub>	z <sub>6</sub>	
I	E <sub>1</sub>	Jedburgh	3,10,[15]	z <sub>29</sub>	–	
II	E <sub>1</sub>		3,10	z <sub>29</sub>	–	
I	E <sub>1</sub>	Everleigh	3,10	z <sub>29</sub>	e,n,x	
II	E <sub>1</sub>		3,10	z <sub>29</sub>	e,n,x	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>1</sub>	Zongo	3,10	Z <sub>35</sub>	1,7	
I	E <sub>1</sub>	Shannon	3,10	Z <sub>35</sub>	1,w	
I	E <sub>1</sub>	Cairina	3,10	Z <sub>35</sub>	Z <sub>6</sub>	
I	E <sub>1</sub>	Macallen	3,10	Z <sub>36</sub>	–	
II	E <sub>1</sub>	Mpila	3,10	Z <sub>38</sub>	Z <sub>42</sub>	
I	E <sub>1</sub>	Bolombo	3,10	Z <sub>38</sub>	[Z <sub>6</sub> ]	
II	E <sub>1</sub>	Winchester	3,10	Z <sub>39</sub>	1,[5],7	
I	E <sub>1</sub>	Adabraka	3,10	Z <sub>4</sub> ,Z <sub>23</sub>	[1,7]	
I	E <sub>1</sub>	Wagadugu	3,10	Z <sub>4</sub> ,Z <sub>23</sub>	Z <sub>6</sub>	
I	E <sub>1</sub>	Florian	3,10,[15]	Z <sub>4</sub> ,Z <sub>24</sub>	–	
II	E <sub>1</sub>		3,10	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	E <sub>1</sub>	Pietersburg	3,10,[15,34]	Z <sub>69</sub>	1,7	
I	E <sub>2</sub>	Clichy	3,15	a	1,5	IP combined Clichy with Goelzau (3,10:a:1,5) to form Goelzau 3,10,[15]:a:1,5. Clichy is now called Goelzau var. O 15+ by IP.
I	E <sub>2</sub>	Rosenthal	3,15	b	1,5	IP combined Rosenthal and unnamed 3,15,34:b:1,5 with Butantan (3,10:b:1,5) to form Butantan 3,10,[15],[15,34]:b:1,5. Rosenthal is now called Butantan var. O 15+ by IP.
I	E <sub>2</sub>	Pankow	3,15	d	1,5	IP combined Pankow with Shangani (3,10:d:1,5) to form Shangani 3,10,15:d:1,5. Pankow is now called Shangani var. O 15+ by IP.
I	E <sub>2</sub>	Eschersheim	3,15	d	e,n,x	IP combined Eschersheim with Souza (3,10:d:e,n,x) to form Souza 3,10,[15]:d:e,n,x. Eschersheim is now called Souza var. O 15+ by IP.
I	E <sub>2</sub>	Goerlitz	3,15	e,h	1,2	IP combined Goerlitz with Vejle (3,10:e,h:1,2) to form Vejle 3,10,15:e,h:1,2. Goerlitz is now called Vejle var. O 15+ by IP.
I	E <sub>2</sub>	Newhaw	3,15	e,h	1,5	IP combined Newhaw and Arkansas (3,15,34:e,h:1,5) with Muenster (3,10:e,h:1,5) to form Muenster 3,10,[15],[15,34]:e,h:1,5. Newhaw is now called Muenster var. O 15+ by IP.
I	E <sub>2</sub>	Newington	3,15	e,h	1,6	IP combined Newington and Minneapolis (3,15,34:e,h:1,6) with Anatum (3,10:e,h:1,6) to form Anatum 3,10,[15],[15,34]:e,h:1,6. Newington is now called Anatum var. O 15+ by IP.
I	E <sub>2</sub>	Selandia	3,15	e,h	1,7	IP combined Selandia with Nyborg (3,10:e,h:1,7) to form Nyborg 3,10,[15]:e,h:1,7. Selandia is now called Nyborg var. O 15+ by IP.
I	E <sub>2</sub>	Cambridge	3,15	e,h	1,w	IP combined Cambridge and Wildwood (3,15,34:e,h:1,w) with Meleagridis (3,10:e,h:1,w) to form Meleagridis 3,10,[15],[15,34]:e,h:1,w. Cambridge is now called Meleagridis var. O 15+ by IP.
I	E <sub>2</sub>	Drypool	3,[15],[15,34]	g,m,s	–	IP combined Drypool 3,15:g,m,s:- and Drypool var. O34+ with Amsterdam (3,10:g,m,s:-) to form Amsterdam 3,10,[15],[15,34]:g,m,s:-. Drypool is now called Amsterdam var. O 15+ or O 15+, 34+ by IP.
I	E <sub>2</sub>	Halmstad	3,15	g,s,t	–	IP combined Halmstad and Canoga (3,15,34:g,s,t:-) with Westhampton (3,10:g,s,t:-) to form Westhampton 3,10,[15],[15,34]:g,s,t:-. Halmstad is now called Westhampton var. O 15+ by IP.
I	E <sub>2</sub>	Nancy	3,15	l,v	1,2	IP combined Nancy with Nchanga (3,10:l,v:1,2) to form Nchanga 3,10,[15]:l,v:1,2. Nancy is now called Nchanga var. O 15+ by IP.

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>2</sub>	Portsmouth	3,15	l,v	1,6	IP combined Portsmouth with London (3,10:l,v:1,6) to form London 3,10,[15]:l,v:1,6. Portsmouth is now called London var. O 15+ by IP.
I	E <sub>2</sub>	Newbrunswick	3,15	l,v	1,7	IP combined Newbrunswick and Menhaden (3,15,34:l,v:1,7) with Give (3,10:l,v:1,7) to form Give 3,10,[15],[15,34]:[d],l,v:[d],1,7. Newbrunswick is now called Give var. O 15+ by IP.
I	E <sub>2</sub>	Kinshasa	3,15	l,z <sub>13</sub>	1,5	IP combined Kinshasa with Uganda (3,10:l,z <sub>13</sub> :1,5) to form Uganda 3,10,[15]:l,z <sub>13</sub> :1,5. Kinshasa is now called Uganda var. O 15+ by IP.
I	E <sub>2</sub>	Lanka	3,15	r	z <sub>6</sub>	IP combined Lanka with Weltevreden (3,10:r:z <sub>6</sub> ) to form Weltevreden 3,10,[15]:r:z <sub>6</sub> . Lanka is now called Weltevreden var. O 15+ by IP.
I	E <sub>2</sub>	Hamilton	3,15	Rz <sub>27</sub>	–	IP combined Hamilton with Goerlitz (3,15:e,h:1,2) and Vejle (3,10:e,h:1,2) to form Vejle 3,10,15:e,h:1,2:Rz <sub>27</sub> . Hamilton is now called Vejle var. Rz <sub>27</sub> +. The name Hamilton has been dropped.
I	E <sub>2</sub>	Tuebingen	3,15	y	1,2	IP combined Tuebingen with Amager (e,10:y:1,2) to form Amager 3,10,[15]:y:1,2. Tuebingen is now called Amager var. O 15+ by IP.
I	E <sub>2</sub>	Binza	3,15	y	1,5	IP combined Binza and Thomasville (3,15,34:y:1,5) with Orion (3,10:y:1,5) to form Orion 3,10,[15],[15,34]:y:1,5. Binza is now called Orion var. O 15+ by IP.
I	E <sub>2</sub>	Tournai	3,15	y	z <sub>6</sub>	IP combined Tournai with Stockholm (3,10:y:z <sub>6</sub> ) to form Stockholm 3,10,[15]:y:z <sub>6</sub> . Tournai is now called Stockholm var. O 15+ by IP.
I	E <sub>2</sub>	Manila	3,15	z <sub>10</sub>	1,5	IP combined Manila and Illinois (3,15,34:z <sub>10</sub> :1,5) with Lexington (3,10:z <sub>10</sub> :1,5) to form Lexington 3,10,[15],[15,34]:z <sub>10</sub> :1,5. Manila is now called Lexington var. O 15+ by IP.
I	E <sub>3</sub>	Khartoum	3,15,34	a	1,7	IP combined Khartoum with Oxford (3,10:a:1,7) to form Oxford 3,10,[15],[15,34]:a:1,7. Khartoum is now called Oxford var. O 15+ by IP. CDC has no 3,15:a:1,7. Khartoum was found by IP with colonies containing O 3,15.
I	E <sub>3</sub>	Arkansas	3,15,34	e,h	1,5	IP combined Arkansas and Newhaw (3,15:e,h:1,5) with Muenster (3,10:e,h:1,5) to form Muenster 3,10,[15],[15,34]:e,h:1,5. Arkansas is now called Muenster var. O 15+, 34+ by IP.
I	E <sub>3</sub>	Minneapolis	3,15,34	e,h	1,6	IP combined Minneapolis and Newington (3,15:e,h:1,6) with Anatum (3,10:e,h:1,6) to form Anatum 3,10,[15],[15,34]:e,h:1,6. Minneapolis is now called Anatum var. O 15+ by IP.
I	E <sub>3</sub>	Wildwood	3,15,34	e,h	1,w	IP combined Wildwood and Cambridge (3,15:e,h:l,w) with Meleagridis (3,10:e,h:l,w) to form Meleagridis 3,10,[15],[15,34]:e,h:l,w. Wildwood is now called Meleagridis var. O 15+, 34+ by IP.
I	E <sub>3</sub>	Drypool	3,15,34	g,m,s	–	IP combined Drypool (3,15:g,m,s:-) and Drypool var. O 34+ with Amsterdam (3,10:g,m,s:-) to form Amsterdam 3,10,[15],[15,34]:g,m,s:-. Drypool is now called Amsterdam var. O 15+ or var. O 15+,34+ by IP.
I	E <sub>3</sub>	Canoga	3,15,34	g,s,t	-	IP combined Canoga and Halmstad (3,15:g,s,t:-) with Westhampton (3,10:g,s,t:-) to form Westhampton 3,10,[15],[15,34]:g,s,t:-. Canoga is now called Westhampton var. O 15+, 34+ by IP.

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>3</sub>	Menhaden	3,15,34	l,v	1,7	IP combined Menhaden with Give (3,10:l,v:1,7) and Newbrunswick (3,15:l,v:1,7) to form Give 3,10,[15],[15,34]:[d],l,v:1,7. Menhaden is now called Give var. O 15+, 34+ by IP.
I	E <sub>3</sub>	Thomasville	3,15,34	y	1,5	IP combined Thomasville and Binza (3,15:y:1,5) with Orion (3,10:y:1,5) to form Orion 3,10,[15],[15,34]:y:1,5. Thomasville is now called Orion var. O 15+ by IP.
I	E <sub>3</sub>	Illinois	3,15,34	z <sub>10</sub>	1,5	IP combined Illinois and Manila (3,15:z <sub>10</sub> :1,5) with Lexington (3,10:z <sub>10</sub> :1,5) to form Lexington 3,10,[15],[15,34]:z <sub>10</sub> :1,5. Illinois is now called Lexington var. O 15+, 34+ by IP.
I	E <sub>4</sub>	Niumi	1,3,19	a	1,5	
I	E <sub>4</sub>	Juba	1,3,19	a	1,7	
I	E <sub>4</sub>	Gwoza	1,3,19	a	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Alkmaar	1,3,19	a	l,w	
I	E <sub>4</sub>	Gnesta	1,3,19	b	1,5	
I	E <sub>4</sub>	Visby	1,3,19	b	1,6	
I	E <sub>4</sub>	Tambacounda	1,3,19	b	e,n,x	
I	E <sub>4</sub>	Kande	1,3,19	b	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Broughton	1,3,19	b	l,w	
I	E <sub>4</sub>	Chittagong	1,3,10,19	b	z <sub>35</sub>	
I	E <sub>4</sub>	Accra	1,3,19	b	z <sub>6</sub>	
I	E <sub>4</sub>	Eastglam	1,3,19	c	1,5	
I	E <sub>4</sub>	Bida	1,3,19	c	1,6	
I	E <sub>4</sub>	Madiago	1,3,19	c	1,7	
I	E <sub>4</sub>	Ahmadi	1,3,19	d	1,5	
I	E <sub>4</sub>	Liverpool	1,3,19	d	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Tilburg	1,3,19	d	l,w	Tilburg may possess H phase Rz <sub>49</sub> .
I	E <sub>4</sub>	Niloese	1,3,19	d	z <sub>6</sub>	
I	E <sub>4</sub>	Vilvoorde	1,3,19	e,h	1,5	
I	E <sub>4</sub>	Hayindogo	1,3,19	e,h	1,6	
I	E <sub>4</sub>	Sanktmarx	1,3,19	e,h	1,7	
I	E <sub>4</sub>	Sao	1,3,19	e,h	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Calabar	1,3,19	e,h	l,w	
I	E <sub>4</sub>	Rideau	1,3,19	f,g	-	
I	E <sub>4</sub>	Bilu	(1),3,10,(19)	f,g,t	1,(2),7	
I	E <sub>4</sub>	Petahtikva	1,3,19	f,g,t	1,7	
I	E <sub>4</sub>	Maiduguri	1,3,19	f,g,t	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Kouka	1,3,19	g,m,[t]	-	
I	E <sub>4</sub>	Dessau	1,3,15,19	g,s,t	-	
I	E <sub>4</sub>	Senftenberg	1,3,19	g,[s],t		Senftenberg may possess H phase Rz <sub>37</sub> or Rz <sub>43</sub> or Rz <sub>45</sub> or Rz <sub>46</sub> . Simsbury (1,3,19:Rz <sub>27</sub> :-) is now considered an H phase Rz <sub>27</sub> of Senftenberg.
I	E <sub>4</sub>	Stratford	1,3,19	i	1,2	
I	E <sub>4</sub>	Chichester	1,3,19	i	1,6	
I	E <sub>4</sub>	Machaga	1,3,19	i	e,n,x	
I	E <sub>4</sub>	Avonmouth	1,3,19	i	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Zuilen	1,3,19	i	l,w	
I	E <sub>4</sub>	Taksony	1,3,19	[i]	z <sub>6</sub>	
I	E <sub>4</sub>	Bethune	1,3,19	k	1,7	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>4</sub>	Ngor	1,3,19	l,v	1,5	
I	E <sub>4</sub>	Parkroyal	1,3,19	l,v	1,7	
I	E <sub>4</sub>	Svedvi	1,3,19	l,v	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Fulda	1,3,19	l,w	1,5	
I	E <sub>4</sub>	Westerstede	1,3,19	l,z <sub>13</sub>	[1,2]	
I	E <sub>4</sub>	Winterthur	1,3,19	l,z <sub>13</sub>	1,6	
I	E <sub>4</sub>	Lokstedt	1,3,19	l,z <sub>13</sub> ,z <sub>28</sub>	1,2	
I	E <sub>4</sub>	Stuivenberg	1,3,19	l,z <sub>13</sub> ,z <sub>28</sub>	1,5	
I	E <sub>4</sub>	Bedford	1,3,19	l,z <sub>13</sub> ,z <sub>28</sub>	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Tomelilla	1,3,19	l,z <sub>28</sub>	1,7	
I	E <sub>4</sub>	Kindia	1,3,19	l,z <sub>28</sub>	e,n,x	
I	E <sub>4</sub>	Cannstatt	1,3,19	m,t	–	
I	E <sub>4</sub>	Yalding	1,3,19	r	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Fareham	1,3,19	r,i	l,w	
I	E <sub>4</sub>	Simsbury	1,3,19	Rz <sub>27</sub>	–	IP combined Simsbury with Senftenberg 1,3,19:g,[s],t:-. Simsbury is now considered an R phase of Senftenberg. The name Simsbury has been dropped.
I	E <sub>4</sub>	Gatineau	1,3,19	y	1,5	
I	E <sub>4</sub>	Thies	1,3,19	y	1,7	
I	E <sub>4</sub>	Cannonhill	1,3, <u>15</u> ,19	y	e,n,x	
I	E <sub>4</sub>	Kinson	1,3,19	y	e,n,x	
I	E <sub>4</sub>	Slade	1,3,19	y	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Krefeld	1,3,19	y	l,w	
I	E <sub>4</sub>	Korlebu	1,3,19	z	1,5	
I	E <sub>4</sub>	Kainji	1,3,19	z	1,6	
I	E <sub>4</sub>	Lerum	1,3,19	z	1,7	
I	E <sub>4</sub>	Schoeneberg	1,3,19	z	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Carno	1,3,19	z	l,w	
I	E <sub>4</sub>	Hongkong	1,3,19	z	z <sub>6</sub>	
I	E <sub>4</sub>	Dallgow	1,3,19	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Llandoff	1,3,19	z <sub>29</sub>	[z <sub>6</sub> ]	
I	E <sub>4</sub>	Ochiogu	1,3,19	z <sub>38</sub>	[e,n,z <sub>15</sub> ]	
I	E <sub>4</sub>	Ilugun	1,3,10,19	z <sub>4</sub> ,z <sub>23</sub>	z <sub>6</sub>	
I	E <sub>4</sub>	Sambre	1,3,19	z <sub>4</sub> ,z <sub>24</sub>	–	
II	F		11	–	1,5	
I	F	Gallen	11	a	1,2	
I	F	Marseille	11	a	1,5	
VI	F		11	a	1,5	
I	F	Toowong	11	a	1,7	
II	F	Montgomery	11	a,[d]	[d]:e,n,z <sub>15</sub>	
I	F	Luciana	11	a	e,n,z <sub>15</sub>	
I	F	Epinay	11	a	l,z <sub>13</sub> ,z <sub>28</sub>	
II	F	Glencairn	11	a	z <sub>6</sub> :z <sub>42</sub>	
I	F	Atento	11	b	1,2	
I	F	Leeuwarden	11	b	1,5	
I	F	Wohlen	11	b	1,6	
VI	F		11	b	1,7	
VI	F	Srinagar	11	b	e,n,x	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	F	Pharr	11	b	e,n,z <sub>15</sub>	
I	F	Chiredzi	11	c	1,5	
I	F	Woodinville	11	c	e,n,x	
II	F		11	c	e,n,z <sub>15</sub>	
I	F	Ati	11	d	1,2	
I	F	Gustavia	11	d	1,5	
I	F	Chandans	11	d	e,n,x:[r]	
I	F	Pennsylvania	11	d	e,n,z <sub>15</sub>	
I	F	Findorff	11	d	z <sub>6</sub>	
I	F	Chingola	11	e,h	1,2	
I	F	Adamstua	11	e,h	1,6	
I	F	Redhill	11	e,h	1,z <sub>13</sub> ,z <sub>28</sub>	
II	F	Grabouw	11	g,[m],s,t	[z <sub>39</sub> ]	
I	F	Missouri	11	g,s,t	–	
IV	F	Mundsborg	11	g,z <sub>51</sub>	–	
I	F	Aberdeen	11	i	1,2	
I	F	Brijbhumi	11	i	1,5	
I	F	Heerlen	11	i	1,6	
I	F	Veneziana	11	i	e,n,x	
I	F	Pretoria	11	k	1,2	
I	F	Abaetetuba	11	k	1,5	
I	F	Sharon	11	k	1,6	
I	F	Colobane	11	k	1,7	
I	F	Kisarawe	11	k	e,n,x,[z <sub>15</sub> ]	
I	F	Mannheim	11	k	1,w	
I	F	Amba	11	k	1,z <sub>13</sub> ,z <sub>28</sub>	
IIIb	F		11	k	z <sub>53</sub>	(Ar. 17:29:25)
I	F	Stendal	11	l,v	1,2	
I	F	Maracaibo	11	l,v	1,5	
I	F	Fann	11	l,v	e,n,x	
I	F	Bullbay	11	l,v	e,n,z <sub>15</sub>	
IIIb	F		11	l,v	z <sub>53</sub>	(Ar. 17:23:25)
IIIb	F		11	l,v	z	(Ar. 17:23:31). May possess H phase Rz <sub>56</sub> (Ar. 38).
I	F	Glidji	11	l,w	1,5	
I	F	Connecticut	11	1,z <sub>13</sub> ,z <sub>28</sub>	1,5	
I	F	Osnabrueck	11	1,z <sub>13</sub> ,z <sub>28</sub>	e,n,x	
II	F	Huila	11	1,z <sub>28</sub>	e,n,x	
I	F	Moers	11	m,t	–	
II	F	Lincoln	11	m,t	e,n,x	
I	F	Senegal	11	r	1,5	
I	F	Rubislaw	11	r	[e,n,x]	
I	F	Volta	11	r	1,z <sub>13</sub> ,z <sub>28</sub>	
I	F	Euston	11	r,i	e,n,x,z <sub>15</sub>	
I	F	Solt	11	y	1,5	
I	F	Jalisco	11	y	1,7	
I	F	Herzliya	11	y	e,n,x	
I	F	Crewe	11	z	1,5	
I	F	Maroua	11	z	1,7	
II	F		11	z	e,n,x	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	F	Soutpan	11	z	Z <sub>39</sub>	
I	F	Nyanza	11	z	Z <sub>6</sub>	
I	F	Wentworth	11	Z <sub>10</sub>	1,2	
I	F	Straengnaes	11	Z <sub>10</sub>	1,5	
I	F	Telhashomer	11	Z <sub>10</sub>	e,n,x	
I	F	Lene	11	Z <sub>38</sub>	–	
IIIa	F		11	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 17:1,2,5:-)
IV	F	Parera	11	Z <sub>4</sub> ,Z <sub>23</sub>	–	
I	F	Remete	11	Z <sub>4</sub> ,Z <sub>23</sub>	1,6	
I	F	Etterbeek	11	Z <sub>4</sub> ,Z <sub>23</sub>	e,n,Z <sub>15</sub>	
I	F	Yehuda	11	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IV	F		11	Z <sub>4</sub> ,Z <sub>32</sub>	–	
I	F	Maastricht	11	Z <sub>41</sub>	1,2	
II	G		13,23	–	1,6	
I	G	Chagoua	<u>1</u> ,13,23	a	1,5	
II	G		<u>1</u> ,13,23	a	1,5	
I	G	Mim	13,22	a	1,6	
II	G		13,22	a	e,n,x	
I	G	Wyldegreen	<u>1</u> ,13,23	a	1,w	
I	G	Marshall	13,22	a	1,Z <sub>13</sub> ,Z <sub>28</sub>	
II	G	Tygerberg	<u>1</u> ,13,23	a	Z <sub>42</sub>	
I	G	Atlanta	13,23	b	–	Atlanta was combined with Mississippi ( <u>1</u> ,13,23:b:1,5). The name Atlanta has been dropped.
I	G	Ibadan	13,22	b	1,5	
I	G	Mississippi	<u>1</u> ,13,23	b	[1,5]	
II	G	Acres	<u>1</u> ,13,23	b	[1,5]:Z <sub>42</sub>	
I	G	Bracknell	13,23	b	1,6	
I	G	Oudwijk	13,22	b	1,6	
I	G	Rottnest	<u>1</u> ,13,22	b	1,7	
I	G	Ullevi	<u>1</u> ,13,23	b	e,n,x	
I	G	Vaertan	13,22	b	e,n,x	
I	G	Bahati	13,22	b	e,n,Z <sub>15</sub>	
I	G	Durham	13,23	b	e,n,Z <sub>15</sub>	
II	G		<u>1</u> ,13,22	b	Z <sub>42</sub>	
I	G	Haouaria	13,22	c	e,n,x,Z <sub>15</sub>	
I	G	Handen	<u>1</u> ,13,23	d	1,2	
I	G	Mishmarhaemek	<u>1</u> ,13,23	d	1,5	
I	G	Friedenau	13,22	d	1,6	
I	G	Wichita	1,13,23	d	1,6	Wichita may possess H phase Rz <sub>37</sub> .
I	G	Grumpensis	<u>1</u> ,13,23	d	1,7	
II	G		13,23	d	e,n,x	
I	G	Diguel	<u>1</u> ,13,22	d	e,n,Z <sub>15</sub>	
I	G	Teitelkebir	13,23	d	e,n,Z <sub>15</sub>	
II	G		<u>1</u> ,13,23	d	e,n,Z <sub>15</sub>	
I	G	Putten	13,23	d	1,w	
I	G	Isuge	13,23	d	Z <sub>6</sub>	
I	G	Tschangu	<u>1</u> ,13,23	e,h	1,5	
I	G	Willemstad	1,13,22	e,h	1,6	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	G	Vridi	<u>1</u> ,13,23	e,h	l,w	
II	G	Epping	<u>1</u> ,13,23	e,n,x	1,[5],7	
I	G	Raus	13,22	f,g	e,n,x	
I	G	Havana	<u>1</u> ,13,23	f,g,[s]	–	Havana may possess H phase RZ <sub>45</sub> or RZ <sub>79</sub> .
IIIa	G		<u>1</u> ,13,23	g,51	–	(Ar. 18:13,14:-)
I	G	Bron	13,22	g,m	[e,n,z15]	
II	G		<u>1</u> ,13,23	g,m,s,t	1,5	
II	G	Luanshya	<u>1</u> ,13,23	g,m,s,t	[e,n,x]	IP combined Luanshya with Kraaifontein ( <u>1</u> ,13,23:g,m,t:[e,n,x]) to form Luanshya <u>1</u> ,13,23:g,m,[s],t:[e,n,x].
II	G		<u>1</u> ,13,23	g,[s],t	z <sub>42</sub>	
II	G	Limbe	<u>1</u> ,13,22	g,m,t	[1,5]	
II	G	Kraaifontein	<u>1</u> ,13,23	g,m,t	[e,n,x]	IP combined Kraaifontein with Luanshya ( <u>1</u> ,13,23:g,m,s,t:[e,n,x]) to form Luanshya <u>1</u> ,13,23:g,m,[s],t:[e,n,x]. The name Kraaifontein has been dropped.
I	G	Agbeni	<u>1</u> ,13,23	g,m,[t]	–	
I	G	Congo	13,23	g,[m],[s],t	–	IP calls Congo 13,23:g,m,s,t:-.
I	G	Newyork	13,22	g,s,t	–	
I	G	Okatie	13,23	g,[s],t	–	
II	G		<u>1</u> ,13,23	g,m,s,t	z <sub>42</sub>	
II	G	Gojenberg	<u>1</u> ,13,23	g,t	1,5	
II	G	Rotterdam	<u>1</u> ,13,22	g,t	1,5	
V	G		<u>1</u> ,13,22	i	–	
I	G	Idikan	<u>1</u> ,13,23	i	1,5	
I	G	Jukestown	13,23	i	e,n,z <sub>15</sub>	
I	G	Kedougou	<u>1</u> ,13,23	i	l,w	
I	G	Marburg	13,23	k	–	
II	G		13,22	k	1,5:z <sub>42</sub>	
II	G		13,23	k	z <sub>41</sub>	
I	G	Lovelace	13,22	l,v	1,5	
IIIb	G		13,22	l,v	1,5,7	(Ar. 18:23:30)
I	G	Borbeck	13,22	l,v	1,6	
I	G	Nanga	<u>1</u> ,13,23	l,v	e,n,z <sub>15</sub>	
II	G		13,23	l,w	e,n,x	
II	G		13,22	l,z <sub>28</sub>	1,5	
II	G		13,23	l,z <sub>28</sub>	1,5	
II	G	Vredelust	<u>1</u> ,13,23	l,z <sub>28</sub>	z <sub>42</sub>	
II	G		13,23	l,z <sub>28</sub>	z <sub>6</sub>	
I	G	Kintambo	<u>1</u> ,13,23	m,t	–	
I	G	Washington	13,22	m,t	–	
II	G	Katesgrove	<u>1</u> ,13,23	m,t	1,5	
II	G	Worcester	1,13,23	m,t	e,n,x	
II	G	Boulders	<u>1</u> ,13,23	m,t	z <sub>42</sub>	
II	G		13,22	m,t	z <sub>42</sub> :z <sub>39</sub>	
V	G		13,22	r	–	
I	G	Adjame	13,23	r	1,6	
I	G	Linton	13,23	r	e,n,z <sub>15</sub>	
I	G	Tanger	<u>1</u> ,13,22	y	1,6	
I	G	Yarrabah	13,23	y	1,7	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	G	Ordonez	<u>1</u> ,13,23	y	1,w	
I	G	Tunis	<u>1</u> ,13,23	y	z <sub>6</sub>	
II	G		13,22	z	–	
II	G	Nachshonim	<u>1</u> ,13,23	z	1,5	
I	G	Farmsen	13,23	z	1,6	
I	G	Poona	<u>1</u> ,13,22	z	1,6	Poona may possess H phase Rz <sub>59</sub> .
I	G	Bristol	13,22	z	1,7	
I	G	Tanzania	<u>1</u> ,13,22	z	e,n,z <sub>15</sub>	
I	G	Worthington	<u>1</u> ,13,23	z	1,w	Worthington may possess H phase Rz <sub>45</sub> .
II	G		<u>1</u> ,13,23	z	z <sub>42</sub>	
I	G	Roodepoort	<u>1</u> ,13,22	z <sub>10</sub>	1,5	
I	G	Demerara	13,23	z <sub>10</sub>	1,w	
II	G		<u>1</u> ,13,22	z <sub>10</sub>	z <sub>6</sub>	
I	G	Agoueve	13,22	z <sub>29</sub>	–	
I	G	Cubana	<u>1</u> ,13,23	z <sub>29</sub>	–	Cubana may possess H phase Rz <sub>37</sub> or Rz <sub>43</sub> .
II	G	Clifton	13,22	z <sub>29</sub>	1,5	
II	G		<u>1</u> ,13,23	z <sub>29</sub>	1,5	
II	G	Goodwood	13,22	z <sub>29</sub>	e,n,x	
II	G		<u>1</u> ,13,23	z <sub>29</sub>	e,n,x	
I	G	Mampong	13,22	z <sub>35</sub>	1,6	
I	G	Anna	13,23	z <sub>35</sub>	e,n,z <sub>15</sub>	
I	G	Nimes	13,22	z <sub>35</sub>	e,n,z <sub>15</sub>	
I	G	Fanti	13,23	z <sub>38</sub>	–	
I	G	Leiden	13,22	z <sub>38</sub>	–	
II	G		<u>1</u> ,13,23	z <sub>39</sub>	1,5,7	
II	G		13,22	z <sub>39</sub>	1,7	
I	G	Ajiobo	13,23	z <sub>4</sub> ,z <sub>23</sub>	–	
IIIa	G		13,22	z <sub>4</sub> ,z <sub>23</sub>	–	(Ar. 18:1,2,5:-)
I	G	Ried	<u>1</u> ,13,22	z <sub>4</sub> ,z <sub>23</sub>	[e,n,z <sub>15</sub> ]	
IIIa	G		13,23	z <sub>4</sub> ,z <sub>23</sub> ,z <sub>32</sub>	–	(Ar. 18:1,6,7:-). CDC would call this 1,6,7,9.
I	G	Romanby	<u>1</u> ,13,23	z <sub>4</sub> ,z <sub>24</sub>	–	
IIIa	G		<u>1</u> ,13,23	z <sub>4</sub> ,z <sub>24</sub>	–	(Ar. 18:1,3,11:-)
II	G	Stevenage	<u>1</u> ,13,23	[z <sub>42</sub> ]	1,[5],7	
I	H	Garba	1,6,14,25	a	1,5	
VI	H		[1],6,14	a	1,5	
VI	H	Ferlac	1,6,14,25	a	e,n,x	
I	H	Banjul	1,6,14,25	a	e,n,z <sub>15</sub>	
I	H	Ndjamena	1,6,14,25	b	1,2	
I	H	Kuntair	1,6,14,25	b	1,5	
I	H	Tucson	[1],6,14,[25]	b	[1,7]	
IIIb	H		(6),14	b	e,n,x,z <sub>15</sub>	(Ar. 7a,7c:43:28)
I	H	Blijdorp	1,6,14,25	c	1,5	
I	H	Kassberg	1,6,14,25	c	1,6	
I	H	Runby	1,6,14,25	c	e,n,x	
I	H	Minna	1,6,14,25	c	1,w	
I	H	Finkenwerder	[1],6,14,[25]	d	1,5	
I	H	Heves	6,14,[24]	d	1,5	
I	H	Woodhull	1,6,14,25	d	1,6	
I	H	Florida	[1],6,14,[25]	d	1,7	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	H	Midway	6,14,24	d	1,7	
I	H	Charity	[1],6,14,[25]	d	e,n,x	
I	H	Lindern	6,14,[24]	d	e,n,x	
I	H	Teko	1,6,14,25	d	e,n,z <sub>15</sub>	
I	H	Encino	1,6,14,25	d	l,z <sub>13</sub> ,z <sub>28</sub>	
I	H	Albuquerque	1,6,14,24	d	z <sub>6</sub>	
I	H	Bahrenfeld	6,14,24	e,h	1,5	
I	H	Onderstepoort	1,6,14,[25]	e,h	1,5	
I	H	Magumeri	1,6,14,25	e,h	1,6	
I	H	Beaudesert	[1],6,14,[25]	e,h	1,7	
I	H	Warragul	[1],6,14,[25]	g,m	–	
I	H	Caracas	[1],6,14,[25]	g,m,s	–	
I	H	Sylvania	[1],6,14,[25]	g,p	–	
I	H	Catanzaro	6,14	g,s,t	–	
I	H	Mampeza	1,6,14,25	i	1,5	
I	H	Buzu	[1],6,14,[25]	i	1,7	
I	H	Schalkwijk	6,14,[24]	i	e,n,..	
I	H	Moussoro	1,6,14,25	i	e,n,z <sub>15</sub>	
I	H	Harburg	[1],6,14,[25]	k	1,5	
II	H		6,14,[24]	k	1,6	
II	H		6,14	k	[e,n,x]	
IIIb	H		(6),14	k	z	(Ar. 7a,7c:29:31)
IIIb	H		(6),14	k	z <sub>53</sub>	(Ar. 7a,7c:29:25)
II	H		1,6,14	k	z <sub>6</sub> :z <sub>42</sub>	
I	H	Boecker	[1],6,14,[25]	l,v	1,7	
I	H	Horsham	1,6,14,[25]	l,v	e,n,x	
IIIb	H		(6),14	l,v	z	(Ar. 7a,7c:23:31)
IIIb	H		(6),14	l,v	z <sub>35</sub>	(Ar. 7a,7c:23:21)
IIIb	H		(6),14	l,v	z <sub>53</sub>	(Ar. 7a,7c:23:25)
I	H	Aflao	1,6,14,25	l,z <sub>28</sub>	e,n,x	
I	H	Kaitaan	1,6,14,25	m,t	–	
II	H	Rooikrantz	1,6,14	m,t	1,5	
II	H	Emmerich	6,14	[m,t]	e,n,x	
IIIb	H		(6),14	r	z	(Ar. 7a,7c:24:31)
I	H	Istoria	1,6,14,25	r,i	1,5	
I	H	Surat	[1],6,14,[25]	[r],[i]	e,n,z <sub>15</sub>	
I	H	Carrau	6,14,[24]	y	1,7	
I	H	Madelia	1,6,14,25	y	1,7	
I	H	Fischerkietz	1,6,14,25	y	e,n,x	
I	H	Mornington	1,6,14,25	y	e,n,z <sub>15</sub>	
I	H	Homosassa	1,6,14,25	z	1,5	
I	H	Kanifing	1,6,14,25	z	1,6	
I	H	Soahanina	6,14,24	z	e,n,x	
I	H	Sundsvall	[1],6,14,[25]	z	e,n,x	
I	H	Royan	1,6,14,25	z	e,n,z <sub>15</sub>	
I	H	Poano	1,6,14,25	z	l,z <sub>13</sub> ,z <sub>28</sub>	
I	H	Nessa	1,6,14,25	z <sub>10</sub>	1,2	
VI	H	Bornheim	1,6,14,25	z <sub>10</sub>	1,(2),7	Bornheim was formerly in Subspecies II.
II	H	Simonstown	1,6,14	z <sub>10</sub>	1,5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIb	H		(6),14	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 7a,7c:27:28)
IIIb	H		(6),14	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 7a,7c:27:25)
II	H	Slangkop	1,6,14	Z <sub>10</sub>	Z <sub>6</sub> :Z <sub>42</sub>	
IIIb	H		(6),14	Z <sub>10</sub>	z:[Z <sub>53</sub> ]	(Ar. 7a,7c:27:31:[25])
I	H	Potosi	6,14	Z <sub>36</sub>	1,5	
I	H	Sara	1,6,14,25	Z <sub>38</sub>	[e,n,x]	
IV	H		6,14	Z <sub>4</sub> ,Z <sub>23</sub>	–	
I	H	Arapahoe	1,6,14	Z <sub>4</sub> ,Z <sub>23</sub>	1,5	
I	H	Bouso	1,6,14,25	Z <sub>4</sub> ,Z <sub>23</sub>	[e,n,Z <sub>15</sub> ]	
I	H	Chichiri	6,14,24	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	H	Uzaramo	1,6,14,25	Z <sub>4</sub> ,Z <sub>24</sub>	–	
II	H		1,6,14	Z <sub>42</sub>	1,6	
IIIb	H		(6),14	Z <sub>52</sub>	e,n,x,Z <sub>15</sub>	(Ar. 7a,7c:26:28)
IIIb	H		(6),14	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 7a,7c:26:21)
I	I	Hannover	16	a	1,2	
I	I	Brazil	16	a	1,5	
I	I	Amunigun	16	a	1,6	
I	I	Nyeko	16	a	1,7	
I	I	Togba	16	a	e,n,x	
I	I	Fischerhuetten	16	a	e,n,Z <sub>15</sub>	
I	I	Heron	16	a	Z <sub>6</sub>	
I	I	Hull	16	b	1,2	
I	I	Wa	16	b	1,5	
I	I	Glasgow	16	b	1,6	
I	I	Hvittingfoss	16	b	e,n,x	
II	I		16	b	e,n,x	
I	I	Sangera	16	b	e,n,Z <sub>15</sub>	
I	I	Vege sack	16	b	1,w	
II	I		16	b	Z <sub>39</sub>	
II	I		16	b	Z <sub>42</sub>	
I	I	Malstatt	16	b	Z <sub>6</sub>	
I	I	Vancouver	16	c	1,5	
I	I	Gafsa	16	c	1,6	
I	I	Shamba	16	c	e,n,x	
I	I	Hithergreen	16	c	e,n,Z <sub>15</sub>	
I	I	Yoruba	16	c	1,w	
I	I	Oldenburg	16	d	1,2	
I	I	Sculcoates	16	d	1,5	
II	I		16	d	1,5	
I	I	Sherbrooke	16	d	1,6	
I	I	Gaminara	16	d	1,7	
I	I	Barranquilla	16	d	e,n,x	
I	I	Nottingham	16	d	e,n,Z <sub>15</sub>	
I	I	Caen	16	d	1,w	
I	I	Barmbek	16	d	Z <sub>6</sub>	
I	I	Malakal	16	e,h	1,2	
I	I	Saboya	16	e,h	1,5	
I	I	Rhydyfelin	16	e,h	e,n,x	
I	I	Weston	16	e,h	Z <sub>6</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	I	Bellville	16	e,n,x	1,(5),7	
II	I	Elsiesrivier	16	[e,n,x]	1,6;z <sub>42</sub>	
I	I	Tees	16	f,g	–	
I	I	Nikolaifleet	16	g,m,s	–	
I	I	Adeoyo	16	g,m,[t]	–	
II	I	Mobeni	16	g,[m],[s],t	[e,n,x]	
II	I		16	g,[m],[s],t	z <sub>42</sub>	
I	I	Cardoner	16	g,s,t	–	
II	I	Merseyside	16	g,t	[1,5]	
I	I	Amina	16	i	1,5	
I	I	Agbara	16	i	1,6	
I	I	Wisbech	16	i	1,7	
I	I	Frankfurt	16	i	e,n,z <sub>15</sub>	
I	I	Pisa	16	i	1,w	
IIIb	I		16	i	z <sub>35</sub>	(Ar. 25:33:21)
I	I	Abobo	16	i	z <sub>6</sub>	
I	I	Szentes	16	k	1,2	
I	I	Nuatja	16	k	e,n,x	
I	I	Orientalis	16	k	e,n,z <sub>15</sub>	
IIIb	I		16	k	z	(Ar. 25:29:31)
IIIb	I		16	k	z <sub>53</sub>	(Ar. 25:29:25)
IIIb	I		16	(k)	z <sub>35</sub>	(Ar. 25:22:21)
IIIb	I		16	l,v	1,5,7	(Ar. 25:23:30)
I	I	Shanghai	16	l,v	1,6	
I	I	Welikade	16	l,v	1,7	
I	I	Salford	16	l,v	e,n,x	
I	I	Burgas	16	l,v	e,n,z <sub>15</sub>	
IIIb	I		16	l,v	z <sub>35</sub>	(Ar. 25:23:21)
IIIb	I		16	l,v	z <sub>53</sub>	(Ar. 25:23:25)
I	I	Losangeles	16	l,v	z <sub>6</sub>	
IIIb	I		16	l,v	z:[z <sub>61</sub> ]	(Ar. 25:23:31:[41])
I	I	Zigong	16	l,w	1,5	
I	I	Westeinde	16	l,w	1,6	
I	I	Brooklyn	16	l,w	e,n,x	
I	I	Lomnava	16	l,w	e,n,z <sub>15</sub>	
II	I	Noordhoek	16	l,w	z <sub>6</sub>	
I	I	Mandera	16	l,z <sub>13</sub>	e,n,z <sub>15</sub>	
I	I	Battle	16	l,z <sub>13</sub> ,z <sub>28</sub>	1,6	
I	I	Ablogame	16	l,z <sub>13</sub> ,z <sub>28</sub>	z <sub>6</sub>	
I	I	Enugu	16	l,[z <sub>13</sub> ],z <sub>28</sub>	[1,5]	
II	I	Sarepta	16	l,z <sub>28</sub>	z <sub>42</sub>	
I	I	Mpouto	16	m,t	–	
II	I		16	m,t	e,n,x	
II	I	Rowburton	16	m,t	[z <sub>42</sub> ]	
I	I	Ivory	16	r	1,6	
I	I	Brunflo	16	r	1,7	
I	I	Annedal	16	r,i	e,n,x	
I	I	Zwickau	16	r,i	e,n,z <sub>15</sub>	
I	I	Rovaniemi	16	r,[i]	1,5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	I	Saphra	16	y	1,5	
I	I	Akuafo	16	y	1,6	
I	I	Kikoma	16	y	e,n,x	
I	I	Avignon	16	y	e,n,z <sub>15</sub>	
I	I	Fortlamy	16	z	1,6	
I	I	Lingwala	16	z	1,7	
II	I	Louwbester	16	z	[e,n,x]	
I	I	Brevik	16	z	e,n,[x],z <sub>15</sub>	
II	I		16	z	z <sub>42</sub>	
I	I	Bouake	16	z	z <sub>6</sub>	
I	I	Badagry	16	z <sub>10</sub>	1,5	
IIIb	I		16	z <sub>10</sub>	1,5,7	(Ar. 25:27:30)
I	I	Lisboa	16	z <sub>10</sub>	1,6	
IIIb	I		16	z <sub>10</sub>	e,n,x,z <sub>15</sub>	(Ar. 25:27:28)
I	I	Redlands	16	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	I	Angouleme	16	z <sub>10</sub>	z <sub>6</sub>	
I	I	Saloniki	16	z <sub>29</sub>	–	
II	I		16	z <sub>29</sub>	1,5	
II	I	Jacksonville	16	z <sub>29</sub>	[e,n,x]	
I	I	Trier	16	z <sub>35</sub>	1,6	
I	I	Dakota	16	z <sub>35</sub>	e,n,z <sub>15</sub>	
IV	I		16	z <sub>36</sub>	–	
I	I	Naware	16	z <sub>38</sub>	–	
I	I	Grancanaria	16	z <sub>39</sub>	[1,6]	Grancanaria can be d-tartrate neg., dulcitol neg., ONPG pos., and anaerogenic.
II	I	Haddon	16	z <sub>4</sub> ,z <sub>23</sub>	–	
IV	I	Ochsenzoll	16	z <sub>4</sub> ,z <sub>23</sub>	–	
I	I	Kibi	16	z <sub>4</sub> ,z <sub>23</sub>	[1,6]	
II	I		16	z <sub>4</sub> ,z <sub>24</sub>	–	
IV	I	Chameleon	16	z <sub>4</sub> ,z <sub>32</sub>	–	
II	I	Woodstock	16	z <sub>42</sub>	1,[5],7	
IIIb	I		16	z <sub>52</sub>	z <sub>35</sub>	(Ar. 25:26:21)
II	I		16	z <sub>6</sub>	1,6	
I	J	Bonames	17	a	1,2	
I	J	Jangwani	17	a	1,5	
I	J	Kinondoni	17	a	e,n,x	
I	J	Kirkee	17	b	1,2	
I	J	Dahra	17	b	1,5	
II	J	Hillbrow	17	b	e,n,x,z <sub>15</sub>	
I	J	Bignona	17	b	e,n,z <sub>15</sub>	
II	J		17	b	z <sub>6</sub>	
I	J	Victoriaborg	17	c	1,6	
II	J	Woerden	17	c	z <sub>39</sub>	
I	J	Berlin	17	d	1,5	
I	J	Niamey	17	d	1,w	
I	J	Jubilee	17	e,h	1,2	
II	J		17	e,n,x,z <sub>15</sub>	1,[5],7	
II	J	Verity	17	e,n,x,z <sub>15</sub>	1,6	
II	J	Bleaton	17	(f),g,t	[e,n,x,z <sub>15</sub> ]	IP has dropped f.

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	J		17	g,m,s,t	–	
I	J	Lowestoft	17	g,s,t	–	
II	J		17	g,t	Z <sub>39</sub>	
I	J	Ahanou	17	i	1,7	
IIIb	J		17	i	Z <sub>35</sub>	(Ar. 12:33:21)
II	J		17	k	–	
I	J	Irenea	17	k	1,5	
I	J	Warri	17	k	1,7	
I	J	Matadi	17	k	e,n,x	
I	J	Zaria	17	k	e,n,Z <sub>15</sub>	
IIIb	J		17	k	z	(Ar. 12:29:32)
I	J	Morotai	17	l,v	1,2	
I	J	Michigan	17	l,v	1,5	
I	J	Lancaster	17	l,v	1,7	
I	J	Carmel	17	l,v	e,n,x	
IIIb	J		17	l,v	e,n,x,Z <sub>15</sub>	(Ar. 12:23:28)
IIIb	J		17	l,v	Z <sub>35</sub>	(Ar. 12:23:21)
I	J	Granlo	17	l,Z <sub>28</sub>	e,n,x	
I	J	Bama	17	m,t	–	
II	J		17	m,t	–	
I	J	Lode	17	r	1,2	
IIIb	J		17	r	z	(Ar. 12:24:31)
II	J		17	y	–	
I	J	Hadejia	17	y	e,n,Z <sub>15</sub>	
I	J	Gori	17	z	1,2	
I	J	Warengo	17	z	1,5	
I	J	Tchamba	17	z	e,n,Z <sub>15</sub>	
II	J	Constantia	17	z	l,w:Z <sub>42</sub>	
I	J	Djibouti	17	Z <sub>10</sub>	e,n,x	
IIIb	J		17	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar.12:27:28). May possess H phase Rz <sub>56</sub> (Ar. 38).
IIIb	J		17	Z <sub>10</sub>	z	(Ar. 12:27:31)
I	J	Kandla	17	Z <sub>29</sub>	–	
IIIa	J		17	Z <sub>29</sub>	–	(Ar. 12:16,17,18:-)
IV	J		17	Z <sub>29</sub>	–	
IIIa	J		17	Z <sub>36</sub>	–	(Ar. 12:17,20:-)
IV	J		17	Z <sub>36</sub>	–	
IIIa	J		17	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 12:1,2,5:- and 12:1,2,6:-)
IIIa	J		17	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 12:1,6,7,9:-)
IIIa	J		17	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 12:1,3,11:-)
IIIa	J		17	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 12:1,6,7:- and 12:1,7,8:-)
I	K	Cotia	18	–	1,6	
I	K	Brazos	6,14,18	a	e,n,Z <sub>15</sub>	
I	K	Fluntern	6,14,18	b	1,5	
I	K	Rawash	6,14,18	c	e,n,x	
I	K	Groenekan	18	d	1,5	
I	K	Usumbura	6,14,18	d	1,7	
I	K	Pontypridd	18	g,m	–	
IIIa	K		18	g,Z <sub>51</sub>	–	(Ar. 7a,7b:13,14:-)
I	K	Memphis	18	k	1,5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIb	K		18	(k)	Z <sub>53</sub>	(Ar. 7a,7b:22:25)
IIIb	K		18	(k)	Z <sub>54</sub>	(Ar. 7a,7b:22:34)
IIIb	K		18	l,v	e,n,x,Z <sub>15</sub>	(Ar. 7a,7b:23:28)
I	K	Orlando	18	l,v	e,n,Z <sub>15</sub>	
IIIb	K		18	l,v	z	(Ar. 7a,7b:23:31)
IIIb	K		18	l,v	Z <sub>53</sub>	(Ar. 7a,7b:23:25)
I	K	Toulon	18	l,w	e,n,Z <sub>15</sub>	
I	K	Langenhorn	18	m,t	–	
II	K		18	m,t	1,5	
IIIb	K		18	r	z	(Ar. 7a,7b:24,31)
II	K		18	y	e,n,x,Z <sub>15</sub>	
I	K	Potengi	18	z	–	
I	K	Leer	18	Z <sub>10</sub>	1,5	
IIIb	K		18	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 7a,7b:27:28)
I	K	Carnac	18	Z <sub>10</sub>	Z <sub>6</sub>	
II	K	Zeist	18	Z <sub>10</sub>	Z <sub>6</sub>	
II	K	Beloha	18	Z <sub>36</sub>	–	
IV	K		18	Z <sub>36</sub> ,Z <sub>38</sub>	–	
I	K	Sinthia	18	Z <sub>38</sub>	–	
II	K		18	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	K		18	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 7a,7b:1,2,5:- and 7a,7b:1,2,6:-)
I	K	Cerro	18	Z <sub>4</sub> ,Z <sub>23</sub>	[1,5]	Cerro was combined with Siegburg (6,14,18:Z <sub>4</sub> ,Z <sub>23</sub> : [1,5]) and called Cerro. Cerro may possess H phase RZ <sub>45</sub> .
I	K	Siegburg	6,14,18	Z <sub>4</sub> ,Z <sub>23</sub>	[1,5]	IP combined Siegburg with Cerro (18:Z <sub>4</sub> ,Z <sub>23</sub> : [1,5]) to form Cerro 6,14,18:Z <sub>4</sub> ,Z <sub>23</sub> : [1,5]. Siegburg is now called Cerro var. O 14+. The name Siegburg has been dropped.
I	K	Aarhus	18	Z <sub>4</sub> ,Z <sub>23</sub>	Z <sub>64</sub>	
I	K	Blukwa	18	Z <sub>4</sub> ,Z <sub>24</sub>	–	
II	K		18	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIa	K	Shomron	18	Z <sub>4</sub> ,Z <sub>32</sub>	–	Shomron was formerly in Subspecies II, but is now combined with <i>Arizona</i> 7a,7b:1,7,8:-. The name Shomron has been dropped.
IIIa	K		18	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 7a,7b:1,6,7:- and 7a,7b:1,7,8:-)
I	K	Delmenhorst	18	Z <sub>71</sub>	–	
I	L	Assen	21	a	[1,5]	
II	L		21	b	1,5	
I	L	Ghana	21	b	1,6	
I	L	Minnesota	21	b	e,n,x	Minnesota may possess H phase RZ <sub>33</sub> or RZ <sub>49</sub> .
I	L	Hydra	21	c	1,6	
I	L	Rhone	21	c	e,n,x	
II	L		21	c	e,n,x	
IIIb	L		21	c	e,n,x,Z <sub>15</sub>	(Ar. 22:32:28)
I	L	Spartel	21	d	1,5	
I	L	Magwa	21	d	e,n,x	
I	L	Madison	21	d	Z <sub>6</sub>	
I	L	Good	21	f,g	e,n,x	
II	L		21	g,[m],[s],t	–	
IIIa	L		21	g,Z <sub>51</sub>	–	(Ar. 22:13,14:-)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IV	L		21	g,Z <sub>51</sub>	–	
I	L	Diourbel	21	i	1,2	
IIIb	L		21	i	1,5,7	(Ar. 22:33:30)
IIIb	L		21	i	e,n,x,Z <sub>15</sub>	(Ar. 22:33:28)
IIIb	L		21	k	e,n,x,Z <sub>15</sub>	(Ar. 22:29:28)
IIIb	L		21	k	z	(Ar. 22:29:31)
IIIb	L		21	l,v	z	(Ar. 22:23:31)
IIIb	L		21	l,v	Z <sub>57</sub>	(Ar. 22:23:40)
I	L	Keve	21	l,w	–	
I	L	Jambur	21	l,Z <sub>28</sub>	e,n,Z <sub>15</sub>	
II	L		21	m,t	–	
I	L	Mountmagnet	21	r	–	
IIIb	L		21	r	z	(Ar. 22:24:31)
I	L	Ibaragi	21	y	1,2	
I	L	Ruiru	21	y	e,n,x	
II	L		21	z	–	CDC does not have this.
IIIb	L		21	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 22:27:28)
IIIb	L		21	Z <sub>10</sub>	z	(Ar. 22:27:31)
IIIb	L		21	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 22:27:25)
II	L	Wandsbek	21	Z <sub>10</sub>	[Z <sub>6</sub> ]	
IIIa	L		21	Z <sub>29</sub>	–	(Ar. 22:16,17,18:-)
I	L	Gambaga	21	Z <sub>35</sub>	e,n,Z <sub>15</sub>	
IV	L		21	Z <sub>36</sub>	–	
I	L	Baguida	21	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	L		21	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 22:1,2,5:- and 22:1,2,6:-)
IV	L	Soesterberg	21	Z <sub>4</sub> ,Z <sub>23</sub>	–	
II	L	Gwaai	21	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIa	L		21	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 22:1,3,11:-)
IV	L		21	Z <sub>4</sub> ,Z <sub>32</sub>	–	
IIIb	L		21	Z <sub>65</sub>	e,n,x,Z <sub>15</sub>	(Ar. 22:32:28)
I	M	Solna	28	a	1,5	
I	M	Dakar	28	a	1,6	
I	M	Bakau	28	a	1,7	
I	M	Seattle	28	a	e,n,x	
II	M		28	a	e,n,x	
I	M	Honelis	28	a	e,n,Z <sub>15</sub>	
I	M	Dibra	28	a	Z <sub>6</sub>	
I	M	Moero	28	b	1,5	
I	M	Ashanti	28	b	1,6	
I	M	Bokanjac	28	b	1,7	
I	M	Soumbédioune	28	b	e,n,x	
II	M		28	b	e,n,x	
I	M	Langford	28	b	e,n,Z <sub>15</sub>	
II	M	Kaltenhausen	28	b	Z <sub>6</sub>	
I	M	Hermannswerder	28	c	1,5	
I	M	Eberswalde	28	c	1,6	
I	M	Halle	28	c	1,7	
I	M	Dresden	28	c	e,n,x	
I	M	Wedding	28	c	e,n,Z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	M	Techimani	28	c	Z <sub>6</sub>	
I	M	Amoutive	28	d	1,5	
I	M	Hatfield	28	d	1,6	
I	M	Mundonobo	28	d	1,7	
I	M	Mocamedes	28	d	e,n,x	
I	M	Patience	28	d	e,n,Z <sub>15</sub>	
I	M	Cullingworth	28	d	l,w	
I	M	Kpeme	28	e,h	1,7	
I	M	Gozo	28	e,h	e,n,Z <sub>15</sub>	
II	M		28	e,n,x	1,7	
I	M	Friedrichsfelde	28	f,g	-	
I	M	Yardley	28	g,m	1,6	
I	M	Abadina	28	g,m	[e,n,Z <sub>15</sub> ]	
I	M	Croft	28	g,m,s	[e,n,Z <sub>15</sub> ]	
II	M		28	g,m,t	e,n,x	
II	M		28	g,m,t	Z <sub>39</sub>	
II	M	Llandudno	28	g,[m],[s],t	1,5	
I	M	Ona	28	g,s,t	-	
II	M		28	g,s,t	e,n,x	
I	M	Doorn	28	i	1,2	
I	M	Cotham	28	i	1,5	
I	M	Volksmarsdorf	28	i	1,6	
I	M	Dieuppeul	28	i	1,7	
I	M	Warnemuende	28	i	e,n,x	
I	M	Kuessel	28	i	e,n,Z <sub>15</sub>	
I	M	Douala	28	i	l,w	
I	M	Guildford	28	k	1,2	
I	M	Ilala	28	k	1,5	
I	M	Adamstown	28	k	1,6	
I	M	Ikeja	28	k	1,7	
I	M	Taunton	28	k	e,n,x	
I	M	Ank	28	k	e,n,Z <sub>15</sub>	
I	M	Leoben	28	l,v	1,5	
I	M	Vitkin	28	l,v	e,n,x	
I	M	Nashua	28	l,v	e,n,Z <sub>15</sub>	
I	M	Ramsey	28	l,w	1,6	
I	M	Catalunia	28	l,Z <sub>13</sub> ,Z <sub>28</sub>	1,5	
I	M	Penilla	28	l,Z <sub>13</sub> ,Z <sub>28</sub>	e,n,Z <sub>15</sub>	
II	M		28	l,Z <sub>28</sub>	1,5	
I	M	Fajara	28	l,Z <sub>28</sub>	e,n,x	
I	M	Morillons	28	m,t	1,6	
II	M		28	m,t	[e,n,x]	
I	M	Vinohrady	28	m,t	[e,n,Z <sub>15</sub> ]	
I	M	Bassadji	28	r	1,6	
I	M	Kibusi	28	r	e,n,x	
II	M	Oevelgoenne	28	r	e,n,Z <sub>15</sub>	
I	M	Fairfield	28	r	l,w	
I	M	Banco	28	r,i	1,7	
I	M	Chicago	28	r,[i]	1,5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	M	Sanktgeorg	28	r,[i]	e,n,z <sub>15</sub>	
I	M	Oskarshamn	28	y	1,2	
I	M	Nima	28	y	1,5	
I	M	Pomona	28	y	1,7	Pomona may possess H phases Rz <sub>60</sub> , Rz <sub>70</sub> or Rz <sub>80</sub> .
I	M	Kitenge	28	y	e,n,x	
I	M	Telaviv	28	y	e,n,z <sub>15</sub>	
I	M	Shomolu	28	y	1,w	
I	M	Selby	28	y	z <sub>6</sub>	
I	M	Vanier	28	z	1,5	
II	M		28	z	1,5	
I	M	Doel	28	z	1,6	
I	M	Ezra	28	z	1,7	
I	M	Brisbane	28	z	e,n,z <sub>15</sub>	
II	M	Ceres	28	z	z <sub>39</sub>	
I	M	Rogy	28	z <sub>10</sub>	1,2	
I	M	Farakan	28	z <sub>10</sub>	1,5	
I	M	Libreville	28	z <sub>10</sub>	1,6	
I	M	Malaysia	28	z <sub>10</sub>	1,7	
I	M	Umbilo	28	z <sub>10</sub>	e,n,x	
I	M	Luckenwalde	28	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	M	Moroto	28	z <sub>10</sub>	1,w	
IIIb	M		28	z <sub>10</sub>	z	(Ar. 35:27:31)
IIIb	M		28	z <sub>10</sub>	z:[z <sub>57</sub> ]	(Ar. 35:27:31:[40])
I	M	Djermaia	28	z <sub>29</sub>	–	
II	M		28	z <sub>29</sub>	1,5	
II	M		28	z <sub>29</sub>	e,n,x	
I	M	Konolfingen	28	z <sub>35</sub>	1,6	
I	M	Babili	28	z <sub>35</sub>	1,7	
I	M	Santander	28	z <sub>35</sub>	e,n,z <sub>15</sub>	
I	M	Aderike	28	z <sub>38</sub>	[e,n,z <sub>15</sub> ]	
I	M	Cannobio	28	z <sub>4</sub> ,z <sub>23</sub>	1,5	
I	M	Teltow	28	z <sub>4</sub> ,z <sub>23</sub>	1,6	
I	M	Babelsberg	28	z <sub>4</sub> ,z <sub>23</sub>	[e,n,z <sub>15</sub> ]	
I	N	Overvecht	30	a	1,2	
I	N	Zehlendorf	30	a	1,5	
I	N	Guarapiranga	30	a	e,n,x	
I	N	Doulassame	30	a	e,n,z <sub>15</sub>	
II	N	Odijk	30	a	z <sub>39</sub>	
I	N	Louga	30	b	1,2	
I	N	Aschersleben	30	b	1,5	
I	N	Urbana	30	b	e,n,x	
I	N	Neudorf	30	b	e,n,z <sub>15</sub>	
II	N		30	b	z <sub>6</sub>	
I	N	Zaire	30	c	1,7	
I	N	Morningside	30	c	e,n,z <sub>15</sub>	
II	N		30	c	z <sub>39</sub>	
I	N	Messina	30	d	1,5	
I	N	Livulu	30	e,h	1,2	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	N	Torhout	30	e,h	1,5	
I	N	Giessen	30	g,m,s	–	
II	N		30	g,m,s	e,n,x	
I	N	Godesberg	30	g,m,[t]	–	
I	N	Sternschanze	30	g,s,t	–	Sternschanze may possess H phase Rz <sub>59</sub> .
II	N	Slatograd	30	g,t	–	
I	N	Wayne	30	g,z <sub>51</sub>	–	
I	N	Landau	30	i	1,2	
I	N	Morehead	30	i	1,5	
I	N	Mjordan	30	i	e,n,z <sub>15</sub>	
I	N	Soerenga	30	i	1,w	
I	N	Hilversum	30	k	1,2	
I	N	Ramatgan	30	k	1,5	
I	N	Aqua	30	k	1,6	
I	N	Angoda	30	k	e,n,x	
II	N		30	k	e,n,x,z <sub>15</sub>	
I	N	Odozi	30	k	e,n,[x],z <sub>15</sub>	
I	N	Ligeo	30	l,v	1,2	
I	N	Donna	30	l,v	1,5	
I	N	Ockenheim	30	l,z <sub>13</sub> ,z <sub>28</sub>	1,6	
I	N	Morocco	30	l,z <sub>13</sub> ,z <sub>28</sub>	e,n,z <sub>15</sub>	
II	N		30	l,z <sub>28</sub>	z <sub>6</sub>	
II	N		30	m,t	–	
I	N	Grandhaven	30	r	1,2	
I	N	Gege	30	r	1,5	
I	N	Matopeni	30	y	1,2	
I	N	Bietri	30	y	1,5	
I	N	Steinplatz	30	y	1,6	
I	N	Baguirmi	30	y	e,n,x	
I	N	Nijmegen	30	y	e,n,z <sub>15</sub>	
I	N	Sada	30	z <sub>10</sub>	1,2	
I	N	Senneville	30	z <sub>10</sub>	1,5	
I	N	Kumasi	30	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	N	Aragua	30	z <sub>29</sub>	–	
I	N	Kokoli	30	z <sub>35</sub>	1,6	
I	N	Wuiti	30	z <sub>35</sub>	e,n,z <sub>15</sub>	
I	N	Ago	30	z <sub>38</sub>	–	
II	N		30	z <sub>39</sub>	1,7	
I	N	Stoneferry	30	z <sub>4</sub> ,z <sub>23</sub>	–	
I	N	Bodjonegoro	30	z <sub>4</sub> ,z <sub>24</sub>	–	
II	N		30	z <sub>6</sub>	1,6	
I	O	Umhlatazana	35	a	e,n,z <sub>15</sub>	
I	O	Tchad	35	b	–	
I	O	Gouloumbo	35	c	1,5	
I	O	Yolo	35	c	[e,n,z <sub>15</sub> ]	
II	O		35	d	1,5	
I	O	Dembe	35	d	1,w	
I	O	Gassi	35	e,h	z <sub>6</sub>	
I	O	Adelaide	35	f,g	–	Adelaide may possess H phase Rz <sub>27</sub> .

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	O	Ealing	35	g,m,s	–	
II	O		35	g,m,s,t	–	
I	O	Ebrie	35	g,m,t	–	
I	O	Anecho	35	g,s,t	–	
I	O	Agodi	35	g,t	–	
II	O		35	g,t	1,5	
II	O		35	g,t	Z <sub>42</sub>	
IIIa	O		35	g,Z <sub>51</sub>	–	(Ar. 20:13,14:-)
IIIb	O		35	i	e,n,x,Z <sub>15</sub>	(Ar. 20:33:28)
I	O	Gambia	35	i	e,n,Z <sub>15</sub>	
I	O	Bandia	35	i	l,w	
IIIb	O		35	i	z	(Ar. 20:33:31)
IIIb	O		35	i	Z <sub>35</sub>	(Ar. 20:33:21)
IIIb	O		35	i	Z <sub>53</sub>	(Ar. 20:33:25)
IIIb	O		35	k	e,n,x,Z <sub>15</sub>	(Ar. 20:29:28)
IIIb	O		35	k	z	(Ar. 20:29:31)
IIIb	O		35	k	Z <sub>53</sub>	(Ar. 20:29:25). May possess H phase RZ <sub>50</sub> (Ar.42).
IIIb	O		35	(k)	z	(Ar. 20:22:31)
IIIb	O		35	(k)	Z <sub>35</sub>	(Ar. 20:22:21)
IIIb	O		35	l,v	1,5,7	(Ar. 20:23:30)
IIIb	O		35	l,v	e,n,x,Z <sub>15</sub>	(Ar. 20:23:28)
IIIb	O		35	l,v	Z <sub>35</sub>	(Ar. 20:23:21)
II	O		35	l,Z <sub>28</sub>	–	
I	O	Monschau	35	m,t	-	
II	O		35	m,t	-	
IIIb	O		35	r	e,n,x,Z <sub>15</sub>	(Ar. 20:24:28)
I	O	Massakory	35	r	l,w	
IIIb	O		35	r	z	(Ar. 20:24:31)
IIIb	O		35	r	Z <sub>35</sub>	(Ar. 20:24:21)
IIIb	O		35	r	Z <sub>61</sub>	(Ar. 20:24:41)
I	O	Camberene	35	Z <sub>10</sub>	1,5	
I	O	Enschede	35	Z <sub>10</sub>	l,w	
IIIb	O		35	Z <sub>10</sub>	Z <sub>35</sub>	(Ar. 20:27:21)
I	O	Ligna	35	Z <sub>10</sub>	Z <sub>6</sub>	
I	O	Widemarsh	35	Z <sub>29</sub>	–	
IIIa	O		35	Z <sub>29</sub>	–	(Ar. 20:16,17,18:-)
II	O	Utbremen	35	Z <sub>29</sub>	e,n,x	
IIIa	O		35	Z <sub>36</sub>	–	(Ar. 20:17,20:-)
I	O	Haga	35	Z <sub>38</sub>	–	
I	O	Alachua	35	Z <sub>4,Z<sub>23</sub></sub>	–	Alachua may possess H phase RZ <sub>37</sub> or RZ <sub>45</sub> .
IIIa	O		35	Z <sub>4,Z<sub>23</sub></sub>	–	(Ar. 20:1,2,6:-)
I	O	Westphalia	35	Z <sub>4,Z<sub>24</sub></sub>	–	
IIIa	O		35	Z <sub>4,Z<sub>32</sub></sub>	–	(Ar. 20:1,7,8:-)
IIIb	O		35	Z <sub>52</sub>	1,5,7	(Ar. 20:26:30)
IIIb	O		35	Z <sub>52</sub>	e,n,x,Z <sub>15</sub>	(Ar. 20:26:28)
IIIb	O		35	Z <sub>52</sub>	z	(Ar. 20:26:31)
IIIb	O		35	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 20:26:21)
II	P		38	b	1,2	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	P	Rittersbach	38	b	e,n,z <sub>15</sub>	
I	P	Sheffield	38	c	1,5	
I	P	Kidderminster	38	c	1,6	
II	P	Carletonville	38	d	[1,5]	
I	P	Thiaroye	38	e,h	1,2	
I	P	Kasenyi	38	e,h	1,5	
I	P	Korovi	38	g,m,[s]	–	
II	P	Foulpointe	38	g,t	–	
IIIa	P		38	g,z <sub>51</sub>	–	(Ar. 16:13,14:-)
IV	P		38	g,z <sub>51</sub>	–	
I	P	Mgulani	38	i	1,2	
I	P	Lansing	38	i	1,5	
IIIb	P		38	i	z	(Ar. 16:33:31)
IIIb	P		38	i	z <sub>53</sub>	(Ar. 16:33:25)
I	P	Echa	38	k	1,2	
I	P	Mango	38	k	1,5	
I	P	Inverness	38	k	1,6	
I	P	Njala	38	k	e,n,x	
IIIb	P		38	k	e,n,x,z <sub>15</sub>	(Ar. 16:29:28)
IIIb	P		38	k	z	(Ar. 16:29:31)
IIIb	P		38	k	z <sub>53</sub>	(Ar. 16:29:25)
IIIb	P		38	(k)	1,5,7	(Ar. 16:22:30)
IIIb	P		38	(k)	z	(Ar. 16:22:31)
IIIb	P		38	(k)	z <sub>35</sub>	(Ar. 16:22:21). May possess H phase RZ <sub>56</sub> (Ar. 38).
IIIb	P		38	(k)	z <sub>54</sub>	(Ar. 16:22:34)
IIIb	P		38	(k)	z <sub>55</sub>	(Ar. 16:22:37)
I	P	Alger	38	l,v	1,2	
I	P	Kimberley	38	l,v	1,5	
I	P	Roan	38	l,v	e,n,x	
IIIb	P		38	l,v	z	(Ar. 16:23:31)
IIIb	P		38	l,v	z <sub>35</sub>	(Ar. 16:23:21)
IIIb	P		38	l,v	z <sub>53</sub> : [z <sub>54</sub> ]	(Ar. 16:23:25:[34])
I	P	Rothenburgsort	38	m,t	–	
I	P	Lindi	38	r	1,5	
IIIb	P		38	r	1,5,7	(Ar. 16:24:30)
I	P	Emmastad	38	r	1,6	
IIIb	P		38	r	e,n,x,z <sub>15</sub>	(Ar. 16:24:28)
IIIb	P		38	r	z:[z <sub>57</sub> ]	(Ar. 16:24:31:[40])
IIIb	P		38	r	z <sub>35</sub>	(Ar. 16:24:21)
I	P	Freetown	38	y	1,5	
I	P	Colombo	38	y	1,6	
I	P	Perth	38	y	e,n,x	
I	P	Stachus	38	z	–	
I	P	Neunkirchen	38	z <sub>10</sub>	–	
IIIb	P		38	z <sub>10</sub>	z	(Ar. 16:27:31)
IIIb	P		38	z <sub>10</sub>	z <sub>53</sub>	(Ar. 16:27:25)
I	P	Klouto	38	z <sub>38</sub>	–	
IIIa	P		38	z <sub>4</sub> ,z <sub>23</sub>	–	(Ar. 16:1,2,6:-)
IV	P		38	z <sub>4</sub> ,z <sub>23</sub>	–	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	P	Yoff	38	Z <sub>4</sub> ,Z <sub>23</sub>	1,2	
I	P	Bangkok	38	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIb	P		38	Z <sub>47</sub>	Z <sub>53</sub>	(Ar. 16:39:25)
IIIb	P		38	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 16:26:21)
IIIb	P		38	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 16:26:25)
IIIb	P		38	Z <sub>53</sub>	–	(Ar. 16:25:-). May possess H phase RZ <sub>50</sub> (Ar. 42) or RZ <sub>76</sub> (Ar. RZ <sub>76</sub> ). CDC does not have monophasic.
IIIa	P		38	Z <sub>61</sub>	–	(Ar. 16:41:-)
IIIb	P		38	Z <sub>61</sub>	Z <sub>53</sub>	(Ar. 16:41:25)
II	Q		39	–	1,7	
II	Q		39	a	Z <sub>39</sub>	
I	Q	Wandsworth	39	b	1,2	
I	Q	Abidjan	39	b	1,w	
II	Q		39	c	e,n,x	
I	Q	Logone	39	d	1,5	
I	Q	Mara	39	e,h	[1,5]	
II	Q		39	e,n,x	1,7	
II	Q		39	g,m,t	–	
I	Q	Hofit	39	i	1,5	
I	Q	Cumberland	39	i	e,n,x	
I	Q	Champaign	39	k	1,5	Champaign may possess H phase RZ <sub>48</sub>
II	Q		39	l,v	1,5	
I	Q	Kokomlemle	39	l,v	e,n,x	
I	Q	Oerlikon	39	l,v	e,n,Z <sub>15</sub>	
II	Q	Mondeor	39	l,Z <sub>28</sub>	e,n,x	
II	Q		39	l,Z <sub>28</sub>	Z <sub>39</sub>	
II	Q		39	m,t	e,n,x	
I	Q	Cook	39	RZ <sub>48</sub>	1,5	IP combined Cook with Champaign (39:k:1,5). The name Cook has been dropped.
I	Q	Anfo	39	y	1,2	
I	Q	Windermere	39	y	1,5	
I	Q	Hegau	39	Z <sub>10</sub>	–	
I	R	Shikmonah	40	a	1,5	
II	R		1,40	a	1,5	
II	R	Springs	40	a	Z <sub>39</sub>	
I	R	Greiz	40	a	Z <sub>6</sub>	
II	R		1,40	a	Z <sub>6</sub>	
II	R		40	b	–	
I	R	Riogrande	40	b	1,5	
I	R	Saugus	40	b	1,7	
I	R	Johannesburg	1,40	b	e,n,x	
I	R	Duval	1,40	b	e,n,Z <sub>15</sub>	
I	R	Benguella	40	b	Z <sub>6</sub>	
II	R	Suarez	1,40	c	e,n,x,Z <sub>15</sub>	
II	R		1,40	c	Z <sub>39</sub>	
II	R	Ottershaw	40	d	–	
I	R	Driffield	1,40	d	1,5	
I	R	Tilene	1,40	e,h	1,2	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	R		<u>1</u> ,40	e,n,x	1,[5],7	
II	R		<u>1</u> ,40	e,n,x,Z <sub>15</sub>	1,6	
I	R	Bijlmer	<u>1</u> ,40	g,m	–	
II	R	Alsterdorf	<u>1</u> ,40	g,m,[s],t	[1,5]	
II	R	Boksborg	40	g,[m],s,[t]	e,n,x	
II	R		<u>1</u> ,40	g,[m],[s],t	Z <sub>42</sub>	
II	R		<u>1</u> ,40	g,t	[e,n,x]	
II	R		<u>1</u> ,40	g,t	1,5	
II	R		<u>1</u> ,40	g,t	e,n,x,Z <sub>15</sub>	
II	R		40	g,t	Z <sub>39</sub>	
IV	R	Seminole	<u>1</u> ,40	g,Z <sub>51</sub>	–	
IIIb	R		40	g,Z <sub>51</sub>	[e,n,x,Z <sub>15</sub> ]	(Ar. 10a,10b:13,14:[28])
IIIb	R		40	i	1,5,7	(Ar. 10a,10b:33:30)
I	R	Goulfey	<u>1</u> ,40	k	1,5	
I	R	Allandale	<u>1</u> ,40	k	1,6	
I	R	Hann	40	k	e,n,x	
II	R	Sunnydale	<u>1</u> ,40	k	e,n,x,Z <sub>15</sub>	
IIIb	R		40	k	Z <sub>53</sub>	(Ar. 10a,10b:29:25)
II	R		40	k	Z <sub>6</sub>	
IIIb	R		40	k	z:Z <sub>57</sub>	(Ar. 10a,10b:29:31:40)
I	R	Millesi	<u>1</u> ,40	l,v	1,2	
I	R	Canary	40	l,v	1,6	
IIIb	R		40	l,v	z	(Ar. 10a,10b,(10c):23:31)
IIIb	R		40	l,v	Z <sub>53</sub>	(Ar. 10a,10b:23:25)
I	R	Overchurch	<u>1</u> ,40	l,w	[1,2]	
I	R	Tiko	40	l,Z <sub>13</sub> ,Z <sub>28</sub>	1,2	
I	R	Bukavu	<u>1</u> ,40	l,Z <sub>28</sub>	1,5	
II	R		<u>1</u> ,40	l,Z <sub>28</sub>	1,5:Z <sub>42</sub>	
I	R	Santhiaba	40	l,Z <sub>28</sub>	1,6	
II	R		<u>1</u> ,40	l,Z <sub>28</sub>	Z <sub>39</sub>	
IV	R		40	m,t	–	
II	R		40	m,t	Z <sub>39</sub>	
II	R		<u>1</u> ,40	m,t	Z <sub>42</sub>	
I	R	Odiene	40	y	1,5	
II	R	Bulawayo	<u>1</u> ,40	z	1,5	
I	R	Casamance	40	z	e,n,x	
II	R		<u>1</u> ,40	z	Z <sub>39</sub>	
II	R		40	z	Z <sub>42</sub>	
I	R	Nowawes	40	z	Z <sub>6</sub>	
II	R		<u>1</u> ,40	z	Z <sub>6</sub>	
IIIb	R		40	Z <sub>10</sub>	Z <sub>35</sub>	(Ar. 10a,10b:27:21)
I	R	Trotha	40	Z <sub>10</sub>	Z <sub>6</sub>	
I	R	Omifisan	40	Z <sub>29</sub>	–	
IIIa	R		40	Z <sub>29</sub>	–	(Ar. 10a,10b:16,18:-)
V	R		<u>1</u> ,40	Z <sub>35</sub>	–	
II	R	Fandran	<u>1</u> ,40	Z <sub>35</sub>	e,n,x,Z <sub>15</sub>	
I	R	Yekepa	<u>1</u> ,40	Z <sub>35</sub>	e,n,Z <sub>15</sub>	
IIIa	R		40	Z <sub>35</sub>	–	(Ar. 10a,10b:17,20:-)
II	R		40	Z <sub>39</sub>	1,5:Z <sub>42</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	R	Grunty	<u>1</u> ,40	Z <sub>39</sub>	1,6	
II	R		40	Z <sub>39</sub>	1,7	
IIIa	R		40	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 10a,10b:1,2,5:-; 10a,10b:1,2,5,6:-; and 10a,10b:1,2,6:-)
IV	R	Sachsenwald	<u>1</u> ,40	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	R		40	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 10a,10b:1,3,11:-)
IV	R	Degania var. subsp. IV	40	Z <sub>4</sub> ,Z <sub>24</sub>	–	
II	R	Degania	40	Z <sub>4</sub> ,Z <sub>24</sub>	[Z <sub>39</sub> ]	
IIIa	R		40	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 10a,10b:1,2,10:-; 10a,10c:1,2,10:-; and 10a,10b:1,7,8:-)
IIIa	R		40	Z <sub>4</sub> ,Z <sub>32</sub>	–	
IV	R	Bern	<u>1</u> ,40	Z <sub>4</sub> ,Z <sub>32</sub>	–	
I	R	Karamoja	40	Z <sub>41</sub>	1,2	
II	R		<u>1</u> ,40	Z <sub>42</sub>	1,6	
II	R		<u>1</u> ,40	[Z <sub>42</sub> ]	1,(5),7	
II	R		<u>1</u> ,40	Z <sub>6</sub>	1,5	
V	R		40	Z <sub>81</sub>	–	H Z <sub>81</sub> was formerly H a in <i>S. bongori</i> .
II	S		41	–	1,6	
I	S	Burundi	41	a	–	
II	S	Vietnam var. subsp. II	41	b	–	
II	S		41	b	[1,5]	
I	S	Vaugirard	41	b	1,6	
VI	S		41	b	1,7	
I	S	Sica	41	b	e,n,Z <sub>15</sub>	
I	S	Vietnam	41	b	Z <sub>6</sub>	
IIIb	S		41	c	e,n,x,Z <sub>15</sub>	(Ar. 13:32:28)
II	S		41	c	[Z <sub>6</sub> ]	
I	S	Egusi	41	d	[1,5]	
II	S	Hennepin	41	d	Z <sub>6</sub>	
II	S		41	g,m,s,t	Z <sub>6</sub>	
II	S	Lethe	41	g,t	–	
IIIa	S		41	g,Z <sub>51</sub>	–	(Ar. 13:13,14:-)
I	S	Samaru	41	i	1,5	
I	S	Verona	41	i	1,6	
I	S	Ferlo	41	k	1,6	
II	S		41	k	1,6	
II	S		41	k	[Z <sub>6</sub> ]	
IIIb	S		41	(k)	[Z <sub>35</sub> ]	(Ar. 13:22:[21])
II	S		41	l,Z <sub>13</sub> ,Z <sub>28</sub>	e,n,x,Z <sub>15</sub>	
I	S	Leatherhead	41	m,t	1,6	
I	S	Lubumbashi	41	r	1,5	
II	S	Dubrovnik	41	z	1,5	
II	S	Negev	41	Z <sub>10</sub>	1,2	
I	S	Leipzig	41	Z <sub>10</sub>	1,5	
I	S	Landala	41	Z <sub>10</sub>	1,6	
I	S	Inpraw	41	Z <sub>10</sub>	e,n,x	
II	S	Lurup	41	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	
II	S	Lichtenberg	41	Z <sub>10</sub>	[Z <sub>6</sub> ]	
I	S	Lodz	41	Z <sub>29</sub>	–	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIa	S		41	Z <sub>29</sub>	–	(Ar. 13:16,17,18:-)
IV	S		41	Z <sub>29</sub>	–	
I	S	Ahoutoue	41	Z <sub>35</sub>	1,6	
IIIa	S		41	Z <sub>36</sub>	–	(Ar. 13:17,20:-)
I	S	Offa	41	Z <sub>38</sub>	–	
IIIa	S		41	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 13:1,2,5:- and 13:1,2,6:-)
IV	S	Waycross var. subsp. IV41	41	Z <sub>4</sub> ,Z <sub>23</sub>	–	
I	S	Waycross	41	Z <sub>4</sub> ,Z <sub>23</sub>	[e,n,Z <sub>15</sub> ]	
IIIa	S		41	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 13:1,6,7,9:-)
IIIa	S		41	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 13:1,3,11:-)
I	S	Ipswich	41	Z <sub>4</sub> ,Z <sub>24</sub>	[1,5]	
IIIa	S		41	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 13:1,6,7:- and 13:1,7,8:-)
IV	S		41	Z <sub>52</sub>	–	
I	T	Faji	1,42	a	e,n,Z <sub>15</sub>	
II	T	Chinovum	42	b	1,5	
I	T	Orbe	42	b	1,6	
II	T	Uphill	42	b	e,n,x,Z <sub>15</sub>	
I	T	Tomegbe	1,42	b	e,n,Z <sub>15</sub>	
I	T	Egusitoo	1,42	b	Z <sub>6</sub>	
II	T		42	b	Z <sub>6</sub>	
I	T	Antwerpen	1,42	c	e,n,Z <sub>15</sub>	
I	T	Kampala	1,42	c	Z <sub>6</sub>	
II	T		42	d	Z <sub>6</sub>	
II	T		42	[e,n,x]	1,6	
II	T	Fremantle	42	(f),g,t	–	
IIIa	T		42	g,Z <sub>51</sub>	–	(Ar. 15:13,14:-)
IV	T		1,42	g,Z <sub>51</sub>	–	
I	T	Maricopa	1,42	g,Z <sub>51</sub>	1,5	
I	T	Borromea	42	i	1,6	
I	T	Kaneshie	1,42	i	1,w	
I	T	Middlesbrough	1,42	i	Z <sub>6</sub>	
IIIb	T		42	k	–	(Ar. 15:29:-)
I	T	Haferbreite	42	k	[1,6]	
IIIb	T		42	k	e,n,x,Z <sub>15</sub>	(Ar. 15:29:28)
IIIb	T		42	k	z	(Ar. 15:29:31)
IIIb	T		42	k	Z <sub>35</sub>	(Ar. 15:29:21)
I	T	Gwale	1,42	k	Z <sub>6</sub>	
IIIb	T		42	(k)	Z <sub>35</sub>	(Ar. 15:22:21)
IIIb	T		42	l,v	1,5,7	(Ar. 15:23:30)
II	T	Portbech	42	l,v	e,n,x,Z <sub>15</sub>	
IIIb	T		42	l,v	e,n,x,Z <sub>15</sub>	(Ar. 15:23:28)
I	T	Coogee	42	l,v	e,n,Z <sub>15</sub>	
IIIb	T		42	l,v	z	(Ar. 15:23:31)
IIIb	T		42	l,v	Z <sub>53</sub>	(Ar. 15:23:25)
II	T		1,42	l,w	e,n,x	
II	T		42	l,[z <sub>13</sub> ],Z <sub>28</sub>	[z <sub>6</sub> ]	
I	T	Waral	1,42	m,t	–	
II	T		42	m,t	[e,n,x,Z <sub>15</sub> ]	
II	T	Nairobi	42	r	–	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIb	T		42	r	–	(Ar. 15:24:-). May possess H phase RZ <sub>50</sub> (Ar. 42).
I	T	Sipane	<u>1</u> ,42	r	e,n,Z <sub>15</sub>	
I	T	Brive	<u>1</u> ,42	r	l,w	
IIIb	T		42	r	z	(Ar. 15:24:31)
IIIb	T		42	r	Z <sub>53</sub>	(Ar. 15:24:25)
I	T	Spalentor	<u>1</u> ,42	y	e,n,Z <sub>15</sub>	
I	T	Harvestehude	<u>1</u> ,42	y	Z <sub>6</sub>	
II	T	Detroit	42	z	1,5	
I	T	Ursenbach	<u>1</u> ,42	z	1,6	
II	T	Rand	42	z	e,n,x,Z <sub>15</sub>	
I	T	Melbourne	42	z	e,n,Z <sub>15</sub>	
II	T	Nuernberg	42	z	Z <sub>6</sub>	
IIIb	T		42	Z <sub>10</sub>	–	(Ar. 15:27:-). May possess H phase RZ <sub>56</sub> (Ar. 38) and RZ <sub>50</sub> (Ar. 42).
II	T		42	Z <sub>10</sub>	1,2	
II	T		42	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	
IIIb	T		42	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 15:27:28)
IIIb	T		42	Z <sub>10</sub>	z	(Ar. 15:27:31)
IIIb	T		42	Z <sub>10</sub>	Z <sub>35</sub>	(Ar. 15:27:21)
I	T	Loenga	<u>1</u> ,42	Z <sub>10</sub>	Z <sub>6</sub>	
II	T		42	Z <sub>10</sub>	Z <sub>6</sub>	
IIIb	T		42	Z <sub>10</sub>	Z <sub>67</sub>	(Ar. 15:27:46)
I	T	Djama	<u>1</u> ,42	Z <sub>29</sub>	[1,5]	
I	T	Kahla	<u>1</u> ,42	Z <sub>35</sub>	1,6	
I	T	Tema	<u>1</u> ,42	Z <sub>35</sub>	Z <sub>6</sub>	
I	T	Weslaco	42	Z <sub>36</sub>	–	
IV	T		42	Z <sub>36</sub>	–	
I	T	Vogan	<u>1</u> ,42	Z <sub>38</sub>	Z <sub>6</sub>	
IIIa	T		42	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 15:1,2,5:- and 15:1,2,6:-)
I	T	Gera	<u>1</u> ,42	Z <sub>4</sub> ,Z <sub>23</sub>	[1,6]	
I	T	Broc	42	Z <sub>4</sub> ,Z <sub>23</sub>	e,n,Z <sub>15</sub>	
I	T	Toricada	<u>1</u> ,42	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIa	T		42	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 15:1,3,11:-)
IV	T		<u>1</u> ,42	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	T	Taset	<u>1</u> ,42	Z <sub>41</sub>	–	
IIIb	T		42	Z <sub>52</sub>	z	(Ar. 15:26:31)
II	T		42	Z <sub>6</sub>	1,6	
I	U	Graz	43	a	1,2	
I	U	Berkeley	43	a	1,5	
II	U		43	a	1,5	
II	U		43	a	Z <sub>6</sub>	
II	U	Kommetje	43	b	Z <sub>42</sub>	
I	U	Montreal	43	c	1,5	
I	U	Orleans	43	d	1,5	
II	U		43	d	e,n,x,Z <sub>15</sub>	
II	U		43	d	Z <sub>39</sub>	
II	U		43	d	Z <sub>42</sub>	
II	U		43	e,n,x,Z <sub>15</sub>	1,(5),7	
II	U		43	e,n,x,Z <sub>15</sub>	1,6	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	U	Milwaukee	43	f,g,[t]	–	
II	U	Mosselbay	43	g,m,[s],t	[z <sub>42</sub> ]	
II	U	Veddel	43	g,t	–	
II	U		43	g,t	1,5	
IIIa	U		43	g,z <sub>51</sub>	–	(Ar. 21:13,14:-)
IV	U		43	g,z <sub>51</sub>	–	
II	U		43	g,z <sub>62</sub>	e,n,x	
I	U	Mbao	43	i	1,2	
I	U	Voulte	43	i	e,n,x	
I	U	Thetford	43	k	1,2	
I	U	Ahuza	43	k	1,5	
IIIb	U		43	k	z	(Ar. 21:29:31)
IIIb	U		43	l,v	z <sub>53</sub> :[Rz <sub>56</sub> ]	(Ar. 21:23:25:[38])
I	U	Sudan	43	l,z <sub>13</sub>	–	
II	U		43	l,z <sub>13</sub> ,z <sub>28</sub>	1,5	
IIIb	U		43	r	e,n,x,z <sub>15</sub>	(Ar. 21:24:28)
IIIb	U		43	r	z	(Ar. 21:24:31)
IIIb	U		43	r	z <sub>53</sub>	(Ar. 21:24:25)
I	U	Farcha	43	y	1,2	
I	U	Kingabwa	43	y	1,5	
I	U	Ogbete	43	z	1,5	
II	U		43	z	1,5	
I	U	Arusha	43	z	e,n,z <sub>15</sub>	
I	U	Adana	43	z <sub>10</sub>	1,5	
I	U	Makiling	43	z <sub>29</sub>	–	
IV	U		43	z <sub>29</sub>	–	
II	U		43	z <sub>29</sub>	e,n,x	
II	U		43	z <sub>29</sub>	z <sub>42</sub>	
I	U	Ahepe	43	z <sub>35</sub>	1,6	
IIIa	U		43	z <sub>36</sub>	–	(Ar. 21:17,20:-)
IV	U	Volksdorf	43	z <sub>36</sub> ,z <sub>38</sub>	–	
I	U	Irigny	43	z <sub>38</sub>	–	
IIIa	U		43	z <sub>4</sub> ,z <sub>23</sub>	–	(Ar. 21:1,2,5:- and 21:1,2,6:-)
IV	U	Houten	43	z <sub>4</sub> ,z <sub>23</sub>	–	
IV	U		43	z <sub>4</sub> ,z <sub>23</sub>	–	
IIIa	U		43	z <sub>4</sub> ,z <sub>24</sub>	–	(Ar. 21:1,3,11:-)
IV	U		43	z <sub>4</sub> ,z <sub>24</sub>	–	
IV	U	Tuindorp	43	z <sub>4</sub> ,z <sub>32</sub>	–	
II	U	Bunnik	43	z <sub>42</sub>	[1,5,7]	
IIIb	U		43	z <sub>52</sub>	z <sub>53</sub>	(Ar. 21:26:25)
I	V	Niakhar	44	a	1,5	
I	V	Tiergarten	44	a	e,n,x	
I	V	Niarembe	44	a	1,w	
I	V	Sedgwick	44	b	e,n,z <sub>15</sub>	
I	V	Madigan	44	c	1,5	
I	V	Quebec	44	c	e,n,z <sub>15</sub>	
I	V	Bobo	44	d	1,5	
I	V	Kermel	44	d	e,n,x	
I	V	Fischerstrasse	44	d	e,n,z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	V	Palamaner	<u>1</u> ,44	d	Z <sub>35</sub>	
II	V		<u>1</u> ,44	e,n,x	1,6	
I	V	Vleuten	44	f,g	–	
I	V	Gamaba	<u>1</u> ,44	g,m,[s]	–	
I	V	Splott	44	g,s,t	–	
II	V		44	g,t	Z <sub>42</sub>	
I	V	Carswell	44	g,Z <sub>51</sub>	–	
IV	V		44	g,Z <sub>51</sub>	–	
I	V	Maritzburg	<u>1</u> ,44	i	e,n,Z <sub>15</sub>	
I	V	Lawra	44	k	e,n,Z <sub>15</sub>	
I	V	Malika	44	l,Z <sub>28</sub>	1,5	
I	V	Muguga	44	m,t	–	
V	V	Camdeni	44	r	–	
I	V	Brefet	44	r	e,n,Z <sub>15</sub>	
I	V	Bolama	44	z	e,n,x	
I	V	Uhlenhorst	44	z	1,w	
I	V	Guinea	<u>1</u> ,44	Z <sub>10</sub>	[1,7]	
I	V	Llobregat	<u>1</u> ,44	Z <sub>10</sub>	e,n,x	
I	V	Zinder	44	Z <sub>29</sub>	–	
IV	V		44	Z <sub>29</sub>	–	
II	V		44	Z <sub>29</sub>	e,n,x:Z <sub>42</sub>	
IV	V		44	Z <sub>36</sub> ,[Z <sub>38</sub> ]	–	
I	V	Koketime	44	Z <sub>38</sub>	–	
V	V		44	Z <sub>39</sub>	–	
II	V	Clovelly	<u>1</u> ,44	Z <sub>39</sub>	[e,n,x,Z <sub>15</sub> ]	
I	V	Kua	44	Z <sub>4</sub> ,Z <sub>23</sub>	–	
II	V		44	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	V		44	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 1,3:1,2,5:- and 1,3:1,2,6:-)
IV	V		44	Z <sub>4</sub> ,Z <sub>23</sub>	–	
I	V	Ploufragan	1,44	Z <sub>4</sub> ,Z <sub>23</sub>	e,n,Z <sub>15</sub>	
IIIa	V		44	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 1,3:1,6,7,9:-)
I	V	Christiansborg	44	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIa	V		44	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 1,3:1,3,11:-)
IV	V		44	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIa	V		44	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 1,3:1,2,10:- and 1,3:1,7,8:-). IP calls Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub> , Ar. 1,2,10.
IV	V	Lohbruegge	44	Z <sub>4</sub> ,Z <sub>32</sub>	–	
VI	W	Vrindaban	45	a	e,n,x	
I	W	Meekatharra	45	a	e,n,Z <sub>15</sub>	
II	W	Ejeda	45	a	Z <sub>10</sub>	
I	W	Riverside	45	b	1,5	
I	W	Fomeco	45	b	e,n,Z <sub>15</sub>	
I	W	Deversoir	45	c	e,n,x	
I	W	Dugbe	45	d	1,6	
I	W	Karachi	45	d	e,n,x	
I	W	Warmesen	45	d	e,n,Z <sub>15</sub>	
I	W	Suellendorf	45	f,g	–	
II	W	Windhoek	45	g,m,s,t	1,5	
II	W	Bremen	45	g,m,s,t	e,n,x	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	W	Tornow	45	g,m,[s],[t]	–	
II	W	Perinet	45	g,m,t	e,n,x,Z <sub>15</sub>	
I	W	Binningen	45	g,s,t	–	
IIIa	W		45	g,Z <sub>51</sub>	–	(Ar. 11:13,14:-)
IV	W		45	g,Z <sub>51</sub>	–	
I	W	Verviers	45	k	1,5	
I	W	Casablanca	45	k	1,7	
I	W	Cairns	45	k	e,n,Z <sub>15</sub>	
I	W	Imo	45	l,v	[e,n,Z <sub>15</sub> ]	
I	W	Apapa	45	m,t	–	
II	W		45	m,t	1,5	
I	W	Kofandoka	45	r	e,n,Z <sub>15</sub>	
II	W		45	z	1,5	
I	W	Yopougon	45	z	e,n,Z <sub>15</sub>	
II	W	Klapmuts	45	z	Z <sub>39</sub>	
I	W	Jodhpur	45	Z <sub>29</sub>	–	
IIIa	W		45	Z <sub>29</sub>	–	(Ar. 11:16,18:-)
II	W		45	Z <sub>29</sub>	1,5	
II	W		45	Z <sub>29</sub>	e,n,x	
II	W		45	Z <sub>29</sub>	Z <sub>42</sub>	
I	W	Lattenkamp	45	Z <sub>35</sub>	1,5	
I	W	Balcones	45	Z <sub>36</sub>	–	
IV	W		45	Z <sub>36</sub> ,Z <sub>38</sub>	–	
IIIa	W		45	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 11:1,2,5:-)
IV	W		45	Z <sub>4</sub> ,Z <sub>23</sub>	–	
I	W	Transvaal	45	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIa	W		45	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 11:1,3,11:-)
IIIa	W		45	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 11:1,7,8:-)
II	X	Bilthoven	47	a	[1,5]	
II	X		47	a	e,n,x,Z <sub>15</sub>	
I	X	Saka	47	b	–	IP combined Saka with Sya (47:b:z <sub>6</sub> ) and called it Sya.
I	X	Wenatchee	47	b	1,2	
II	X	Phoenix	47	b	1,5	
II	X	Khami	47	b	[e,n,x,Z <sub>15</sub> ]	
I	X	Sya	47	b	z <sub>6</sub>	
II	X		47	b	z <sub>6</sub>	
IIIb	X		47	c	1,5,7	(Ar. 28:32:30)
I	X	Kodjovi	47	c	[1,6]	Kodjovi may possess H phase Rz <sub>78</sub> .
IIIb	X		47	c	e,n,x,Z <sub>15</sub> : [Z <sub>57</sub> ]	(Ar. 23:32:28 and 28:32:28:[40])
IIIb	X		47	c	z	(Ar. 28:32:31)
IIIb	X		47	c	Z <sub>35</sub>	(Ar. 28:32:21)
II	X		47	d	e,n,x,Z <sub>15</sub>	
I	X	Stellingen	47	d	[e,n,x]	
II	X	Quimbamba	47	d	Z <sub>39</sub>	
II	X		47	e,n,x,Z <sub>15</sub>	1,6	
I	X	Sljeme	1,47	f,g	–	
I	X	Anie	47	(g),m,t	–	IP combined Anie with Mesbit (47:m,t:[e,n,Z <sub>15</sub> ]) and called it Mesbit.

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	X	Luke	1,47	g,m	–	
II	X		47	g,t	e,n,x	
IIIa	X		47	g,z <sub>51</sub>	–	(Ar. 28:13,14)
IIIb	X		47	i	e,n,x,z <sub>15</sub>	(Ar. 23:33:28). May possess H phase RZ <sub>50</sub> (Ar. 42).
I	X	Bergen	47	i	e,n,z <sub>15</sub>	
IIIb	X		47	i	z	(Ar. 28:33:31)
IIIb	X		47	i	z <sub>35</sub>	(Ar. 23:33:21 and 28:33:21)
IIIb	X		47	i	z <sub>53</sub> :[z <sub>57</sub> ]	(Ar. 23:33:25 and 28:33:25:[40])
I	X	Staoueli	47	k	1,2	
I	X	Bootle	47	k	1,5	
IIIb	X		47	k	1,5,7	(Ar. 28:29:30)
I	X	Dahomey	47	k	1,6	Dahomey may possess H phase RZ <sub>58</sub> .
IIIb	X		47	k	e,n,x,z <sub>15</sub>	(Ar. 28:29:28)
I	X	Lyon	47	k	e,n,z <sub>15</sub>	
IIIb	X		47	k	z	(Ar. 28:29:31)
IIIb	X		47	k	z <sub>35</sub>	(Ar. 23:29:21)
IIIb	X		47	k	z <sub>53</sub>	(Ar. 23:29:25)
IV	X		47	l,v	–	
IIIb	X		47	l,v	1,5,(7)	(Ar. 23:23:30). May possess H phase RZ <sub>50</sub> (Ar. 42).
IIIa	X		47	l,v	e,n,x,z <sub>15</sub>	(Ar. 28:23:28)
IIIb	X		47	l,v	z	(Ar. 23:23:31)
IIIb	X		47	l,v	z <sub>35</sub>	(Ar. 28:23:21)
IIIb	X		47	l,v	z <sub>53</sub>	(Ar. 28:23:25)
IIIb	X		47	l,v	z <sub>57</sub>	(Ar. 28:23:40)
I	X	Teshie	1,47	l,z <sub>13</sub> ,z <sub>28</sub>	e,n,z <sub>15</sub>	
I	X	Mesbit	47	m,t	[e,n,z <sub>15</sub> ]	
IIIa	X		47	r	–	(Ar. 23:24:-). CDC does not have this.
I	X	Dapango	47	r	1,2	
IIIb	X		47	r	1,5,7	(Ar. 23:24:30)
IIIb	X		47	r	z	(Ar. 23:24:31)
IIIb	X		47	r	z <sub>35</sub>	(Ar. 23:24:21 and 28:24:21)
IIIb	X		47	r	z <sub>53</sub>	(Ar. 23:24:25). May possess H phase RZ <sub>74</sub> (Ar. RZ <sub>74</sub> )
IIIb	X		47	r	z <sub>53</sub> :RZ <sub>50</sub> :z <sub>60</sub>	(Ar. 28:24:25:42:44). Not in IP book.
IIIb	X		47	r	z <sub>53</sub> :[z <sub>60</sub> ]	(Ar. 23:24:25:[44]). May possess H phase RZ <sub>70</sub> and RZ <sub>72</sub> (Ar. RZ <sub>70</sub> or RZ <sub>72</sub> ).
I	X	Moualine	47	y	1,6	
I	X	Blitta	47	y	e,n,x	
I	X	Mountpleasant	47	z	1,5	
I	X	Kaolack	47	z	1,6	
II	X		47	z	e,n,x,z <sub>15</sub>	
II	X	Chersina	47	z	z <sub>6</sub>	
IIIb	X		47	z <sub>10</sub>	1,5,7	(Ar. 28:27:30)
IIIb	X		47	z <sub>10</sub>	z	(Ar. 28:27:31)
IIIb	X		47	z <sub>10</sub>	z <sub>35</sub>	(Ar. 28:27:21)
I	X	Ekpoui	47	z <sub>29</sub>	–	
IIIa	X		47	z <sub>29</sub>	–	(Ar. 28:16,18:-)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	X		47	Z <sub>29</sub>	e,n,x,Z <sub>15</sub>	
I	X	Bingerville	47	Z <sub>35</sub>	e,n,Z <sub>15</sub>	
IV	X		47	Z <sub>36</sub>	–	
I	X	Alexanderplatz	47	Z <sub>38</sub>	–	
IIIa	X		47	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 28:1,2,5:-)
I	X	Tabligbo	47	Z <sub>4</sub> ,Z <sub>23</sub>	[e,n,Z <sub>15</sub> ]	
I	X	Binche	47	Z <sub>4</sub> ,Z <sub>23</sub>	l,w	
I	X	Bere	47	Z <sub>4</sub> ,Z <sub>23</sub>	[Z <sub>6</sub> ]	Bere may possess H phase Rz45.
I	X	Tamberma	47	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	X	Quinhon	47	Z <sub>44</sub>	–	
IIIb	X		47	Z <sub>52</sub>	1,5,7	(Ar. 28:26:30)
IIIb	X		47	Z <sub>52</sub>	e,n,x,Z <sub>15</sub>	(Ar. 28:26:28)
IIIb	X		47	Z <sub>52</sub>	z	(Ar. 28:26:31)
IIIb	X		47	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 28:26:21)
II	X		47	Z <sub>6</sub>	1,6	
I	Y	Hisingen	48	a	1,5,7	
IIIb	Y		48	a	Z <sub>35</sub>	(Ar. 5:35:21). Not in 1992 IP, but is in Bergey.
II	Y		48	a	Z <sub>39</sub>	
II	Y		48	a	Z <sub>6</sub>	
V	Y		48	b	–	
II	Y		48	b	[Z <sub>6</sub> ]	
IIIb	Y		48	c	z	(Ar. 29:32:31)
II	Y	Etosha	48	d	1,11	Etosha was not considered a new serotype by Kauffmann and is not used.
II	Y		48	d	1,2	
II	Y	Hagenbeck	48	d	Z <sub>6</sub>	
I	Y	Fitzroy	48	e,h	1,5	
II	Y	Hammonia	48	e,n,x,Z <sub>15</sub>	Z <sub>6</sub>	
II	Y	Erlangen	48	g,m,t	–	
IIIa	Y		48	g,Z <sub>51</sub>	–	(Ar. 5:13,14:-)
IV	Y	Marina	48	g,Z <sub>51</sub>	–	
IIIb	Y	Sydney	48	i	z	Sydney was formerly in subspecies II, but it is now combined with <i>Arizona</i> 5:33:31. The name Sydney has been dropped.
IIIb	Y		48	i	Z <sub>35</sub> : [Z <sub>57</sub> ]	(Ar. 29:33:21:[40])
IIIb	Y		48	i	Z <sub>53</sub>	(Ar. 5:33:25)
IIIb	Y		48	i	Z <sub>61</sub>	(Ar. 5,29:33:41)
IIIb	Y		48	i	z: [Z <sub>72</sub> ]	(Ar. 5,29:33:31: [Z <sub>72</sub> ]). CDC does not have Z <sub>72</sub> strain.
IIIb	Y		48	k	1,5,(7)	(Ar. 5:29:30)
II	Y		48	k	e,n,x,Z <sub>15</sub>	
IIIb	Y		48	k	e,n,x,Z <sub>15</sub>	(Ar. 5:29:28)
I	Y	Dahlem	48	k	e,n,Z <sub>15</sub>	
IIIb	Y		48	k	z	(Ar. 5,29:29:31)
IIIb	Y		48	k	Z <sub>35</sub> : [RZ <sub>75</sub> ]	(Ar. [5:29:21:RZ <sub>75</sub> ]). CDC does not have RZ <sub>75</sub> .
IIIb	Y		48	k	Z <sub>53</sub>	(Ar. 5,29:29:25)
IIIb	Y		48	(k)	Z <sub>53</sub>	(Ar. 5:22:25 and Ar. 5,29:22:25). Called 5:22:25 by IP.
II	Y	Sakaraha	48	[k]	Z <sub>39</sub>	
I	Y	Australia	48	l,v	1,5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIb	Y		48	l,v	1,5,(7)	(Ar. 5:23:30). May possess H phase Rz <sub>47</sub> or Rz <sub>50</sub> (Ar. 39 or 42).
IIIb	Y		48	l,v	z	(Ar. 5:29:23:31)
IIIb	Y		48	r	e,n,x,Z <sub>15</sub>	(Ar. 5:24:28)
IIIb	Y		48	r	z	(Ar. 5:29:24:31)
I	Y	Toucra	48	z	1,5	Toucra may possess H phase Rz <sub>58</sub> .
II	Y		48	z	1,5	
VI	Y		48	Z <sub>10</sub>	1,5	
II	Y	Ngozi	48	Z <sub>10</sub>	[1,5]	
I	Y	Isaszeg	48	Z <sub>10</sub>	e,n,x	
IIIb	Y		48	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 5:27:28)
IIIb	Y		48	Z <sub>10</sub>	z	(Ar. 5:29:27:31)
II	Y		48	Z <sub>29</sub>	-	
IIIa	Y		48	Z <sub>29</sub>	-	(Ar. 5:16,18). This is not in IP book, but is on Rohde's list.
IV	Y		48	Z <sub>29</sub>	-	
V	Y	Bongor	48	Z <sub>35</sub>	-	
IIIb	Y		48	Z <sub>35</sub>	Z <sub>52</sub>	(Ar. 5:21:26)
IIIa	Y		48	Z <sub>36</sub>	-	(Ar. 5:29:17,20:-)
IV	Y		48	Z <sub>36</sub> ,[Z <sub>38</sub> ]	-	
V	Y		48	Z <sub>39</sub>	-	
IIIa	Y		48	Z <sub>4</sub> ,Z <sub>23</sub>	-	(Ar. 5:1,2,5:-; 5:1,2,5,6:-; and 5:1,6:-)
IV	Y		48	Z <sub>4</sub> ,Z <sub>23</sub>	-	
IIIa	Y		48	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	-	(Ar. 5:1,6,7,9:-). IP calls this 5:1,6,7:-.
I	Y	Djakarta	48	Z <sub>4</sub> ,Z <sub>24</sub>	-	
IIIa	Y		48	Z <sub>4</sub> ,Z <sub>24</sub>	-	(Ar. 5:1,3,11:-)
IIIa	Y		48	Z <sub>4</sub> ,Z <sub>32</sub>	-	(Ar. 5:1,6,7:-; 5:1,7,8:-; and Ar. 5:1,2,10:-). IP calls Z <sub>4</sub> ,Z <sub>32</sub> , Ar. 1,7,8; and would call Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub> , Ar. 1,2,10.
IV	Y		48	Z <sub>4</sub> ,Z <sub>32</sub>	-	
V	Y	Balboa	48	Z <sub>41</sub>	-	
IIIb	Y		48	Z <sub>52</sub>	e,n,x,Z <sub>15</sub>	(Ar. 29:26:28)
IIIb	Y		48	Z <sub>52</sub>	z	(Ar. 5:26:31)
V	Y		48	Z <sub>65</sub>	-	
V	Y		48	Z <sub>81</sub>	-	
IV	Z		50	a	-	
IV	Z		50	b	-	
I	Z	Rochdale	50	b	e,n,x	
II	Z		50	b	Z <sub>6</sub>	
IV	Z		50	d	-	
I	Z	Hemingford	50	d	1,5	Hemingford may possess H phase Rz <sub>82</sub> .
II	Z	Krugersdorp	50	e,n,x	1,7	
II	Z	Namib	50	g,[m],s,t	[1,5]	
IV	Z	Wassenaar	50	g,Z <sub>51</sub>	-	
II	Z		50	g,Z <sub>62</sub>	e,n,x	
IIIb	Z		50	i	1,5,7	(Ar. 9a,9c:33:30)
IIIb	Z		50	i	e,n,x,Z <sub>15</sub>	(Ar. 9a,9c:33:28)
IIIb	Z		50	i	z	(Ar. 9a,9c:33:31)
IIIb	Z		50	k	1,5,7	(Ar. 9a,9c:29:30)
IIIb	Z		50	k	e,n,x,Z <sub>15</sub>	(Ar. 9a,9c:29:28)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	Z		50	k	e,n,x:Z <sub>42</sub>	
IIIb	Z		50	k	z	(Ar. 9a,9b:29:31 and 9a,9c:29:31). Ar. 9a,9b may possess H phase Rz <sub>50</sub> (Ar. 42).
IIIb	Z		50	k	Z <sub>35</sub>	(Ar. 9a,9b:29:21)
IIIb	Z		50	k	Z <sub>53</sub>	(Ar. 9a,9b:29:25 and 9a,9c:29:25). IP and Rohde only list the 9a,9c.
II	Z	Seaforth	50	k	Z <sub>6</sub>	
IIIb	Z		50	(k)	z	(Ar. 9a,9b:22:31)
IIIb	Z		50	(k)	Z <sub>35</sub>	(Ar. 9a,9b:22:21)
I	Z	Fass	50	l,v	1,2	
IIIb	Z		50	l,v	e,n,x,Z <sub>15</sub>	(Ar. 9a,9b:23:28)
IIIb	Z		50	l,v	z	(Ar. 9a,9b:23:31 and 9a,9c:23:31). IP only lists 9a,9c.
IIIb	Z		50	l,v	Z <sub>35</sub>	(Ar. 9a,9c:23:21)
II	Z		50	l,w	e,n,x,Z <sub>15</sub> :Z <sub>42</sub>	
II	Z		50	l,Z <sub>28</sub>	Z <sub>42</sub>	
II	Z	Atra	50	m,t	Z <sub>6</sub> :Z <sub>42</sub>	
IIIb	Z		50	r	1,5,(7)	(Ar. 9a,9b:24:30)
IIIb	Z		50	r	e,n,x,Z <sub>15</sub>	(Ar. 9a,9c:24:28)
IIIb	Z		50	r	z	(Ar. 9a,9b:24:31 and 9a,9c:24:31).
IIIb	Z		50	r	Z <sub>35</sub>	(Ar. 9a,9b:24:21). May possess H phase Rz <sub>58</sub> (Ar. Rz58). This is not in IP book, but is on Rohde's list.
IIIb	Z		50	r	Z <sub>53</sub>	(Ar. 9a,9b:24:25). May possess H phase Rz <sub>50</sub> (Ar. 42). This is not in IP book, but is on Rohde's list.
I	Z	Dougi	50	y	1,6	
II	Z	Greenside	50	z	e,n,x	
IIIb	Z		50	Z <sub>10</sub>	z	(Ar. 9a,9c:27:31). May possess H phase Rz <sub>56</sub> (Ar. 38).
IIIb	Z		50	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 9a,9c:27:25)
II	Z	Hooggraven	50	Z <sub>10</sub>	Z <sub>6</sub> :Z <sub>42</sub>	
I	Z	Ivorycoast	50	Z <sub>29</sub>	–	
IIIa	Z		50	Z <sub>29</sub>	–	(Ar. 9a,9b:16,18:-)
IIIa	Z		50	Z <sub>36</sub>	–	(Ar. 9a,9b:17,20:-)
IIIa	Z		50	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 9a,9b:1,2,5:- and 9a,9b:1,2,6:-)
IV	Z	Flint	50	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	Z		50	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 9a,9b:1,2,10:-). Called 9a,9b:1,6,7:- by IP and Rohde.
IIIa	Z		50	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 9a,9b:1,3,11:-)
IV	Z		50	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIa	Z		50	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 9a,9b:1,2,10; 9a,9b:1,6,7:-; and 9a,9b:1,7,8:-). 9a,9b:1,2,10:- and 9a,9b:1,7,8:- used by IP and Rohde.
IV	Z	Bonaire	50	Z <sub>4</sub> ,Z <sub>32</sub>	–	
II	Z	Faure	50	Z <sub>42</sub>	1,7	
IIIb	Z		50	Z <sub>52</sub>	1,5,7	(Ar. 9a,9b:26:30 and 9a,9c:26:30)
IIIb	Z		50	Z <sub>52</sub>	z	(Ar. 9a,9b:26:31 and 9a,9c:26:31)
IIIb	Z		50	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 9a,9b:26:21 and 9a,9c:26:21)
IIIb	Z		50	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 9a,9b:26:25 and 9a,9c:26:25)
II	51	Roggeveld	51	–	1,7	
I	51	Tione	51	a	e,n,x	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IV	51		51	b	–	
I	51	Karaya	51	b	1,5	
II	51		51	c	–	
I	51	Gokul	1,51	d	[1,5]	
I	51	Meskin	51	e,h	1,2	
II	51		51	g,s,t	e,n,x	
IIIa	51		51	g,z <sub>51</sub>	–	(Ar. 1,2:13,14:-)
I	51	Kabete	51	i	1,5	
I	51	Dan	51	k	e,n,z <sub>15</sub>	
IIIb	51		51	k	z <sub>35</sub>	(Ar. 1,2:29:21)
I	51	Harcourt	51	l,v	1,2	
I	51	Overschie	51	l,v	1,5	
I	51	Dadzie	51	l,v	e,n,x	
IIIb	51		51	l,v	z	(Ar. 1,2:23:31)
I	51	Moundou	51	l,z <sub>28</sub>	1,5	
II	51	Askraal	51	l,z <sub>28</sub>	[z <sub>6</sub> ]	
I	51	Lutetia	51	r,i	l,z <sub>13</sub> ,z <sub>28</sub>	
I	51	Antsalova	51	z	1,5	
I	51	Treforest	1,51	z	1,6	
I	51	Lechler	51	z	e,n,z <sub>15</sub>	
I	51	Bergues	51	z <sub>10</sub>	1,5	
II	51		51	z <sub>29</sub>	e,n,x,z <sub>15</sub>	
IIIa	51		51	z <sub>4</sub> ,z <sub>23</sub>	–	(Ar. 1,2:1,2,5:- and 1,2:1,2,6:-)
IV	51	Harmelen	51	z <sub>4</sub> ,z <sub>23</sub>	–	
IIIa	51		51	z <sub>4</sub> ,z <sub>24</sub>	–	(Ar. 1,2:1,3,11:-)
IIIa	51		51	z <sub>4</sub> ,z <sub>32</sub>	–	(Ar. 1,2:1,7,8:-)
I	52	Uithof	52	a	1,5	
I	52	Ord	52	a	e,n,z <sub>15</sub>	
I	52	Molesey	52	b	1,5	
I	51	Flottbek	52	b	[e,n,x]	
II	52		52	c	k	
IIIb	52		52	c	k	(Ar. 31:32:29). This is not in IP book, but is on Rohde's list.
I	52	Utrecht	52	d	1,5	
II	52		52	d	e,n,x,z <sub>15</sub>	
II	52		52	d	z <sub>39</sub>	CDC does not have this.
I	52	Butare	52	e,h	1,6	
I	52	Derkle	52	e,h	1,7	
I	52	Saintemarie	52	g,t	–	
II	52		52	g,t	–	
I	52	Bordeaux	52	k	1,5	
IIIb	52		52	k	z <sub>35</sub>	(Ar. 31:29:21)
IIIb	52		52	k	z <sub>53</sub>	(Ar. 31:29:25)
IIIb	52		52	(k)	z <sub>35</sub>	(Ar. 31:22:21)
IIIb	52		52	l,v	z <sub>53</sub>	(Ar. 31:23:25)
II	52		52	z	z <sub>39</sub>	
II	52	Wilhemstrasse	52	z <sub>44</sub>	1,5	IP combined Wilhemstrasse with Lobatsi (52:z44:1,5,7). The name Wilhemstrasse has been dropped.

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	52	Lobatsi	52	Z <sub>44</sub>	1,5,7	
IIIb	52		52	Z <sub>52</sub>	z	(Ar. 31:26:31)
II	53		53	c	1,5	
II	53		53	d	1,5	
II	53		1,53	d	Z <sub>39</sub>	
II	53		53	d	Z <sub>42</sub>	
IIIa	53		53	g,Z <sub>51</sub>	–	(Ar. 1,4:13,14:-)
IV	53		1,53	g,Z <sub>51</sub>	–	
IIIb	53		53	i	z	(Ar. 1,4:33:31)
IIIb	53		53	k	e,n,x,Z <sub>15</sub>	(Ar. 1,4:29:28)
IIIb	53		53	k	z	(Ar. 1,4:29:31)
IIIb	53		53	(k)	z	(Ar. 1,4:22:31)
IIIb	53		53	(k)	Z <sub>35</sub>	(Ar. 1,4:22:21)
IIIb	53		53	l,v	e,n,x,Z <sub>15</sub>	(Ar. 1,4:23:28)
IIIb	53		53	l,v	Z <sub>35</sub>	(Ar. 1,4:23:21)
II	53		53	l,Z <sub>28</sub>	e,n,x	
II	53	Midhurst	53	l,Z <sub>28</sub>	Z <sub>39</sub>	
II	53		53	l,Z <sub>28</sub>	Z <sub>6</sub>	
IIIb	53		53	r	z	(Ar. 1,4:24:31)
IIIb	53		53	r	Z <sub>35</sub>	(Ar. 1,4:24:21)
IIIb	53		53	r	Z <sub>68</sub>	(Ar. 1,4:24:47). This was formerly called Z <sub>56</sub> (Ar. 38), but was changed to Z <sub>68</sub> (Ar. 47).
II	53		53	z	1,5	
IIIb	53		53	z	1,5,(7)	(Ar. 1,4:30:31)
II	53		53	z	Z <sub>6</sub>	
IIIb	53		53	Z <sub>10</sub>	Z <sub>35</sub>	(Ar. 1,4:27:21)
IIIa	53		53	Z <sub>29</sub>	–	(Ar. 1,4:16,18:-)
IV	53	Bockenheim	1,53	Z <sub>36</sub> ,Z <sub>38</sub>	–	
IIIa	53		53	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 1,4:1,2,5:- and 1,4:1,2,6:-)
IV	53		53	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	53		53	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 1,4:1,6,7:- and 1,4:1,6,7,9:-)
II	53	Humber	53	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIa	53		53	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 1,4:1,3,11:-)
IIIa	53		53	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 1,4:1,6,7:-). IP combined this with 53:z <sub>4</sub> ,z <sub>23</sub> ,z <sub>32</sub> :- (Ar. 1,4:1,6,7,9:-).
IIIb	53		53	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 1,4:26:21)
IIIb	53		53	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 1,4:26:25)
I	54	Tonev	21,54	b	e,n,x	
I	54	Winnipeg	54	e,h	1,5	
I	54	Rossleben	54	e,h	1,6	
I	54	Borreze	54	f,g,s	–	
I	54	Uccle	3,54	g,s,t	–	
I	54	Poeseldorf	8,20,54	i	Z <sub>6</sub>	
I	54	Ochsenwerder	6,7,54	k	1,5	
I	54	Newholland	4,12,54	m,t	–	
I	54	Czernyring	54	r	1,5	
I	54	Steinwerder	3,15,54	y	1,5	
I	54	Canton	54	Z <sub>10</sub>	e,n,x	
I	54	Barry	54	Z <sub>10</sub>	e,n,Z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	54	Yerba	54	Z <sub>4</sub> ,Z <sub>23</sub>	–	
II	55	Tranoroa	55	k	Z <sub>39</sub>	
II	56	Artis	56	b	–	
II	56		56	d	–	
II	56		56	e,n,x	1,7	
II	56		56	l,v	Z <sub>39</sub>	
II	56		56	l,Z <sub>28</sub>	–	
II	56		56	z	Z <sub>6</sub>	
II	56		56	Z <sub>10</sub>	e,n,x	
IIIa	56		56	Z <sub>29</sub>	–	(Ar. 14:16,18:-)
IIIa	56		56	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 14:1,2,5:- and 14:1,2,6:-)
IIIa	56		56	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	-	(Ar. 14:1,6,7,9:-)
II	57		57	a	Z <sub>42</sub>	
I	57	Antonio	57	a	Z <sub>6</sub>	
I	57	Maryland	57	b	1,7	
I	57	Batonrouge	57	b	e,n,Z <sub>15</sub>	
IIIb	57		57	c	z:[Z <sub>60</sub> ]	(Ar. 34:32:31:[44])
II	57		57	d	1,5	
II	57		57	g,[m],s,t	Z <sub>42</sub>	
II	57		57	g,t	-	
IIIb	57		57	i	e,n,x,Z <sub>15</sub>	(Ar. 34:33:28)
IIIb	57		57	i	z	(Ar. 34:33:31)
IIIb	57		57	k	e,n,x,Z <sub>15</sub>	(Ar. 34:29:28). CDC does not have this and not on Rohde's list.
IIIb	57		57	Z <sub>10</sub>	z	(Ar. 34:27:31)
II	57	Locarno	57	Z <sub>29</sub>	Z <sub>42</sub>	
II	57	Manombo	57	Z <sub>39</sub>	e,n,x,Z <sub>15</sub>	
IV	57		57	Z <sub>4</sub> ,Z <sub>23</sub>	-	
II	57	Tokai	57	Z <sub>42</sub>	1,6:Z <sub>53</sub>	
II	58		58	a	Z <sub>6</sub>	
II	58		58	b	1,5	
II	58		58	c	Z <sub>6</sub>	
II	58		58	d	Z <sub>6</sub>	
IIIb	58		58	i	e,n,x,Z <sub>15</sub>	(Ar. 1,33:33:28)
IIIb	58		58	k	z	(Ar. 1,33:29:31)
IIIb	58		58	l,v	e,n,x,Z <sub>15</sub>	(Ar. 1,33:23:28)
IIIb	58		58	l,v	Z <sub>35</sub>	(Ar. 1,33:23:21)
II	58	Basel	58	l,Z <sub>13</sub> ,Z <sub>28</sub>	1,5	
II	58		58	l,Z <sub>13</sub> ,Z <sub>28</sub>	Z <sub>6</sub>	
IIIb	58		58	r	e,n,x,Z <sub>15</sub>	(Ar. 1,33:24:28)
IIIb	58		58	r	z	(Ar. 1,33:24:31)
IIIb	58		58	r	Z <sub>53</sub>	(Ar. 1,33:24:25). May possess H phase RZ <sub>47</sub> (Ar. 39) or RZ <sub>57</sub> (Ar. 40) or RZ <sub>70</sub> (Ar. RZ <sub>70</sub> ).
II	58		58	Z <sub>10</sub>	1,6	
IIIb	58		58	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 1,33:27:28)
IIIb	58		58	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 1,33:27:25). May possess H phase RZ <sub>50</sub> (Ar. 42).
II	58		58	Z <sub>10</sub>	Z <sub>6</sub>	
II	58		58	Z <sub>39</sub>	e,n,x,Z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIb	58		58	Z <sub>52</sub>	z	(Ar. 1,33:26:31)
IIIb	58		58	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 1,33:26:21)
II	58		58	Z <sub>6</sub>	1,6	
IIIb	59		59	c	e,n,x,Z <sub>15</sub>	(Ar. 19:32:28)
IIIb	59		59	i	e,n,x,Z <sub>15</sub>	(Ar. 19:33:28)
IIIb	59		59	i	z	(Ar. 19:33:31)
IIIb	59		59	i	Z <sub>35</sub>	(Ar. 19:33:21)
II	59	Betioky	59	k	(z)	
IIIb	59		59	k	Z <sub>53</sub>	(Ar. 19:29:25)
IIIb	59		59	(k)	e,n,x,Z <sub>15</sub>	(Ar. 19:22:28)
IIIb	59		59	(k)	z	(Ar. 19:22:31)
IIIb	59		59	(k)	Z <sub>35</sub>	(Ar. 19:22:21)
IIIb	59		59	l,v	z	(Ar. 19:23:31)
IIIb	59		59	l,v	Z <sub>53</sub>	(Ar. 19:23:25)
IIIb	59		59	r	Z <sub>35</sub>	(Ar. 19:24:21)
II	59		1,59	z	Z <sub>6</sub>	
IIIb	59		59	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 19:27:25)
IIIb	59		59	Z <sub>10</sub>	Z <sub>57</sub>	(Ar. 19:27:40)
IIIa	59		59	Z <sub>29</sub>	–	(Ar. 19:16,18:-)
IIIa	59		59	Z <sub>36</sub>	–	(Ar. 19:17,20:-)
IIIa	59		59	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 19:1,2,5:- and 19:1,2,6:-)
IIIb	59		59	Z <sub>52</sub>	[Z <sub>53</sub> ]	(Ar. 19:26:[25])
II	60		60	b	[1,16]	
II	60	Setubal	60	g,m,t	Z <sub>6</sub>	
IIIb	60		60	i	–	(Ar. 24:33:-). May possess H phase Rz <sub>50</sub> (Ar. 42).
IIIb	60		60	i	e,n,x,Z <sub>15</sub>	(Ar. 24:33:28)
IIIb	60		60	i	Z <sub>35</sub>	(Ar. 24:33:21)
IIIb	60		60	k	z	(Ar. 24:29:31)
IIIb	60		60	k	Z <sub>35</sub>	(Ar. 24:29:21)
IIIb	60		60	(k)	Z <sub>53</sub>	(Ar. 24:22:25)
IIIb	60		60	l,v	z	(Ar. 24:23:31)
IIIb	60		60	r	e,n,x,Z <sub>15</sub>	(Ar. 24:24:28)
IIIb	60		60	r	z	(Ar. 24:24:31)
IIIb	60		60	r	Z <sub>35</sub>	(Ar. 24:24:21)
IIIb	60		60	r	Z <sub>53</sub>	(Ar. 24:24:25)
II	60	Luton	60	z	e,n,x	
IIIb	60		60	Z <sub>10</sub>	z	(Ar. 24:27:31)
IIIb	60		60	Z <sub>10</sub>	Z <sub>35</sub>	(Ar. 24:27:21)
IIIb	60		60	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 24:27:25)
II	60		60	Z <sub>29</sub>	e,n,x	
V	60		60	Z <sub>41</sub>	–	
IIIb	60		60	Z <sub>52</sub>	1,5,[7]	(Ar. 24:26:30)
IIIb	60		60	Z <sub>52</sub>	z	(Ar. 24:26:31)
IIIb	60		60	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 24:26:21)
IIIb	60		60	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 24:26:25)
IIIb	61		61	c	1,5,(7)	(Ar. 26:32:30)
IIIb	61		61	c	Z <sub>35</sub>	(Ar. 26:32:21)
IIIb	61		61	i	e,n,x,Z <sub>15</sub>	(Ar. 26:33:28)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIb	61	Eilbeck	61	i	z	Eilbeck was formerly in subspecies II, but is now combined with <i>Arizona</i> 26:33:31. The name Eilbeck has been dropped.
IIIb	61		61	i	z	(Ar. 26:33:31)
IIIb	61		61	i	Z <sub>35</sub>	(Ar. 26:33:21)
IIIb	61		61	i	Z <sub>53</sub>	(Ar. 26:33:25)
IIIb	61		61	k	1,5,(7)	(Ar. 26:29:30)
IIIb	61		61	k	Z <sub>35</sub>	(Ar. 26:29:21). CDC does not have this.
IIIb	61		61	(k)	Z <sub>53</sub>	(Ar. 26:22:25)
IIIb	61		61	l,v	1,5,7:[Z <sub>57</sub> ]	(Ar. 26:23:30:[40])
IIIb	61		61	l,v	z	(Ar. 26:23:31)
IIIb	61		61	l,v	Z <sub>35</sub>	(Ar. 26:23:21)
IIIb	61		61	r	1,5,7	(Ar. 26:24:30)
IIIb	61		61	r	z	(Ar. 26:24:31)
IIIb	61		61	r	Z <sub>35</sub>	(Ar. 26:24:21)
IIIb	61		61	r	Z <sub>53</sub>	(Ar. 26:24:25). May possess H phase RZ <sub>47</sub> (Ar. 39).
IIIb	61		61	Z <sub>10</sub>	Z <sub>35</sub>	(Ar. 26:27:21)
V	61		61	Z <sub>35</sub>	-	
IIIb	61		61	Z <sub>52</sub>	1,5,7	(Ar. 26:26:30)
IIIb	61		61	Z <sub>52</sub>	z	(Ar. 26:26:31)
IIIb	61		61	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 26:26:21)
IIIb	61		61	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 26:26:25)
IIIa	62		62	g,Z <sub>51</sub>	-	(Ar. 6:13,14:-)
IIIa	62		62	Z <sub>29</sub>	-	(Ar. 6:17,18:-)
IIIa	62		62	Z <sub>36</sub>	-	(Ar. 6:17,20:-)
IIIa	62		62	Z <sub>4</sub> ,Z <sub>23</sub>	-	(Ar. 6:1,2,5:-)
IIIa	62		62	Z <sub>4</sub> ,Z <sub>32</sub>	-	(Ar. 6:1,7,8:-)
IIIa	63		63	g,Z <sub>51</sub>	-	(Ar. 8:13,14:-)
IIIa	63		63	Z <sub>36</sub>	-	(Ar. 8:17,20:-)
IIIa	63		63	Z <sub>4</sub> ,Z <sub>23</sub>	-	(Ar. 8:1,2,5:-)
IIIa	63		63	Z <sub>4</sub> ,Z <sub>32</sub>	-	(Ar. 8:1,7,8:-)
II	65		65	-	1,6	
IIIb	65		65	c	1,5,7	(Ar. 30:32:30)
IIIb	65		65	c	z	(Ar. 30:32:31)
IIIb	65		65	c	Z <sub>53</sub>	(Ar. 30:32:25)
II	65		65	g,t	-	
IIIb	65		65	i	e,n,x,Z <sub>15</sub>	(Ar. 30:33:28)
IIIb	65		65	(k)	z	(Ar. 30:22:31)
IIIb	65		65	(k)	Z <sub>35</sub>	(Ar. 30:22:21)
IIIb	65		65	(k)	Z <sub>53</sub>	(Ar. 30:22:25)
IIIb	65		65	l,v	e,n,x,Z <sub>15</sub>	(Ar. 30:23:28)
IIIb	65		65	l,v	z	(Ar. 30:23:31)
IIIb	65		65	l,v	Z <sub>35</sub>	(Ar. 30:23:21)
IIIb	65		65	l,v	Z <sub>53</sub>	(Ar. 30:23:25)
IIIb	65		65	r	Z <sub>35</sub>	(Ar. 30:24:21)
IIIb	65		65	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 30:27:28)
IIIb	65		65	Z <sub>10</sub>	z	(Ar. 30:27:31)
IIIb	65		65	Z <sub>52</sub>	e,n,x,Z <sub>15</sub>	(Ar. 30:26:28)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIb	65		65	Z <sub>52</sub>	z	(Ar. 30:26:31)
IIIb	65		65	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 30:26:21)
IIIb	65		65	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 30:26:25)
V	66	Maregrosso	66	Z <sub>35</sub>	–	
V	66		66	Z <sub>39</sub>	–	
V	66	Brookfield	66	Z <sub>41</sub>	–	
V	66	Malawi	66	Z <sub>65</sub>	–	
V	66		66	Z <sub>81</sub>	–	
I	67	Crossness	67	r	1,2	

## Appendix B

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	K	Aarhus	18	Z <sub>4</sub> ,Z <sub>23</sub>	Z <sub>64</sub>	
I	C <sub>2</sub>	Aba	6,8	i	e,n,Z <sub>15</sub>	
I	M	Abadina	28	g,m	[e,n,Z <sub>15</sub> ]	
I	F	Abaetetuba	11	k	1,5	
I	F	Aberdeen	11	i	1,2	
I	Q	Abidjan	39	b	l,w	
I	I	Ablogame	16	l,Z <sub>13</sub> ,Z <sub>28</sub>	Z <sub>6</sub>	
I	I	Abobo	16	i	Z <sub>6</sub>	
I	B	Abony	1,4,5,12	b	e,n,x	IP combined Sladun (1,4,12,27:b:e,n,x) with Abony to form Abony 1,4,[5],12,27:b:e,n,x (gelatin neg.).
I	B	Abortusbovis	1,4,12,27	b	e,n,x	Gelatin pos., mucate pos.1-4 days. Abortusbovis was combined with Abony (1,4,5,12:b:e,n,x), gelatin neg. The name Abortusbovis was dropped.
I	B	Abortuscanis	4,5,12	b	Rz <sub>5</sub>	Abortuscanis was combined with Paratyphi B (1,4,[5],12:b:1,2) in 1938 and the name Abortuscanis was dropped.
I	B	Abortusequi	4,12	–	e,n,x	
I	B	Abortusovis	4,12	c	1,6	
I	E <sub>4</sub>	Accra	1,3,19	b	z <sub>6</sub>	
I	D <sub>2</sub>	Ackwepe	9,46	l,w	–	
II	G	Acres	1,13,23	b	[1,5]:z <sub>42</sub>	
I	E <sub>1</sub>	Adabraka	3,10	Z <sub>4</sub> ,Z <sub>23</sub>	[1,7]	
I	M	Adamstown	28	k	1,6	
I	F	Adamstua	11	e,h	1,6	
I	U	Adana	43	Z <sub>10</sub>	1,5	
I	O	Adelaide	35	f,g	–	Adelaide may possess H phase Rz <sub>27</sub> .
I	I	Adeoyo	16	g,m,[t]	–	
I	M	Aderike	28	Z <sub>38</sub>	[e,n,Z <sub>15</sub> ]	
I	C <sub>1</sub>	Adime	6,7	b	1,6	
I	G	Adjame	13,23	r	1,6	
I	C <sub>1</sub>	Aequatoria	6,7	Z <sub>4</sub> ,Z <sub>23</sub>	[e,n,Z <sub>15</sub> ]	
I	C <sub>2</sub>	Aesch	6,8	Z <sub>60</sub>	1,2	
I	H	Aflao	1,6,14,25	l,Z <sub>28</sub>	e,n,x	
I	B	Africana	4,12	r,i	l,w	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	C <sub>1</sub>	Afula	6,7	f,g,t	e,n,x	
I	B	Agama	4,12	i	1,6	
I	I	Agbara	16	i	1,6	
I	G	Agbeni	1,13,23	g,m,[t]	–	
I	E <sub>1</sub>	Agege	3,10	c	e,n,z <sub>15</sub>	
I	N	Ago	30	z <sub>38</sub>	–	
I	O	Agodi	35	g,t	–	
I	B	Agona	1,4,[5],12	f,g,s	[1,2]	
I	G	Agoueve	13,22	z <sub>29</sub>	–	
I	J	Ahanou	17	i	1,7	
I	U	Ahepe	43	z <sub>35</sub>	1,6	
I	E <sub>4</sub>	Ahmadi	1,3,19	d	1,5	
I	S	Ahoutoue	41	z <sub>35</sub>	1,6	
I	U	Ahuza	43	k	1,5	
I	G	Ajiobo	13,23	z <sub>4</sub> ,z <sub>23</sub>	–	
I	C <sub>2</sub>	Akanji	6,8	r	1,7	
I	I	Akuafo	16	y	1,6	
I	D <sub>1</sub>	Alabama	9,12	c	e,n,z <sub>15</sub>	
I	O	Alachua	35	z <sub>4</sub> ,z <sub>23</sub>	–	Alachua may possess H phase Rz <sub>37</sub> or Rz <sub>45</sub> .
I	C <sub>2</sub>	Alagbon	8	y	1,7	
I	C <sub>1</sub>	Alamo	6,7	g,z <sub>51</sub>	1,5	
I	C <sub>2</sub>	Albany	8,20	z <sub>4</sub> ,z <sub>24</sub>	–	Albany may possess H phase Rz <sub>45</sub> .
I	B	Albert	4,12	z <sub>10</sub>	e,n,x	
I	H	Albuquerque	1,6,14,24	d	z <sub>6</sub>	
II	E <sub>1</sub>	Alexander	3,10	z	1,5	
I	X	Alexanderplatz	47	z <sub>38</sub>	–	
I	C <sub>2</sub>	Alexanderpolder	8	c	l,w	
I	E <sub>1</sub>	Alfort	3,10	f,g	e,n,x	
I	P	Alger	38	l,v	1,2	
I	E <sub>4</sub>	Alkmaar	1,3,19	a	l,w	
I	R	Allandale	1,40	k	1,6	
I	E <sub>1</sub>	Allerton	3,10	b	1,6	
I	C <sub>2</sub>	Alminko	8,20	g,s,t	–	
II	R	Alsterdorf	1,40	g,m,[s],t	[1,5]	
I	B	Altendorf	4,12	c	1,7	IP combined Womba (4,12,27:c:1,7) with Altendorf to form Altendorf 4,12,27:c:1,7.
I	C <sub>2</sub>	Altona	8,20	r,[i]	z <sub>6</sub>	Pikine (8,20:r:z <sub>6</sub> ) was combined with Altona and called Altona.
I	E <sub>1</sub>	Amager	3,10	y	1,2	IP combined Tuebingen (3,15:y:1,2) with Amager to form Amager 3,10,[15]:y:1,2. Amager may possess H phase Rz <sub>45</sub> .
I	F	Amba	11	k	l,z <sub>13</sub> ,z <sub>28</sub>	
I	C <sub>1</sub>	Amersfoort	6,7	d	e,n,x	IP combined Omderman (6,7,14:d:e,n,x) with Amersfoort to form Amersfoort 6,7,14:d:e,n,x.
I	C <sub>2</sub>	Amherstiana	8	l,(v)	1,6	
I	I	Amina	16	i	1,5	
I	E <sub>1</sub>	Aminatu	3,10	a	1,2	
I	E <sub>1</sub>	Amounderness	3,10	i	1,5	
I	M	Amoutive	28	d	1,5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>1</sub>	Amsterdam	3,10	g,m,s	–	IP combined Drypool (3,15:g,m,s) and Drypool var. O 34+ with Amsterdam to form Amsterdam 3,10,[15],[15,34]:g,m,s:-.
I	I	Amunigun	16	a	1,6	
I	E <sub>1</sub>	Anatum	3,10	e,h	1,6	IP combined Newington (3,15:e,h:1,6) and Minneapolis (3,15,34:e,h:1,6) with Anatum to form Anatum 3,10,[15],[15,34]:e,h:1,6.
I	E <sub>1</sub>	Anderlecht	3,10	c	1,w	
I	O	Anecho	35	g,s,t	–	
I	Q	Anfo	39	y	1,2	
I	C <sub>2</sub>	Angers	8,20	Z <sub>35</sub>	Z <sub>6</sub>	
I	N	Angoda	30	k	e,n,x	
II	D <sub>1</sub>	Angola	1,9,12	z	Z <sub>6</sub>	
I	I	Angouleme	16	Z <sub>10</sub>	Z <sub>6</sub>	
I	X	Anie	47	(g),m,t	–	IP combined Anie with Mesbit (47:m,t:[e,n,Z <sub>15</sub> ]) and called it Mesbit.
I	M	Ank	28	k	e,n,Z <sub>15</sub>	
I	G	Anna	13,23	Z <sub>35</sub>	e,n,Z <sub>15</sub>	
I	I	Annedal	16	r,i	e,n,x	
I	D <sub>1</sub>	Antarctica	9,12	g,Z <sub>63</sub>	–	
I	57	Antonio	57	a	Z <sub>6</sub>	
I	51	Antsalova	51	z	1,5	
I	T	Antwerpen	1,42	c	e,n,Z <sub>15</sub>	
I	W	Apapa	45	m,t	–	
I	C <sub>2</sub>	Apeyeme	8,20	Z <sub>38</sub>	–	
I	N	Aqua	30	k	1,6	
I	N	Aragua	30	Z <sub>29</sub>	–	
I	H	Arapahoe	1,6,14	Z <sub>4</sub> ,Z <sub>23</sub>	1,5	
I	C <sub>1</sub>	Ardwick	6,7,14	f,g	–	IP combined Ardwick with Rissen (6,7:f,g:-) to form Rissen 6,7,14:f,g:-. Ardwick is now called Rissen var. O 14+ by IP.
I	B	Arechavaleta	4,[5],12	a	[1,7]	
IV	C <sub>1</sub>	Argentina	6,7	Z <sub>36</sub>	–	
I	E <sub>3</sub>	Arkansas	3,15,34	e,h	1,5	IP combined Arkansas and Newhaw (3,15:e,h:1,5) with Muenster (3,10:e,h:1,5) to form Muenster 3,10,[15],[15,34]:e,h:1,5. Arkansas is now called Muenster var. O 15+, 34+ by IP.
II	56	Artis	56	b	–	
I	U	Arusha	43	z	e,n,Z <sub>15</sub>	
I	N	Aschersleben	30	b	1,5	
I	M	Ashanti	28	b	1,6	
II	51	Askraal	51	1,Z <sub>28</sub>	[Z <sub>6</sub> ]	
I	L	Assen	21	a	[1,5]	
I	E <sub>1</sub>	Assinie	3,10	1,w	Z <sub>6</sub>	Assinie may possess H phase RZ <sub>45</sub> .
I	E <sub>1</sub>	Asylanta	3,10	c	1,2	
I	C <sub>2</sub>	Atakpame	8,20	e,h	1,7	
I	F	Atento	11	b	1,2	
I	C <sub>1</sub>	Athinai	6,7	i	e,n,Z <sub>15</sub>	
I	F	Ati	11	d	1,2	
I	G	Atlanta	13,23	b	–	Atlanta was combined with Mississippi (1,13,23:b:1,5). The name Atlanta has been dropped.

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	Z	Atra	50	m,t	Z <sub>6</sub> :Z <sub>42</sub>	
I	C <sub>1</sub>	Augustenborg	6,7, <u>14</u>	i	1,2	
I	C <sub>1</sub>	Austin	6,7	a	1,7	
I	Y	Australia	48	l,v	1,5	
I	I	Avignon	16	y	e,n,Z <sub>15</sub>	
I	E <sub>4</sub>	Avonmouth	1,3,19	i	e,n,Z <sub>15</sub>	
I	B	Ayinde	<u>1,4,12,27</u>	d	Z <sub>6</sub>	
I	B	Ayton	<u>1,4,12,27</u>	l,w	Z <sub>6</sub>	
I	B	Azteca	4, <u>5</u> ,12, <u>27</u>	l,v	1,5	
I	M	Babelsberg	28	Z <sub>4</sub> ,Z <sub>23</sub>	[e,n,Z <sub>15</sub> ]	
I	M	Babili	28	Z <sub>35</sub>	1,7	
II	C <sub>1</sub>	Bacongo	6,7	Z <sub>36</sub>	Z <sub>42</sub>	
I	I	Badagry	16	Z <sub>10</sub>	1,5	
I	L	Baguida	21	Z <sub>4</sub> ,Z <sub>23</sub>	–	
I	N	Baguirmi	30	y	e,n,x	
I	G	Bahati	13,22	b	e,n,Z <sub>15</sub>	
I	H	Bahrenfeld	6,14,24	e,h	1,5	
I	C <sub>1</sub>	Baiboukoum	6,7	k	1,7	
I	D <sub>2</sub>	Baildon	9,46	a	e,n,x	
I	M	Bakau	28	a	1,7	
V	Y	Balboa	48	Z <sub>41</sub>	–	
I	W	Balcones	45	Z <sub>36</sub>	–	
I	B	Ball	<u>1,4,12,27</u>	y	e,n,x	IP combined Ball with Ruki (4,5,12:y:e,n,x) and Dalat (4,5,27:y:e,n,x) to form Ball <u>1,4,[5],12,27</u> :y:e,n,x.
I	J	Bama	17	m,t	–	
I	D <sub>2</sub>	Bamboye	9,46	b	l,w	
I	D <sub>2</sub>	Bambylor	9,46	z	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Banalia	6,8	b	Z <sub>6</sub>	
I	B	Banana	<u>1,4,[5],12</u>	m,t	[1,5]	Banana may possess H phase Rz <sub>45</sub> .
I	M	Banco	28	r,i	1,7	
I	O	Bandia	35	i	l,w	
I	P	Bangkok	38	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	D <sub>1</sub>	Bangui	9,12	d	e,n,Z <sub>15</sub>	
I	H	Banjul	1,6,14,25	a	e,n,Z <sub>15</sub>	
II	C <sub>2</sub>	Baragwanath	6,8	m,t	1,5	
I	C <sub>2</sub>	Bardo	8	e,h	1,2	
I	C <sub>1</sub>	Bareilly	6,7, <u>14</u>	y	1,5	
I	C <sub>2</sub>	Bargny	8, <u>20</u>	i	1,5	
I	I	Barmbek	16	d	Z <sub>6</sub>	
I	I	Barranquilla	16	d	e,n,x	
I	54	Barry	54	Z <sub>10</sub>	e,n,Z <sub>15</sub>	
II	58	Basel	58	1,Z <sub>13</sub> ,Z <sub>28</sub>	1,5	
I	D	Basingstoke	9,46	Z <sub>35</sub>	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Bassa	6,8	m,t	–	
I	M	Bassadji	28	r	1,6	
I	57	Batonrouge	57	b	e,n,Z <sub>15</sub>	
I	I	Battle	16	1,Z <sub>13</sub> ,Z <sub>28</sub>	1,6	
I	C <sub>2</sub>	Bazenheid	8, <u>20</u>	Z <sub>10</sub>	1,2	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	C <sub>2</sub>	Be	8,20	a	[z <sub>6</sub> ]	
I	H	Beaudesert	[1],6,14,[25]	e,h	1,7	
II	B	Bechuana	1,4,12,27	g,[m],t	[1,5]	
I	E <sub>4</sub>	Bedford	1,3,19	1,z <sub>13</sub> ,z <sub>28</sub>	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Belem	6,8	c	e,n,x	
I	C <sub>2</sub>	Belfast	6,8	c	1,7	
I	C <sub>2</sub>	Bellevue	8	z <sub>4</sub> ,z <sub>23</sub>	1,7	
I	I	Bellville	16	e,n,x	1,(5),7	
II	K	Beloha	18	z <sub>36</sub>	–	
I	E <sub>1</sub>	Benfica	3,10	b	e,n,x	
I	R	Benguella	40	b	z <sub>6</sub>	
I	D <sub>2</sub>	Benin	9,46	y	1,7	
I	C <sub>2</sub>	Benue	6,8	y	1,w	
I	X	Bere	47	z <sub>4</sub> ,z <sub>23</sub>	[z <sub>6</sub> ]	Bere may possess H phase Rz <sub>45</sub> .
I	D <sub>2</sub>	Bergedorf	9,46	e,h	1,2	
I	X	Bergen	47	i	e,n,z <sub>15</sub>	
I	51	Bergues	51	z <sub>10</sub>	1,5	
I	U	Berkeley	43	a	1,5	
I	J	Berlin	17	d	1,5	
IV	R	Bern	1,40	z <sub>4</sub> ,z <sub>32</sub>	–	
I	D <sub>1</sub>	Berta	1,9,12	[f],g,t	–	
I	E <sub>1</sub>	Bessi	3,10	i	e,n,x	
I	E <sub>4</sub>	Bethune	1,3,19	k	1,7	
II	59	Betioky	59	k	(z)	
I	E <sub>1</sub>	Biafra	3,10	z <sub>10</sub>	z <sub>6</sub>	
I	E <sub>4</sub>	Bida	1,3,19	c	1,6	
I	N	Bietri	30	y	1,5	
I	J	Bignona	17	b	e,n,z <sub>15</sub>	
I	R	Bijlmer	1,40	g,m	–	
II	X	Bilthoven	47	a	[1,5]	
I	E <sub>4</sub>	Bilu	(1),3,10,(19)	f,g,t	1,(2),7	
I	X	Binche	47	z <sub>4</sub> ,z <sub>23</sub>	1,w	
I	X	Bingerville	47	z <sub>35</sub>	e,n,z <sub>15</sub>	
I	W	Binningen	45	g,s,t	–	
I	E <sub>2</sub>	Binza	3,15	y	1,5	IP combined Binza and Thomasville (3,15,34:y:1,5) with Orion (3,10:y:1,5) to form Orion 3,10,[15],[15,34]:y:1,5. Binza is now called Orion var. O 15+ by IP.
I	C <sub>1</sub>	Birkenhead	6,7	c	1,6	
I	E <sub>1</sub>	Birmingham	3,10	d	1,w	
I	B	Bispebjerg	1,4,[5],12	a	e,n,x	
I	B	Bissau	4,12	c	e,n,x	
II	D <sub>1</sub>	Blankenese	1,9,12	b	z <sub>6</sub>	
II	J	Bleaton	17	(f),g,t	[e,n,x,z <sub>15</sub> ]	IP has dropped f.
I	D <sub>1</sub>	Blegdam	9,12	g,m,q	–	
I	H	Blijdorp	1,6,14,25	c	1,5	
I	X	Blitta	47	y	e,n,x	
I	C <sub>2</sub>	Blockley	6,8	k	1,5	Blockley may possess H phase Rz <sub>58</sub>
II	C <sub>1</sub>	Bloemfontein	6,7	b	[e,n,x]:z <sub>42</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>1</sub>	Bloomsbury	3,10	g,t	1,5	
I	K	Blukwa	18	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	V	Bobo	44	d	1,5	
I	B	Bochum	4,[5],12	r	1,w	
IV	53	Bockenheim	1,53	Z <sub>36</sub> ,Z <sub>38</sub>	–	
I	N	Bodjonegoro	30	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	H	Boecker	[1],6,14,[25]	l,v	1,7	
I	M	Bokanjac	28	b	1,7	
II	R	Boksburg	40	g,[m],s,[t]	e,n,x	
I	V	Bolama	44	z	e,n,x	
I	E <sub>1</sub>	Bolombo	3,10	Z <sub>38</sub>	[Z <sub>6</sub> ]	
I	E <sub>1</sub>	Bolton	3,10	y	e,n,Z <sub>15</sub>	
IV	Z	Bonaire	50	Z <sub>4</sub> ,Z <sub>32</sub>	–	
I	J	Bonames	17	a	1,2	
I	C <sub>2</sub>	Bonariensis	6,8	i	e,n,x	
V	Y	Bongor	48	Z <sub>35</sub>	–	
I	C <sub>1</sub>	Bonn	6,7	l,v	e,n,x	
I	X	Bootle	47	k	1,5	
I	G	Borbeck	13,22	l,v	1,6	
I	52	Bordeaux	52	k	1,5	
VI	H	Bornheim	1,6,14,25	Z <sub>10</sub>	1,(2),7	Bornheim was formerly in Subspecies II.
I	C <sub>1</sub>	Bornum	6,7, <u>14</u>	Z <sub>38</sub>	–	IP combined Bornum with Lille (6,7:z <sub>38</sub> :-) to form Lille 6,7, <u>14</u> :z <sub>38</sub> :-. Bornum is now called Lille var. O 14+ by IP.
I	54	Borreze	54	f,g,s	–	
I	T	Borromea	42	i	1,6	
I	I	Bouake	16	z	Z <sub>6</sub>	
II	G	Boulders	<u>1</u> ,13,23	m,t	Z <sub>42</sub>	
I	D <sub>1</sub>	Bournemouth	9,12	e,h	1,2	
I	H	Bouso	1,6,14,25	Z <sub>4</sub> ,Z <sub>23</sub>	[e,n,Z <sub>15</sub> ]	
I	C <sub>2</sub>	Bovismorbificans	6,8	r,[i]	1,5	
I	G	Bracknell	13,23	b	1,6	
I	B	Bradford	4,12, <u>27</u>	r	1,5	
I	C <sub>1</sub>	Braenderup	6,7, <u>14</u>	e,h	e,n,Z <sub>15</sub>	
I	B	Brancaster	<u>1</u> ,4,12, <u>27</u>	Z <sub>29</sub>	–	
I	B	Brandenburg	<u>1</u> ,4,[5],12, <u>27</u>	l,v	e,n,Z <sub>15</sub>	
I	I	Brazil	16	a	1,5	
I	K	Brazos	<u>6</u> ,14,18	a	e,n,Z <sub>15</sub>	
I	C <sub>1</sub>	Brazzaville	6,7	b	1,2	
I	C <sub>2</sub>	Breda	6,8	Z <sub>4</sub> ,Z <sub>23</sub>	e,n,x	
I	B	Bredeney	<u>1</u> ,4,12, <u>27</u>	l,v	1,7	Bredeney may possess H phase Rl <sub>Z<sub>40</sub></sub> instead of l,v.
I	V	Brefet	44	r	e,n,Z <sub>15</sub>	
II	W	Bremen	45	g,m,s,t	e,n,x	
I	C <sub>2</sub>	Breukelen	6,8	l,Z <sub>13</sub> ,[Z <sub>28</sub> ]	e,n,Z <sub>15</sub>	
I	I	Brevik	16	z	e,n,[x],Z <sub>15</sub>	
I	B	Brezany	<u>1</u> ,4,12, <u>27</u>	d	1,6	
I	F	Brijbhumi	11	i	1,5	
I	C <sub>2</sub>	Brikama	8, <u>20</u>	r,[i]	1,w	
I	M	Brisbane	28	z	e,n,Z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	G	Bristol	13,22	z	1,7	
I	T	Brive	<u>1</u> ,42	r	1,w	
I	T	Broc	42	Z <sub>4</sub> ,Z <sub>23</sub>	e,n,Z <sub>15</sub>	
I	G	Bron	13,22	g,m	[e,n,Z <sub>15</sub> ]	
I	C <sub>2</sub>	Bronx	6,8	c	1,6	
V	66	Brookfield	66	Z <sub>41</sub>	–	
I	I	Brooklyn	16	l,w	e,n,x	
I	E <sub>4</sub>	Broughton	1,3,19	b	1,w	
I		Broxbourne				Combined with Wien. The name Broxbourne has been dropped.
I	C <sub>1</sub>	Bruck	6,7	z	1,w	
I	C <sub>2</sub>	Brunei	8, <u>20</u>	y	1,5	
I	I	Brunflo	16	r	1,7	
I	C <sub>2</sub>	Bsilla	6,8	r	1,2	
I	B	Budapest	<u>1</u> ,4,12, <u>27</u>	g,t	–	
I	R	Bukavu	<u>1</u> ,40	l,Z <sub>28</sub>	1,5	
I	C <sub>2</sub>	Bukuru	6,8	b	1,w	
II	R	Bulawayo	<u>1</u> ,40	z	1,5	
I	C <sub>2</sub>	Bulgaria	6,8	y	1,6	
I	F	Bullbay	11	l,v	e,n,Z <sub>15</sub>	
I	C <sub>1</sub>	Bulovka	6,7	Z <sub>44</sub>	–	
II	U	Bunnik	43	Z <sub>42</sub>	[1,5,7]	
I	I	Burgas	16	l,v	e,n,Z <sub>15</sub>	
I	S	Burundi	41	a	–	
I	B	Bury	4,12, <u>27</u>	c	Z <sub>6</sub>	
I	C <sub>1</sub>	Businga	6,7	z	e,n,Z <sub>15</sub>	
I	E <sub>1</sub>	Butantan	3,10	b	1,5	IP combined Rosenthal (3,15:b:1,5) and unnamed 3,15,34:b:1,5 with Butantan to form Butantan 3,10,[15],[15,34]:b:1,5.
I	52	Butare	52	e,h	1,6	
I	H	Buzu	[1],6,14,[25]	i	1,7	
I	I	Caen	16	d	1,w	
I	E <sub>1</sub>	Cairina	3,10	Z <sub>35</sub>	Z <sub>6</sub>	
I	W	Cairns	45	k	e,n,Z <sub>15</sub>	
I	B	Cairo	<u>1</u> ,4,12, <u>27</u>	d	1,2	IP combined Cairo with Stanley (4,5,12:d:1,2) to form Stanley <u>1</u> ,4,[5],12, <u>27</u> :d:1,2. The name Cairo has been dropped.
I	E <sub>4</sub>	Calabar	1,3,19	e,h	1,w	
II	B	Caledon	<u>1</u> ,4,12, <u>27</u>	g,[m],[s],t	e,n,x	
I	B	California	4,12	g,m,t	–	
II	C <sub>1</sub>	Calvinia	6,7	a	Z <sub>42</sub>	
I	O	Camberene	35	Z <sub>10</sub>	1,5	
I	D <sub>1</sub>	Camberwell	9,12	r	1,7	
I	E <sub>2</sub>	Cambridge	3,15	e,h	1,w	IP combined Cambridge and Wildwood (3,15,34:e,h:l,w) with Meleagridis (3,10:e,h:l,w) to form Meleagridis 3,10,[15],[15,34]:e,h:l,w. Cambridge is now called Meleagridis var. O 15+ by IP.
V	V	Camdeni	44	r	–	
I	D <sub>1</sub>	Campinense	9,12	r	e,n,Z <sub>15</sub>	
I	B	Canada	4,12, <u>27</u>	b	1,6	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	R	Canary	40	l,v	1,6	
II	D <sub>1</sub>	Canastel	9,12	z <sub>29</sub>	[1,5]	
I	M	Cannobio	28	z <sub>4</sub> ,z <sub>23</sub>	1,5	
I	E <sub>4</sub>	Cannonhill	1,3,15,19	y	e,n,x	
I	E <sub>4</sub>	Cannstatt	1,3,19	m,t	–	
I	E <sub>3</sub>	Canoga	3,15,34	g,s,t	–	IP combined Canoga and Halmstad (3,15:g,s,t-) with Westhampton (3,10:g,s,t-) to form Westhampton 3,10,[15],[15,34]:g,s,t-.- Canoga is now called Westhampton var. O 15+, 34+ by IP.
I	54	Canton	54	z <sub>10</sub>	e,n,x	
II	C <sub>1</sub>	Cape	6,7	z <sub>6</sub>	1,7	
I	H	Caracas	[1],6,14,[25]	g,m,s	–	
I	C <sub>1</sub>	Cardiff	6,7	k	R1,10	IP combined Cardiff that contains H phase R1,10 with Thompson (6,7,14:k:1,5).
I	I	Cardoner	16	g,s,t	–	
II	P	Carletonville	38	d	[1,5]	
I	J	Carmel	17	l,v	e,n,x	
I	K	Carnac	18	z <sub>10</sub>	z <sub>6</sub>	
I	E <sub>4</sub>	Carno	1,3,19	z	1,w	
I	H	Carrau	6,14,[24]	y	1,7	
I	V	Carswell	44	g,z <sub>51</sub>	–	
I	W	Casablanca	45	k	1,7	
I	R	Casamance	40	z	e,n,x	
I	M	Catalunia	28	l,z <sub>13</sub> ,z <sub>28</sub>	1,5	
I	H	Catanzaro	6,14	g,s,t	–	
I	C <sub>1</sub>	Cayar	6,7	z	e,n,x	
II	M	Ceres	28	z	z <sub>39</sub>	
I	K	Cerro	18	z <sub>4</sub> ,z <sub>23</sub>	[1,5]	Cerro was combined with Siegburg (6,14,18:z <sub>4</sub> ,z <sub>23</sub> : [1,5]) and called Cerro. Cerro may possess H phase Rz <sub>45</sub> .
I	D <sub>2</sub>	Ceyco	9,46	k	z <sub>35</sub>	
I	G	Chagoua	1,13,23	a	1,5	
I	C <sub>2</sub>	Chailey	6,8	z <sub>4</sub> ,z <sub>23</sub>	[e,n,z <sub>15</sub> ]	
IV	I	Chameleon	16	z <sub>4</sub> ,z <sub>32</sub>	–	
I	Q	Champaign	39	k	1,5	Champaign may possess H phase Rz <sub>48</sub>
I	F	Chandans	11	d	e,n,x:[r]	
I	H	Charity	[1],6,14,[25]	d	e,n,x	
I	C <sub>2</sub>	Charlottenburg	6,8	k	e,n,z <sub>15</sub>	
I	D <sub>2</sub>	Cheltenham	9,46	b	1,5	
II	X	Chersina	47	z	z <sub>6</sub>	
I	B	Chester	1,4,[5],12	e,h	e,n,x	
I	M	Chicago	28	r,[i]	1,5	
I	E <sub>4</sub>	Chichester	1,3,19	i	1,6	
I	H	Chichiri	6,14,24	z <sub>4</sub> ,z <sub>24</sub>	–	
I	C <sub>1</sub>	Chile	6,7	z	1,2	
I	C <sub>2</sub>	Chincol	6,8	g,m,[s]	[e,n,x]	
I	F	Chingola	11	e,h	1,2	
II	T	Chinovum	42	b	1,5	
I	F	Chiredzi	11	c	1,5	
I	E <sub>4</sub>	Chittagong	1,3,10,19	b	z <sub>35</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	C <sub>1</sub>	Choleraesuis	6,7	c	1,5	H <sub>2</sub> S negative
I	C <sub>1</sub>	Choleraesuis var. Kunzendorf		6,7	c	[1,5] H <sub>2</sub> S positive
I	C <sub>2</sub>	Chomedey	8,20	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	V	Christiansborg	44	z <sub>4</sub> ,z <sub>24</sub>	–	
II	E <sub>1</sub>	Chudleigh	3,10	e,n,x	1,7	
I	B	Clackamas	4,12	l,v,[z13]	1,6	IP does not get z <sub>13</sub>
I	D <sub>1</sub>	Claibornei	1,9,12	k	1,5	
I	E <sub>1</sub>	Clerkenwell	3,10	z	1,w	
I	C <sub>2</sub>	Cleveland	6,8	z <sub>10</sub>	1,7	
I	E <sub>2</sub>	Clichy	3,15	a	1,5	IP combined Clichy with Goelzau (3,10:a:1,5) to form Goelzau 3,10,[15]:a:1,5. Clichy is now called Goelzau var. O 15+ by IP.
II	G	Clifton	13,22	z <sub>29</sub>	1,5	
II	V	Clovelly	1,44	z <sub>39</sub>	[e,n,x,z <sub>15</sub> ]	
I	D <sub>2</sub>	Cochin	9,46	k	1,5	
I	C <sub>2</sub>	Cocody	8,20	r,i	e,n,z <sub>15</sub>	
I	B	Coeln	4,[5],12	y	1,2	
I	C <sub>1</sub>	Coleypark	6,7,14	a	1,w	
I	C <sub>1</sub>	Colindale	6,7	r	1,7	
I	F	Colobane	11	k	1,7	
I	P	Colombo	38	y	1,6	
I	C <sub>1</sub>	Colorado	6,7	l,w	1,5	
I	C <sub>1</sub>	Concord	6,7	l,v	1,2	
I	G	Congo	13,23	g,[m],[s],t	–	IP calls Congo 13,23:g,m,s,t:-.
I	F	Connecticut	11	l,z <sub>13</sub> ,z <sub>28</sub>	1,5	
II	J	Constantia	17	z	l,w:z <sub>42</sub>	
I	T	Coogee	42	l,v	e,n,z <sub>15</sub>	
I	Q	Cook	39	Rz <sub>48</sub>	1,5	IP combined Cook with Champaign (39:k:1,5). The name Cook has been dropped.
I	E <sub>1</sub>	Coquilhatville	3,10	z <sub>10</sub>	1,7	
I	C <sub>1</sub>	Coromandel	6,7	l,v	z <sub>35</sub>	
I	C <sub>2</sub>	Corvallis	8,20	z <sub>4</sub> ,z <sub>23</sub>	[z <sub>6</sub> ]	
I	M	Cotham	28	i	1,5	
I	K	Cotia	18	–	1,6	
I	C <sub>2</sub>	Cremieu	6,8	e,h	1,[6]	
I	F	Crewe	11	z	1,5	
I	M	Croft	28	g,m,s	[e,n,z <sub>15</sub> ]	
I	67	Crossness	67	r	1,2	
I	G	Cubana	1,13,23	z <sub>29</sub>	–	Cubana may possess H phase Rz <sub>37</sub> or Rz <sub>43</sub> .
I	E <sub>1</sub>	Cuckmere	3,10	i	1,2	
I	M	Cullingworth	28	d	1,w	
I	Q	Cumberland	39	i	e,n,x	
I	C <sub>2</sub>	Curacao	6,8	a	1,6	
I	C <sub>2</sub>	Cyprus	6,8	i	1,w	
I	54	Czernyring	54	r	1,5	
I	C <sub>2</sub>	Daarle	6,8	y	e,n,x	
I	C <sub>2</sub>	Dabou	8,20	z <sub>4</sub> ,z <sub>23</sub>	l,w	
I	51	Dadzie	51	l,v	e,n,x	
I	Y	Dahlem	48	k	e,n,z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	X	Dahomey	47	k	1,6	Dahomey may possess H phase Rz <sub>58</sub> .
I	J	Dahra	17	b	1,5	
I	M	Dakar	28	a	1,6	
I	I	Dakota	16	Z <sub>35</sub>	e,n,Z <sub>15</sub>	
I	B	Dalat	4,5,27	y	e,n,x	Dalat was combined with Ball. The name Dalat has been dropped.
I	E <sub>4</sub>	Dallgow	1,3,19	Z <sub>10</sub>	e,n,Z <sub>15</sub>	
I	C <sub>1</sub>	Damman	6,7	a	Z <sub>6</sub>	
I	51	Dan	51	k	e,n,Z <sub>15</sub>	
I	X	Dapango	47	r	1,2	
II	D <sub>1</sub>	Daressalaam	1,9,12	l,w	e,n,x	
I	C <sub>2</sub>	Daula	8,20	z	Z <sub>6</sub>	
I	C <sub>1</sub>	Daytona	6,7	k	1,6	
I	C <sub>1</sub>	Decatur	6,7	c	1,5	IP has dropped Decatur and calls it dulcitol positive, mucate positive variant of Choleraesuis.
I	D <sub>2</sub>	Deckstein	9,46	r	1,7	
II	R	Degania	40	Z <sub>4</sub> ,Z <sub>24</sub>	[Z <sub>39</sub> ]	
IV	R	Degania var. subsp. IV	40	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	K	Delmenhorst	18	Z <sub>71</sub>	–	
I	O	Dembe	35	d	1,w	
I	G	Demerara	13,23	Z <sub>10</sub>	1,w	
I	C <sub>1</sub>	Denver	6,7	a	e,n,Z <sub>15</sub>	
I	B	Derby	1,4,[5],12	f,g	[1,2]	
I	52	Derkle	52	e,h	1,7	
I	E <sub>4</sub>	Dessau	1,3,15,19	g,s,t	–	
II	T	Detroit	42	z	1,5	
I	W	Deversoir	45	c	e,n,x	
I	M	Dibra	28	a	Z <sub>6</sub>	
I	M	Dieuppeul	28	i	1,7	
I	G	Diguel	1,13,22	d	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Diogoye	8,20	Z <sub>41</sub>	Z <sub>6</sub>	
I	L	Diourbel	21	i	1,2	
I	Y	Djakarta	48	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	T	Djama	1,42	Z <sub>29</sub>	[1,5]	
I	C <sub>2</sub>	Djelfa	8	b	1,2	
I	M	Djermaia	28	Z <sub>29</sub>	–	
I	J	Djibouti	17	Z <sub>10</sub>	e,n,x	
I	C <sub>1</sub>	Djugu	6,7	Z <sub>10</sub>	e,n,x	
I	D <sub>2</sub>	Doba	9,46	a	e,n,Z <sub>15</sub>	
I	M	Doel	28	z	1,6	
I	C <sub>2</sub>	Doncaster	6,8	a	1,5	
I	N	Donna	30	l,v	1,5	
I	M	Doorn	28	i	1,2	
I	M	Douala	28	i	1,w	
I	Z	Dougi	50	y	1,6	
I	N	Doulassame	30	a	e,n,Z <sub>15</sub>	
I	M	Dresden	28	c	e,n,x	
I	R	Driffield	1,40	d	1,5	
I	B	Drogana	1,4,12,27	r,(i)	e,n,Z <sub>15</sub>	IP calls Drogana r,i.

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>2</sub>	Drypool	3,[15],[15,34]	g,m,s	–	IP combined Drypool 3, <u>15</u> :g,m,s:- and Drypool var. O34+ with Amsterdam (3,10:g,m,s:-) to form Amsterdam 3,10,[15],[15,34]:g,m,s:-. Drypool is now called Amsterdam var. O 15+ or O 15+, 34+ by IP.
I	D <sub>1</sub>	Dublin	<u>1,9,12</u> ,[Vi]	g,p	–	
II	S	Dubrovnik	41	z	1,5	
I	C <sub>2</sub>	Duesseldorf	6,8	z <sub>4</sub> ,Z <sub>24</sub>	–	
I	W	Dugbe	45	d	1,6	
I	B	Duisburg	<u>1,4,12,27</u>	d	e,n,Z <sub>15</sub>	
II	D <sub>2</sub>	Duivenhoks	9,46	g,[m],[s],t	[e,n,x]	
I	E <sub>1</sub>	Dumfries	3,10	r,[i]	1,6	
I	C <sub>2</sub>	Dunkwa	6,8	d	1,7	
I	D <sub>1</sub>	Durban	9,12	a	e,n,Z <sub>15</sub>	
II	B	Durbanville	<u>1,4,12,27</u>	[Z <sub>39</sub> ]	1,[5],7	
I	G	Durham	13,23	b	e,n,Z <sub>15</sub>	
I	R	Duval	1,40	b	e,n,Z <sub>15</sub>	
I	O	Ealing	35	g,m,s	–	
I	D <sub>1</sub>	Eastbourne	<u>1,9,12</u>	e,h	1,5	
I	E <sub>4</sub>	Eastglam	1,3,19	c	1,5	
I	M	Eberswalde	28	c	1,6	
I	C <sub>2</sub>	Eboko	6,8	b	1,7	
I	O	Ebrie	35	g,m,t	–	
I	P	Echa	38	k	1,2	
I	C <sub>1</sub>	Edinburg	6,7	b	1,5	
I	C <sub>2</sub>	Edmonton	6,8	l,v	e,n,Z <sub>15</sub>	
I	S	Egusi	41	d	[1,5]	
I	T	Egusitoo	<u>1,42</u>	b	z <sub>6</sub>	
IIIb	61	Eilbeck	61	i	z	Eilbeck was formerly in subspecies II, but is now combined with <i>Arizona</i> 26:33:31. The name Eilbeck has been dropped.
I	C <sub>1</sub>	Eimsbuettel	6,7, <u>14</u>	d	l,w	IP combined Eimsbuettel with Livingstone (6,7:d:l,w) to form Livingstone 6,7, <u>14</u> :d:l,w. Eimbuettel is now called Livingstone var. O 14+ by IP.
I	C <sub>1</sub>	Eingedi	6,7	f,g,t	1,2,7	
II	W	Ejeda	45	a	Z <sub>10</sub>	
I	B	Eko	4,12	e,h	1,6	
I	D <sub>2</sub>	Ekotedo	9,46	z <sub>4</sub> ,Z <sub>23</sub>	–	
I	X	Ekpoui	47	Z <sub>29</sub>	–	
I	E <sub>1</sub>	Elisabethville	3,10	r	1,7	
I	D <sub>1</sub>	Elokate	9,12	c	1,7	
I	D <sub>1</sub>	Elomrane	<u>1,9,12</u>	z38	–	
II	I	Elsiesrivier	16	[e,n,x]	1,6:Z <sub>42</sub>	
I	C <sub>2</sub>	Emek	8, <u>20</u>	g,m,s	–	
I	P	Emmastad	38	r	1,6	
II	H	Emmerich	6,14	[m,t]	e,n,x	
	I	H	Encino	1,6,14,25	d	1,Z <sub>13</sub> ,Z <sub>28</sub>
I	O	Enschede	35	Z <sub>10</sub>	l,w	
I	B	Entebbe	<u>1,4,12,27</u>	z	z <sub>6</sub>	
I	D <sub>1</sub>	Enteritidis	<u>1,9,12</u>	[f],g,m,[p],[t]	[1,7]	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	I	Enugu	16	1,[z13],z <sub>28</sub>	[1,5]	
I	E <sub>1</sub>	Epicrates	3,10	b	1,w	
I	F	Epinay	11	a	1,z <sub>13</sub> ,z <sub>28</sub>	
I	B	Eppendorf	<u>1,4,12,27</u>	d	1,5	
II	G	Epping	<u>1,13,23</u>	e,n,x	1,[5],7	
II	Y	Erlangen	48	g,m,t	–	
I	C <sub>1</sub>	Escanaba	6,7	k	e,n,z <sub>15</sub>	
I	D <sub>1</sub>	Eschberg	9,12	d	1,7	
I	E <sub>2</sub>	Eschersheim	3, <u>15</u>	d	e,n,x	IP combined Eschersheim with Souza (3,10:d:e,n,x) to form Souza 3,10,[15]:d:e,n,x. Eschersheim is now called Souza var. O 15+ by IP.
I	C <sub>1</sub>	Eschweiler	6,7	z <sub>10</sub>	1,6	
I	B	Essen	4,12	g,m	–	
II	Y	Etosha	48	d	1,11	Etosha was not considered a new serotype by Kauffmann and is not used.
I	F	Etterbeek	11	z <sub>4</sub> ,z <sub>23</sub>	e,n,z <sub>15</sub>	
I	F	Euston	11	r,i	e,n,x,z <sub>15</sub>	
I	E <sub>1</sub>	Everleigh	3,10	z <sub>29</sub>	e,n,x	
I	M	Ezra	28	z	1,7	
I	M	Fairfield	28	r	1,w	
I	M	Fajara	28	1,z <sub>28</sub>	e,n,x	
I	T	Faji	<u>1,42</u>	a	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Falkensee	3,10	i	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Fallowfield	3,10	1,z <sub>13</sub> ,z <sub>28</sub>	e,n,z <sub>15</sub>	
II	R	Fandran	<u>1,40</u>	z <sub>35</sub>	e,n,x,z <sub>15</sub>	
I	F	Fann	11	l,v	e,n,x	
I	G	Fanti	13,23	z <sub>38</sub>	–	
I	M	Farakan	28	z <sub>10</sub>	1,5	
I	U	Farcha	43	y	1,2	
I	E <sub>4</sub>	Fareham	1,3,19	r,i	1,w	
I	G	Farmsen	13,23	z	1,6	
I	B	Farsta	4,12	i	e,n,x	
I	Z	Fass	50	l,v	1,2	
II	Z	Faure	50	z <sub>42</sub>	1,7	
I	C <sub>2</sub>	Fayed	6,8	1,w	1,2	
VI	H	Ferlac	1,6,14,25	a	e,n,x	
I	S	Ferlo	41	k	1,6	
I	C <sub>2</sub>	Ferruch	8	e,h	1,5	
I	B	Finaghy	4,12	y	1,6	
II	E <sub>1</sub>	Finchley	3,10	z	e,n,x	
I	F	Findorff	11	d	z <sub>6</sub>	
I	H	Finkenwerder	[1],6,14,[25]	d	1,5	
I	I	Fischerhuetten	16	a	e,n,z <sub>15</sub>	
I	H	Fischerkietz	1,6,14,25	y	e,n,x	
I	V	Fischerstrasse	44	d	e,n,z <sub>15</sub>	
I	Y	Fitzroy	48	e,h	1,5	
IV	Z	Flint	50	z <sub>4</sub> ,z <sub>23</sub>	–	
I	E <sub>1</sub>	Florian	3,10,[15]	z <sub>4</sub> ,z <sub>24</sub>	–	
I	H	Florida	[1],6,14,[25]	d	1,7	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	51	Flottbek	52	b	[e,n,x]	
I	K	Fluntern	6,14,18	b	1,5	
I	W	Fomeco	45	b	e,n,z <sub>15</sub>	
I	I	Fortlamy	16	z	1,6	
I	B	Fortune	<u>1,4,12,27</u>	z <sub>10</sub>	z <sub>6</sub>	
II	P	Foulpointe	38	g,t	–	
I	D <sub>1</sub>	Franken	<u>1,9,12</u>	z <sub>60</sub>	z <sub>67</sub>	
I	I	Frankfurt	16	i	e,n,z <sub>15</sub>	
I	P	Freetown	38	y	1,5	
I	E <sub>1</sub>	Freiburg	3,10	l,z <sub>13</sub>	1,2	
II	T	Fremantle	42	(f),g,t	–	
I	D <sub>2</sub>	Fresno	9,46	z <sub>38</sub>	–	
I	G	Friedenau	13,22	d	1,6	
I	M	Friedrichsfelde	28	f,g	–	
I	D <sub>1</sub>	Frintrop	1,9,12	b	1,5	
I	E <sub>1</sub>	Fufu	3,10	z	1,5	
II	E <sub>1</sub>	Fuhlsbuettel	3,10	l,v	z <sub>6</sub>	
I	E <sub>4</sub>	Fulda	1,3,19	l,w	1,5	
I	B	Fulica	4,[5],12	a	1,5	
I	B	Fyris	4,[5],12	l,v	1,2	
I	C <sub>1</sub>	Gabon	6,7	l,w	1,2	
I	I	Gafsa	16	c	1,6	
I	C <sub>1</sub>	Galiema	6,7, <u>14</u>	k	1,2	
I	E <sub>1</sub>	Galil	3,10	a	e,n,z <sub>15</sub>	
I	F	Gallen	11	a	1,2	
I	D <sub>1</sub>	Gallinarum	<u>1,9,12</u>	–	–	Gallinarum must be identified biochemically.
I	V	Gamaba	<u>1,44</u>	g,m,[s]	–	
I	L	Gambaga	21	z <sub>35</sub>	e,n,z <sub>15</sub>	
I	O	Gambia	35	i	e,n,z <sub>15</sub>	
I	I	Gaminara	16	d	1,7	
I	H	Garba	1,6,14,25	a	1,5	
I	C <sub>1</sub>	Garoli	6,7	i	1,6	
I	O	Gassi	35	e,h	z <sub>6</sub>	
I	D <sub>2</sub>	Gateshead	9,46	g,s,t	–	
I	E <sub>4</sub>	Gatineau	1,3,19	y	1,5	
I	C <sub>1</sub>	Gatow	6,7	y	1,7	
I	C <sub>2</sub>	Gatuni	6,8	b	e,n,x	
I	E <sub>1</sub>	Gbadago	3,10,[ <u>15</u> ]	c	1,5	
I	C <sub>1</sub>	Gdansk	6,7	l,v	z <sub>6</sub>	IP combined Gelsenkirchen (6,7, <u>14</u> :l,v:z <sub>6</sub> ) with Gdansk to form Gdansk 6,7, <u>14</u> :l,v:z <sub>6</sub> .
I	N	Gege	30	r	1,5	
I	C <sub>1</sub>	Gelsenkirchen	6,7, <u>14</u>	l,v	z <sub>6</sub>	IP combined Gelsenkirchen with Gdansk (6,7:l,v:z <sub>6</sub> ) to form Gdansk 6,7, <u>14</u> :l,v:z <sub>6</sub> . Gelsenkirchen is now called Gdansk var. O 14+ by IP.
I	C <sub>1</sub>	Georgia	6,7	b	e,n,z <sub>15</sub>	
I	T	Gera	<u>1,42</u>	z <sub>4</sub> ,z <sub>23</sub>	[1,6]	
I	D <sub>2</sub>	Geraldton	9,46	l,v	1,6	
II	C <sub>2</sub>	Germiston	6,8	m,t	e,n,x	
I	L	Ghana	21	b	1,6	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	N	Giessen	30	g,m,s	–	
II	C <sub>1</sub>	Gilbert	6,7	z <sub>39</sub>	1,5,7	
I	E <sub>1</sub>	Give	3,10	l,v	1,7	IP combined Newbrunswick (3,15:l,v:1,7) and Menhaden (3,15,34:l,v:1,7) with Give to form Give 3,10,[15],[15,34]:[d],l,v:[d],1,7. Give may possess H phase d; Rl,z <sub>40</sub> ; or Rz <sub>77</sub> .
I	C <sub>2</sub>	Giza	8,20	y	1,2	
I	I	Glasgow	16	b	1,6	
II	F	Glencairn	11	a	z <sub>6</sub> ;z <sub>42</sub>	
I	F	Glidji	11	l,w	1,5	
I	C <sub>2</sub>	Glostrup	6,8	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	B	Gloucester	1,4,12,27	i	1,w	
I	E <sub>4</sub>	Gnesta	1,3,19	b	1,5	
I	N	Godesberg	30	g,m,[t]	–	
I	E <sub>1</sub>	Goelzau	3,10	a	1,5	IP combined Clichy (3,15:a:1,5) with Goelzau to form Goelzau 3,10,[15]:a:1,5.
I	E <sub>2</sub>	Goerlitz	3,15	e,h	1,2	IP combined Goerlitz with Vejle (3,10:e,h:1,2) to form Vejle 3,10,15:e,h:1,2. Goerlitz is now called Vejle var. O 15+ by IP.
I	D <sub>1</sub>	Goeteborg	9,12	c	1,5	
I	D <sub>1</sub>	Goettingen	9,12	l,v	e,n,z <sub>15</sub>	
II	G	Gojenberg	1,13,23	g,t	1,5	
I	51	Gokul	1,51	d	[1,5]	
I	C <sub>2</sub>	Goldcoast	6,8	r	1,w	
I	C <sub>1</sub>	Goma	6,7	z <sub>4</sub> ;z <sub>23</sub>	z <sub>6</sub>	
I	C <sub>1</sub>	Gombe	6,7,14	d	e,n,z <sub>15</sub>	
I	M	Good	21	f,g	e,n,x	
II	G	Goodwood	13,22	z <sub>29</sub>	e,n,x	
I	J	Gori	17	z	1,2	
I	R	Goulfey	1,40	k	1,5	
I	O	Gouloumbo	35	c	1,5	
I	D <sub>1</sub>	Goverdhan	9,12	k	1,6	
I	M	Gozo	28	e,h	e,n,z <sub>15</sub>	
II	F	Grabouw	11	g,[m],s,t	[z <sub>39</sub> ]	
I	C <sub>1</sub>	Grampian	6,7	r	1,w	
I	I	Grancanaria	16	z <sub>39</sub>	[1,6]	Grancanaria can be d-tartrate neg., dulcitol neg., ONPG pos., and anaerogenic.
I	N	Grandhaven	30	r	1,2	
I	J	Granlo	17	l,z <sub>28</sub>	e,n,x	
I	U	Graz	43	a	1,2	
II	Z	Greenside	50	z	e,n,x	
I	R	Greiz	40	a	z <sub>6</sub>	
I	K	Groenekan	18	d	1,5	
		Group A	1,2,12			
		Group B	4,12; 1,4,5,12; or 1,4,12,27			
		Group C1	6,7,[Vi] or 6,7,14			
		Group C2	6,8			IP combined C <sub>2</sub> and C <sub>3</sub> .
		Group C3	8; or 8,20			IP combined with C <sub>2</sub> .
		Group D1	1,9,12			
		Group D2	9,46			

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
		Group D3	1,9,12,46,27			
		Group E1	3,10			
		Group E2	3,15			IP combined E <sub>2</sub> and E <sub>3</sub> with E <sub>1</sub> .
		Group E3	3,15,34			IP combined E <sub>2</sub> and E <sub>3</sub> with E <sub>1</sub> .
		Group E4	1,3,19			
		Group F	11			
		Group G	13,22 or 13,23			
		Group H	6,14; 6,14,24; or 1,6,14,25			
		Group I	16			
		Group J	17			
		Group K	18			
		Group L	21			
		Group M	28			
		Group N	30			
		Group O	35			
		Group P	38			
		Group Q	39			
		Group R	40			
		Group S	41			
		Group T	42			
		Group U	43			
		Group V	44			
		Group W	45			
		Group X	47			
		Group Y	48			
		Group Z	50			
		Group 51	51			
		Group 52	52			
		Group 53	53			
		Group 54	54			
		Group 55	55			
		Group 56	56			
		Group 57	57			
		Group 58	58			
		Group 59	59			
		Group 60	60			
		Group 61	61			
		Group 62	62			
		Group 63	63			
		Group 65	65			
		Group 66	66			
		Group 67	67			
I	G	Grumpensis	1,13,23	d	1,7	
II	R	Grunty	1,40	z <sub>39</sub>	1,6	
I	N	Guarapiranga	30	a	e,n,x	
I	D <sub>2</sub>	Guerin	9,46	e,h	z <sub>6</sub>	
I	M	Guildford	28	k	1,2	
I	V	Guinea	1,44	z <sub>10</sub>	[1,7]	
I	F	Gustavia	11	d	1,5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	L	Gwaai	21	z <sub>4</sub> ,z <sub>24</sub>	–	
I	T	Gwale	<u>1</u> ,42	k	z <sub>6</sub>	
I	E <sub>4</sub>	Gwoza	1,3,19	a	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Haardt	8	k	1,5	
II	D <sub>2</sub>	Haarlem	9,46	z	e,n,x	
I	C <sub>2</sub>	Hadar	6,8	z <sub>10</sub>	e,n,x	
II	I	Haddon	16	z <sub>4</sub> ,z <sub>23</sub>	–	
I	J	Hadejia	17	y	e,n,z <sub>15</sub>	
I	B	Haduna	4,12	l,z <sub>13</sub> ,[z <sub>28</sub> ]	1,6	
I	C <sub>1</sub>	Haelsingborg	6,7	m,p,t,[u]	–	
I	T	Haferbreite	42	k	[1,6]	
I	O	Haga	35	z <sub>38</sub>	–	
II	Y	Hagenbeck	48	d	z <sub>6</sub>	
I	B	Haifa	<u>1</u> ,4,[5],12	z <sub>10</sub>	1,2	
I	M	Halle	28	c	1,7	
I	B	Hallfold	<u>1</u> ,4,12,27	c	l,w	
I	E <sub>2</sub>	Halmstad	3, <u>15</u>	g,s,t	–	IP combined Halmstad and Canoga (3, <u>15</u> ,34:g,s,t:-) with Westhampton (3,10:g,s,t:-) to form Westhampton 3,10,[ <u>15</u> ],[ <u>15</u> ,34]:g,s,t:-. Halmstad is now called Westhampton var. O 15+ by IP.
II	D <sub>1</sub>	Hamburg	<u>1</u> ,9,12	g,t	–	IP combined Hamburg with Manica ( <u>1</u> ,9,12:g,m,s,t:z <sub>42</sub> ) and Muizenberg (9,12:g,m,s,t:1,5) to form S. II <u>1</u> ,9,12:g,m,[s],t:[1,5,7]:[z <sub>42</sub> ].
I	E <sub>2</sub>	Hamilton	3, <u>15</u>	Rz <sub>27</sub>	–	IP combined Hamilton with Goerlitz (3, <u>15</u> :e,h:1,2) and Vejle (3,10:e,h:1,2) to form Vejle 3,10, <u>15</u> :e,h:1,2:Rz <sub>27</sub> . Hamilton is now called Vejle var. Rz <sub>27</sub> +. The name Hamilton has been dropped.
II	Y	Hammonia	48	e,n,x,z <sub>15</sub>	z <sub>6</sub>	
I	G	Handen	<u>1</u> ,13,23	d	1,2	
I	R	Hann	40	k	e,n,x	
I	I	Hannover	16	a	1,2	
I	G	Haouaria	13,22	c	e,n,x,z <sub>15</sub>	
I	H	Harburg	[1],6,14,[25]	k	1,5	
I	51	Harcourt	51	l,v	1,2	
I	E <sub>1</sub>	Harleystreet	3,10	z	1,6	
IV	51	Harmelen	51	z <sub>4</sub> ,z <sub>23</sub>	–	
I	E <sub>1</sub>	Harrisonburg	3,10,[ <u>15</u> ],[ <u>15</u> ,34]	z <sub>10</sub>	1,6	
I	C <sub>1</sub>	Hartford	6,7	y	e,n,x	Hartford may possess H phase Rz <sub>50</sub> or Rz <sub>67</sub> .
I	T	Harvestehude	<u>1</u> ,42	y	z <sub>6</sub>	
I	M	Hatfield	28	d	1,6	
I	B	Hato	4,[5],12	g,m,s	–	
I	G	Havana	<u>1</u> ,13,23	f,g,[s]	–	Havana may possess H phase Rz <sub>45</sub> or Rz <sub>79</sub> .
I	E <sub>4</sub>	Hayindogo	1,3,19	e,h	1,6	
I	F	Heerlen	11	i	1,6	
I	Q	Hegau	39	z <sub>10</sub>	–	
I	B	Heidelberg	<u>1</u> ,4,[5],12	r	1,2	
II	C <sub>1</sub>	Heilbron	6,7	l,z <sub>28</sub>	1,5:[z <sub>42</sub> ]	
II	B	Helsinki	<u>1</u> ,4,12	z <sub>29</sub>	[e,n,x]	
I	Z	Hemingford	50	d	1,5	Hemingford may possess H phase Rz <sub>82</sub> .

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	S	Hennepin	41	d	Z <sub>6</sub>	
I	M	Hermannswerder	28	c	1,5	
I	I	Heron	16	a	Z <sub>6</sub>	
I	C <sub>2</sub>	Herston	6,8	d	e,n,Z <sub>15</sub>	
I	F	Herzliya	11	y	e,n,x	
I	B	Hessarek	4,12,27	a	1,5	
I	H	Heves	6,14,[24]	d	1,5	
I	C <sub>2</sub>	Hidalgo	6,8	r,i	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Hiduddify	6,8	1,Z <sub>13</sub> ,Z <sub>28</sub>	1,5	
II	J	Hillbrow	17	b	e,n,x,Z <sub>15</sub>	
I	D <sub>2</sub>	Hillegersberg	9,46	Z <sub>35</sub>	1,5	
I	D <sub>2</sub>	Hillingdon	9,46	g,m	–	
I	C <sub>1</sub>	Hillsborough	6,7	Z <sub>41</sub>	l,w	
I	N	Hilversum	30	k	1,2	
I	C <sub>2</sub>	Hindmarsh	8,20	r	1,5	
I	Y	Hisingen	48	a	1,5,7	
I	C <sub>1</sub>	Hissar	6,7,14	c	1,2	
I	I	Hithergreen	16	c	e,n,Z <sub>15</sub>	
I	E <sub>1</sub>	Hoboken	3,10	i	l,w	
I	Q	Hofit	39	i	1,5	
I	E <sub>1</sub>	Hoghton	3,10	1,Z <sub>13</sub> ,Z <sub>28</sub>	Z <sub>6</sub>	
I	C <sub>2</sub>	Holcomb	6,8	l,v	e,n,x	
I	H	Homosassa	1,6,14,25	z	1,5	
I	M	Honelis	28	a	e,n,Z <sub>15</sub>	
I	E <sub>4</sub>	Hongkong	1,3,19	z	Z <sub>6</sub>	
II	Z	Hooggraven	50	Z <sub>10</sub>	Z <sub>6</sub> :Z <sub>42</sub>	
I	H	Horsham	1,6,14,[25]	l,v	e,n,x	
IV	U	Houten	43	Z <sub>4</sub> ,Z <sub>23</sub>	–	
I	E <sub>1</sub>	Huddinge	3,10	z	1,7	
II	D <sub>1</sub>	Hueningen	9,12	z	Z <sub>39</sub>	
I	B	Huettwilen	1,4,12	a	l,w	
II	F	Huila	11	1,Z <sub>28</sub>	e,n,x	
I	I	Hull	16	b	1,2	
II	53	Humber	53	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	E <sub>1</sub>	Huvudsta	3,10	b	1,7	
I	I	Hvittingfoss	16	b	e,n,x	
I	L	Hydra	21	c	1,6	
I	G	Ibadan	13,22	b	1,5	
I	L	Ibaragi	21	y	1,2	
I	G	Idikan	1,13,23	i	1,5	
I	E <sub>1</sub>	Ikayi	3,10,[15]	c	1,6	Ikayi Var. O 15+ was described after E <sub>1</sub> and E <sub>2</sub> were combined.
I	M	Ikeja	28	k	1,7	
I	M	Ilala	28	k	1,5	
I	E <sub>3</sub>	Illinois	3,15,34	Z <sub>10</sub>	1,5	IP combined Illinois and Manila (3,15:z <sub>10</sub> :1,5) with Lexington (3,10:z <sub>10</sub> :1,5) to form Lexington 3,10,[15],[15,34]:z <sub>10</sub> :1,5. Illinois is now called Lexington var. O 15+, 34+ by IP.
I	E <sub>4</sub>	Ilugun	1,3,10,19	Z <sub>4</sub> ,Z <sub>23</sub>	Z <sub>6</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	W	Imo	45	l,v	[e,n,z <sub>15</sub> ]	
I	C <sub>2</sub>	Inchpark	6,8	y	1,7	
I	D <sub>2</sub>	India	9,46	l,v	1,5	
I	B	Indiana	<u>1,4,12</u>	z	1,7	
I	C <sub>1</sub>	Infantis	6,7, <u>14</u>	r	1,5	Infantis may possess H phase Rz <sub>49</sub> .
I	C <sub>1</sub>	Inganda	6,7	z <sub>10</sub>	1,5	
I	D <sub>2</sub>	Inglis	9,46	z <sub>10</sub>	e,n,x	
I	S	Inpraw	41	z <sub>10</sub>	e,n,x	
I	P	Inverness	38	k	1,6	
I	D <sub>1</sub>	Ipeko	9,12	c	1,6	
I	S	Ipswich	41	z <sub>4</sub> ,z <sub>24</sub>	[1,5]	
I	D <sub>2</sub>	Irchel	9,46	y	e,n,x	
I	J	Irenea	17	k	1,5	
I	U	Irigny	43	z <sub>38</sub>	–	
I	C <sub>1</sub>	Irumu	6,7	l,v	1,5	
I	C <sub>1</sub>	Isangi	6,7, <u>14</u>	d	1,5	
I	Y	Isaszeg	48	z <sub>10</sub>	e,n,x	
II	E <sub>1</sub>	Islington	3,10	g,t	–	
I	D <sub>1</sub>	Israel	9,12	e,h	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Istanbul	8	z <sub>10</sub>	e,n,x	
I	H	Istoria	1,6,14,25	r,i	1,5	
I	G	Isuge	13,23	d	z <sub>6</sub>	
I	D <sub>1</sub>	Italiana	9,12	l,v	R1,11	IP combined Italiana that contains H phase R1,11 with Panama (1,9,12:l,v:1,5). The name Italiana has been dropped.
I	D <sub>1</sub>	Itami	9,12	l,z <sub>13</sub>	1,5	
I	B	Ituri	<u>1,4,12</u>	z <sub>10</sub>	1,5	
I	D <sub>2</sub>	Itutaba	9,46	c	z <sub>6</sub>	
I	I	Ivory	16	r	1,6	
I	Z	Ivorycoast	50	z <sub>29</sub>	–	
II	I	Jacksonville	16	z <sub>29</sub>	[e,n,x]	
I	D <sub>1</sub>	Jaffna	<u>1,9,12</u>	d	z <sub>35</sub>	
I	B	Jaja	4,12, <u>27</u>	z <sub>4</sub> ,z <sub>23</sub>	–	IP combined Jaja with Stanleyville (1,4,[5],12:z <sub>4</sub> ,z <sub>23</sub> :1,5) to form Stanleyville 1,4,[5],12, <u>27</u> :z <sub>4</sub> ,z <sub>23</sub> :1,5. Jaja is now called Stanleyville var. O 27+. The name Jaja has been dropped.
I	F	Jalisco	11	y	1,7	
I	D <sub>1</sub>	Jamaica	9,12	r	1,5	
I	L	Jambur	21	l,z <sub>28</sub>	e,n,z <sub>15</sub>	
I	J	Jangwani	17	a	1,5	
I	B	Java	<u>1,4,5,12</u>	b	[1,2], (tartrate +)	IP calls Java, Paratyphi B var. Java. Java is often monophasic in the U.S. May possess H phase Rz <sub>33</sub> .
I	D <sub>1</sub>	Javiana	<u>1,9,12</u>	l,z <sub>28</sub>	1,5	
I	E <sub>1</sub>	Jedburgh	3,10, <u>15</u>	z <sub>29</sub>	–	
I	B	Jericho	<u>1,4,12,27</u>	c	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Jerusalem	6,7, <u>14</u>	z <sub>10</sub>	l,w	
I	E <sub>1</sub>	Joal	3,10	l,z <sub>28</sub>	1,7	
I	W	Jodhpur	45	z <sub>29</sub>	–	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	B	Joenkoeping	4,5,12	g,s,t	–	IP combined Joenkoeping with Kingston (1,4,12,27:g,s,t:-) to form Kingston 1,4,[5],12,27:g,s,t:-. The name Joenkoeping has been dropped.
I	R	Johannesburg	1,40	b	e,n,x	
I	B	Jos	1,4,12,27	y	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Juba	1,3,19	a	1,7	
I	J	Jubilee	17	e,h	1,2	
I	G	Jukestown	13,23	i	e,n,z <sub>15</sub>	
I	B	Kaapstad	4,12	e,h	1,7	
I	51	Kabete	51	i	1,5	
I	C <sub>1</sub>	Kaduna	6,7,14	c	e,n,z <sub>15</sub>	
I	T	Kahla	1,42	z <sub>35</sub>	1,6	
I	E <sub>4</sub>	Kainji	1,3,19	z	1,6	
I	H	Kaitaan	1,6,14,25	m,t	–	
I	B	Kalamu	4,[5],12	z <sub>4</sub> ,z <sub>24</sub>	[1,5]	
I	E <sub>1</sub>	Kalina	3,10	b	1,2	
I	C <sub>2</sub>	Kallo	6,8	k	1,2	
II	M	Kaltenhausen	28	b	z <sub>6</sub>	
I	C <sub>2</sub>	Kalumburu	6,8	z	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Kambole	6,7	d	1,[2],7	
I	B	Kamoru	4,12,27	y	z <sub>6</sub>	
I	T	Kampala	1,42	c	z <sub>6</sub>	
I	E <sub>4</sub>	Kande	1,3,19	b	e,n,z <sub>15</sub>	
I	J	Kandla	17	z <sub>29</sub>	–	
I	T	Kaneshie	1,42	i	l,w	
I	H	Kanifing	1,6,14,25	z	1,6	
I	B	Kano	1,4,12,27	l,z <sub>13</sub> ,z <sub>28</sub>	e,n,x	
I	X	Kaolack	47	z	1,6	
I	D <sub>1</sub>	Kapemba	9,12	l,v	1,7	
I	W	Karachi	45	d	e,n,x	
I	R	Karamoja	40	z <sub>41</sub>	1,2	
I	51	Karaya	51	b	1,5	
I	P	Kasenyi	38	e,h	1,5	
I	H	Kassberg	1,6,14,25	c	1,6	
II	G	Katesgrove	1,13,23	m,t	1,5	
I	G	Kedougou	1,13,23	i	l,w	
I	C <sub>2</sub>	Kentucky	8,20	i	z <sub>6</sub>	
I	C <sub>1</sub>	Kenya	6,7	l,z <sub>13</sub>	e,n,x	
I	V	Kermel	44	d	e,n,x	
I	L	Keve	21	l,w	–	
II	X	Khami	47	b	[e,n,x,z <sub>15</sub> ]	
I	E <sub>3</sub>	Khartoum	3,15,34	a	1,7	IP combined Khartoum with Oxford (3,10:a:1,7) to form Oxford 3,10,[15],[15,34]:a:1,7. Khartoum is now called Oxford var. O 15+ by IP. CDC has no 3,15:a:1,7. Khartoum was found by IP with colonies containing O 3,15.
I	B	Kiambu	4,12	z	1,5	
I	I	Kibi	16	z <sub>4</sub> ,z <sub>23</sub>	[1,6]	
I	M	Kibusi	28	r	e,n,x	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	P	Kidderminster	38	c	1,6	
I	A	Kiel	<u>1,2</u> ,12	g,p	–	
I	I	Kikoma	16	y	e,n,x	
II	B	Kilwa	4,12	l,w	e,n,x	
I	P	Kimberley	38	l,v	1,5	
I	D <sub>1</sub>	Kimpese	9,12	z	1,6	
I	B	Kimuenza	<u>1,4,12,27</u>	l,v	e,n,x	
I	E <sub>4</sub>	Kindia	1,3,19	l,z <sub>28</sub>	e,n,x	
I	U	Kingabwa	43	y	1,5	
I	B	Kingston	<u>1,4,12,27</u>	g,s,t	[1,2]	IP combined Kingston with Joenkoeping (4,5,12:g,s,t:-) to form Kingston <u>1,4,[5],12,27</u> :g,s,t:[1,2]. Kingston may possess H phase Rz <sub>27</sub> or Rz <sub>43</sub> .
I	J	Kinondoni	17	a	e,n,x	
I	E <sub>2</sub>	Kinshasa	3,15	l,z <sub>13</sub>	1,5	IP combined Kinshasa with Uganda (3,10:l,z <sub>13</sub> :1,5) to form Uganda 3,10, <u>[15]</u> :l,z <sub>13</sub> :1,5. Kinshasa is now called Uganda var. O 15+ by IP.
I	E <sub>4</sub>	Kinson	1,3,19	y	e,n,x	
I	G	Kintambo	<u>1,13,23</u>	m,t	–	
I	J	Kirkee	17	b	1,2	
I	B	Kisangani	<u>1,4,[5],12</u>	a	1,2	
I	F	Kisarawe	11	k	e,n,x,[z <sub>15</sub> ]	
I	C <sub>1</sub>	Kisii	6,7	d	1,2	
I	M	Kitenge	28	y	e,n,x	
I	C <sub>1</sub>	Kivu	6,7	d	1,6	
II	W	Klapmuts	45	z	z <sub>39</sub>	
I	P	Klouto	38	z <sub>38</sub>	–	
II	B	Kluetjenfelde	4,12	d	e,n,x	
I	X	Kodjovi	47	c	[1,6]	Kodjovi may possess H phase Rz <sub>78</sub> .
I	B	Koenigstuhl	<u>1,4,[5],12</u>	z	e,n,z <sub>15</sub>	
I	A	Koessen	2,12	l,v	1,5	
I	W	Kofandoka	45	r	e,n,z <sub>15</sub>	
I	V	Koketime	44	z <sub>38</sub>	–	
I	N	Kokoli	30	z <sub>35</sub>	1,6	
I	Q	Kokomlele	39	l,v	e,n,x	
I	D <sub>2</sub>	Kolar	9,46	b	z <sub>35</sub>	
I	C <sub>2</sub>	Kolda	8, <u>20</u>	z <sub>35</sub>	1,2	
II	U	Kommetje	43	b	z <sub>42</sub>	
I	M	Konolfingen	28	z <sub>35</sub>	1,6	
I	C <sub>2</sub>	Konstanz	8	b	e,n,x	
I	C <sub>2</sub>	Korbol	8, <u>20</u>	b	1,5	
I	E <sub>4</sub>	Korlebu	1,3,19	z	1,5	
I	P	Korovi	38	g,m,[s]	–	
I	C <sub>1</sub>	Kortrijk	6,7	l,v	1,7	
I	C <sub>2</sub>	Kottbus	6,8	e,h	1,5	
I	C <sub>1</sub>	Kotte	6,7	b	z <sub>35</sub>	
I	D <sub>1</sub>	Kotu	9,12	l,z <sub>28</sub>	1,6	
I	E <sub>4</sub>	Kouka	1,3,19	g,m,[t]	–	
I	C <sub>1</sub>	Koumra	6,7	b	1,7	
I	M	Kpeme	28	e,h	1,7	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	G	Kraaifontein	<u>1</u> ,13,23	g,m,t	[e,n,x]	IP combined Kraaifontein with Luanshya ( <u>1</u> ,13,23:g,m,s,t:[e,n,x]) to form Luanshya <u>1</u> ,13,23:g,m,[s],t:[e,n,x]. The name Kraaifontein has been dropped.
IV	C <sub>1</sub>	Kralendyk	6,7	z <sub>4</sub> ,z <sub>24</sub>	–	
I	C <sub>2</sub>	Kralingen	8,20	y	z <sub>6</sub>	
I	E <sub>4</sub>	Krefeld	1,3,19	y	l,w	
I	E <sub>1</sub>	Kristianstad	3,10	z <sub>10</sub>	e,n,z <sub>15</sub>	
II	Z	Krugersdorp	50	e,n,x	1,7	
I	V	Kua	44	z <sub>4</sub> ,z <sub>23</sub>	–	
I	B	Kubacha	<u>1</u> ,4,12,27	l,z <sub>13</sub> ,z <sub>28</sub>	1,7	
I	M	Kuessel	28	i	e,n,z <sub>15</sub>	
II	D <sub>1</sub>	Kuilsrivier	<u>1</u> ,9,12	g,m,s,t	e,n,x	
I	N	Kumasi	30	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	B	Kunduchi	<u>1</u> ,4,[5],12,27	l,[z13],[z28]	[1,2]	
I	H	Kuntair	1,6,14,25	b	1,5	
I	C <sub>2</sub>	Kuru	6,8	z	l,w	
I	C <sub>2</sub>	Labadi	8,20	d	z <sub>6</sub>	
I	B	Lagos	<u>1</u> ,4,[5],12	i	1,5	
I	E <sub>1</sub>	Lamberhurst	3,10	e,h	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Lamin	3,10	l,z <sub>28</sub>	e,n,x	
I	J	Lancaster	17	l,v	1,7	
I	S	Landala	41	z <sub>10</sub>	1,6	
I	N	Landau	30	i	1,2	
I	E <sub>1</sub>	Landwasser	3,10	z	z <sub>6</sub>	
I	K	Langenhorn	18	m,t	–	
I	E <sub>1</sub>	Langensalza	3,10	y	l,w	
I	M	Langford	28	b	e,n,z <sub>15</sub>	
I	E <sub>2</sub>	Lanka	3,15	r	z <sub>6</sub>	IP combined Lanka with Weltevreden (3,10:r:z <sub>6</sub> ) to form Weltevreden 3,10,[15]:r:z <sub>6</sub> . Lanka is now called Weltevreden var. O 15+ by IP.
I	P	Lansing	38	i	1,5	
I	C <sub>1</sub>	Larochelle	6,7	e,h	1,2	
I	W	Lattenkamp	45	z <sub>35</sub>	1,5	
I	D <sub>1</sub>	Lawndale	<u>1</u> ,9,12	z	1,5	
I	V	Lawra	44	k	e,n,z <sub>15</sub>	
I	S	Leatherhead	41	m,t	1,6	
I	51	Lechler	51	z	e,n,z <sub>15</sub>	
I	K	Leer	18	z <sub>10</sub>	1,5	
I	F	Leeuwarden	11	b	1,5	
I	B	Legon	<u>1</u> ,4,12,27	c	1,5	
I	G	Leiden	13,22	z <sub>38</sub>	–	
I	S	Leipzig	41	z <sub>10</sub>	1,5	
I	C <sub>2</sub>	Leith	6,8	a	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Lekke	3,10	d	1,6	
I	F	Lene	11	z <sub>38</sub>	–	
I	M	Leoben	28	l,v	1,5	
I	C <sub>1</sub>	Leopoldville	6,7,14	b	z <sub>6</sub>	
I	E <sub>4</sub>	Lerum	1,3,19	z	1,7	
II	S	Lethe	41	g,t	–	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>1</sub>	Lexington	3,10	z <sub>10</sub>	1,5	IP combined Lexington with Manila (3,15:z <sub>10</sub> :1,5) and Illinois (3,15,34:z <sub>10</sub> :1,5) to form Lexington 3,10,[15],[15,34]:z <sub>10</sub> :1,5. Lexington may possess H phase Rz <sub>49</sub> .
I	C <sub>2</sub>	Lezennes	6,8	z <sub>4</sub> ,z <sub>23</sub>	1,7	
I	M	Libreville	28	z <sub>10</sub>	1,6	
II	S	Lichtenberg	41	z <sub>10</sub>	[z <sub>6</sub> ]	
I	N	Ligeo	30	l,v	1,2	
I	O	Ligna	35	z <sub>10</sub>	z <sub>6</sub>	
I	C <sub>1</sub>	Lika	6,7	i	1,7	
I	C <sub>1</sub>	Lille	6,7	z <sub>38</sub>	–	IP combined Lille with Bornum (6,7,14:z <sub>38</sub> :–) to form Lille 6,7,14:z <sub>38</sub> :–.
II	G	Limbe	1,13,22	g,m,t	[1,5]	
I	B	Limete	1,4,12,27	b	1,5	
II	F	Lincoln	11	m,t	e,n,x	
I	C <sub>2</sub>	Lindenburg	6,8	i	1,2	
I	H	Lindern	6,14,[24]	d	e,n,x	
I	P	Lindi	38	r	1,5	
II	D <sub>1</sub>	Lindrick	9,12	e,n,x	1,[5],7	
I	D <sub>2</sub>	Linguere	9,46	b	z <sub>6</sub>	
I	I	Lingwala	16	z	1,7	
I	G	Linton	13,23	r	e,n,z <sub>15</sub>	
I	I	Lisboa	16	z <sub>10</sub>	1,6	
I	D <sub>2</sub>	Lishabi	9,46	z <sub>10</sub>	1,7	
I	C <sub>2</sub>	Litchfield	6,8	l,v	1,2	
I	E <sub>4</sub>	Liverpool	1,3,19	d	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Livingstone	6,7	d	1,w	IP combined Eimsbuettel (6,7,14:d:l,w) with Livingstone to form Livingstone 6,7,14:d:l,w.
I	N	Livulu	30	e,h	1,2	
I	B	Ljubljana	4,12,27	k	e,n,x	
I	E <sub>4</sub>	Llandoff	1,3,19	z <sub>29</sub>	[z <sub>6</sub> ]	
II	M	Llandudno	28	g,[m],[s],t	1,5	
I	V	Llobregat	1,44	z <sub>10</sub>	e,n,x	
I	C <sub>2</sub>	Loanda	6,8	l,v	1,5	
II	52	Lobatsi	52	z <sub>44</sub>	1,5,7	
II	57	Locarno	57	z <sub>29</sub>	z <sub>42</sub>	
I	C <sub>1</sub>	Lockleaze	6,7,14	b	e,n,x	
I	J	Lode	17	r	1,2	
I	S	Lodz	41	z <sub>29</sub>	–	
I	T	Loenga	1,42	z <sub>10</sub>	z <sub>6</sub>	
I	Q	Logone	39	d	1,5	
IV	V	Lohbruegge	44	z <sub>4</sub> ,z <sub>32</sub>	–	
I	E <sub>4</sub>	Lokstedt	1,3,19	1,z <sub>13</sub> ,z <sub>28</sub>	1,2	
I	D <sub>1</sub>	Lomalinda	1,9,12	a	e,n,x	
I	D <sub>1</sub>	Lome	9,12	r	z <sub>6</sub>	
I	C <sub>1</sub>	Lomita	6,7	e,h	1,5	
I	I	Lomnava	16	l,w	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	London	3,10	l,v	1,6	IP combined London with Portsmouth (3,15:l,v:1,6) to form London 3,10,[15]:l,v:1,6.
I	I	Losangeles	16	l,v	z <sub>6</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	B	Loubomo	4,12	z	1,6	
I	N	Louga	30	b	1,2	
I	D <sub>2</sub>	Louisiana	9,46	Z <sub>10</sub>	Z <sub>6</sub>	
II	I	Louwbester	16	z	[e,n,x]	
I	G	Lovelace	13,22	l,v	1,5	
I	J	Lowestoft	17	g,s,t	–	
II	G	Luanshya	<u>1,13,23</u>	g,m,s,t	[e,n,x]	IP combined Luanshya with Kraifontein ( <u>1,13,23</u> :g,m,t:[e,n,x]) to form Luanshya <u>1,13,23</u> :g,m,[s],t:[e,n,x].
I	S	Lubumbashi	41	r	1,5	
I	F	Luciana	11	a	e,n,Z <sub>15</sub>	
I	M	Luckenwalde	28	Z <sub>10</sub>	e,n,Z <sub>15</sub>	
I	X	Luke	1,47	g,m	–	
II	D <sub>2</sub>	Lundby	9,46	b	e,n,x	
II	S	Lurup	41	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	
I	51	Lutetia	51	r,i	1,Z <sub>13</sub> ,Z <sub>28</sub>	
II	60	Luton	60	z	e,n,x	
I	X	Lyon	47	k	e,n,Z <sub>15</sub>	
II	D <sub>2</sub>	Maarssen	9,46	Z <sub>4</sub> ,Z <sub>24</sub>	Z <sub>39</sub> :Z <sub>42</sub>	
I	F	Maastricht	11	Z <sub>41</sub>	1,2	
I	E <sub>1</sub>	Macallen	3,10	Z <sub>36</sub>	–	
I	D <sub>2</sub>	Macclesfield	9,46	g,m,s,t	1,(2),7	
I	E <sub>4</sub>	Machaga	1,3,19	i	e,n,x	
I	H	Madelia	1,6,14,25	y	1,7	
I	E <sub>4</sub>	Madiago	1,3,19	c	1,7	
I	V	Madigan	44	c	1,5	
I	L	Madison	21	d	Z <sub>6</sub>	
I	E <sub>1</sub>	Madjorio	3,10	d	e,n,Z <sub>15</sub>	
I	B	Madras	4,[5],12	m,t	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Magherafelt	8, <u>20</u>	i	1,w	
I	H	Magumeri	1,6,14,25	e,h	1,6	
I	L	Magwa	21	d	e,n,x	
I	D <sub>2</sub>	Mahina	9,46	Z <sub>10</sub>	e,n,Z <sub>15</sub>	
I	E <sub>4</sub>	Maiduguri	1,3,19	f,g,t	e,n,Z <sub>15</sub>	
I	U	Makiling	43	Z <sub>29</sub>	–	
I	C <sub>1</sub>	Makiso	6,7	1,Z <sub>13</sub> ,Z <sub>28</sub>	Z <sub>6</sub>	
II	B	Makoma	<u>1,4,[5],12,27</u>	a	[e,n,x]	
II	B	Makumira	<u>1,4,12,27</u>	e,n,x	1,[5],7	
I	I	Malakal	16	e,h	1,2	
V	66	Malawi	66	Z <sub>65</sub>	–	
I	M	Malaysia	28	Z <sub>10</sub>	1,7	
I	V	Malika	44	1,Z <sub>28</sub>	1,5	
I	C <sub>2</sub>	Malmoe	6,8	i	1,7	
I	I	Malstatt	16	b	Z <sub>6</sub>	
I	H	Mampeza	1,6,14,25	i	1,5	
I	G	Mampong	13,22	Z <sub>35</sub>	1,6	
I	C <sub>2</sub>	Manchester	6,8	l,v	1,7	
I	I	Mandera	16	1,Z <sub>13</sub>	e,n,Z <sub>15</sub>	
I	P	Mango	38	k	1,5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	C <sub>2</sub>	Manhattan	6,8	d	1,5	
II	D <sub>1</sub>	Manica	<u>1,9,12</u>	g,m,s,t	z <sub>42</sub>	IP combined Manica with Hamburg ( <u>1,9,12:g,t:-</u> ) and Muizenberg (9,12:g,m,s,t:1,5) to form S. II <u>1,9,12:g,[m],[s],t:[1,5,7]:[z<sub>42</sub>]</u> .
I	E <sub>2</sub>	Manila	3, <u>15</u>	z <sub>10</sub>	1,5	IP combined Manila and Illinois (3, <u>15,34</u> :z <sub>10</sub> :1,5) with Lexington (3,10:z <sub>10</sub> :1,5) to form Lexington 3,10, <u>[15],[15,34]</u> :z <sub>10</sub> :1,5. Manila is now called Lexington var. O 15+ by IP.
I	F	Mannheim	11	k	1,w	
II	57	Manombo	57	z <sub>39</sub>	e,n,x,z <sub>15</sub>	
I	C <sub>2</sub>	Mapo	6,8	z <sub>10</sub>	1,5	
I	Q	Mara	39	e,h	[1,5]	
I	F	Maracaibo	11	l,v	1,5	
I	G	Marburg	13,23	k	-	
V	66	Maregrosso	66	z <sub>35</sub>	-	
I	T	Maricopa	<u>1,42</u>	g,z <sub>51</sub>	1,5	
I	E <sub>1</sub>	Marienthal	3,10	k	e,n,z <sub>15</sub>	
IV	Y	Marina	48	g,z <sub>51</sub>	-	
I	V	Maritzburg	<u>1,44</u>	i	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Maron	3,10	d	z <sub>35</sub>	
I	F	Maroua	11	z	1,7	
I	F	Marseille	11	a	1,5	
I	G	Marshall	13,22	a	l,z <sub>13</sub> ,z <sub>28</sub>	
I	57	Maryland	57	b	1,7	
I	D <sub>2</sub>	Marylebone	9,46	k	1,2	
I	E <sub>1</sub>	Masembe	3,10	a	e,n,x	Masembe may possess H phase Rz <sub>5</sub> .
I	B	Maska	<u>1,4,12,27</u>	z <sub>41</sub>	e,n,z <sub>15</sub>	
I	O	Massakory	35	r	1,w	
I	B	Massenya	<u>1,4,12,27</u>	k	1,5	
I	J	Matadi	17	k	e,n,x	
I	D <sub>2</sub>	Mathura	9,46	i	e,n,z <sub>15</sub>	
I	N	Matopeni	30	y	1,2	
II	E <sub>1</sub>	Matroosfontein	3,10	a	e,n,x	
I	D <sub>2</sub>	Mayday	9,46	y	z <sub>6</sub>	
I	C <sub>1</sub>	Mbandaka	6,7, <u>14</u>	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	U	Mbao	43	i	1,2	
I	W	Meekatharra	45	a	e,n,z <sub>15</sub>	
I	T	Melbourne	42	z	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Meleagridis	3,10	e,h	1,w	IP combined Meleagridis with Cambridge (3, <u>15</u> :e,h:l,w) and Wildwood (3, <u>15,34</u> :e,h:l,w) to form Meleagridis 3,10, <u>[15],[15,34]</u> :e,h:l,w.
I	K	Memphis	18	k	1,5	
I	C <sub>1</sub>	Menden	6,7	z <sub>10</sub>	1,2	
I	D <sub>1</sub>	Mendoza	9,12	l,v	1,2	
I	E <sub>3</sub>	Menhaden	3, <u>15,34</u>	l,v	1,7	IP combined Menhaden with Give (3,10:l,v:1,7) and Newbrunswick (3, <u>15</u> :l,v:1,7) to form Give 3,10, <u>[15],[15,34]</u> :l,v:1,7. Menhaden is now called Give var. O 15+, 34+ by IP.
I	C <sub>1</sub>	Menston	6,7	g,s,[t]	[1,6]	
II	I	Merseyside	16	g,t	[1,5]	
I	X	Mesbit	47	m,t	[e,n,z <sub>15</sub> ]	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	51	Meskin	51	e,h	1,2	
I	N	Messina	30	d	1,5	
I	C <sub>2</sub>	Mexicana	6,8	d	1,2	Mexicana was combined with Muenchen. The name Mexicana has been dropped.
I	P	Mgulani	38	i	1,2	
I	D <sub>1</sub>	Miami	<u>1</u> ,9,12	a	1,5	Miami must be differentiated from Sendai with biochemical tests. Miami is pos. for H <sub>2</sub> S, citrate, and tartrate; Sendai is neg.
I	J	Michigan	17	l,v	1,5	
I	T	Middlesbrough	<u>1</u> ,42	i	z <sub>6</sub>	
II	53	Midhurst	53	l,z <sub>28</sub>	z <sub>39</sub>	
I	H	Midway	6,14,24	d	1,7	
I	C <sub>1</sub>	Mikawasima	6,7, <u>14</u>	y	e,n,z <sub>15</sub>	Mikawasima may possess H phase Rz <sub>47</sub> or Rz <sub>50</sub> .
I	R	Millesi	<u>1</u> ,40	l,v	1,2	
I	U	Milwaukee	43	f,g,[t]	-	
I	G	Mim	13,22	a	1,6	
I	H	Minna	1,6,14,25	c	1,w	
I	E <sub>3</sub>	Minneapolis	3, <u>15</u> , <u>34</u>	e,h	1,6	IP combined Minneapolis and Newington (3, <u>15</u> :e,h:1,6) with Anatum (3,10:e,h:1,6) to form Anatum 3,10,[ <u>15</u> ],[ <u>15</u> , <u>34</u> ]:e,h:1,6. Minneapolis is now called Anatum var. O 15+ by IP.
I	L	Minnesota	21	b	e,n,x	Minnesota may possess H phase Rz <sub>33</sub> or Rz <sub>49</sub> .
I	G	Mishmarhaemek	<u>1</u> ,13,23	d	1,5	
I	C <sub>1</sub>	Mission	6,7	d	1,5	Mission was combined with Isangi 6,7, <u>14</u> :d:1,5. The name Mission has been dropped.
I	G	Mississippi	<u>1</u> ,13,23	b	[1,5]	
I	F	Missouri	11	g,s,t	-	
I	D <sub>1</sub>	Miyazaki	9,12	l,z <sub>13</sub>	1,7	
II	D <sub>1</sub>	Mjimwema	<u>1</u> ,9,12	b	e,n,x	
I	N	Mjordan	30	i	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Mkamba	6,7	l,v	1,6	
II	I	Mobeni	16	g,[m],[s],t	[e,n,x]	
I	M	Mocamedes	28	d	e,n,x	
I	M	Moero	28	b	1,5	
I	F	Moers	11	m,t	-	
I	E <sub>1</sub>	Mokola	3,10	y	1,7	
I	C <sub>2</sub>	Molade	8, <u>20</u>	z <sub>10</sub>	z <sub>6</sub>	
I	52	Molesey	52	b	1,5	
II	Q	Mondeor	39	l,z <sub>28</sub>	e,n,x	
I	B	Mono	4,12	l,w	1,5	
I	B	Mons	<u>1</u> ,4,12, <u>27</u>	d	l,w	
I	O	Monschau	35	m,t	-	
I	C <sub>1</sub>	Montevideo	6,7, <u>14</u>	g,m,[p],s	[1,2,7]	
II	F	Montgomery	11	a,[d]	[d]:e,n,z <sub>15</sub>	
I	U	Montreal	43	c	1,5	
I	N	Morehead	30	i	1,5	
I	M	Morillons	28	m,t	1,6	
I	N	Morningside	30	c	e,n,z <sub>15</sub>	
I	H	Mornington	1,6,14,25	y	e,n,z <sub>15</sub>	
I	N	Morocco	30	l,z <sub>13</sub> ,z <sub>28</sub>	e,n,z <sub>15</sub>	
I	J	Morotai	17	l,v	1,2	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	M	Moroto	28	z <sub>10</sub>	l,w	
I	D <sub>1</sub>	Moscow	9,12	g,q	–	
II	U	Mosselbay	43	g,m,[s],t	[z <sub>42</sub> ]	
I	X	Moualine	47	y	1,6	
I	51	Moundou	51	l,z <sub>28</sub>	1,5	
I	L	Mountmagnet	21	r	–	
I	X	Mountpleasant	47	z	1,5	
I	H	Moussoro	1,6,14,25	i	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Mowanjum	6,8	z	1,5	
II	E <sub>1</sub>	Mpila	3,10	z <sub>38</sub>	z <sub>42</sub>	
I	I	Mpouto	16	m,t	–	
I	C <sub>2</sub>	Muenchen	6,8	d	1,2	
I	E <sub>1</sub>	Muenster	3,10	e,h	1,5	IP combined Muenster with Newhaw (3,15:e,h:1,5) and Arkansas (3,15,34:e,h:1,5) to form Muenster 3,10,[15],[15,34]:e,h:1,5.
I	V	Muguga	44	m,t	–	
II	D <sub>1</sub>	Muizenberg	9,12	g,m,s,t	1,5	IP combined Muizenberg with Hamburg (1,9,12:g,t:-) and Manica (1,9,12:g,m,s,t:z <sub>42</sub> ) to form S. II 1,9,12:g,[m],[s],t:[1,5,7]:[z <sub>42</sub> ].
I	M	Mundonobo	28	d	1,7	
IV	F	Mundsborg	11	g,z <sub>51</sub>	–	
I	B	Mura	1,4,12	z <sub>10</sub>	l,w	
II	G	Nachshonim	1,13,23	z	1,5	
I	D <sub>1</sub>	Naestved	1,9,12	g,p,s	–	
I	C <sub>2</sub>	Nagoya	6,8	b	1,5	
II	T	Nairobi	42	r	–	
I	B	Nakuru	1,4,12,27	a	z <sub>6</sub>	
II	Z	Namib	50	g,[m],s,t	[1,5]	
I	C <sub>1</sub>	Namibia	6,7	c	e,n,x	
I	E <sub>2</sub>	Nancy	3,15	l,v	1,2	IP combined Nancy with Nchanga (3,10:l,v:1,2) to form Nchanga 3,10,[15]:l,v:1,2. Nancy is now called Nchanga var. O 15+ by IP.
I	C <sub>2</sub>	Nanergou	6,8	g,s,t	–	
I	G	Nanga	1,13,23	l,v	e,n,z <sub>15</sub>	
I	D <sub>2</sub>	Nantes	9,46	y	l,w	
I	D <sub>1</sub>	Napoli	1,9,12	l,z <sub>13</sub>	e,n,x	
I	C <sub>2</sub>	Narashino	6,8	a	e,n,x	
I	M	Nashua	28	l,v	e,n,z <sub>15</sub>	
I	D <sub>1</sub>	Natal	9,12	z <sub>4</sub> ,z <sub>24</sub>	–	
I	I	Naware	16	z <sub>38</sub>	–	
I	E <sub>1</sub>	Nchanga	3,10	l,v	1,2	IP combined Nchanga with Nancy (3,15:l,v:1,2) to form Nchanga 3,10,[15]:l,v:1,2.
I	H	Ndjamena	1,6,14,25	b	1,2	
I	D <sub>1</sub>	Ndolo	1,9,12	d	1,5	
II	D <sub>1</sub>	Neasden	9,12	g,s,t	e,n,x	
I	B	Neftenbach	4,12	z	e,n,x	
II	S	Negev	41	z <sub>10</sub>	1,2	
I	H	Nessa	1,6,14,25	z <sub>10</sub>	1,2	
I	C <sub>1</sub>	Nessziona	6,7	l,z <sub>13</sub>	1,5	
I	N	Neudorf	30	b	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Neukoelin	6,7	l,z <sub>13</sub> ,[z <sub>28</sub> ]	e,n,z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	B	Neumuenster	<u>1,4,12,27</u>	k	1,6	
I	P	Neunkirchen	38	z <sub>10</sub>	-	
I	E <sub>2</sub>	Newbrunswick	3, <u>15</u>	l,v	1,7	IP combined Newbrunswick and Menhaden (3, <u>15,34</u> :l,v:1,7) with Give (3,10:l,v:1,7) to form Give 3,10,[ <u>15</u> ],[ <u>15,34</u> ]:[d],l,v:[d],1,7. Newbrunswick is now called Give var. O 15+ by IP.
I	E <sub>2</sub>	Newhaw	3, <u>15</u>	e,h	1,5	IP combined Newhaw and Arkansas (3, <u>15,34</u> :e,h:1,5) with Muenster (3,10:e,h:1,5) to form Muenster 3,10,[ <u>15</u> ],[ <u>15,34</u> ]:e,h:1,5. Newhaw is now called Muenster var. O 15+ by IP.
I	54	Newholland	4,12,54	m,t	-	
I	E <sub>2</sub>	Newington	<u>3,15</u>	e,h	1,6	IP combined Newington and Minneapolis (3, <u>15,34</u> :e,h:1,6) with Anatum (3,10:e,h:1,6) to form Anatum 3,10,[ <u>15</u> ],[ <u>15,34</u> ]:e,h:1,6. Newington is now called Anatum var. O 15+ by IP.
I	E <sub>1</sub>	Newlands	3,10,[ <u>15,34</u> ]	e,h	e,n,x	
I	D <sub>1</sub>	Newmexico	9,12	g,z <sub>51</sub>	1,5	
I	C <sub>2</sub>	Newport	6,8, <u>20</u>	e,h	1,2	Newport may possess H phase Rz <sub>50</sub> or Rz <sub>58</sub> or Rz <sub>78</sub> or R1,12
I	C <sub>2</sub>	Newport var. Puerto Rico	6,8	-	1,2	
I	E <sub>1</sub>	Newrochelle	3,10	k	1,w	
I	G	Newyork	13,22	g,s,t	-	
I	D <sub>2</sub>	Ngaparou	9,46	z <sub>4</sub> ,z <sub>24</sub>	-	
I	C <sub>1</sub>	Ngili	6,7	z <sub>10</sub>	1,7	
I	E <sub>4</sub>	Ngor	1,3,19	l,v	1,5	
II	Y	Ngozi	48	z <sub>10</sub>	[1,5]	
I	V	Niakhar	44	a	1,5	
I	J	Niamey	17	d	1,w	
I	V	Niarembe	44	a	1,w	
I	C <sub>1</sub>	Nienstedten	6,7, <u>14</u>	b	[1,w]	Nienstedten was combined with Nissii (6,7, <u>14</u> :b:-) and called Nienstedten; then IP combined Nienstedten with Ohio (6,7:b:l,w) to form Ohio (6,7, <u>14</u> :b:[1,w]). Nienstedten is now called Ohio var.O 14+ by IP.
I	C <sub>1</sub>	Nieukerk	6,7, <u>14</u>	d	z <sub>6</sub>	
I	C <sub>1</sub>	Nigeria	6,7	r	1,6	
I	N	Nijmegen	30	y	e,n,z <sub>15</sub>	
I	I	Nikolaifleet	16	g,m,s	-	
I	E <sub>4</sub>	Niloese	1,3,19	d	z <sub>6</sub>	
I	M	Nima	28	y	1,5	
I	G	Nimes	13,22	z <sub>35</sub>	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Nissii	6,7, <u>14</u>	b	-	Nissii was combined with Nienstedten (6,7, <u>14</u> :b:l,w) as a monophasic variant of Nienstedten. Nienstedten is now called a variant of Ohio by IP. The name Nissii has been dropped.
I	A	Nitra	2,12	g,m	-	
I	E <sub>4</sub>	Niumi	1,3,19	a	1,5	
I	P	Njala	38	k	e,n,x	
I	C <sub>1</sub>	Nola	6,7	e,h	1,7	
II	I	Noordhoek	16	l,w	z <sub>6</sub>	
II	B	Nordenham	<u>1,4,12,27</u>	z	e,n,x	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	C <sub>2</sub>	Nordufer	6,8	a	1,7	
I	C <sub>1</sub>	Norton	6,7	i	1,w	
I	C <sub>1</sub>	Norwich	6,7	e,h	1,6	
I	I	Nottingham	16	d	e,n,z <sub>15</sub>	
I	R	Nowawes	40	z	z <sub>6</sub>	
I	C <sub>2</sub>	Noya	8	r	1,7	
I	I	Nuatja	16	k	e,n,x	
II	T	Nuernberg	42	z	z <sub>6</sub>	
I	F	Nyanza	11	z	z <sub>6</sub>	
I	E <sub>1</sub>	Nyborg	3,10	e,h	1,7	IP combined Selandia (3,15:e,h:1,5) with Nyborg to form Nyborg 3,10,[15]:e,h:1,7.
I	I	Nyeko	16	a	1,7	
I	C <sub>1</sub>	Oakey	6,7	m,t	z <sub>64</sub>	
I	C <sub>1</sub>	Oakland	6,7	z	1,6,[7]	
I	C <sub>1</sub>	Obogu	6,7	z <sub>4</sub> ,z <sub>23</sub>	1,5	
I	E <sub>4</sub>	Ochiogu	1,3,19	z <sub>38</sub>	[e,n,z <sub>15</sub> ]	
I	54	Ochsenwerder	6,7,54	k	1,5	
IV	I	Ochsenzoll	16	z <sub>4</sub> ,z <sub>23</sub>	–	
I	N	Ockenheim	30	1,z <sub>13</sub> ,z <sub>28</sub>	1,6	
I	R	Odiene	40	y	1,5	
II	N	Odijk	30	a	z <sub>39</sub>	
I	N	Odozi	30	k	e,n,[x],z <sub>15</sub>	
I	Q	Oerlikon	39	l,v	e,n,z <sub>15</sub>	
II	M	Oevelgoenne	28	r	e,n,z <sub>15</sub>	
I	S	Offa	41	z <sub>38</sub>	–	
I	U	Ogbete	43	z	1,5	
I	C <sub>1</sub>	Ohio	6,7	b	1,w	IP combined Nienstedten (6,7,14:b:[1,w]) with Ohio to form Ohio 6,7,14:b:[1,w]. Ohio may possess H phase Rz <sub>59</sub> .
I	E <sub>1</sub>	Ohlstedt	3,10	y	e,n,x	
I	G	Okatie	13,23	g,[s],t	–	
I	E <sub>1</sub>	Okefoko	3,10	c	z <sub>6</sub>	
I	E <sub>1</sub>	Okerara	3,10	z <sub>10</sub>	1,2	
I	I	Oldenburg	16	d	1,2	
I	D <sub>2</sub>	Olten	9,46	d	e,n,z <sub>15</sub>	
I	C <sub>1</sub>	Omderman	6,7,14	d	e,n,x	IP combined Omderman with Amersfoort (6,7:d:e,n,x) to form Amersfoort 6,7,14:d:e,n,x. Omderman is now called Amersfoort var. O 14+ by IP.
I	R	Omifisan	40	z <sub>29</sub>	–	
I	C <sub>1</sub>	Omuna	6,7	z <sub>10</sub>	z <sub>35</sub>	
I	M	Ona	28	g,s,t	–	
I	D <sub>1</sub>	Onarimon	1,9,12	b	1,2	
I	H	Onderstepoort	1,6,14,[25]	e,h	1,5	
I	E <sub>1</sub>	Onireke	3,10	d	1,7	
I	D <sub>2</sub>	Ontario	9,46	d	1,5	
I	C <sub>1</sub>	Oranienburg	6,7	m,t	–	IP combined Theilallee (6,7,14:m,t:-) with Oranienburg to form Oranienburg 6,7,14:m,t:-. Oranienburg may possess H phase Rz <sub>57</sub> .
I	T	Orbe	42	b	1,6	
I	52	Ord	52	a	e,n,z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	G	Ordonez	<u>1</u> ,13,23	y	1,w	
I	I	Orientalis	16	k	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Orion	3,10	y	1,5	IP combined Binza (3, <u>15</u> :y:1,5) and Thomasville (3, <u>15</u> ,34:y:1,5) with Orion to form Orion 3,10,[ <u>15</u> ],[ <u>15</u> ,34]:y:1,5.
I	C <sub>1</sub>	Oritamerin	6,7	i	1,5	
I	K	Orlando	18	l,v	e,n,z <sub>15</sub>	
I	U	Orleans	43	d	1,5	
BBBI	D <sub>1</sub>	Os	9,12	a	1,6	
I	M	Oskarshamn	28	y	1,2	
I	C <sub>1</sub>	Oslo	6,7, <u>14</u>	a	e,n,x	
I	F	Osnabrueck	11	l,z <sub>13</sub> ,z <sub>28</sub>	e,n,x	
I	C <sub>1</sub>	Othmarschen	6,7, <u>14</u>	g,m,[t]	–	
I	D <sub>1</sub>	Ottawa	<u>1</u> ,9,12	z <sub>41</sub>	1,5	
II	R	Ottershaw	40	d	–	
I	D <sub>2</sub>	Ouakam	9,46	z <sub>29</sub>	–	
I	G	Oudwijk	13,22	b	1,6	
I	R	Overchurch	<u>1</u> ,40	l,w	[1,2]	
I	51	Overschie	51	l,v	1,5	
I	N	Overvecht	30	a	1,2	
I	E <sub>1</sub>	Oxford	3,10	a	1,7	IP combined Khartoum (3, <u>15</u> ,34:a:1,7) with Oxford to form Oxford 3,10,[ <u>15</u> ],[ <u>15</u> ,34]:a:1,7.
I	C <sub>1</sub>	Oyonnax	6,7	y	1,6	
II	C <sub>1</sub>	Oysterbeds	6,7	z	z <sub>42</sub>	
I	C <sub>2</sub>	Pakistan	8	l,v	1,2	
I	V	Palamaner	<u>1</u> ,44	d	z <sub>35</sub>	
I	C <sub>1</sub>	Palime	6,7	z <sub>35</sub>	e,n,z <sub>15</sub>	
I	D <sub>1</sub>	Panama	<u>1</u> ,9,12	l,v	1,5	Panama may possess H phase R1,11
I	E <sub>2</sub>	Pankow	3, <u>15</u>	d	1,5	IP combined Pankow with Shangani (3,10:d:1,5) to form Shangani 3,10, <u>15</u> :d:1,5. Pankow is now called Shangani var. O 15+ by IP.
I	C <sub>1</sub>	Papuana	6,7	r	e,n,z <sub>15</sub>	
I	A	Paratyphi A	<u>1</u> ,2,12	a	[1,5]	
I	B	Paratyphi B	1,4,[5],12	[b]	[1,2]	Paratyphi B is tartrate neg.; Paratyphi B var. Java (CDC calls this S. ser. Java) is often monophasic ( <u>1</u> ,4,5,12:b:-) and is tartrate pos. Paratyphi B and Java may possess H phase RZ <sub>33</sub> .
I	C <sub>1</sub>	Paratyphi C	6,7,[Vi]	c	1,5	
IV	C <sub>2</sub>	Parera	11	z <sub>4</sub> ,z <sub>23</sub>	–	
I	C <sub>2</sub>	Paris	8, <u>20</u>	z <sub>10</sub>	1,5	
I	E <sub>4</sub>	Parkroyal	1,3,19	l,v	1,7	
II	E <sub>1</sub>	Parow	3,10,[ <u>15</u> ]	g,m,s,t	–	
I	B	Pasing	4,12	z <sub>35</sub>	1,5	
I	M	Patience	28	d	e,n,z <sub>15</sub>	
I	D <sub>1</sub>	Penarth	9,12	z <sub>35</sub>	z <sub>6</sub>	
I	M	Penilla	28	l,z <sub>13</sub> ,z <sub>28</sub>	e,n,z <sub>15</sub>	
I	F	Pennsylvania	11	d	e,n,z <sub>15</sub>	
I	D <sub>1</sub>	Pensacola	<u>1</u> ,9,12	m,t	[1,2]	
II	W	Perinet	45	g,m,t	e,n,x,z <sub>15</sub>	
I	P	Perth	38	y	e,n,x	
I	E <sub>4</sub>	Petahtikva	1,3,19	f,g,t	1,7	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	C <sub>2</sub>	Phaliron	8	z	e,n,Z <sub>15</sub>	
I	F	Pharr	11	b	e,n,Z <sub>15</sub>	
II	X	Phoenix	47	b	1,5	
I	E <sub>1</sub>	Pietersburg	3,10,[15,34]	Z <sub>69</sub>	1,7	
I	C <sub>2</sub>	Pikine	8,20	r	Z <sub>6</sub>	Pikine was combined with Altona (8,20:r,[i]:z <sub>6</sub> ). The name Pikine has been dropped.
I	I	Pisa	16	i	1,w	
I	C <sub>1</sub>	Planckendael	6,7	Z <sub>4</sub> ,Z <sub>23</sub>	1,6	
I	V	Ploufragan	1,44	Z <sub>4</sub> ,Z <sub>23</sub>	e,n,Z <sub>15</sub>	
I	D <sub>2</sub>	Plymouth	9,46	d	Z <sub>6</sub>	
I	H	Poano	1,6,14,25	z	1,Z <sub>13</sub> ,Z <sub>28</sub>	
I	54	Poeseldorf	8,20,54	i	Z <sub>6</sub>	
I	C <sub>1</sub>	Poitiers	6,7	z	1,5	
I	M	Pomona	28	y	1,7	Pomona may possess H phases RZ <sub>60</sub> , RZ <sub>70</sub> or RZ <sub>80</sub> .
I	K	Pontypridd	18	g,m	–	
I	G	Poona	1,13,22	z	1,6	Poona may possess H phase Rz59.
I	C <sub>2</sub>	Portanigra	8,20	d	1,7	
II	T	Portbech	42	l,v	e,n,x,Z <sub>15</sub>	
I	D <sub>1</sub>	Portland	9,12	Z <sub>10</sub>	1,5	
I	E <sub>2</sub>	Portsmouth	3,15	l,v	1,6	IP combined Portsmouth with London (3,10:l,v:1,6) to form London 3,10,[15]:l,v:1,6. Portsmouth is now called London var. O 15+ by IP.
I	K	Potengi	18	z	–	
I	H	Potosi	6,14	Z <sub>36</sub>	1,5	
I	C <sub>1</sub>	Potsdam	6,7,14	l,v	e,n,Z <sub>15</sub>	
I	D <sub>2</sub>	Potto	9,46	i	Z <sub>6</sub>	
I	D <sub>1</sub>	Powell	9,12	y	1,7	
I	C <sub>2</sub>	Praha	6,8	y	e,n,Z <sub>15</sub>	
I	E <sub>1</sub>	Pramiso	3,10	c	1,7	
I	C <sub>2</sub>	Presov	6,8	b	e,n,Z <sub>15</sub>	
I	B	Preston	1,4,12	z	1,w	
I	F	Pretoria	11	k	1,2	
I	D <sub>1</sub>	Pullorum	1,9,12	–	–	IP combined Pullorum with Gallinarum (1,9,12:-:-). They must be identified biochemically.
I	G	Putten	13,23	d	1,w	
I	V	Quebec	44	c	e,n,Z <sub>15</sub>	
I	D <sub>2</sub>	Quentin	9,46	d	1,6	
II	X	Quimbamba	47	d	Z <sub>39</sub>	
I	X	Quinhon	47	Z <sub>44</sub>	–	
I	C <sub>2</sub>	Quiniela	6,8	c	e,n,Z <sub>15</sub>	
I	N	Ramatgan	30	k	1,5	
I	M	Ramsey	28	l,w	1,6	
II	T	Rand	42	z	e,n,x,Z <sub>15</sub>	
I	G	Raus	13,22	f,g	e,n,x	
I	K	Rawash	6,14,18	c	e,n,x	
I	B	Reading	1,4,[5],12	e,h	[1,5]	
I	C <sub>2</sub>	Rechovot	8,20	e,h	Z <sub>6</sub>	
I	C <sub>1</sub>	Redba	6,7	Z <sub>10</sub>	Z <sub>6</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	F	Redhill	11	e,h	1,z <sub>13</sub> ,z <sub>28</sub>	
I	I	Redlands	16	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	E <sub>1</sub>	Regent	3,10	f,g,[s]	[1,6]	
I	B	Reinickendorf	4,12	1,z <sub>28</sub>	e,n,x	
I	F	Remete	11	z <sub>4</sub> ,z <sub>23</sub>	1,6	
I	C <sub>2</sub>	Remiremont	8,20	z <sub>10</sub>	1,w	
I	B	Remo	1,4,12,27	r	1,7	
I	C <sub>2</sub>	Reubeuss	8,20	g,m,t	–	
II	D <sub>1</sub>	Rhodesiense	9,12	d	e,n,x	
I	L	Rhone	21	c	e,n,x	
I	I	Rhydyfelin	16	e,h	e,n,x	
I	C <sub>1</sub>	Richmond	6,7	y	1,2	
I	E <sub>4</sub>	Rideau	1,3,19	f,g	–	
I	D <sub>1</sub>	Ridge	9,12	c	z <sub>6</sub>	
I	G	Ried	1,13,22	z <sub>4</sub> ,z <sub>23</sub>	[e,n,z <sub>15</sub> ]	
I	C <sub>1</sub>	Riggil	6,7	g,t	–	
I	R	Riogrande	40	b	1,5	
I	C <sub>1</sub>	Rissen	6,7	f,g	–	IP combined Ardwick (6,7,14:f:g:-) with Rissen for form Rissen 6,7,14:f:g:-.
I	P	Rittersbach	38	b	e,n,z <sub>15</sub>	
I	W	Riverside	45	b	1,5	
I	P	Roan	38	l,v	e,n,x	
I	Z	Rochdale	50	b	e,n,x	
II	51	Roggeveld	51	–	1,7	
I	M	Rogy	28	z <sub>10</sub>	1,2	
I	G	Romanby	1,13,23	z <sub>4</sub> ,z <sub>24</sub>	–	
I	G	Roodepoort	1,13,22	z <sub>10</sub>	1,5	
II	H	Rooikrantz	1,6,14	m,t	1,5	
I	E <sub>2</sub>	Rosenthal	3,15	b	1,5	IP combined Rosenthal and unnamed 3,15,34:b:1,5 with Butantan (3,10:b:1,5) to form Butantan 3,10,[15],[15,34]:b:1,5. Rosenthal is now called Butantan var. O 15+ by IP.
I	54	Rossleben	54	e,h	1,6	
I	D <sub>1</sub>	Rostock	1,9,12	g,p,u	–	
IV	C <sub>1</sub>	Roterberg	6,7	z <sub>4</sub> ,z <sub>23</sub>	–	
I	P	Rothenburgsort	38	m,t	–	
II	G	Rotterdam	1,13,22	g,t	1,5	
I	G	Rottnest	1,13,22	b	1,7	
I	I	Rovaniemi	16	r,[i]	1,5	
II	I	Rowburton	16	m,t	[z <sub>42</sub> ]	
I	H	Royan	1,6,14,25	z	e,n,z <sub>15</sub>	
I	D <sub>1</sub>	Ruanda	9,12	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	F	Rubislaw	11	r	[e,n,x]	
I	L	Ruiru	21	y	e,n,x	
I	B	Ruki	4,5,12	y	e,n,x	IP combined Ruki with Ball (1,4,12,27:y:e,n,x) and Dalat (4,5,27:y:e,n,x) to form Ball 1,4,[5],12,27:y:e,n,x. The name Ruki has been dropped.
I	C <sub>1</sub>	Rumford	6,7	z <sub>38</sub>	1,2	
I	H	Runby	1,6,14,25	c	e,n,x	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>1</sub>	Rutgers	3,10	R <sub>1</sub> ,Z <sub>40</sub>	1,7	Rutgers has been dropped from the scheme and the H phase R <sub>1</sub> ,Z <sub>40</sub> is now considered an R phase of Give.
I	E <sub>1</sub>	Ruzizi	3,10	l,v	e,n,Z <sub>15</sub>	
I	D <sub>1</sub>	Saarbruecken	<u>1</u> ,9,12	a	1,7	
I	I	Saboya	16	e,h	1,5	
IV	R	Sachsenwald	<u>1</u> ,40	Z <sub>4</sub> ,Z <sub>23</sub>	–	
I	N	Sada	30	Z <sub>10</sub>	1,2	
I	52	Saintemarie	52	g,t	–	
I	B	Saintpaul	<u>1</u> ,4,[5],12	e,h	1,2	
I	X	Saka	47	b	–	IP combined Saka with Sya (47:b:z <sub>6</sub> ) and called it Sya.
II	Y	Sakaraha	48	[k]	Z <sub>39</sub>	
I	I	Salford	16	l,v	e,n,x	
I	B	Salinatis	4,12	d,e,h	d,e,n,Z <sub>15</sub>	IP states that Salinatis was combined with Duisburg (1,4,12,27:d:e,n,Z <sub>15</sub> ). This is incorrect; IP should have stated that it was combined with Sandiego (4,[5],12:e,h:e,n,Z <sub>15</sub> ), because Salinatis loses the d and becomes Sandiego.
I	I	Saloniki	16	Z <sub>29</sub>	–	
I	S	Samaru	41	i	1,5	
I	E <sub>4</sub>	Sambre	1,3,19	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	B	Sandiego	4,[5],12	[e,h]	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Sadow	6,8	f,g	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Sanga	8	b	1,7	
I	D <sub>2</sub>	Sangalkam	9,46	m,t	–	
I	I	Sangera	16	b	e,n,Z <sub>15</sub>	
I	C <sub>1</sub>	Sanjuan	6,7	a	1,5	
I	M	Sanktgeorg	28	r,[i]	e,n,Z <sub>15</sub>	
I	E <sub>4</sub>	Sanktmarx	1,3,19	e,h	1,7	
I	M	Santander	28	Z <sub>35</sub>	e,n,Z <sub>15</sub>	
I	R	Santhiaba	40	l,Z <sub>28</sub>	1,6	
I	C <sub>2</sub>	Santiago	8, <u>20</u>	c	e,n,x	
I	E <sub>4</sub>	Sao	1,3,19	e,h	e,n,Z <sub>15</sub>	
I	I	Saphra	16	y	1,5	
I	H	Sara	1,6,14,25	Z <sub>38</sub>	[e,n,x]	
I	B	Sarajane	<u>1</u> ,4,[5],12, <u>27</u>	d	e,n,x	
II	I	Sarepta	16	l,Z <sub>28</sub>	Z <sub>42</sub>	
I	R	Saugus	40	b	1,7	
I	H	Schalkwijk	6,14,[24]	i	e,n,..	
I	B	Schleissheim	4,12, <u>27</u>	b	–	
I	E <sub>4</sub>	Schoeneberg	1,3,19	z	e,n,Z <sub>15</sub>	
I	C <sub>1</sub>	Schwabach	6,7	c	1,7	
I	B	Schwarzengrund	<u>1</u> ,4,12, <u>27</u>	d	1,7	
I	C <sub>2</sub>	Schwerin	6,8	k	e,n,x	
I	I	Sculcoates	16	d	1,5	
II	Z	Seaforth	50	k	Z <sub>6</sub>	
I	M	Seattle	28	a	e,n,x	
I	V	Sedgwick	44	b	e,n,Z <sub>15</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>1</sub>	Seegefeld	3,10	r,[i]	1,2	
I	E <sub>1</sub>	Sekondi	3,10	e,h	z <sub>6</sub>	
I	E <sub>2</sub>	Selandia	3, <u>15</u>	e,h	1,7	IP combined Selandia with Nyborg (3,10:e,h:1,7) to form Nyborg 3,10,[15]:e,h:1,7. Selandia is now called Nyborg var. O 15+ by IP.
I	M	Selby	28	y	z <sub>6</sub>	
IV	R	Seminole	<u>1</u> ,40	g,z <sub>51</sub>	–	
I	D <sub>1</sub>	Sendai	<u>1</u> ,9,12	a	1,5	Sendai must be differentiated from Miami with biochemical tests. Sendai is neg. for H <sub>2</sub> S, citrate, and tartrate; Miami is pos.
I	F	Senegal	11	r	1,5	
I	E <sub>4</sub>	Senftenberg	1,3,19	g,[s],t		Senftenberg may possess H phase Rz <sub>37</sub> or Rz <sub>43</sub> or Rz <sub>45</sub> or Rz <sub>46</sub> . Simsbury (1,3,19:Rz <sub>27</sub> :-) is now considered an H phase Rz <sub>27</sub> of Senftenberg.
I	N	Senneville	30	z <sub>10</sub>	1,5	
I	D <sub>1</sub>	Seremban	9,12	i	1,5	
I	E <sub>1</sub>	Serrekunda	3,10	k	1,7	
II	60	Setubal	60	g,m,t	z <sub>6</sub>	
I	I	Shamba	16	c	e,n,x	
I	E <sub>1</sub>	Shangani	3,10	d	1,5	IP combined Pankow (3, <u>15</u> :d:1,5) with Shangani to form Shangani 3,10,[15]:d:1,5.
BBI	I	Shanghai	16	l,v	1,6	
I	E <sub>1</sub>	Shannon	3,10	z <sub>35</sub>	1,w	
I	F	Sharon	11	k	1,6	
I	P	Sheffield	38	c	1,5	
I	I	Sherbrooke	16	d	1,6	
I	R	Shikmonah	40	a	1,5	
I	C <sub>2</sub>	Shiplay	8, <u>20</u>	b	e,n,z <sub>15</sub>	
I	M	Shomolu	28	y	1,w	
IIIa	K	Shomron	18	z <sub>4</sub> ,z <sub>32</sub>	–	Shomron was formerly in Subspecies II, but is now combined with <i>Arizona 7a,7b:1,7,8</i> :-. The name Shomron has been dropped.
I	D <sub>2</sub>	Shoreditch	9,46	r	e,n,z <sub>15</sub>	
I	B	Shubra	4,[5],12	z	1,2	
I	S	Sica	41	b	e,n,z <sub>15</sub>	
I	K	Siegburg	<u>6</u> , <u>14</u> ,18	z <sub>4</sub> ,z <sub>23</sub>	[1,5]	IP combined Siegburg with Cerro (18:z <sub>4</sub> ,z <sub>23</sub> :[1,5]) to form Cerro <u>6</u> , <u>14</u> ,18:z <sub>4</sub> ,z <sub>23</sub> :[1,5]. Siegburg is now called Cerro var. O 14+. The name Siegburg has been dropped.
I	E <sub>1</sub>	Simi	3,10	r	e,n,z <sub>15</sub>	
II	H	Simonstown	1,6,14	z <sub>10</sub>	1,5	
I	E <sub>4</sub>	Simsbury	1,3,19	Rz <sub>27</sub>	–	IP combined Simsbury with Senftenberg 1,3,19:g,[s],t:-. Simsbury is now considered an R phase of Senftenberg. The name Simsbury has been dropped.
I	E <sub>1</sub>	Sinchew	3,10	l,v	z <sub>35</sub>	
I	C <sub>2</sub>	Sindelfingen	8, <u>20</u>	y	1,w	
I	C <sub>1</sub>	Singapore	6,7	k	e,n,x	
I	E <sub>1</sub>	Sinstorf	3,10	l,v	1,5	
I	K	Sinthia	18	z <sub>38</sub>	–	
I	T	Sipane	<u>1</u> ,42	r	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Skansen	6,8	b	1,2	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>4</sub>	Slade	1,3,19	y	e,n,z <sub>15</sub>	
I	B	Sladun	<u>1,4,12,27</u>	b	e,n,x	IP combined Sladun with Abony (1,4,5,12:b:e,n,x) to form Abony 1,4,[5],12,27:b:e,n,x. Sladun is now called Abony var. O 27+. The name Sladun has been
II	H	Slangkop	<u>1,6,14</u>	z <sub>10</sub>	z <sub>6</sub> ;z <sub>42</sub>	
II	N	Slatograd	30	g,t	–	
I	X	Sljeme	1,47	f,g	–	
I	B	Sloterdijk	<u>1,4,12,27</u>	z <sub>35</sub>	z <sub>6</sub>	
I	H	Soahanina	6,14,24	z	e,n,x	
I	N	Soerenga	30	i	l,w	
IV	L	Soesterberg	21	z <sub>4</sub> ,z <sub>23</sub>	–	
II	B	Sofia	<u>1,4,12,27</u>	b	[e,n,x]	
I	D <sub>2</sub>	Sokode	9,46	r	z <sub>6</sub>	
I	M	Solna	28	a	1,5	
I	F	Solt	11	y	1,5	
I	C <sub>1</sub>	Somone	6,7	z <sub>4</sub> ,z <sub>24</sub>	–	
I	M	Soumbedioune	28	b	e,n,x	
I	B	Southampton	<u>1,4,12,27</u>	r	z <sub>6</sub>	
I	E <sub>1</sub>	Southbank	3,10, <u>15</u> ,34	m,t	[1,6]	
II	F	Soutpan	11	z	z <sub>39</sub>	
I	E <sub>1</sub>	Souza	3,10	d	e,n,x	IP combined Eschersheim (3, <u>15</u> :d:e,n,x) with Souza to form Souza 3,10,[ <u>15</u> ]:d:e,n,x.
I	T	Spalentor	<u>1,42</u>	y	e,n,z <sub>15</sub>	
I	L	Spartel	21	d	1,5	
I	V	Splott	44	g,s,t	–	
II	R	Springs	40	a	z <sub>39</sub>	
VI	F	Srinagar	11	b	e,n,x	
I	P	Stachus	38	z	–	
I	B	Stanley	4,5,12	d	1,2	IP combined Cairo (1,4,12,27:d:1,2) with Stanley to form Stanley <u>1,4</u> ,[5],12,27:d:1,2
I	B	Stanleyville	<u>1,4</u> ,[5],12	z <sub>4</sub> ,z <sub>23</sub>	[1,2]	IP combined Jaja (4,12,27:z <sub>4</sub> ,z <sub>23</sub> :-) with Stanleyville to form Stanleyville <u>1,4</u> ,[5],12,27:z <sub>4</sub> ,z <sub>23</sub> : [1,2].
I	X	Staoueli	47	k	1,2	
I	N	Steinplatz	30	y	1,6	
I	54	Steinwerder	3, <u>15</u> ,54	y	1,5	
II	D <sub>1</sub>	Stellenbosch	<u>1,9,12</u>	z	1,7	
I	X	Stellingen	47	d	[e,n,x]	
I	F	Stendal	11	l,v	1,2	
I	N	Sternschanze	30	g,s,t	–	Sternschanze may possess H phase Rz <sub>59</sub> .
I	C <sub>2</sub>	Sterrenbos	6,8	d	e,n,x	
II	G	Stevenage	<u>1,13,23</u>	[z <sub>42</sub> ]	1,[5],7	
II	E <sub>1</sub>	Stikland	3,10	m,t	e,n,x	
I	E <sub>1</sub>	Stockholm	3,10	y	z <sub>6</sub>	IP combined Tournai (3, <u>15</u> :y:z <sub>6</sub> ) with Stockholm to form Stockholm 3,10,[ <u>15</u> ]:y:z <sub>6</sub> .
I	N	Stoneferry	30	z <sub>4</sub> ,z <sub>23</sub>	–	
I	E <sub>1</sub>	Stormont	3,10	d	1,2	
I	C <sub>2</sub>	Stourbridge	6,8	b	1,6	
I	F	Straengnaes	11	z <sub>10</sub>	1,5	
I	D <sub>2</sub>	Strasbourg	9,46	d	1,7	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	E <sub>4</sub>	Stratford	1,3,19	i	1,2	
I	C <sub>1</sub>	Strathcona	6,7	1,z <sub>13</sub> ,z <sub>28</sub>	1,7	
I	E <sub>4</sub>	Stuivenberg	1,3,19	1,z <sub>13</sub> ,z <sub>28</sub>	1,5	
I	C <sub>1</sub>	Stuttgart	6,7,14	i	z <sub>6</sub>	
II	R	Suarez	1,40	c	e,n,x,z <sub>15</sub>	
I	E <sub>1</sub>	Suberu	3,10	g,m	–	
I	U	Sudan	43	1,z <sub>13</sub>	–	
II	D <sub>1</sub>	Suederelbe	1,9,12	b	z <sub>39</sub>	
I	W	Suelldorf	45	f,g	–	
II	C <sub>1</sub>	Sullivan	6,7	z <sub>42</sub>	1,7	
I	H	Sundsvall	[1],6,14,[25]	z	e,n,x	
I	C <sub>2</sub>	Sunncove	8	y	e,n,x	
II	R	Sunnydale	1,40	k	e,n,x,z <sub>15</sub>	
I	H	Surat	[1],6,14,[25]	[r],[i]	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Svedvi	1,3,19	l,v	e,n,z <sub>15</sub>	
I	X	Sya	47	b	z <sub>6</sub>	
IIIb	Y	Sydney	48	i	z	Sydney was formerly in subspecies II, but it is now combined with <i>Arizona</i> 5:33:31. The name Sydney has been dropped.
I	H	Sylvania	[1],6,14,[25]	g,p	–	
I	I	Szentes	16	k	1,2	
I	X	Tabligbo	47	z <sub>4</sub> ,z <sub>23</sub>	[e,n,z <sub>15</sub> ]	
I	C <sub>2</sub>	Tado	8,20	c	z <sub>6</sub>	
II	E <sub>1</sub>	Tafelbaai	3,10	z	z <sub>39</sub>	
I	B	Tafo	1,4,12,27	z <sub>35</sub>	1,7	
I	C <sub>2</sub>	Takoradi	6,8	i	1,5	
I	E <sub>4</sub>	Taksony	1,3,19	[i]	z <sub>6</sub>	
I	C <sub>2</sub>	Tallahassee	6,8	z <sub>4</sub> ,z <sub>32</sub>	–	
I	C <sub>2</sub>	Tamale	8,20	z <sub>29</sub>	[e,n,z <sub>15</sub> ]	
I	E <sub>4</sub>	Tambacounda	1,3,19	b	e,n,x	
I	X	Tamberma	47	z <sub>4</sub> ,z <sub>24</sub>	–	
I	C <sub>1</sub>	Tamilnadu	6,7	z <sub>41</sub>	z <sub>35</sub>	
I	C <sub>1</sub>	Tampico	6,7	z <sub>36</sub>	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Tananarive	6,8	y	1,5	
I	G	Tanger	1,13,22	y	1,6	
I	G	Tanzania	1,13,22	z	e,n,z <sub>15</sub>	
I	D <sub>1</sub>	Tarshyne	9,12	d	1,6	
I	T	Taset	1,42	z <sub>41</sub>	–	
I	M	Taunton	28	k	e,n,x	
I	O	Tchad	35	b	–	
I	J	Tchamba	17	z	e,n,z <sub>15</sub>	
I	M	Techimani	28	c	z <sub>6</sub>	
I	B	Teddington	1,4,12,27	y	1,7	
I	I	Tees	16	f,g	–	
I	B	Tejas	4,12	z <sub>36</sub>	–	
I	H	Teko	1,6,14,25	d	e,n,z <sub>15</sub>	
I	M	Telaviv	28	y	e,n,z <sub>15</sub>	
I	G	Teitelkebir	13,23	d	e,n,z <sub>15</sub>	
I	F	Telhashomer	11	z <sub>10</sub>	e,n,x	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	M	Teltow	28	z <sub>4</sub> ,z <sub>23</sub>	1,6	
I	T	Tema	<u>1</u> ,42	z <sub>35</sub>	z <sub>6</sub>	
I	C <sub>1</sub>	Tennessee	6,7, <u>14</u>	z <sub>29</sub>	[1,2,7]	
I	B	Tennyson	4,5,12	g,z <sub>51</sub>	e,n,z <sub>15</sub>	
I	X	Teshie	1,47	l,z <sub>13</sub> ,z <sub>28</sub>	e,n,z <sub>15</sub>	
I	B	Texas	4,[5],12	k	e,n,z <sub>15</sub>	
I	B	Thayngen	<u>1</u> ,4,12, <u>27</u>	z <sub>41</sub>	1,(2),5	
I	U	Thetford	43	k	1,2	
I	P	Thiaroye	38	e,h	1,2	
I	C <sub>1</sub>	Thielallee	6,7, <u>14</u>	m,t	-	IP combined Thielallee with Oranienburg (6,7:m,t:-) to form Oranienburg 6,7, <u>14</u> :m,t:-. Thielallee is now called Oranienburg var. O 14+ by IP.
I	E <sub>4</sub>	Thies	1,3,19	y	1,7	
I	E <sub>3</sub>	Thomasville	3, <u>15</u> , <u>34</u>	y	1,5	IP combined Thomasville and Binza (3, <u>15</u> :y:1,5) with Orion (3,10:y:1,5) to form Orion 3,10,[ <u>15</u> ],[ <u>15</u> , <u>34</u> ]:y:1,5. Thomasville is now called Orion var. O 15+ by IP.
I	C <sub>1</sub>	Thompson	6,7, <u>14</u>	[k]	[1,5]	IP combined Cardiff that contains H phase R1,10 (6,7:k:R1,10) with Thompson.
I	C <sub>1</sub>	Tienba	6,7	z <sub>35</sub>	1,6	
I	V	Tiergarten	44	a	e,n,x	
I	R	Tiko	40	l,z <sub>13</sub> ,z <sub>28</sub>	1,2	
I	E <sub>4</sub>	Tilburg	1,3,19	d	1,w	Tilburg may possess H phase Rz <sub>49</sub> .
I	R	Tilene	<u>1</u> ,40	e,h	1,2	
I	B	Tinda	<u>1</u> ,4,5, <u>27</u>	a	e,n,z <sub>15</sub>	
I	51	Tione	51	a	e,n,x	
I	I	Togba	16	a	e,n,x	
I	B	Togo	4,12	l,w	1,6	
II	57	Tokai	57	z <sub>42</sub>	1,6:z <sub>53</sub>	
I	B	Tokoin	4,12	z <sub>10</sub>	e,n,z <sub>15</sub>	
I	T	Tomegbe	<u>1</u> ,42	b	e,n,z <sub>15</sub>	
I	E <sub>4</sub>	Tomelilla	1,3,19	l,z <sub>28</sub>	1,7	
I	54	Tonev	21,54	b	e,n,x	
I	F	Toowong	11	a	1,7	
I	N	Torhout	30	e,h	1,5	
I	T	Toricada	<u>1</u> ,42	z <sub>4</sub> ,z <sub>24</sub>	-	
I	W	Tornow	45	g,m,[s],[t]	-	
I	D <sub>2</sub>	Toronto	9,46	l,v	e,n,x;[z <sub>44</sub> ]	
II	C <sub>1</sub>	Tosamanga	6,7	z	1,5	
I	Y	Toucra	48	z	1,5	Toucra may possess H phase Rz <sub>58</sub> .
I	K	Toulon	18	l,w	e,n,z <sub>15</sub>	
I	C <sub>2</sub>	Tounouma	8, <u>20</u>	b	z <sub>6</sub>	
I	E <sub>2</sub>	Tournai	3, <u>15</u>	y	z <sub>6</sub>	IP combined Tournai with Stockholm (3,10:y:z <sub>6</sub> ) to form Stockholm 3,10,[ <u>15</u> ]:y:z <sub>6</sub> . Tournai is now called Stockholm var. O 15+ by IP.
I	B	Trachau	4,12, <u>27</u>	y	1,5	
II	55	Tranoroa	55	k	z <sub>39</sub>	
I	W	Transvaal	45	z <sub>4</sub> ,z <sub>24</sub>	-	
I	B	Travis	4,[5],12	g,z <sub>51</sub>	1,7	
I	51	Treforest	<u>1</u> ,51	z	1,6	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	D1	Treguier	<u>1</u> ,9,12	Z <sub>10</sub>	Z <sub>6</sub>	
I	I	Trier	16	Z <sub>35</sub>	1,6	
I	D <sub>2</sub>	Trimdon	9,46	Z <sub>35</sub>	Z <sub>6</sub>	
I	B	Tripoli	<u>1</u> ,4,12, <u>27</u>	b	Z <sub>6</sub>	
I	R	Trotha	40	Z <sub>10</sub>	Z <sub>6</sub>	
I	E <sub>1</sub>	Truro	3,10	i	1,7	
I	G	Tschangu	<u>1</u> ,13,23	e,h	1,5	
I	B	Tsevie	4,12	i	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Tshiongwe	6,8	e,h	e,n,Z <sub>15</sub>	
I	H	Tucson	[1],6,14,[25]	b	[1,7]	
I	B	Tudu	4,12	Z <sub>10</sub>	1,6	
I	E <sub>2</sub>	Tuebingen	3, <u>15</u>	y	1,2	IP combined Tuebingen with Amager (e,10:y:1,2) to form Amager 3,10,[ <u>15</u> ]:y:1,2. Tuebingen is now called Amager var. O 15+ by IP.
IV	U	Tuindorp	43	Z <sub>4</sub> ,Z <sub>32</sub>	–	
II	C <sub>2</sub>	Tular	6,8	a	Z <sub>52</sub>	
I	B	Tumodi	<u>1</u> ,4,12	i	Z <sub>6</sub>	
I	G	Tunis	<u>1</u> ,13,23	y	Z <sub>6</sub>	
II	G	Tygerberg	<u>1</u> ,13,23	a	Z <sub>42</sub>	
I	D <sub>1</sub>	Typhi	9,12,[Vi]	d	–	Typhi may possess H phase R <sub>j</sub> or R <sub>Z<sub>66</sub></sub> .
I	B	Typhimurium	<u>1</u> ,4,5,12	i	1,2,[7]	
I	B	Typhimurium var. Copenhagen	<u>1</u> ,4,12	i	1,2	
I	C <sub>1</sub>	Typhisuis	6,7	c	1,5	Typhisuis is a bioserotype found in pigs. It is like Choleraesuis except tartrate negative.
I	B	Tyresoe	4,12	1,[Z <sub>13</sub> ],Z <sub>28</sub>	1,5	
I	54	Uccle	3,54	g,s,t	–	
I	E <sub>1</sub>	Uganda	3,10, <u>15</u>	1,Z <sub>13</sub>	1,5	
I	E <sub>1</sub>	Ughelli	3,10	r	1,5	
I	V	Uhlenhorst	44	z	1,w	
I	52	Uithof	52	a	1,5	
I	G	Ullevi	<u>1</u> ,13,23	b	e,n,x	
I	M	Umbilo	28	Z <sub>10</sub>	e,n,x	
I	C <sub>1</sub>	Umhlali	6,7	a	1,6	
I	O	Umhlatazana	35	a	e,n,Z <sub>15</sub>	
I	C <sub>2</sub>	Uno	6,8	Z <sub>29</sub>	[e,n,Z <sub>15</sub> ]	
II	T	Uphill	42	b	e,n,x,Z <sub>15</sub>	
I	B	Uppsala	4,12, <u>27</u>	b	1,7	
I	N	Urbana	30	b	e,n,x	
I	T	Ursenbach	<u>1</u> ,42	z	1,6	
I	K	Usumbura	<u>6</u> ,14,18	d	1,7	
I	C <sub>2</sub>	Utah	6,8	c	1,5	
II	O	Utbremen	35	Z <sub>29</sub>	e,n,x	
I	52	Utrecht	52	d	1,5	
I	H	Uzaramo	1,6,14,25	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	G	Vaertan	13,22	b	e,n,x	
I	C <sub>2</sub>	Valdosta	6,8	a	1,2	
I	I	Vancouver	16	c	1,5	
I	M	Vanier	28	z	1,5	
I	S	Vaugirard	41	b	1,6	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	U	Veddel	43	g,t	–	
I	I	Veogesack	16	b	1,w	
I	E <sub>1</sub>	Vejle	3,10	e,h	1,2	IP combined Goerlitz (3,15:e,h:1,2) with Vejle to form Vejle 3,10,[15]:e,h:1,2.
I	B	Vellore	<u>1,4,12,27</u>	z <sub>10</sub>	z <sub>35</sub>	
I	F	Veneziana	11	i	e,n,x	
II	J	Verity	17	e,n,x,z <sub>15</sub>	1,6	
I	S	Verona	41	i	1,6	
I	W	Verviers	45	k	1,5	
I	D <sub>1</sub>	Victoria	<u>1,9,12</u>	l,w	1,5	
I	J	Victoriaborg	17	c	1,6	
I	S	Vietnam	41	b	z <sub>6</sub>	
II	S	Vietnam var. subsp. II	41	b	–	
I	E <sub>4</sub>	Vilvoorde	1,3,19	e,h	1,5	
I	M	Vinohrady	28	m,t	[e,n,z <sub>15</sub> ]	
I	C <sub>1</sub>	Virchow	6,7	r	1,2	
I	C <sub>2</sub>	Virginia	8	d	[1,2]	
I	E <sub>4</sub>	Visby	1,3,19	b	1,6	
I	M	Vitkin	28	l,v	e,n,x	
I	V	Vleuten	44	f,g	–	
I	T	Vogan	<u>1,42</u>	z <sub>38</sub>	z <sub>6</sub>	
IV	U	Volksdorf	43	z <sub>36</sub> ,z <sub>38</sub>	–	
I	M	Volksmarsdorf	28	i	1,6	
I	F	Volta	11	r	l,z <sub>13</sub> ,z <sub>28</sub>	
I	B	Vom	<u>1,4,12,27</u>	l,[z <sub>13</sub> ],[z <sub>28</sub> ]	e,n,z <sub>15</sub>	
I	U	Voulte	43	i	e,n,x	
II	G	Vredelust	<u>1,13,23</u>	l,z <sub>28</sub>	z <sub>42</sub>	
I	G	Vridi	<u>1,13,23</u>	e,h	1,w	
VI	W	Vrindaban	45	a	e,n,x	
I	B	Vuadens	<u>4,12,27</u>	z <sub>4</sub> ,z <sub>23</sub>	z <sub>6</sub>	
I	I	Wa	16	b	1,5	
I	D <sub>2</sub>	Waedenswil	9,46	e,h	1,5	
I	E <sub>1</sub>	Wagadugu	3,10	z <sub>4</sub> ,z <sub>23</sub>	z <sub>6</sub>	
I	B	Wagenia	<u>1,4,12,27</u>	b	e,n,z <sub>15</sub>	
II	L	Wandsbek	21	z <sub>10</sub>	[z <sub>6</sub> ]	
I	Q	Wandsworth	39	b	1,2	
I	D <sub>1</sub>	Wangata	<u>1,9,12</u>	z <sub>4</sub> ,z <sub>23</sub>	[1,7]	
I	T	Waral	<u>1,42</u>	m,t	–	
I	J	Warengo	17	z	1,5	
I	W	Warmesen	45	d	e,n,z <sub>15</sub>	
I	M	Warnemuende	28	i	e,n,x	
I	C <sub>2</sub>	Warnow	6,8	i	1,6	
I	H	Warragul	[1],6,14,[25]	g,m	–	
I	J	Warri	17	k	1,7	
I	G	Washington	13,22	m,t	–	
IV	Z	Wassenaar	50	g,z <sub>51</sub>	–	
I	S	Waycross	41	z <sub>4</sub> ,z <sub>23</sub>	[e,n,z <sub>15</sub> ]	
IV	S	Waycross var. subsp. IV	41	z <sub>4</sub> ,z <sub>23</sub>	–	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	N	Wayne	30	g,z <sub>51</sub>	–	
I	M	Wedding	28	c	e,n,z <sub>15</sub>	
I	I	Welikade	16	l,v	1,7	
I	E <sub>1</sub>	Weltevreden	3,10	r	z <sub>6</sub>	IP combined Lanka 3,15:r:z <sub>6</sub> with Weltevreden to form Weltevreden 3,10,[15]:r:z <sub>6</sub> .
I	X	Wenatchee	47	b	1,2	
I	F	Wentworth	11	z <sub>10</sub>	1,2	
I	D <sub>2</sub>	Wernigerode	9,46	f,g	–	
I	T	Weslaco	42	z <sub>36</sub>	–	
I	D <sub>1</sub>	Westafrica	9,12	e,h	1,7	
I	I	Westeinde	16	l,w	1,6	
I	E <sub>4</sub>	Westerstede	1,3,19	l,z <sub>13</sub>	[1,2]	
I	E <sub>1</sub>	Westhampton	3,10	g,s,t	–	IP combined Halmstad (3,15:g,s,t:-) and Canoga (3,15,34:g,s,t:-) with Westhampton to form Westhampton 3,10,[15],[15,34]:g,s,t:-. Westhampton may possess H phase R <sub>Z37</sub> or R <sub>Z43</sub> or R <sub>Z45</sub> .
I	E <sub>1</sub>	Westminster	3,10,[15]	b	z <sub>35</sub>	CDC has no 3,10:b:z <sub>35</sub> .
I	I	Weston	16	e,h	z <sub>6</sub>	
II	E <sub>1</sub>	Westpark	3,10	l,z <sub>28</sub>	e,n,x	
I	O	Westphalia	35	z <sub>4</sub> ,z <sub>24</sub>	–	
I	E <sub>1</sub>	Weybridge	3,10	d	z <sub>6</sub>	
I	G	Wichita	1,13,23	d	1,6	Wichita may possess H phase R <sub>Z37</sub> .
I	O	Widemarsh	35	z <sub>29</sub>	–	
I	B	Wien	1,4,12,27	b	1,w	
I	C <sub>1</sub>	Wil	6,7	d	l,z <sub>13</sub> ,z <sub>28</sub>	
I	E <sub>3</sub>	Wildwood	3,15,34	e,h	1,w	IP combined Wildwood and Cambridge (3,15:e,h:l,w) with Meleagridis (3,10:e,h:l,w) to form Meleagridis 3,10,[15],[15,34]:e,h:l,w. Wildwood is now called Meleagridis var. O 15+, 34+ by IP.
I	B	Wilhelmsburg	1,4,[5],12,27	z <sub>38</sub>	[e,n,z <sub>15</sub> ]	
II	52	Wilhemstrasse	52	z <sub>44</sub>	1,5	IP combined Wilhemstrasse with Lobatsi (52:z <sub>44</sub> :1,5,7). The name Wilhemstrasse has been dropped.
I	G	Willemstad	1,13,22	e,h	1,6	
I	E <sub>1</sub>	Wilmington	3,10	b	z <sub>6</sub>	
I	E <sub>1</sub>	Wimborne	3,10	k	1,2	
II	E <sub>1</sub>	Winchester	3,10	z <sub>39</sub>	1,[5],7	
I	Q	Windermere	39	y	1,5	
II	W	Windhoek	45	g,m,s,t	1,5	
I	C <sub>2</sub>	Wingrove	6,8	c	1,2	
I	B	Winneba	4,12	r	1,6	
I	54	Winnipeg	54	e,h	1,5	
I	C <sub>1</sub>	Winston	6,7	m,t	1,6	
I	E <sub>4</sub>	Winterthur	1,3,19	l,z <sub>13</sub>	1,6	
I	C <sub>2</sub>	Wippra	6,8	z <sub>10</sub>	z <sub>6</sub>	
I	I	Wisbech	16	i	1,7	
II	J	Woerden	17	c	z <sub>39</sub>	
I	F	Wohlen	11	b	1,6	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
I	B	Womba	4,12, <u>27</u>	c	1,7	IP combined Womba with Altendorf (4,12:c:1,7) to form Altendorf 4,12,27:c:1,7. Womba is now called Altendorf var. O 27+. The name Womba has been dropped.
I	H	Woodhull	1,6,14,25	d	1,6	
I	F	Woodinville	11	c	e,n,x	
II	I	Woodstock	16	Z <sub>42</sub>	1,[5],7	
I	D <sub>2</sub>	Worb	9,46	b	e,n,x	
II	G	Worcester	1,13,23	m,t	e,n,x	
I	G	Worthington	<u>1</u> ,13,23	z	1,w	Worthington may possess H phase Rz <sub>45</sub> .
I	N	Wuiti	30	Z <sub>35</sub>	e,n,Z <sub>15</sub>	
I	D <sub>2</sub>	Wuppertal	9,46	Z <sub>41</sub>	–	
I	G	Wyldegreen	<u>1</u> ,13,23	a	1,w	
II	D <sub>1</sub>	Wynberg	<u>1</u> ,9,12	Z <sub>39</sub>	1,7	
I	E <sub>1</sub>	Yaba	3,10,[ <u>15</u> ]	b	e,n,Z <sub>15</sub>	
I	E <sub>4</sub>	Yalding	1,3,19	r	e,n,Z <sub>15</sub>	
I	B	Yaounde	<u>1</u> ,4,12, <u>27</u>	Z <sub>35</sub>	e,n,Z <sub>15</sub>	
I	M	Yardley	28	g,m	1,6	
I	C <sub>2</sub>	Yarm	6,8	Z <sub>35</sub>	1,2	
I	G	Yarrabah	13,23	y	1,7	
I	E <sub>1</sub>	Yeerongpilly	3,10	i	z <sub>6</sub>	
I	F	Yehuda	11	Z <sub>4</sub> ,Z <sub>24</sub>	–	
I	R	Yekepa	<u>1</u> ,40	Z <sub>35</sub>	e,n,Z <sub>15</sub>	
I	54	Yerba	54	Z <sub>4</sub> ,Z <sub>23</sub>	–	
I	P	Yoff	38	Z <sub>4</sub> ,Z <sub>23</sub>	1,2	
I	C <sub>2</sub>	Yokoe	8, <u>20</u>	m,t	–	
I	O	Yolo	35	c	[e,n,Z <sub>15</sub> ]	
I	W	Yopougon	45	z	e,n,Z <sub>15</sub>	
I	I	Yoruba	16	c	1,w	
I	C <sub>2</sub>	Yovokome	8, <u>20</u>	d	1,5	
I	E <sub>1</sub>	Yundum	3,10	k	e,n,x	
I	D <sub>2</sub>	Zadar	9,46	b	1,6	
I	D <sub>1</sub>	Zaiman	9,12	l,v	e,n,x	
I	N	Zaire	30	c	1,7	
I	E <sub>1</sub>	Zanzibar	3,10,[ <u>15</u> ]	k	1,5	
I	J	Zaria	17	k	e,n,Z <sub>15</sub>	
I	D <sub>1</sub>	Zega	9,12	d	z <sub>6</sub>	
I	N	Zehlendorf	30	a	1,5	
II	K	Zeist	18	Z <sub>10</sub>	z <sub>6</sub>	
I	C <sub>2</sub>	Zerifin	6,8	Z <sub>10</sub>	1,2	
I	I	Zigong	16	l,w	1,5	
I	V	Zinder	44	Z <sub>29</sub>	–	
I	E <sub>1</sub>	Zongo	3,10	Z <sub>35</sub>	1,7	
II	D <sub>3</sub>	Zuerich	1,9,12,46,27	c	Z <sub>39</sub>	
I	E <sub>4</sub>	Zuilen	1,3,19	i	1,w	
I	I	Zwickau	16	r,i	e,n,Z <sub>15</sub>	
II	B		<u>1</u> ,4,[5],12, <u>27</u>	a	e,n,x	
II	B		<u>1</u> ,4,12, <u>27</u>	a	Z <sub>39</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	B		4,12	b	1,5	
II	B		4,12	e,n,x	1,2,7	
II	B		4,[5],12	f,g,t	Z <sub>6</sub> ,Z <sub>42</sub>	
II	B		4,12	(f),g	–	Not in IP book
II	B		4,12	g,m,t	Z <sub>39</sub>	
II	B		4,12	g,m,t	–	IP calls this monophasic var. of Bechuana.
II	B		4,12	g,Z <sub>62</sub>	–	
II	B		4,12,27	i	Z <sub>35</sub>	
II	B		1,4,12,27	k	1,6	
II	B		1,4,12,27	l,v	e,n,x	
II	B		1,4,12,27	l,v	Z <sub>39</sub>	
II	B		4,12	l,Z <sub>28</sub>	–	
II	B		1,4,12,27	z	1,5	
II	B		4,12	z	1,7	
II	B		4,12	z	Z <sub>39</sub>	
II	B		4,12	–	1,6	
II	C <sub>1</sub>		6,7,14	a	1,5	
II	C <sub>1</sub>		6,7	a	Z <sub>6</sub>	
II	C <sub>1</sub>		6,7	b	Z <sub>39</sub>	
II	C <sub>1</sub>		6,7	d	Z <sub>42</sub>	
II	C <sub>1</sub>		6,7	g,m,[s],t	e,n,x	
II	C <sub>1</sub>		6,7	(g),m,[s],t	[1,5]	
II	C <sub>1</sub>		6,7	g,[m],s,t	[Z <sub>42</sub> ]	
II	C <sub>1</sub>		6,7	g,t	e,n,x:Z <sub>42</sub>	
IV	C <sub>1</sub>		6,7	g,Z <sub>51</sub>	–	
II	C <sub>1</sub>		6,7	k	[z6]	
IIIa	C <sub>1</sub>		6,7	(k)	z:[Z <sub>55</sub> ]	(Ar. 27:22:31:[37])
IIIb	C <sub>1</sub>		6,7	l,v	Z <sub>53</sub>	(Ar. 27:23:25)
II	C <sub>1</sub>		6,7	l,w	1,5,7	
II	C <sub>1</sub>		6,7	l,w	Z <sub>42</sub>	
II	C <sub>1</sub>		6,7	l,Z <sub>28</sub>	e,n,x	
II	C <sub>1</sub>		6,7	l,Z <sub>28</sub>	Z <sub>6</sub>	
II	C <sub>1</sub>		6,7	m,t	–	
II	C <sub>1</sub>		6,7	z	Z <sub>6</sub>	
II	C <sub>1</sub>		6,7	z	Z <sub>39</sub>	
II	C <sub>1</sub>		6,7	Z <sub>4</sub> ,Z <sub>24</sub>	Z <sub>42</sub>	
II	C <sub>1</sub>		6,7	Z <sub>10</sub>	Z <sub>35</sub>	
II	C <sub>1</sub>		6,7	Z <sub>29</sub>	–	
VI	C <sub>1</sub>		6,7	Z <sub>41</sub>	1,7	
II	C <sub>1</sub>		6,7	Z <sub>42</sub>	e,n,x:1,6	
II	C <sub>1</sub>		6,7	–	1,6	
II	C <sub>2</sub>		6,8	a	e,n,x	
II	C <sub>2</sub>		6,8	a	Z <sub>39</sub>	
II	C <sub>2</sub>		6,8	b	1,5	
II	C <sub>2</sub>		6,8	d	Z <sub>6</sub> :Z <sub>42</sub>	
II	C <sub>2</sub>		6,8	f,g,m,t	[e,n,x]	
II	C <sub>2</sub>		6,8	g,m,t	1,7	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	C <sub>2</sub>		6,8	l,v	e,n,x	
II	C <sub>2</sub>		6,8	l,w	z <sub>6</sub> :z <sub>42</sub>	
II	C <sub>2</sub>		6,8	l,z <sub>28</sub>	e,n,x	
II	C <sub>2</sub>		6,8	y	1,6:z <sub>42</sub>	
II	C <sub>2</sub>		6,8	z	1,5	
II	C <sub>2</sub>		6,8	z <sub>29</sub>	1,5	
II	C <sub>2</sub>		8	z <sub>29</sub>	e,n,x:z <sub>42</sub>	
II	C <sub>2</sub>		6,8	z <sub>29</sub>	e,n,x	
II	C <sub>2</sub>		6,8	–	1,5,7	
II	D <sub>1</sub>		9,12	a	1,5	
II	D <sub>1</sub>		<u>1</u> ,9,12	a	e,n,x	
II	D <sub>1</sub>		9,12	a	z <sub>39</sub>	
II	D <sub>1</sub>		<u>1</u> ,9,12	a	z <sub>42</sub>	
II	D <sub>1</sub>		9,12	d	z <sub>39</sub>	
II	D <sub>1</sub>		9,12	e,n,x	1,6	
II	D <sub>1</sub>		<u>1</u> ,9,12	g,m,[s],t	[1,5,7]:[z <sub>42</sub> ]	
II	D <sub>1</sub>		<u>1</u> ,9,12	g,z <sub>62</sub>	–	
II	D <sub>1</sub>		9,12	l,v	e,n,x	
II	D <sub>1</sub>		9,12	l,v	z <sub>39</sub>	
II	D <sub>1</sub>		9,12	l,z <sub>28</sub>	1,5:[z <sub>42</sub> ]	
II	D <sub>1</sub>		9,12	l,z <sub>28</sub>	e,n,x	
II	D <sub>1</sub>		<u>1</u> ,9,12	m,t	1,5	
II	D <sub>1</sub>		<u>1</u> ,9,12	m,t	z <sub>39</sub>	
II	D <sub>1</sub>		9,12	m,t	–	
II	D <sub>1</sub>		<u>1</u> ,9,12	z <sub>29</sub>	e,n,x	
II	D <sub>1</sub>		<u>1</u> ,9,12	z <sub>42</sub>	1,[5],7	
II	D <sub>2</sub>		9,46	e,n,x	1,5,7	
II	D <sub>2</sub>		9,46	g,z <sub>62</sub>	–	
II	D <sub>2</sub>		9,46	m,t	e,n,x	
II	D <sub>2</sub>		9,46	z	1,5	
II	D <sub>2</sub>		9,46	z <sub>10</sub>	z <sub>39</sub>	
II	D <sub>2</sub>		9,46	z <sub>10</sub>	z <sub>6</sub>	
II	D <sub>2</sub>		9,46	z <sub>39</sub>	1,7	
II	D <sub>3</sub>		9,12,46,27	g,t	e,n,x	
II	D <sub>3</sub>		1,9,12,46,27	l,z <sub>13</sub> ,z <sub>28</sub>	z <sub>39</sub>	
II	D <sub>3</sub>		1,9,12,46,27	y	z <sub>39</sub>	
II	D <sub>3</sub>		1,9,12,46,27	z <sub>10</sub>	1,5	
II	D <sub>3</sub>		1,9,12,46,27	z <sub>10</sub>	e,n,x	
II	D <sub>3</sub>		1,9,12,46,27	z <sub>10</sub>	z <sub>39</sub>	
II	D <sub>3</sub>		1,9,12,46,27	z <sub>4</sub> ,z <sub>24</sub>	[1,5]	
II	E <sub>1</sub>		3,10	a	z <sub>39</sub>	
II	E <sub>1</sub>		3,10	b	e,n,x	
II	E <sub>1</sub>		3,10	b	z <sub>39</sub>	
II	E <sub>1</sub>		3,10	d	e,n,x	
II	E <sub>1</sub>		3,10	l,v	e,n,x	
II	E <sub>1</sub>		3,10	l,z <sub>28</sub>	1,5	
II	E <sub>1</sub>		3,10	l,z <sub>28</sub>	z <sub>39</sub>	
II	E <sub>1</sub>		3,10	m,t	1,5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	E <sub>1</sub>		3,10	Z <sub>4</sub> ,Z <sub>24</sub>	–	
II	E <sub>1</sub>		3,10	Z <sub>29</sub>	e,n,x	
II	E <sub>1</sub>		3,10	Z <sub>29</sub>	–	
VI	F		11	a	1,5	
VI	F		11	b	1,7	
II	F		11	c	e,n,Z <sub>15</sub>	
IIIb	F		11	k	z53	(Ar. 17:29:25)
IIIb	F		11	l,v	z	(Ar. 17:23:31). May possess H phase Rz <sub>56</sub> (Ar. 38).
IIIb	F		11	l,v	z <sub>53</sub>	(Ar. 17:23:25)
II	F		11	z	e,n,x	
IIIa	F		11	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 17:1,2,5:-)
IV	F		11	Z <sub>4</sub> ,Z <sub>32</sub>	–	
II	F		11	–	1,5	
II	G		<u>1</u> ,13,23	a	1,5	
II	G		13,22	a	e,n,x	
II	G		<u>1</u> ,13,22	b	Z <sub>42</sub>	
II	G		13,23	d	e,n,x	
II	G		<u>1</u> ,13,23	d	e,n,Z <sub>15</sub>	
II	G		<u>1</u> ,13,23	g,m,s,t	1,5	
II	G		<u>1</u> ,13,23	g,m,s,t	Z <sub>42</sub>	
II	G		<u>1</u> ,13,23	g,[s],t	Z <sub>42</sub>	
IIIa	G		<u>1</u> ,13,23	g,51	–	(Ar. 18:13,14:-)
V	G		<u>1</u> ,13,22	i	–	
II	G		13,22	k	1,5:Z <sub>42</sub>	
II	G		13,23	k	Z <sub>41</sub>	
IIIb	G		13,22	l,v	1,5,7	(Ar. 18:23:30)
II	G		13,23	l,w	e,n,x	
II	G		13,22	l,Z <sub>28</sub>	1,5	
II	G		13,23	l,Z <sub>28</sub>	1,5	
II	G		13,23	l,Z <sub>28</sub>	Z <sub>6</sub>	
II	G		13,22	m,t	Z <sub>42</sub> :Z <sub>39</sub>	
V	G		13,22	r	–	
II	G		13,22	z	–	
II	G		<u>1</u> ,13,23	z	Z <sub>42</sub>	
IIIa	G		13,22	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 18:1,2,5:-)
IIIa	G		13,23	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 18:1,6,7:-). CDC would call this 1,6,7,9.
IIIa	G		<u>1</u> ,13,23	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 18:1,3,11:-)
II	G		<u>1</u> ,13,22	Z <sub>10</sub>	Z <sub>6</sub>	
II	G		<u>1</u> ,13,23	Z <sub>29</sub>	1,5	
II	G		<u>1</u> ,13,23	Z <sub>29</sub>	e,n,x	
II	G		<u>1</u> ,13,23	Z <sub>39</sub>	1,5,7	
II	G		13,22	Z <sub>39</sub>	1,7	
II	G		13,23	–	1,6	
VI	H		[1],6,14	a	1,5	
IIIb	H		(6),14	b	e,n,x,Z <sub>15</sub>	(Ar. 7a,7c:43:28)
II	H		6,14,[24]	k	1,6	
II	H		6,14	k	[e,n,x]	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIb	H		(6),14	k	z	(Ar. 7a,7c:29:31)
II	H		1,6,14	k	Z <sub>6</sub> :Z <sub>42</sub>	
IIIb	H		(6),14	k	Z <sub>53</sub>	(Ar. 7a,7c:29:25)
IIIb	H		(6),14	l,v	z	(Ar. 7a,7c:23:31)
IIIb	H		(6),14	l,v	Z <sub>35</sub>	(Ar. 7a,7c:23:21)
IIIb	H		(6),14	l,v	Z <sub>53</sub>	(Ar. 7a,7c:23:25)
IIIb	H		(6),14	r	z	(Ar. 7a,7c:24:31)
IV	H		6,14	Z <sub>4</sub> ,Z <sub>23</sub>	-	
IIIb	H		(6),14	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 7a,7c:27:28)
IIIb	H		(6),14	Z <sub>10</sub>	z:[Z <sub>53</sub> ]	(Ar. 7a,7c:27:31:[25])
IIIb	H		(6),14	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 7a,7c:27:25)
II	H		1,6,14	Z <sub>42</sub>	1,6	
IIIb	H		(6),14	Z <sub>52</sub>	e,n,x,Z <sub>15</sub>	(Ar. 7a,7c:26:28)
IIIb	H		(6),14	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 7a,7c:26:21)
II	I		16	b	e,n,x	
II	I		16	b	Z <sub>39</sub>	
II	I		16	b	Z <sub>42</sub>	
II	I		16	d	1,5	
II	I		16	g,[m],[s],t	Z <sub>42</sub>	
IIIb	I		16	i	Z <sub>35</sub>	(Ar. 25:33:21)
IIIb	I		16	k	z	(Ar. 25:29:31)
IIIb	I		16	k	Z <sub>53</sub>	(Ar. 25:29:25)
IIIb	I		16	(k)	Z <sub>35</sub>	(Ar. 25:22:21)
IIIb	I		16	l,v	1,5,7	(Ar. 25:23:30)
IIIb	I		16	l,v	z:[Z <sub>61</sub> ]	(Ar. 25:23:31:[41])
IIIb	I		16	l,v	Z <sub>35</sub>	(Ar. 25:23:21)
IIIb	I		16	l,v	Z <sub>53</sub>	(Ar. 25:23:25)
II	I		16	m,t	e,n,x	
II	I		16	z	Z <sub>42</sub>	
II	I		16	Z <sub>4</sub> ,Z <sub>24</sub>	-	
II	I		16	Z <sub>6</sub>	1,6	
IIIb	I		16	Z <sub>10</sub>	1,5,7	(Ar. 25:27:30)
IIIb	I		16	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 25:27:28)
II	I		16	Z <sub>29</sub>	1,5	
IV	I		16	Z <sub>36</sub>	-	
IIIb	I		16	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 25:26:21)
II	J		17	b	Z <sub>6</sub>	
II	J		17	e,n,x,Z <sub>15</sub>	1,[5],7	
II	J		17	g,m,s,t	-	
II	J		17	g,t	Z <sub>39</sub>	
IIIb	J		17	i	Z <sub>35</sub>	(Ar. 12:33:21)
IIIb	J		17	k	z	(Ar. 12:29:32)
II	J		17	k	-	
IIIb	J		17	l,v	e,n,x,Z <sub>15</sub>	(Ar. 12:23:28)
IIIb	J		17	l,v	Z <sub>35</sub>	(Ar. 12:23:21)
II	J		17	m,t	-	
IIIb	J		17	r	z	(Ar. 12:24:31)
II	J		17	y	-	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIa	J		17	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 12:1,2,5:- and 12:1,2,6:-)
IIIa	J		17	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 12:1,6,7,9:-)
IIIa	J		17	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 12:1,3,11:-)
IIIa	J		17	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 12:1,6,7:- and 12:1,7,8:-)
IIIb	J		17	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar.12:27:28). May possess H phase RZ <sub>56</sub> (Ar. 38).
IIIb	J		17	Z <sub>10</sub>	z	(Ar. 12:27:31)
IIIa	J		17	Z <sub>29</sub>	–	(Ar. 12:16,17,18:-)
IV	J		17	Z <sub>29</sub>	–	
IIIa	J		17	Z <sub>36</sub>	–	(Ar. 12:17,20:-)
IV	J		17	Z <sub>36</sub>	–	
IIIa	K		18	g,Z <sub>51</sub>	–	(Ar. 7a,7b:13,14:-)
IIIb	K		18	(k)	Z <sub>53</sub>	(Ar. 7a,7b:22:25)
IIIb	K		18	(k)	Z <sub>54</sub>	(Ar. 7a,7b:22:34)
IIIb	K		18	l,v	e,n,x,Z <sub>15</sub>	(Ar. 7a,7b:23:28)
IIIb	K		18	l,v	z	(Ar. 7a,7b:23:31)
IIIb	K		18	l,v	Z <sub>53</sub>	(Ar. 7a,7b:23:25)
II	K		18	m,t	1,5	
IIIb	K		18	r	z	(Ar. 7a,7b:24,31)
II	K		18	y	e,n,x,Z <sub>15</sub>	
II	K		18	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	K		18	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 7a,7b:1,2,5:- and 7a,7b:1,2,6:-)
II	K		18	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIa	K		18	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 7a,7b:1,6,7:- and 7a,7b:1,7,8:-)
IIIb	K		18	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 7a,7b:27:28)
IV	K		18	z36,z38	–	
II	L		21	b	1,5	
II	L		21	c	e,n,x	
IIIb	L		21	c	e,n,x,Z <sub>15</sub>	(Ar. 22:32:28)
II	L		21	g,[m],[s],t	–	
IIIa	L		21	g,Z <sub>51</sub>	–	(Ar. 22:13,14:-)
IV	L		21	g,Z <sub>51</sub>	–	
IIIb	L		21	i	1,5,7	(Ar. 22:33:30)
IIIb	L		21	i	e,n,x,Z <sub>15</sub>	(Ar. 22:33:28)
IIIb	L		21	k	e,n,x,Z <sub>15</sub>	(Ar. 22:29:28)
IIIb	L		21	k	z	(Ar. 22:29:31)
IIIb	L		21	l,v	z	(Ar. 22:23:31)
IIIb	L		21	l,v	Z <sub>57</sub>	(Ar. 22:23:40)
II	L		21	m,t	–	
IIIb	L		21	r	z	(Ar. 22:24:31)
II	L		21	z	–	CDC does not have this.
IIIa	L		21	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 22:1,2,5:- and 22:1,2,6:-)
IIIa	L		21	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 22:1,3,11:-)
IV	L		21	Z <sub>4</sub> ,Z <sub>32</sub>	–	
IIIb	L		21	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 22:27:28)
IIIb	L		21	Z <sub>10</sub>	z	(Ar. 22:27:31)
IIIb	L		21	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 22:27:25)
IIIa	L		21	Z <sub>29</sub>	–	(Ar. 22:16,17,18:-)
IV	L		21	Z <sub>36</sub>	–	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIb	L		21	Z <sub>65</sub>	e,n,x,Z <sub>15</sub>	(Ar. 22:32:28)
II	M		28	a	e,n,x	
II	M		28	b	e,n,x	
II	M		28	e,n,x	1,7	
II	M		28	g,m,t	e,n,x	
II	M		28	g,m,t	Z <sub>39</sub>	
II	M		28	g,s,t	e,n,x	
II	M		28	l,Z <sub>28</sub>	1,5	
II	M		28	m,t	[e,n,x]	
II	M		28	z	1,5	
IIIb	M		28	Z <sub>10</sub>	z	(Ar. 35:27:31)
IIIb	M		28	Z <sub>10</sub>	z:[Z <sub>57</sub> ]	(Ar. 35:27:31:[40])
II	M		28	Z <sub>29</sub>	1,5	
II	M		28	Z <sub>29</sub>	e,n,x	
II	N		30	b	Z <sub>6</sub>	
II	N		30	c	Z <sub>39</sub>	
II	N		30	g,m,s	e,n,x	
II	N		30	k	e,n,x,Z <sub>15</sub>	
II	N		30	l,Z <sub>28</sub>	Z <sub>6</sub>	
II	N		30	m,t	–	
II	N		30	Z <sub>6</sub>	1,6	
II	N		30	Z <sub>39</sub>	1,7	
II	O		35	d	1,5	
II	O		35	g,m,s,t	–	
II	O		35	g,t	1,5	
II	O		35	g,t	Z <sub>42</sub>	
IIIa	O		35	g,Z <sub>51</sub>	–	(Ar. 20:13,14:-)
IIIb	O		35	i	e,n,x,Z <sub>15</sub>	(Ar. 20:33:28)
IIIb	O		35	i	z	(Ar. 20:33:31)
IIIb	O		35	i	Z <sub>35</sub>	(Ar. 20:33:21)
IIIb	O		35	i	Z <sub>53</sub>	(Ar. 20:33:25)
IIIb	O		35	k	e,n,x,Z <sub>15</sub>	(Ar. 20:29:28)
IIIb	O		35	k	z	(Ar. 20:29:31)
IIIb	O		35	k	Z <sub>53</sub>	(Ar. 20:29:25). May possess H phase R <sub>Z50</sub> (Ar.42).
IIIb	O		35	(k)	z	(Ar. 20:22:31)
IIIb	O		35	(k)	Z <sub>35</sub>	(Ar. 20:22:21)
IIIb	O		35	l,v	1,5,7	(Ar. 20:23:30)
IIIb	O		35	l,v	e,n,x,Z <sub>15</sub>	(Ar. 20:23:28)
IIIb	O		35	l,v	Z <sub>35</sub>	(Ar. 20:23:21)
II	O		35	l,Z <sub>28</sub>	–	
II	O		35	m,t	–	
IIIb	O		35	r	e,n,x,Z <sub>15</sub>	(Ar. 20:24:28)
IIIb	O		35	r	z	(Ar. 20:24:31)
IIIb	O		35	r	Z <sub>35</sub>	(Ar. 20:24:21)
IIIb	O		35	r	Z <sub>61</sub>	(Ar. 20:24:41)
IIIa	O		35	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 20:1,2,6:-)
IIIa	O		35	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 20:1,7,8:-)
IIIb	O		35	Z <sub>10</sub>	Z <sub>35</sub>	(Ar. 20:27:21)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIa	O		35	Z <sub>29</sub>	–	(Ar. 20:16,17,18:-)
IIIa	O		35	Z <sub>36</sub>	–	(Ar. 20:17,20:-)
IIIb	O		35	Z <sub>52</sub>	1,5,7	(Ar. 20:26:30)
IIIb	O		35	Z <sub>52</sub>	e,n,x,Z <sub>15</sub>	(Ar. 20:26:28)
IIIb	O		35	Z <sub>52</sub>	z	(Ar. 20:26:31)
IIIb	O		35	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 20:26:21)
II	P		38	b	1,2	
IIIa	P		38	g,Z <sub>51</sub>	–	(Ar. 16:13,14:-)
IV	P		38	g,Z <sub>51</sub>	–	
IIIb	P		38	i	z	(Ar. 16:33:31)
IIIb	P		38	i	Z <sub>53</sub>	(Ar. 16:33:25)
IIIb	P		38	k	e,n,x,Z <sub>15</sub>	(Ar. 16:29:28)
IIIb	P		38	k	z	(Ar. 16:29:31)
IIIb	P		38	k	Z <sub>53</sub>	(Ar. 16:29:25)
IIIb	P		38	(k)	1,5,7	(Ar. 16:22:30)
IIIb	P		38	(k)	z	(Ar. 16:22:31)
IIIb	P		38	(k)	Z <sub>35</sub>	(Ar. 16:22:21). May possess H phase RZ <sub>56</sub> (Ar. 38).
IIIb	P		38	(k)	Z <sub>54</sub>	(Ar. 16:22:34)
IIIb	P		38	(k)	Z <sub>55</sub>	(Ar. 16:22:37)
IIIb	P		38	l,v	z	(Ar. 16:23:31)
IIIb	P		38	l,v	Z <sub>35</sub>	(Ar. 16:23:21)
IIIb	P		38	l,v	Z <sub>53</sub> : [Z <sub>54</sub> ]	(Ar. 16:23:25:[34])
IIIb	P		38	r	1,5,7	(Ar. 16:24:30)
IIIb	P		38	r	e,n,x,Z <sub>15</sub>	(Ar. 16:24:28)
IIIb	P		38	r	z: [Z <sub>57</sub> ]	(Ar. 16:24:31:[40])
IIIb	P		38	r	Z <sub>35</sub>	(Ar. 16:24:21)
IIIa	P		38	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 16:1,2,6:-)
IV	P		38	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIb	P		38	Z <sub>10</sub>	z	(Ar. 16:27:31)
IIIb	P		38	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 16:27:25)
IIIb	P		38	Z <sub>47</sub>	Z <sub>53</sub>	(Ar. 16:39:25)
IIIb	P		38	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 16:26:21)
IIIb	P		38	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 16:26:25)
IIIb	P		38	Z <sub>53</sub>	–	(Ar. 16:25:-). May possess H phase RZ <sub>50</sub> (Ar. 42) or RZ <sub>76</sub> (Ar. RZ <sub>76</sub> ). CDC does not have monophasic.
IIIa	P		38	Z <sub>61</sub>	–	(Ar. 16:41:-)
IIIb	P		38	Z <sub>61</sub>	Z <sub>53</sub>	(Ar. 16:41:25)
II	Q		39	a	Z <sub>39</sub>	
II	Q		39	c	e,n,x	
II	Q		39	e,n,x	1,7	
II	Q		39	g,m,t	–	
II	Q		39	l,v	1,5	
II	Q		39	l,Z <sub>28</sub>	Z <sub>39</sub>	
II	Q		39	m,t	e,n,x	
II	Q		39	–	1,7	
II	R		1,40	a	1,5	
II	R		1,40	a	Z <sub>6</sub>	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	R		40	b	–	
II	R		<u>1</u> ,40	c	Z <sub>39</sub>	
II	R		<u>1</u> ,40	e,n,x	1,[5],7	
II	R		<u>1</u> ,40	e,n,x,Z <sub>15</sub>	1,6	
II	R		<u>1</u> ,40	g,t	1,5	
II	R		<u>1</u> ,40	g,t	[e,n,x]	
II	R		<u>1</u> ,40	g,t	e,n,x,Z <sub>15</sub>	
II	R		40	g,t	Z <sub>39</sub>	
II	R		<u>1</u> ,40	g,[m],[s],t	Z <sub>42</sub>	
IIIb	R		40	g,Z <sub>51</sub>	[e,n,x,Z <sub>15</sub> ]	(Ar. 10a,10b:13,14:[28])
IIIb	R		40	i	1,5,7	(Ar. 10a,10b:33:30)
IIIb	R		40	k	z:Z <sub>57</sub>	(Ar. 10a,10b:29:31:40)
II	R		40	k	Z <sub>6</sub>	
IIIb	R		40	k	Z <sub>53</sub>	(Ar. 10a,10b:29:25)
IIIb	R		40	l,v	z	(Ar. 10a,10b,(10c):23:31)
IIIb	R		40	l,v	Z <sub>53</sub>	(Ar. 10a,10b:23:25)
II	R		<u>1</u> ,40	l,Z <sub>28</sub>	1,5:Z <sub>42</sub>	
II	R		<u>1</u> ,40	l,Z <sub>28</sub>	Z <sub>39</sub>	
II	R		40	m,t	Z <sub>39</sub>	
II	R		<u>1</u> ,40	m,t	Z <sub>42</sub>	
IV	R		40	m,t	–	
II	R		<u>1</u> ,40	z	Z <sub>6</sub>	
II	R		<u>1</u> ,40	z	Z <sub>39</sub>	
II	R		40	z	Z <sub>42</sub>	
IIIa	R		40	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 10a,10b:1,2,5:-; 10a,10b:1,2,5,6:-; and 10a,10b:1,2,6:-)
IIIa	R		40	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 10a,10b:1,3,11:-)
IV	R		40	Z <sub>4</sub> ,Z <sub>24</sub>	–	Also called Degania var. subsp. IV.
IIIa	R		40	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 10a,10b:1,2,10:-; 10a,10c:1,2,10:-; and 10a,10b:1,7,8:-)
IIIa	R		40	Z <sub>4</sub> ,Z <sub>32</sub>	–	
II	R		<u>1</u> ,40	Z <sub>6</sub>	1,5	
IIIb	R		40	Z <sub>10</sub>	Z <sub>35</sub>	(Ar. 10a,10b:27:21)
IIIa	R		40	Z <sub>29</sub>	–	(Ar. 10a,10b:16,18:-)
V	R		<u>1</u> ,40	Z <sub>35</sub>	–	
IIIa	R		40	Z <sub>35</sub>	–	(Ar. 10a,10b:17,20:-)
II	R		40	Z <sub>39</sub>	1,5:Z <sub>42</sub>	
II	R		40	Z <sub>39</sub>	1,7	
II	R		<u>1</u> ,40	Z <sub>42</sub>	1,6	
II	R		<u>1</u> ,40	[Z <sub>42</sub> ]	1,(5),7	
V	R		40	Z <sub>81</sub>	–	H Z <sub>81</sub> was formerly H a in <i>S. bongori</i> .
II	S		41	b	[1,5]	
VI	S		41	b	1,7	
IIIb	S		41	c	e,n,x,Z <sub>15</sub>	(Ar. 13:32:28)
II	S		41	c	[z6]	
II	S		41	g,m,s,t	Z <sub>6</sub>	
IIIa	S		41	g,Z <sub>51</sub>	–	(Ar. 13:13,14:-)
II	S		41	k	1,6	
II	S		41	k	[z6]	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIb	S		41	(k)	[z <sub>35</sub> ]	(Ar. 13:22:[21])
II	S		41	l,z <sub>13</sub> ,z <sub>28</sub>	e,n,x,z <sub>15</sub>	
IIIa	S		41	z <sub>4</sub> ,z <sub>23</sub>	–	(Ar. 13:1,2,5:- and 13:1,2,6:-)
IV	S		41	z <sub>4</sub> ,z <sub>23</sub>	–	Also called Waycross var. subsp. IV.
IIIa	S		41	z <sub>4</sub> ,z <sub>23</sub> ,z <sub>32</sub>	–	(Ar. 13:1,6,7,9:-)
IIIa	S		41	z <sub>4</sub> ,z <sub>24</sub>	–	(Ar. 13:1,3,11:-)
IIIa	S		41	z <sub>4</sub> ,z <sub>32</sub>	–	(Ar. 13:1,6,7:- and 13:1,7,8:-)
IIIa	S		41	z <sub>29</sub>	–	(Ar. 13:16,17,18:-)
IV	S		41	z <sub>29</sub>	–	
IIIa	S		41	z <sub>36</sub>	–	(Ar. 13:17,20:-)
IV	S		41	z <sub>52</sub>	–	
II	S		41	–	1,6	
II	T		42	b	z <sub>6</sub>	
II	T		42	d	z <sub>6</sub>	
II	T		42	[e,n,x]	1,6	
IIIa	T		42	g,z <sub>51</sub>	–	(Ar. 15:13,14:-)
IV	T		1,42	g,z <sub>51</sub>	–	
IIIb	T		42	k	–	(Ar. 15:29:-)
IIIb	T		42	k	e,n,x,z <sub>15</sub>	(Ar. 15:29:28)
IIIb	T		42	k	z	(Ar. 15:29:31)
IIIb	T		42	k	z <sub>35</sub>	(Ar. 15:29:21)
IIIb	T		42	(k)	z <sub>35</sub>	(Ar. 15:22:21)
IIIb	T		42	l,v	1,5,7	(Ar. 15:23:30)
IIIb	T		42	l,v	e,n,x,z <sub>15</sub>	(Ar. 15:23:28)
IIIb	T		42	l,v	z	(Ar. 15:23:31)
IIIb	T		42	l,v	z <sub>53</sub>	(Ar. 15:23:25)
II	T		1,42	l,w	e,n,x	
II	T		42	l,[z <sub>13</sub> ],z <sub>28</sub>	[z <sub>6</sub> ]	
II	T		42	m,t	[e,n,x,z <sub>15</sub> ]	
IIIb	T		42	r	–	(Ar. 15:24:-). May possess H phase Rz <sub>50</sub> (Ar. 42).
IIIb	T		42	r	z	(Ar. 15:24:31)
IIIb	T		42	r	z <sub>53</sub>	(Ar. 15:24:25)
IIIa	T		42	z <sub>4</sub> ,z <sub>23</sub>	–	(Ar. 15:1,2,5:- and 15:1,2,6:-)
IIIa	T		42	z <sub>4</sub> ,z <sub>24</sub>	–	(Ar. 15:1,3,11:-)
IV	T		1,42	z <sub>4</sub> ,z <sub>24</sub>	–	
II	T		42	z <sub>6</sub>	1,6	
IIIb	T		42	z <sub>10</sub>	–	(Ar. 15:27:-). May possess H phase Rz <sub>56</sub> (Ar. 38) and Rz <sub>50</sub> (Ar. 42).
II	T		42	z <sub>10</sub>	1,2	
II	T		42	z <sub>10</sub>	e,n,x,z <sub>15</sub>	
IIIb	T		42	z <sub>10</sub>	e,n,x,z <sub>15</sub>	(Ar. 15:27:28)
IIIb	T		42	z <sub>10</sub>	z	(Ar. 15:27:31)
II	T		42	z <sub>10</sub>	z <sub>6</sub>	
IIIb	T		42	z <sub>10</sub>	z <sub>35</sub>	(Ar. 15:27:21)
IIIb	T		42	z <sub>10</sub>	z <sub>67</sub>	(Ar. 15:27:46)
IV	T		42	z <sub>36</sub>	–	
IIIb	T		42	z <sub>52</sub>	z	(Ar. 15:26:31)
II	U		43	a	1,5	

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	U		43	a	Z <sub>6</sub>	
II	U		43	d	e,n,x,Z <sub>15</sub>	
II	U		43	d	Z <sub>39</sub>	
II	U		43	d	Z <sub>42</sub>	
II	U		43	e,n,x,Z <sub>15</sub>	1,(5),7	
II	U		43	e,n,x,Z <sub>15</sub>	1,6	
II	U		43	g,t	1,5	
IIIa	U		43	g,Z <sub>51</sub>	–	(Ar. 21:13,14:-)
IV	U		43	g,Z <sub>51</sub>	–	
II	U		43	g,Z <sub>62</sub>	e,n,x	
IIIb	U		43	k	z	(Ar. 21:29:31)
IIIb	U		43	l,v	z <sub>53</sub> :[Rz <sub>56</sub> ]	(Ar. 21:23:25:[38])
II	U		43	l,Z <sub>13</sub> ,Z <sub>28</sub>	1,5	
IIIb	U		43	r	e,n,x,Z <sub>15</sub>	(Ar. 21:24:28)
IIIb	U		43	r	z	(Ar. 21:24:31)
IIIb	U		43	r	z <sub>53</sub>	(Ar. 21:24:25)
II	U		43	z	1,5	
IV	U		43	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	U		43	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 21:1,2,5:- and 21:1,2,6:-)
IIIa	U		43	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 21:1,3,11:-)
IV	U		43	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IV	U		43	Z <sub>29</sub>	–	
II	U		43	Z <sub>29</sub>	e,n,x	
II	U		43	Z <sub>29</sub>	Z <sub>42</sub>	
IIIa	U		43	Z <sub>36</sub>	–	(Ar. 21:17,20:-)
IIIb	U		43	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 21:26:25)
II	V		44	e,n,x	1,6	
II	V		44	g,t	Z <sub>42</sub>	
IV	V		44	g,Z <sub>51</sub>	–	
II	V		44	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	V		44	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 1,3:1,2,5:- and 1,3:1,2,6:-)
IV	V		44	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	V		44	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 1,3:1,6,7,9:-)
IIIa	V		44	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 1,3:1,3,11:-)
IV	V		44	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIa	V		44	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 1,3:1,2,10:- and 1,3:1,7,8:-). IP calls Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub> , Ar. 1,2,10.
IV	V		44	Z <sub>29</sub>	–	
II	V		44	Z <sub>29</sub>	e,n,x:Z <sub>42</sub>	
IV	V		44	Z <sub>36</sub> ,[Z <sub>38</sub> ]	–	
V	V		44	Z <sub>39</sub>	–	
IIIa	W		45	g,Z <sub>51</sub>	–	(Ar. 11:13,14:-)
IV	W		45	g,Z <sub>51</sub>	–	
II	W		45	m,t	1,5	
II	W		45	z	1,5	
IIIa	W		45	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 11:1,2,5:-)
IV	W		45	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	W		45	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 11:1,3,11:-)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIa	W		45	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 11:1,7,8:-)
IIIa	W		45	Z <sub>29</sub>	–	(Ar. 11:16,18:-)
II	W		45	Z <sub>29</sub>	1,5	
II	W		45	Z <sub>29</sub>	e,n,x	
II	W		45	Z <sub>29</sub>	Z <sub>42</sub>	
IV	W		45	Z <sub>36</sub> ,Z <sub>38</sub>	–	
II	X		47	a	e,n,x,Z <sub>15</sub>	
II	X		47	b	Z <sub>6</sub>	
IIIb	X		47	c	1,5,7	(Ar. 28:32:30)
IIIb	X		47	c	e,n,x,Z <sub>15</sub> : [Z <sub>57</sub> ]	(Ar. 23:32:28 and 28:32:28:[40])
IIIb	X		47	c	z	(Ar. 28:32:31)
IIIb	X		47	c	Z <sub>35</sub>	(Ar. 28:32:21)
II	X		47	d	e,n,x,Z <sub>15</sub>	
II	X		47	e,n,x,Z <sub>15</sub>	1,6	
II	X		47	g,t	e,n,x	
IIIa	X		47	g,Z <sub>51</sub>	–	(Ar. 28:13,14)
IIIb	X		47	i	e,n,x,Z <sub>15</sub>	(Ar. 23:33:28). May possess H phase RZ <sub>50</sub> (Ar. 42).
IIIb	X		47	i	z	(Ar. 28:33:31)
IIIb	X		47	i	Z <sub>35</sub>	(Ar. 23:33:21 and 28:33:21)
IIIb	X		47	i	Z <sub>53</sub> : [Z <sub>57</sub> ]	(Ar. 23:33:25 and 28:33:25:[40])
IIIb	X		47	k	1,5,7	(Ar. 28:29:30)
IIIb	X		47	k	e,n,x,Z <sub>15</sub>	(Ar. 28:29:28)
IIIb	X		47	k	z	(Ar. 28:29:31)
IIIb	X		47	k	Z <sub>35</sub>	(Ar. 23:29:21)
IIIb	X		47	k	Z <sub>53</sub>	(Ar. 23:29:25)
IV	X		47	l,v	–	
IIIb	X		47	l,v	1,5,(7)	(Ar. 23:23:30). May possess H phase RZ <sub>50</sub> (Ar. 42).
IIIa	X		47	l,v	e,n,x,Z <sub>15</sub>	(Ar. 28:23:28)
IIIb	X		47	l,v	z	(Ar. 23:23:31)
IIIb	X		47	l,v	Z <sub>35</sub>	(Ar. 28:23:21)
IIIb	X		47	l,v	Z <sub>53</sub>	(Ar. 28:23:25)
IIIb	X		47	l,v	Z <sub>57</sub>	(Ar. 28:23:40)
IIIa	X		47	r	–	(Ar. 23:24:-). CDC does not have this.
IIIb	X		47	r	1,5,7	(Ar. 23:24:30)
IIIb	X		47	r	z	(Ar. 23:24:31)
IIIb	X		47	r	Z <sub>35</sub>	(Ar. 23:24:21 and 28:24:21)
IIIb	X		47	r	Z <sub>53</sub>	(Ar. 23:24:25). May possess H phase RZ <sub>74</sub> (Ar. RZ74)
IIIb	X		47	r	Z <sub>53</sub> : [Z <sub>60</sub> ]	(Ar. 23:24:25:[44]). May possess H phase RZ <sub>70</sub> and RZ <sub>72</sub> (Ar. RZ <sub>70</sub> or RZ <sub>72</sub> ).
IIIb	X		47	r	Z <sub>53</sub> :RZ <sub>50</sub> :Z <sub>60</sub>	(Ar. 28:24:25:42:44). Not in IP book.
II	X		47	z	e,n,x,Z <sub>15</sub>	
IIIa	X		47	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 28:1,2,5:-)
II	X		47	Z <sub>6</sub>	1,6	
IIIb	X		47	Z <sub>10</sub>	1,5,7	(Ar. 28:27:30)
IIIb	X		47	Z <sub>10</sub>	z	(Ar. 28:27:31)
IIIb	X		47	Z <sub>10</sub>	Z <sub>35</sub>	(Ar. 28:27:21)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIa	X		47	Z <sub>29</sub>	–	(Ar. 28:16,18:-)
II	X		47	Z <sub>29</sub>	e,n,x,Z <sub>15</sub>	
IV	X		47	Z <sub>36</sub>	–	
IIIb	X		47	Z <sub>52</sub>	1,5,7	(Ar. 28:26:30)
IIIb	X		47	Z <sub>52</sub>	e,n,x,Z <sub>15</sub>	(Ar. 28:26:28)
IIIb	X		47	Z <sub>52</sub>	z	(Ar. 28:26:31)
IIIb	X		47	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 28:26:21)
II	Y		48	a	Z <sub>6</sub>	
IIIb	Y		48	a	Z <sub>35</sub>	(Ar. 5:35:21). Not in 1992 IP, but is in Bergey.
II	Y		48	a	Z <sub>39</sub>	
V	Y		48	b	–	
II	Y		48	b	[Z <sub>6</sub> ]	
IIIb	Y		48	c	z	(Ar. 29:32:31)
II	Y		48	d	1,2	
IIIa	Y		48	g,Z <sub>51</sub>	–	(Ar. 5:13,14:-)
IIIb	Y		48	i	z:[Z <sub>72</sub> ]	(Ar. 5,29:33:31:[Z <sub>72</sub> ]). CDC does not have Z <sub>72</sub> strain.
IIIb	Y		48	i	Z <sub>35</sub> :[Z <sub>57</sub> ]	(Ar. 29:33:21:[40])
IIIb	Y		48	i	Z <sub>53</sub>	(Ar. 5:33:25)
IIIb	Y		48	i	Z <sub>61</sub>	(Ar. 5,29:33:41)
IIIb	Y		48	k	1,5,(7)	(Ar. 5:29:30)
II	Y		48	k	e,n,x,Z <sub>15</sub>	
IIIb	Y		48	k	e,n,x,Z <sub>15</sub>	(Ar. 5:29:28)
IIIb	Y		48	k	z	(Ar. 5,29:29:31)
IIIb	Y		48	k	Z <sub>35</sub> :[RZ <sub>75</sub> ]	(Ar. [5:29:21:RZ <sub>75</sub> ]). CDC does not have RZ <sub>75</sub> .
IIIb	Y		48	k	Z <sub>53</sub>	(Ar. 5,29:29:25)
IIIb	Y		48	(k)	Z <sub>53</sub>	(Ar. 5:22:25 and Ar. 5,29:22:25). Called 5:22:25 by IP.
IIIb	Y		48	l,v	1,5,(7)	(Ar. 5:23:30). May possess H phase RZ <sub>47</sub> or RZ <sub>50</sub> (Ar. 39 or 42).
IIIb	Y		48	l,v	z	(Ar. 5,29:23:31)
IIIb	Y		48	r	e,n,x,Z <sub>15</sub>	(Ar. 5:24:28)
IIIb	Y		48	r	z	(Ar. 5,29:24:31)
II	Y		48	z	1,5	
IIIa	Y		48	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 5:1,2,5:-; 5:1,2,5,6:-; and 5:1,6:-)
IV	Y		48	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	Y		48	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 5:1,6,7,9:-). IP calls this 5:1,6,7:-.
IIIa	Y		48	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 5:1,3,11:-)
IIIa	Y		48	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 5:1,6,7:-; 5:1,7,8:-; and Ar. 5:1,2,10:-). IP calls Z <sub>4</sub> ,Z <sub>32</sub> , Ar. 1,7,8; and would call Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub> , Ar. 1,2,10.
IV	Y		48	Z <sub>4</sub> ,Z <sub>32</sub>	–	
VI	Y		48	Z <sub>10</sub>	1,5	
IIIb	Y		48	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 5:27:28)
IIIb	Y		48	Z <sub>10</sub>	z	(Ar. 5,29:27:31)
II	Y		48	Z <sub>29</sub>	–	
IIIa	Y		48	Z <sub>29</sub>	–	(Ar. 5:16,18). This is not in IP book, but is on Rohde's list.
IV	Y		48	Z <sub>29</sub>	–	
IIIb	Y		48	Z <sub>35</sub>	Z <sub>52</sub>	(Ar. 5:21:26)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIa	Y		48	Z <sub>36</sub>	–	(Ar. 5,29:17,20:-)
IV	Y		48	Z <sub>36</sub> , [Z <sub>38</sub> ]	–	
V	Y		48	Z <sub>39</sub>	-	
IIIb	Y		48	Z <sub>52</sub>	e,n,x,z15	(Ar. 29:26:28)
IIIb	Y		48	Z <sub>52</sub>	z	(Ar. 5:26:31)
V	Y		48	Z <sub>65</sub>	-	
V	Y		48	Z <sub>81</sub>	-	
IV	Z		50	a	-	
IV	Z		50	b	-	
II	Z		50	b	Z <sub>6</sub>	
IV	Z		50	d	-	
II	Z		50	g,Z <sub>62</sub>	e,n,x	
IIIb	Z		50	i	1,5,7	(Ar. 9a,9c:33:30)
IIIb	Z		50	i	e,n,x,Z <sub>15</sub>	(Ar. 9a,9c:33:28)
IIIb	Z		50	i	z	(Ar. 9a,9c:33:31)
IIIb	Z		50	k	1,5,7	(Ar. 9a,9c:29:30)
IIIb	Z		50	k	e,n,x,Z <sub>15</sub>	(Ar. 9a,9c:29:28)
II	Z		50	k	e,n,x:Z <sub>42</sub>	
IIIb	Z		50	k	z	(Ar. 9a,9b:29:31 and 9a,9c:29:31). Ar. 9a,9b may possess H phase RZ <sub>50</sub> (Ar. 42).
IIIb	Z		50	k	Z <sub>35</sub>	(Ar. 9a,9b:29:21)
IIIb	Z		50	k	Z <sub>53</sub>	(Ar. 9a,9b:29:25 and 9a,9c:29:25). IP and Rohde only list the 9a,9c.
IIIb	Z		50	(k)	z	(Ar. 9a,9b:22:31)
IIIb	Z		50	(k)	Z <sub>35</sub>	(Ar. 9a,9b:22:21)
IIIb	Z		50	l,v	e,n,x,Z <sub>15</sub>	(Ar. 9a,9b:23:28)
IIIb	Z		50	l,v	z	(Ar. 9a,9b:23:31 and 9a,9c:23:31). IP only lists 9a,9c.
IIIb	Z		50	l,v	Z <sub>35</sub>	(Ar. 9a,9c:23:21)
II	Z		50	l,w	e,n,x,Z <sub>15</sub> :Z <sub>42</sub>	
II	Z		50	l,Z <sub>28</sub>	Z <sub>42</sub>	
IIIb	Z		50	r	1,5,(7)	(Ar. 9a,9b:24:30)
IIIb	Z		50	r	e,n,x,Z <sub>15</sub>	(Ar. 9a,9c:24:28)
IIIb	Z		50	r	z	(Ar. 9a,9b:24:31 and 9a,9c:24:31).
IIIb	Z		50	r	Z <sub>35</sub>	(Ar. 9a,9b:24:21). May possess H phase RZ <sub>58</sub> (Ar. RZ58). This is not in IP book, but is on Rohde's list.
IIIb	Z		50	r	Z <sub>53</sub>	(Ar. 9a,9b:24:25). May possess H phase RZ <sub>50</sub> (Ar. 42). This is not in IP book, but is on Rohde's list.
IIIa	Z		50	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 9a,9b:1,2,5:- and 9a,9b:1,2,6:-)
IIIa	Z		50	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 9a,9b:1,2,10:-). Called 9a,9b:1,6,7:- by IP and Rohde.
IIIa	Z		50	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 9a,9b:1,3,11:-)
IV	Z		50	Z <sub>4</sub> ,Z <sub>24</sub>	–	
IIIa	Z		50	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 9a,9b:1,2,10; 9a,9b:1,6,7:-; and 9a,9b:1,7,8:-). 9a,9b:1,2,10:- and 9a,9b:1,7,8:- used by IP and Rohde.
IIIb	Z		50	Z <sub>10</sub>	z	(Ar. 9a,9c:27:31). May possess H phase RZ <sub>56</sub> (Ar. 38).
IIIb	Z		50	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 9a,9c:27:25)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIa	Z		50	Z <sub>29</sub>	–	(Ar. 9a,9b:16,18:-)
IIIa	Z		50	Z <sub>36</sub>	–	(Ar. 9a,9b:17,20:-)
IIIb	Z		50	Z <sub>52</sub>	1,5,7	(Ar. 9a,9b:26:30 and 9a,9c:26:30)
IIIb	Z		50	Z <sub>52</sub>	Z	(Ar. 9a,9b:26:31 and 9a,9c:26:31)
IIIb	Z		50	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 9a,9b:26:21 and 9a,9c:26:21)
IIIb	Z		50	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 9a,9b:26:25 and 9a,9c:26:25)
IV	51		51	b	–	
II	51		51	c	–	
II	51		51	g,s,t	e,n,x	
IIIa	51		51	g,Z <sub>51</sub>	–	(Ar. 1,2:13,14:-)
IIIb	51		51	k	Z <sub>35</sub>	(Ar. 1,2:29:21)
IIIb	51		51	l,v	Z	(Ar. 1,2:23:31)
IIIa	51		51	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 1,2:1,2,5:- and 1,2:1,2,6:-)
IIIa	51		51	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 1,2:1,3,11:-)
IIIa	51		51	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 1,2:1,7,8:-)
II	51		51	Z <sub>29</sub>	e,n,x,Z <sub>15</sub>	
II	52		52	c	k	
IIIb	52		52	c	k	(Ar. 31:32:29). This is not in IP book, but is on Rohde's list.
II	52		52	d	e,n,x,Z <sub>15</sub>	
II	52		52	d	Z <sub>39</sub>	CDC does not have this.
II	52		52	g,t	–	
IIIb	52		52	k	Z <sub>35</sub>	(Ar. 31:29:21)
IIIb	52		52	k	Z <sub>53</sub>	(Ar. 31:29:25)
IIIb	52		52	(k)	Z <sub>35</sub>	(Ar. 31:22:21)
IIIb	52		52	l,v	Z <sub>53</sub>	(Ar. 31:23:25)
II	52		52	Z	Z <sub>39</sub>	
IIIb	52		52	Z <sub>52</sub>	Z	(Ar. 31:26:31)
II	53		53	c	1,5	
II	53		53	d	1,5	
II	53		1,53	d	Z <sub>39</sub>	
II	53		53	d	Z <sub>42</sub>	
IIIa	53		53	g,Z <sub>51</sub>	–	(Ar. 1,4:13,14:-)
IV	53		1,53	g,Z <sub>51</sub>	–	
IIIb	53		53	i	Z	(Ar. 1,4:33:31)
IIIb	53		53	k	e,n,x,Z <sub>15</sub>	(Ar. 1,4:29:28)
IIIb	53		53	k	Z	(Ar. 1,4:29:31)
IIIb	53		53	(k)	Z	(Ar. 1,4:22:31)
IIIb	53		53	(k)	Z <sub>35</sub>	(Ar. 1,4:22:21)
IIIb	53		53	l,v	e,n,x,Z <sub>15</sub>	(Ar. 1,4:23:28)
IIIb	53		53	l,v	Z <sub>35</sub>	(Ar. 1,4:23:21)
II	53		53	l,Z <sub>28</sub>	e,n,x	
II	53		53	l,Z <sub>28</sub>	Z <sub>6</sub>	
IIIb	53		53	r	Z	(Ar. 1,4:24:31)
IIIb	53		53	r	Z <sub>35</sub>	(Ar. 1,4:24:21)
IIIb	53		53	r	Z <sub>68</sub>	(Ar. 1,4:24:47). This was formerly called Z <sub>56</sub> (Ar. 38), but was changed to Z <sub>68</sub> (Ar. 47).
II	53		53	Z	1,5	
IIIb	53		53	Z	1,5,(7)	(Ar. 1,4:30:31)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	53		53	z	Z <sub>6</sub>	
IIIa	53		53	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 1,4:1,2,5:- and 1,4:1,2,6:-)
IV	53		53	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIa	53		53	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 1,4:1,6,7,9:-)
IIIa	53		53	Z <sub>4</sub> ,Z <sub>24</sub>	–	(Ar. 1,4:1,3,11:-)
IIIa	53		53	Z <sub>4</sub> ,Z <sub>32</sub>	–	(Ar. 1,4:1,6,7:-). IP combined this with 53:Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub> :- (Ar. 1,4:1,6,7,9:-).
IIIb	53		53	Z <sub>10</sub>	Z <sub>35</sub>	(Ar. 1,4:27:21)
IIIa	53		53	Z <sub>29</sub>	–	(Ar. 1,4:16,18:-)
IIIb	53		53	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 1,4:26:21)
IIIb	53		53	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 1,4:26:25)
II	56		56	d	–	
II	56		56	e,n,x	1,7	
II	56		56	l,v	Z <sub>39</sub>	
II	56		56	l,Z <sub>28</sub>	–	
II	56		56	z	Z <sub>6</sub>	
IIIa	56		56	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 14:1,2,5:- and 14:1,2,6:-)
IIIa	56		56	Z <sub>4</sub> ,Z <sub>23</sub> ,Z <sub>32</sub>	–	(Ar. 14:1,6,7,9:-)
II	56		56	Z <sub>10</sub>	e,n,x	
IIIa	56		56	Z <sub>29</sub>	–	(Ar. 14:16,18:-)
II	57		57	a	Z <sub>42</sub>	
IIIb	57		57	c	z:[Z <sub>60</sub> ]	(Ar. 34:32:31:[44])
II	57		57	d	1,5	
II	57		57	g,[m],s,t	Z <sub>42</sub>	
II	57		57	g,t	–	
IIIb	57		57	i	e,n,x,Z <sub>15</sub>	(Ar. 34:33:28)
IIIb	57		57	i	z	(Ar. 34:33:31)
IIIb	57		57	k	e,n,x,Z <sub>15</sub>	(Ar. 34:29:28). CDC does not have this and not on Rohde's list.
IV	57		57	Z <sub>4</sub> ,Z <sub>23</sub>	–	
IIIb	57		57	Z <sub>10</sub>	z	(Ar. 34:27:31)
II	58		58	a	Z <sub>6</sub>	
II	58		58	b	1,5	
II	58		58	c	Z <sub>6</sub>	
II	58		58	d	Z <sub>6</sub>	
IIIb	58		58	i	e,n,x,Z <sub>15</sub>	(Ar. 1,33:33:28)
IIIb	58		58	k	z	(Ar. 1,33:29:31)
IIIb	58		58	l,v	e,n,x,Z <sub>15</sub>	(Ar. 1,33:23:28)
IIIb	58		58	l,v	Z <sub>35</sub>	(Ar. 1,33:23:21)
II	58		58	l,Z <sub>13</sub> ,Z <sub>28</sub>	Z <sub>6</sub>	
IIIb	58		58	r	e,n,x,Z <sub>15</sub>	(Ar. 1,33:24:28)
IIIb	58		58	r	z	(Ar. 1,33:24:31)
IIIb	58		58	r	Z <sub>53</sub>	(Ar. 1,33:24:25). May possess H phase RZ <sub>47</sub> (Ar. 39) or RZ <sub>57</sub> (Ar. 40) or RZ <sub>70</sub> (Ar. RZ <sub>70</sub> ).
II	58		58	Z <sub>6</sub>	1,6	
II	58		58	Z <sub>10</sub>	1,6	
IIIb	58		58	Z <sub>10</sub>	e,n,x,Z <sub>15</sub>	(Ar. 1,33:27:28)
II	58		58	Z <sub>10</sub>	Z <sub>6</sub>	
IIIb	58		58	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 1,33:27:25). May possess H phase RZ <sub>50</sub> (Ar. 42).

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
II	58		58	Z <sub>39</sub>	e,n,x,Z <sub>15</sub>	
IIIb	58		58	Z <sub>52</sub>	z	(Ar. 1,33:26:31)
IIIb	58		58	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 1,33:26:21)
IIIb	59		59	c	e,n,x,Z <sub>15</sub>	(Ar. 19:32:28)
IIIb	59		59	i	e,n,x,Z <sub>15</sub>	(Ar. 19:33:28)
IIIb	59		59	i	z	(Ar. 19:33:31)
IIIb	59		59	i	Z <sub>35</sub>	(Ar. 19:33:21)
IIIb	59		59	k	Z <sub>53</sub>	(Ar. 19:29:25)
IIIb	59		59	(k)	e,n,x,Z <sub>15</sub>	(Ar. 19:22:28)
IIIb	59		59	(k)	z	(Ar. 19:22:31)
IIIb	59		59	(k)	Z <sub>35</sub>	(Ar. 19:22:21)
IIIb	59		59	l,v	z	(Ar. 19:23:31)
IIIb	59		59	l,v	Z <sub>53</sub>	(Ar. 19:23:25)
IIIb	59		59	r	Z <sub>35</sub>	(Ar. 19:24:21)
II	59		1,59	z	Z <sub>6</sub>	
IIIa	59		59	Z <sub>4</sub> ,Z <sub>23</sub>	–	(Ar. 19:1,2,5:- and 19:1,2,6:-)
IIIb	59		59	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 19:27:25)
IIIb	59		59	Z <sub>10</sub>	Z <sub>57</sub>	(Ar. 19:27:40)
IIIa	59		59	Z <sub>29</sub>	–	(Ar. 19:16,18:-)
IIIa	59		59	Z <sub>36</sub>	–	(Ar. 19:17,20:-)
IIIb	59		59	Z <sub>52</sub>	[Z <sub>53</sub> ]	(Ar. 19:26:[25])
II	60		60	b	[1,16]	
IIIb	60		60	i	–	(Ar. 24:33:-). May possess H phase Rz <sub>50</sub> (Ar. 42).
IIIb	60		60	i	e,n,x,Z <sub>15</sub>	(Ar. 24:33:28)
IIIb	60		60	i	Z <sub>35</sub>	(Ar. 24:33:21)
IIIb	60		60	k	z	(Ar. 24:29:31)
IIIb	60		60	k	Z <sub>35</sub>	(Ar. 24:29:21)
IIIb	60		60	(k)	Z <sub>53</sub>	(Ar. 24:22:25)
IIIb	60		60	l,v	z	(Ar. 24:23:31)
IIIb	60		60	r	e,n,x,Z <sub>15</sub>	(Ar. 24:24:28)
IIIb	60		60	r	z	(Ar. 24:24:31)
IIIb	60		60	r	Z <sub>35</sub>	(Ar. 24:24:21)
IIIb	60		60	r	Z <sub>53</sub>	(Ar. 24:24:25)
IIIb	60		60	Z <sub>10</sub>	z	(Ar. 24:27:31)
IIIb	60		60	Z <sub>10</sub>	Z <sub>35</sub>	(Ar. 24:27:21)
IIIb	60		60	Z <sub>10</sub>	Z <sub>53</sub>	(Ar. 24:27:25)
II	60		60	Z <sub>29</sub>	e,n,x	
V	60		60	Z <sub>41</sub>	–	
IIIb	60		60	Z <sub>52</sub>	1,5,[7]	(Ar. 24:26:30)
IIIb	60		60	Z <sub>52</sub>	z	(Ar. 24:26:31)
IIIb	60		60	Z <sub>52</sub>	Z <sub>35</sub>	(Ar. 24:26:21)
IIIb	60		60	Z <sub>52</sub>	Z <sub>53</sub>	(Ar. 24:26:25)
IIIb	61		61	c	1,5,(7)	(Ar. 26:32:30)
IIIb	61		61	c	Z <sub>35</sub>	(Ar. 26:32:21)
IIIb	61		61	i	e,n,x,Z <sub>15</sub>	(Ar. 26:33:28)
IIIb	61		61	i	z	(Ar. 26:33:31)
IIIb	61		61	i	Z <sub>35</sub>	(Ar. 26:33:21)
IIIb	61		61	i	Z <sub>53</sub>	(Ar. 26:33:25)

SUBSPECIES	O ANTIGEN GROUP	SEROTYPE	O ANTIGENS	PHASE 1	PHASE 2	NOTE
IIIb	61		61	k	1,5,(7)	(Ar. 26:29:30)
IIIb	61		61	k	z <sub>35</sub>	(Ar. 26:29:21). CDC does not have this.
IIIb	61		61	(k)	z <sub>53</sub>	(Ar. 26:22:25)
IIIb	61		61	l,v	1,5,7:[z <sub>57</sub> ]	(Ar. 26:23:30:[40])
IIIb	61		61	l,v	z	(Ar. 26:23:31)
IIIb	61		61	l,v	z <sub>35</sub>	(Ar. 26:23:21)
IIIb	61		61	r	1,5,7	(Ar. 26:24:30)
IIIb	61		61	r	z	(Ar. 26:24:31)
IIIb	61		61	r	z <sub>35</sub>	(Ar. 26:24:21)
IIIb	61		61	r	z <sub>53</sub>	(Ar. 26:24:25). May possess H phase Rz <sub>47</sub> (Ar. 39).
IIIb	61		61	z <sub>10</sub>	z <sub>35</sub>	(Ar. 26:27:21)
V	61		61	z <sub>35</sub>	-	
IIIb	61		61	z <sub>52</sub>	1,5,7	(Ar. 26:26:30)
IIIb	61		61	z <sub>52</sub>	z	(Ar. 26:26:31)
IIIb	61		61	z <sub>52</sub>	z <sub>35</sub>	(Ar. 26:26:21)
IIIb	61		61	z <sub>52</sub>	z <sub>53</sub>	(Ar. 26:26:25)
IIIa	62		62	g,z <sub>51</sub>	-	(Ar. 6:13,14:-)
IIIa	62		62	z <sub>4</sub> ,z <sub>23</sub>	-	(Ar. 6:1,2,5:-)
IIIa	62		62	z <sub>4</sub> ,z <sub>32</sub>	-	(Ar. 6:1,7,8:-)
IIIa	62		62	z <sub>29</sub>	-	(Ar. 6:17,18:-)
IIIa	62		62	z <sub>36</sub>	-	(Ar. 6:17,20:-)
IIIa	63		63	g,z <sub>51</sub>	-	(Ar. 8:13,14:-)
IIIa	63		63	z <sub>4</sub> ,z <sub>23</sub>	-	(Ar. 8:1,2,5:-)
IIIa	63		63	z <sub>4</sub> ,z <sub>32</sub>	-	(Ar. 8:1,7,8:-)
IIIa	63		63	z <sub>36</sub>	-	(Ar. 8:17,20:-)
IIIb	65		65	c	1,5,7	(Ar. 30:32:30)
IIIb	65		65	c	z	(Ar. 30:32:31)
IIIb	65		65	c	z <sub>53</sub>	(Ar. 30:32:25)
II	65		65	g,t	-	
IIIb	65		65	i	e,n,x,z <sub>15</sub>	(Ar. 30:33:28)
IIIb	65		65	(k)	z	(Ar. 30:22:31)
IIIb	65		65	(k)	z <sub>35</sub>	(Ar. 30:22:21)
IIIb	65		65	(k)	z <sub>53</sub>	(Ar. 30:22:25)
IIIb	65		65	l,v	e,n,x,z <sub>15</sub>	(Ar. 30:23:28)
IIIb	65		65	l,v	z	(Ar. 30:23:31)
IIIb	65		65	l,v	z <sub>35</sub>	(Ar. 30:23:21)
IIIb	65		65	l,v	z <sub>53</sub>	(Ar. 30:23:25)
IIIb	65		65	r	z <sub>35</sub>	(Ar. 30:24:21)
IIIb	65		65	z <sub>10</sub>	e,n,x,z <sub>15</sub>	(Ar. 30:27:28)
IIIb	65		65	z <sub>10</sub>	z	(Ar. 30:27:31)
IIIb	65		65	z <sub>52</sub>	e,n,x,z <sub>15</sub>	(Ar. 30:26:28)
IIIb	65		65	z <sub>52</sub>	z	(Ar. 30:26:31)
IIIb	65		65	z <sub>52</sub>	z <sub>35</sub>	(Ar. 30:26:21)
IIIb	65		65	z <sub>52</sub>	z <sub>53</sub>	(Ar. 30:26:25)
II	65		65	-	1,6	
V	66		66	z <sub>39</sub>	-	
V	66		66	z <sub>81</sub>	-	

## Bacto® Shigella Antiserum

### Shigella Antiserum Poly Group A · Shigella Antiserum Poly Group A<sub>1</sub> · Shigella Antiserum Poly Group B · Shigella Antiserum Poly Group C · Shigella Antiserum Poly Group C<sub>1</sub> · Shigella Antiserum Poly Group C<sub>2</sub> · Shigella Antiserum Poly Group D · Alkalescens-Dispar Antiserum Poly

#### Intended Use

Bacto Shigella Antiserum are used for identifying *Shigella* species by the slide agglutination test. Bacto Alkalescens-Dispar Antiserum Poly is used for identifying the Alkalescens-Dispar Group of microorganisms by the slide agglutination test.

#### Summary and Explanation

*Shigella* species cause the human diarrheal disease shigellosis (classic bacillary dysentery). The range of illness is from mild diarrhea to severe dysentery characterized by abdominal cramps and frequent passage of bloody, mucoid stools. While the disease is usually self-limiting, it can

be life threatening to the young, the elderly and malnourished persons. *Shigella* species are carried primarily in humans and are not generally distributed in nature. While transmission is usually direct person-to-person and through contaminated water supplies, foodborne outbreaks do occur.

The genus *Shigella* belongs to the family Enterobacteriaceae. *Shigella* species are facultatively anaerobic, gram-negative bacilli that typically are oxidase negative, lactose negative, H<sub>2</sub>S negative and non-gas producing. *Shigella* and *Escherichia* are genetically related. Certain strains of *E. coli* may resemble *Shigella* biochemically because both can be lactose-negative, nonmotile or non-gas-producing. These anaerogenic, nonmotile types have historically been called the Alkalescens-Dispar group and are presently classified as *E. coli*.

Serological testing with polyvalent and group specific antisera should be used to confirm the identification of isolates that are morphologically and biochemically identified as *Shigella* species. *Shigella* species are nonmotile, so serological identification is based on somatic (“O”) antigens. However, some strains have envelope antigens that prevent agglutination in somatic antisera. Heating the suspension at 100°C for 15-60 minutes destroys these interfering antigens. The four named species or serotypes of *Shigella* are *S. dysenteriae* (10 serovars), *S. flexneri* (six serovars), *S. boydii* (15 serovars) and *S. sonnei*. For a complete and current explanation of the classification of *Shigella*, consult appropriate references.<sup>1</sup>

The Alkalescens-Dispar group of microorganisms is currently recognized as anaerogenic, nonmotile biotypes of *E. coli*. Consult appropriate references for biochemical tests specific for differentiating these strains from *Shigella*.<sup>1-6</sup>

#### User Quality Control

##### Identity Specifications

###### Shigella Poly Antiserum Group A-D

Lyophilized appearance: Light gold to amber, button to powdered cake.

Rehydrated appearance: Light gold to amber, clear liquid.

###### Alkalescens-Dispar Antiserum Poly

Lyophilized appearance: Light gold to amber, button to powdered cake.

Rehydrated appearance: Light gold to amber, clear liquid.

##### Culture Response

Rehydrate Shigella Antiserum Poly Groups A-D and Alkalescens-Dispar Antiserum Poly per label directions. Perform the slide agglutination test using appropriate QC Antigens Shigella Group A-D or Alkalescens-Dispar.

SHIGELLA ANTISERUM	QC ANTIGEN	REACTION
Poly Group A	Shigella Group A	3+
Poly Group A <sub>1</sub>	Shigella Group A <sub>1</sub>	3+
Poly Group B	Shigella Group B	3+
Poly Group C	Shigella Group C	3+
Poly Group C <sub>1</sub>	Shigella Group C <sub>1</sub>	3+
Poly Group C <sub>2</sub>	Shigella Group C <sub>2</sub>	3+
Poly Group D	Shigella Group D	3+
Alkalescens-Dispar Antiserum Poly	Alkalescens-Dispar Group 1	3+

#### Principles of the Procedure

Identification of *Shigella* species includes the isolation of the microorganism, biochemical identification and serological confirmation. Serological confirmation involves the reaction in which the microorganism (antigen) reacts with its corresponding antibody. This *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction is rapid, does not dissociate (high avidity) and binds strongly (high affinity).

Because a microorganism (antigen) may agglutinate with antibodies produced in response to other species, heterologous reactions are possible. These are characterized as weak in strength or slow in formation.

Such unexpected and, perhaps, unpredictable reactions may lead to some confusion in serological identification. Therefore, a positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism. Homologous reactions occur rapidly and are strong. Heterologous reactions form slowly and are weak.

## Reagents

Shigella Antisera Poly and Alkalescens-Dispar Antiserum Poly are lyophilized, polyclonal rabbit antisera containing approximately 0.04% Thimerosal as a preservative.

Shigella Antisera Poly are absorbed when necessary to render each lot of serum as specific as practical. Antisera are absorbed to a certain point without reducing homologous reactions to an unsatisfactory level. They have been absorbed inter- and intra-specifically except that *Shigella* antisera are not prepared from or tested for:

- S. dysenteriae* provisional serotypes,
- S. flexneri* X and Y variants, or
- Alkalescens-Dispar Groups other than types 1-4.

ANTISERUM	REACTS WITH
Shigella Antiserum Poly Group A	<i>S. dysenteriae</i> types 1-7
Shigella Antiserum Poly Group A <sub>1</sub>	<i>S. dysenteriae</i> types 8ab, 8ac, 9, 10
Shigella Antiserum Poly Group B	<i>S. flexneri</i> types 1-6
Shigella Antiserum Poly Group C	<i>S. boydii</i> types 1-7
Shigella Antiserum Poly Group C <sub>1</sub>	<i>S. boydii</i> types 8-11
Shigella Antiserum Poly Group C <sub>2</sub>	<i>S. boydii</i> types 12-15
Shigella Antiserum Poly Group D	<i>S. sonnei</i> I and II
Alkalescens-Dispar Antiserum Poly	Alkalescens-Dispar Groups 1,2,3 and 4

When rehydrated and used as described, each vial of Shigella Antisera Poly and Alkalescens-Dispar Antiserum Poly contains sufficient reagent for 60 slide tests.

## Precautions

- For In Vitro Diagnostic Use.
- Shigella Antiserum Poly Group A**  
**Shigella Antiserum Poly Group A<sub>1</sub>**  
**Shigella Antiserum Poly Group B**  
**Shigella Antiserum Poly Group C**  
**Shigella Antiserum Poly Group C<sub>1</sub>**  
**Shigella Antiserum Poly Group C<sub>2</sub>**  
**Shigella Antiserum Poly Group D**  
**Alkalescens-Dispar Antiserum Poly**  
The Packaging of This Product Contains Dry Natural Rubber.
- Follow proper established laboratory procedure in handling and disposing of infectious materials.

## Storage

Store lyophilized and rehydrated Shigella Antisera Poly and Alkalescens-Dispar Antiserum Poly at 2-8°C. Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

## Expiration Date

The expiration date applies to the product in its intact container when stored as directed. Do not use a product if it fails to meet specifications for identity and performance.

Lyophilized Shigella Antisera Poly and Alkalescens-Dispar Poly are stable through the expiration date on the label when stored as described.

Rehydrated Shigella Antisera Poly and Alkalescens-Dispar Antiserum Poly that are cloudy or have a precipitate at any time during the period of use should be discarded.

## Procedure

### Materials Provided

Shigella Antisera Poly Group A  
Shigella Antisera Poly Group A<sub>1</sub>  
Shigella Antisera Poly Group B  
Shigella Antisera Poly Group C  
Shigella Antisera Poly Group C<sub>1</sub>  
Shigella Antisera Poly Group C<sub>2</sub>  
Shigella Antisera Poly Group D  
Alkalescens-Dispar Antiserum Poly

### Materials Required But Not Provided

Agglutination slides with 1 inch squares  
Applicator sticks  
Waterbath, boiling  
Sterile 0.85% NaCl solution  
QC Antigen Shigella Groups  
QC Antigen Alkalescens-Dispar Group 1

### Reagent Preparation

Equilibrate all materials to room temperature before performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergents.

**Shigella Antisera Poly and Alkalescens-Dispar Antiserum Poly:** To rehydrate, add 3 ml of sterile 0.85% NaCl solution and rotate gently to completely dissolve the contents. The rehydrated antiserum is considered a 1:2 working dilution. Subsequent dilutions are based on this as a starting dilution.

### Specimen Collection and Preparation

From clinical specimens, *Shigella* can be recovered on selective differential media such as Hektoen Enteric Agar or XLD Agar. For specific recommendations, consult appropriate references.<sup>2,3,4</sup> Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as a *Shigella* species. After these criteria are met, serological identification can be performed.

*Shigella* can be recovered from various types of foods when samples are processed to recover injured microorganisms and to prevent overgrowth of competing microorganisms. Consult appropriate references for recommended procedures when testing food samples.<sup>5,6</sup> After following an established protocol, determine that a pure culture of the microorganism has been obtained. Biochemical test reactions

should be consistent with the identification of the organism as a *Shigella* species. After these criteria are met, serological identification can be performed.

### Test Procedure

Use this procedure to test the isolate with each selected Shigella Antisera Poly or Alkalescens-Dispar Antiserum Poly.

1. **Shigella Antiserum:** Dispense 1 drop (35 µl) of the antiserum to be tested on an agglutination slide.
2. **Negative control:** Dispense 1 drop of sterile 0.85% NaCl solution on an agglutination slide.
3. **Test isolate:** Transfer a loopful of growth of the test organism to the drops of antisera and NaCl solution and mix thoroughly.
4. **Positive control:** Dispense 1 drop of the Shigella Antiserum to be tested on an agglutination slide. Add 1 drop of the appropriate QC Antigen.
5. Mix each reaction area with a separate applicator stick and rock for 1 minute. Read for agglutination.

### Results

1. Read and record results as follows:
  - 4+ 100% agglutination; background is clear to slightly hazy.
  - 3+ 75% agglutination; background is slightly cloudy.
  - 2+ 50% agglutination; background is moderately cloudy.
  - 1+ 25% agglutination; background is cloudy.
  - No agglutination.
2. **Positive control:** Should show 3+ or greater agglutination.
3. **Negative control:** Should show no agglutination. If autoagglutination occurs, tests results cannot be reported. To test for autoagglutination, transfer the isolate to selective medium.
4. **Test isolates:** 3+ or greater agglutination within 1-2 minutes is a positive result.
5. If no agglutination occurs or agglutination is weak, follow this procedure to remove blocking envelope antigens:
  - Prepare a dense suspension of the isolate from an agar medium in 3-5 ml of sterile 0.85% NaCl solution.
  - Heat in a boiling waterbath for 30-60 minutes and cool. The suspension should not show precipitation after heating. If this occurs, select another colony for testing.
  - Centrifuge at 1,000 rpm for 10-15 minutes.
  - Aspirate and discard the supernatant.
  - Resuspend the sediment in 0.5 ml sterile 0.85% NaCl solution.
  - Use a drop of the suspension and perform the slide agglutination test as outlined above.
6. A partial (less than 3+) or delayed agglutination reaction should be considered negative.
7. If test results for either the positive control or negative control are not as described, the test is invalid and results cannot be reported.

### Limitations of the Procedure

1. Correct interpretation of serological reactions depends on culture purity, morphological characteristics and biochemical reactions that are consistent with identification of the microorganism as a *Shigella* species.
2. Serological methods alone cannot identify the isolate as a *Shigella* species.
3. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent making a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False-positive reactions may occur.
4. Rough culture isolates do occur and will agglutinate spontaneously, causing agglutination of the negative control (autoagglutination). Smooth colonies must be selected and tested in serological procedures.
5. Shigella Antisera Poly and Alkalescens-Dispar Poly have been tested using cultures taken directly from agar media. These antisera have not been tested using antigen suspensions in NaCl solution or other diluents. If the user applies variations in the recommended steps, each lot of antiserum must be tested with known control cultures to verify expected reactions under the modified procedure.

### References

1. **Ewing, W.H. (ed.).** 1986. Edwards and Ewing's identification of Enterobacteriaceae, 4th ed. Elsevier Science Publishing Co., Inc., New York, NY.
2. **Gray, L. D.** 1995. *Escherichia, Salmonella, Shigella and Yersinia*, p. 450-456. In P. R. Murray, E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
3. **Baron, E. J., L. R. Peterson, and S. M. Finegold.** 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis, MO.
4. **Pezzlo, M. (ed.).** 1992. Aerobic bacteriology, p. 1.0.1-1.20.47. In H. D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
5. **Andrews, W. H., G. A. June, and P. S. Sherrod.** 1995. Shigella, p. 6.01-6.06. In FDA Bacteriological Analytical Manual, 8th ed. AOAC International, Gaithersburg, MD.
6. **Vanderzant, C., and D. F. Splittstoesser (eds.).** 1992. Compendium of methods for the microbiological examination of foods, 3rd edition. American Public Health Association, Washington, D.C.

### Packaging

Shigella Antiserum Poly Group A	3 ml	2834-47
Shigella Antiserum Poly Group A <sub>1</sub>	3 ml	2776-47
Shigella Antiserum Poly Group B	3 ml	2835-47
Shigella Antiserum Poly Group C	3 ml	2836-47
Shigella Antiserum Poly Group C <sub>1</sub>	3 ml	2777-47
Shigella Antiserum Poly Group C <sub>2</sub>	3 ml	2778-47
Shigella Antiserum Poly Group D	3 ml	2837-47
Alkalescens-Dispar Antiserum Poly	3 ml	2838-47

# Bacto® Streptococcus Antigens and Antisera

## Streptococcus Antiserum Group A · Streptococcus Antiserum Group B · Streptococcus Antigen Group A · Streptococcus Antigen Group B

### Intended Use

Bacto Streptococcus Antisera are used in the serological grouping of Group A and Group B streptococci by the capillary tube precipitin technique.

Bacto Streptococcus Antigens are used in the quality control testing of Bacto Streptococcus Antisera Groups A and B.

### Summary and Explanation

Streptococci are gram-positive cocci that are facultative anaerobes. They are catalase negative and may be alpha-, beta- or non-hemolytic.

*Streptococcus pyogenes* (Group A) is the most common cause of bacterial pharyngitis in children. Symptoms include fever, pharyngeal erythema and edema, tonsillar exudate and enlarged cervical lymph nodes. Physical findings alone cannot distinguish between Group A streptococcal pharyngitis and pharyngitis caused by other agents such as viruses or mycoplasma. Other infections caused by Group A streptococci include scarlet fever, impetigo and skin infections that range from mild to severe with toxic shock symptoms and tissue necrosis.

*Streptococcus agalactiae* (Group B streptococci) causes neonatal sepsis and meningitis. Other infections in children and adults include bacteremia, endocarditis and pneumonia.

Identification of Group A and Group B streptococci includes isolation of the microorganism and biochemical and serological identification. Serological identification involves the reaction in which the microorganism (antigen) reacts with its corresponding antibody. This *in vitro* reaction produces fine particles called precipitation.

### Principles of the Procedure

Beta-hemolytic streptococci Group A and Group B have carbohydrate group-specific antigens that can be extracted. Streptococcus Antisera and Streptococcus Antigens are used together in the capillary precipitin test to serologically identify the microorganisms.

### Reagents

**Streptococcus Antisera Groups A and B** are lyophilized, polyclonal rabbit antisera containing approximately 0.02% Thimerosal as a preservative. When rehydrated and used as described, each vial of Streptococcus Antisera contains sufficient reagent for 50 precipitin tests.

**Streptococcus Antigens Groups A and B** are ready-to-use cellular extracts of *S. pyogenes* and *S. agalactiae*, respectively, containing Thimerosal as a preservative. When used as described, each vial of Streptococcus Antigens contains sufficient reagent for 50 precipitin tests.

### Precautions

1. For In Vitro Diagnostic Use.
2. **Streptococcus Antiserum Group A**  
**Streptococcus Antiserum Group B**  
**Streptococcus Antigen Group A**  
**Streptococcus Antigen Group B**  
The Packaging of This Product Contains Dry Natural Rubber.
3. Follow proper established laboratory procedures in handling and disposing of infectious materials.
4. Streptococcus Antigens are not intended for use in the immunization of humans or animals.

### Storage

Store lyophilized and rehydrated Streptococcus Antisera at 2-8°C.

Store Streptococcus Antigens at 2-8°C.

Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.

### User Quality Control

#### Identity Specifications

##### Streptococcus Antisera Groups A and B

Lyophilized Appearance: Light gold to amber, button to powdered cake.

Rehydrated Appearance: Light gold to amber, clear liquid.

##### Streptococcus Antigens Groups A and B

Solution Appearance: Colorless to light yellow, clear liquid.

#### Quality Control Results

Rehydrate Streptococcus Antisera per label directions. Perform the capillary tube precipitin technique using appropriate Streptococcus Antigens.

ANTISERUM	ANTIGEN	REACTION
Streptococcus Antiserum Group A	Streptococcus Antigen Group A	Positive
Streptococcus Antiserum Group A	Streptococcus Antigen Group B	Negative
Streptococcus Antiserum Group B	Streptococcus Antigen Group B	Positive
Streptococcus Antiserum Group B	Streptococcus Antigen Group A	Negative

## Expiration Date

The expiration date applies to the product in its intact container when stored as directed. Do not use a product if it fails to meet specifications. Rehydrated Streptococcus Antiserum that is cloudy or has a precipitate anytime during use should be discarded.

Antigens must be smooth uniform suspensions. Examine antigen vials for precipitation before use. Suspensions with precipitation are not usable and should be discarded.

## Procedure

### Materials Provided

Streptococcus Antiserum Group A  
Streptococcus Antiserum Group B  
Streptococcus Antigen Group A  
Streptococcus Antigen Group B

### Materials Required but not Provided

Capillary tubes  
Sterile distilled or deionized water  
Plasticine block

### Reagent Preparation

Equilibrate all materials to room temperature before performing the tests. Ensure that glassware and pipettes are clean and free of residues such as detergent.

**Streptococcus Antisera:** To rehydrate, add 1 ml sterile distilled or deionized water and rotate gently to completely dissolve the contents. Streptococcus Antigens are ready to use.

### Specimen Collection and Preparation

Group A and Group B streptococci can be recovered on routine culture media such as sheep blood agar. For specific recommendations on isolation and presumptive identification, consult appropriate references.<sup>1,2</sup> Determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as a Group A or Group B *Streptococcus*. After these criteria are met, serological identification can be performed.

**Test antigen extract:** To prepare, extract the carbohydrate, group-specific antigen from a pure culture of the microorganism by the Lancefield hot HCl, autoclave, enzyme or other such method. For specific information on these methods, consult appropriate references.<sup>1,2</sup>

### Test Procedure

Add the antiserum to the capillary tube first so that it will be layered above the extract.

1. **Streptococcus Antiserum:** Dip a capillary tube into the antiserum and allow a column of 2-3 cm to rise into the tube.
2. Holding the forefinger on the top end of the capillary tube, remove the tube from the antiserum vial. Clean the tip with a lint-free tissue to remove excess antiserum. Do not allow air into the tube. If this occurs, discard the tube and begin again.
3. **Test antigen extract:** Dip the capillary tube into the prepared extract until the antiserum and the antigen come in contact with

each other. If an air bubble separates them, discard the tube and repeat steps 1-3.

4. Remove the tube from the extract and invert slightly to allow the column to move to the center of the tube.
5. Wipe excess fluid from the tube and insert in a plasticine block, antiserum end upward. Wipe the capillary tube so that it is free of fingerprints or any material that might interfere with a clear reading.
6. **Positive control:** Repeat steps 1-5, using (in step 3) a Streptococcus Antigen (Group A or B) that is *homologous* to the antiserum used in step 1.
7. **Negative control:** Repeat steps 1-5, using (in step 3) a Streptococcus Antigen (Group A or B) that is *not homologous* to the antiserum used in step 1.
8. Incubate all capillary tubes at  $22 \pm 2^\circ\text{C}$  for 5 minutes. Examine for the formation of a white precipitate at the interface of the antiserum and the antigen. Observe at 5 minute intervals for up to 30 minutes.

### Results

1. A strongly positive reaction develops within 5 minutes, a weaker reaction develops within 30 minutes.
2. Disregard any precipitate that appears after 30 minutes.
3. Precipitation in a tube indicates that the test antigen extract is homologous to the Streptococcus Antiserum Group A or Group B used.
4. Observe at 5 minute intervals within the 30 minute period because the precipitate may dissolve (prozone phenomenon).

### Limitations of the Procedure

1. Correct interpretation of serological reactions depends on culture purity as well as on morphological characteristics and biochemical reactions that are consistent with identification of the microorganism as *S. pyogenes* or *S. agalactiae*.
2. Serological methods alone cannot identify the isolate as *S. pyogenes* or *S. agalactiae*.

### References

1. Murray, P. R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.). 1995. Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
2. Isenberg, H. D. (ed.). 1992. Clinical microbiology procedures handbook, vol. 2. American Society for Microbiology, Washington, D.C.

### Packaging

Streptococcus Antiserum Group A	1 ml	2672-50
Streptococcus Antiserum Group B	1 ml	2741-50
Streptococcus Antigen Group A	1 ml	2978-50
Streptococcus Antigen Group B	1 ml	2979-50

# Bacto® USR Antigen

## Bacto USR Test Control Serum Set

### Intended Use

Bacto USR Antigen is nontreponemal antigen used in the Unheated Serum Reagin (USR) Test.<sup>1</sup>

Bacto USR Test Control Serum Set is standardized human sera used for controlling the USR Test.

### Summary and Explanation

*Treponema pallidum* is the causative agent of syphilis. Syphilis is a chronic infection with clinical manifestations that occur in distinct stages. Specific laboratory tests are recommended for the detection of each stage of the disease.

During the primary stage, treponemes present in the characteristic lesion, a chancre, are detectable by dark-field microscopy<sup>2</sup> or by the Direct Fluorescent Antibody Test for *T. pallidum* (DFA-TP). During the secondary stage, most serological tests for syphilis are reactive and treponemes may be found in the lesions by using dark-field microscopy. The latent period, which is asymptomatic, may last for years. Serological tests are usually reactive in the early latent period, but the reactivity of

non-treponemal tests decreases during the late latent period. Symptoms of the tertiary or late stage of syphilis may occur 10-20 years after initial infection. Approximately 71% of patients in the tertiary stage of syphilis have reactive non-treponemal tests.<sup>3,4</sup> In the tertiary stage, treponemal tests will usually be reactive and are the only basis for diagnosis. The lesions in tertiary syphilis will have few treponemes. Neurosyphilis and late cardiovascular syphilis are complications of tertiary syphilis.

Since the clinical manifestations of syphilis can be confused with other infectious or noninfectious conditions, proper diagnosis must include microscopic examination of lesion material and serological results.<sup>3</sup>

The USR Antigen is a nontreponemal antigen composed of cardioliipin, cholesterol and lecithin. The antigen is a modification of the VDRL Antigen emulsion in a USR suspending solution. The antigen suspension contains choline chloride, which enhances the reactivity of reagin in unheated serum.<sup>5</sup> The USR Test is suitable for qualitative as well as quantitative determinations.

Nontreponemal tests measure reagin, an antibody-like substance that can be detected in syphilitic serum. Reagin is also occasionally found in the serum of persons with other acute or chronic diseases. Reactive nontreponemal tests aid in the diagnosis of latent subclinical syphilis and are effective tools for detecting cases in epidemiological investigations. Nontreponemal tests are superior to treponemal tests for following the response to therapy.<sup>3</sup>

Nontreponemal antigen tests are not entirely specific for syphilis, nor do they have satisfactory sensitivity in all stages of syphilis. Whenever the results of a nontreponemal antigen test disagree with the clinical impression, a treponemal antigen test such as the Fluorescent Treponemal Antibody-Absorption (FTA-ABS)<sup>2,3</sup> should be performed. Nontreponemal tests such as the USR, RPR and VDRL tests are used to screen patient serum. Treponemal tests such as the FTA-ABS are used for confirmation.

The likelihood of obtaining a reactive USR Test result in various stages of untreated syphilis has been reported as follows:<sup>3</sup>

STAGES OF UNTREATED SYPHILIS	% REACTIVE
Primary	80
Secondary	100
Latent	95

### Principles of the Procedure

In the USR Test procedure, the patient's unheated serum is mixed with a buffered saline suspension of USR Antigen containing cardioliipin, lecithin and cholesterol. The combination of reagin and USR Antigen forms microscopic clumping called flocculation.

### Reagents

**USR Antigen** is 0.03% cardioliipin and 0.9% cholesterol dissolved in absolute alcohol with sufficient lecithin (approximately 0.2%) to produce standard reactivity. The antigen is suspended in a solution

### User Quality Control

#### Identity Specifications

##### USR Antigen

Appearance: Milky white, opaque suspension after gentle mixing.

##### Nontreponemal Antigen Reactive Serum

Lyophilized Appearance: White to cream colored, button to powdered cake.

Rehydrated Appearance: Light gold to light amber, clear to slightly opalescent.

##### USR Weakly Reactive Serum

Lyophilized Appearance: White to cream colored, button to powdered cake.

Rehydrated Appearance: Light gold to light amber, clear to slightly opalescent.

##### Nontreponemal Antigen Nonreactive Serum

Lyophilized Appearance: White to cream colored, button to powdered cake.

Rehydrated Appearance: Light gold to light amber, clear to slightly opalescent.

#### Performance Response

Rehydrate the sera contained in the USR Test Control Serum Set per label directions. Perform the USR Test according to the Test Procedure. Each serum in the USR Test Control Serum Set should yield appropriate reactions when tested with the USR Antigen.

Use the USR Antigen suspension only if it produces the expected reactivity with the control sera.

containing EDTA, choline chloride and phosphate with 0.2% Thimerosal as a preservative.<sup>7,8</sup>

**USR Test Control Serum Set** contains 3 ml, each, of the following lyophilized human sera: Nontreponemal Antigen Reactive Serum, USR Weakly Reactive Serum and Nontreponemal Antigen Nonreactive Serum. These reagents are standardized to provide reactive, weakly reactive and nonreactive readings, respectively, when tested according to the USR Test procedure.

## Precautions

1. For In Vitro Diagnostic Use.
2. **WARNING! POTENTIAL BIOHAZARDOUS REAGENTS.** Each donor unit used in preparation of USR Antigen and USR Test Control Serum Set was tested by an FDA approved method for the presence of the antibody to human immunodeficiency virus (HIV) and for hepatitis B surface antigen and found negative (were not repeatedly reactive).

Because no test method can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen.<sup>9</sup>

3. **USR Test Control Serum Set**  
The Packaging of This Product Contains Dry Natural Rubber.
4. Observe universal blood and body fluid precautions in handling and disposing of specimens.<sup>10,11</sup>
5. Follow proper established laboratory procedure in handling and disposing of infectious materials.

## Storage

Store USR Antigen at 2-8°C. If the original 3 ml quantity exceeds what is needed for one testing period, transfer the remainder from the first day's use to one or more aliquot vials and store at 2-8°C.

Store the lyophilized sera in the USR Test Control Serum Set at 2-8°C. Store the rehydrated control sera at 2-8°C or divide into aliquots sufficient for one day of testing and store at -20°C. Do not thaw and refreeze.

## Expiration Date

The expiration date applies to the product in its intact container when stored as directed. Do not use a product if it fails to meet specifications for identity and performance.

## Procedure

### Materials Provided

USR Antigen  
USR Test Control Serum Set

### Materials Required But Not Provided

0.9% saline  
Nondisposable glass syringe, 1-2 cc  
Nondisposable calibrated 18-gauge needles without bevel  
Micropipettor, 50 µl  
Pipettes, serological, graduated to tip:  
1.0 ml, graduated in 1/100 ml  
5.0 ml, graduated in 1/10 ml  
10.0 ml, graduated in 1/10 ml

Slides, 2 x 3 inches with paraffin or ceramic rings approximately 14 mm in diameter and high enough to prevent spillage during rotation.

Slide holder for 2 x 3 inch slides

Mechanical rotator, adjustable to 180 ± 2 rpm circumscribing a circle 19 mm in diameter on a horizontal plane.

Light microscope with 10X ocular and 10X objective

Sterile distilled or deionized water

Absolute alcohol

Acetone

Timer

## Reagent Preparation

USR Antigen is ready to use.

Equilibrate all materials to room temperature (23-29°C) before performing the tests. Ensure that all glassware and pipettes are clean and free of detergent residues.

**USR Test Control Serum Set:** To rehydrate the control sera, add 3 ml sterile distilled or deionized water and rotate gently to completely dissolve the contents.

## Specimen Collection and Preparation

Collect a blood specimen by aseptic venipuncture into a clean, dry tube without anticoagulant. After the specimen has clotted, centrifuge the specimen at 1,500-2,000 rpm for five minutes to obtain test serum. Store serum specimens at room temperature for no longer than 4 hours; for prolonged storage, keep at 2-8°C for up to 5 days or maintain below -20°C. Serum specimens must be clear, free of hemolysis and show no visible evidence of bacterial contamination, such as turbidity or particulate matter. Refer to appropriate references for more information on collection of specimens.<sup>1,3,12</sup>

## Test Procedure

### Preparation of Specific Glassware

Syringes with needles:

1. Prerinse with tap water.
2. Soak and hand wash thoroughly in a glassware detergent solution.
3. Rinse with tap water 6-8 times.
4. Rinse with unused distilled or deionized water.
5. Rinse with absolute alcohol.
6. Rinse with acetone.
7. Air dry until the acetone odor is completely eliminated.
8. Remove needles from syringes for storage.

Ceramic-ringed slides:

1. Prerinse with tap water.
2. Wash with a glassware detergent solution. Avoid prolonged soaking of ceramic-ringed slides in detergent solution because the ceramic rings will become brittle and flake off.
3. Rinse with tap water 3-4 times.
4. Rinse with unused distilled or deionized water.
5. Wipe dry with a clean lint-free cloth. If cleaned slides do not allow serum to spread evenly within the inner surface of the circle, proceed as follows.
6. Scrub the slides with a nonscratching cleanser.
7. Rinse, dry and polish with a clean, lint-free cloth.

**Testing the Accuracy of the Antigen Suspension Needle**

1. The accuracy of the test depends on the amount of antigen suspension used. Check the calibration of the needle periodically to ensure delivery of the correct volume of USR Antigen suspension.
2. For the qualitative and quantitative tests on serum, dispense the antigen suspension from a syringe fitted with an 18-gauge needle without bevel that will deliver  $45 \pm 1$  drops (22 ml) of antigen suspension per ml when held vertically.
3. Place the needle on a 1-2 ml syringe. Fill the syringe with 1 ml of USR Antigen suspension. Holding the syringe in a vertical position, count the number of drops delivered in 1.0 ml. The needle is correctly calibrated if  $45 \pm 1$  drops are delivered in 1.0 ml.
4. Adjust or replace the needle if it does not meet this specification. Repeat calibration on the new or adjusted needle.

**Testing and Storing the USR Antigen Suspension**

1. Store the antigen suspension at 2-8°C.
2. For daily use, withdraw a sufficient amount of suspension for 1 day's testing and store the remainder at 2-8°C. Antigen suspensions must be at room temperature (23-29°C) before use.
2. Test antigen suspension reactivity with the Reactive, Weakly Reactive and Nonreactive control sera. Use the antigen suspension only if it produces the expected reactivity with the control sera.
3. After each day of use, clean the dispensing needle, bottle and syringe by rinsing with water, alcohol and acetone, in that order. Remove the needle from the syringe after cleaning.

**USR Qualitative Slide Test on Serum<sup>1,13</sup>**

For reliable and reproducible test results, the USR Antigen suspension, controls and test specimens must be at 23-29°C when tests are performed.

1. Pipette 50 µl of unheated serum into one ring of a paraffin- or ceramic-ringed slide using a safety pipetting device. Do not use a glass slide with concavities, wells or glass rings. Spread the serum with a circular motion of the pipette tip so that the serum covers the entire inner surface of the paraffin or ceramic ring. Include control sera when performing the test.
2. Gently resuspend the USR Antigen and withdraw the desired quantity with a syringe and needle.
3. Hold the syringe and needle containing the USR Antigen suspension in a vertical position. Dispense several drops to clear the needle of air. Add exactly 1 free-falling drop (22 µl) of antigen suspension to each circle containing serum. Do not allow the needle to touch the serum.
4. Place the slide on the mechanical rotator. Rotate the slide for 4 minutes at  $180 \pm 2$  rpm. If the environment is dry, cover the slides with a moist, humidifying cover during rotation to prevent excessive evaporation.
5. Immediately after rotating the slide, remove the slide from the rotator and read the test results microscopically, using a 10X ocular and a 10X objective.

**Results – Qualitative Test**

1. Read and record results as follows:  
Medium to large clumps - Reactive (R)  
Small clumps - Weakly reactive (WR)  
No clumping or very slight roughness - Nonreactive (N)

2. Verify that the control sera results are as expected. If reactions are not as expected, the test is invalid and results cannot be reported.
3. Perform a quantitative test on all serum specimens that produce Reactive, Weakly Reactive or “rough” Nonreactive results, since prozone reactions are occasionally encountered.

**USR Quantitative Test<sup>1,13</sup> on Serum**

1. To quantitate serum samples to an endpoint titer, prepare serum dilutions on the slide at 1:1, 1:2, 1:4 and 1:8, as follows.
2. Dispense 50 µl of 0.9% saline in circles 2-4. Do not spread the saline.
3. Dispense 50 µl of serum in circles 1 and 2.
4. Mix the saline and the serum in circle 2 by drawing the mixture up and down in the pipette 8 times. Mix gently to prevent bubbles.
5. Transfer 50 µl from circle 2 (1:2) to circle 3 and mix.
6. Transfer 50 µl from circle 3 (1:4) to circle 4 (1:8), mix, and then discard 50 µl from circle 4.
7. Holding the syringe and needle containing the USR Antigen suspension in a vertical position, dispense several drops to clear the needle of air. Then, add exactly 1 free-falling drop (22 µl) of antigen suspension to each circle containing serum. Do not allow the needle to touch the serum.
8. Place the slide on the mechanical rotator. Rotate the slide for 4 minutes at  $180 \pm 2$  rpm. If the environment is dry, cover the slides with a moist, humidifying cover during rotation to prevent excessive evaporation.
9. Immediately after rotating the slide, remove the slide from the rotator and read the test results microscopically using a 10X ocular and a 10X objective.
10. If the highest dilution tested (1:8) is reactive, prepare a 1:8 dilution of the test specimen by adding 0.1 ml of serum to 0.7 ml of 0.9% saline. Mix thoroughly. Retest as in steps 1-9, above.

**Results – Quantitative Test**

Report the titer as the highest dilution that produces a Reactive result.

**Table 1.** Sample quantitative USR Test results.

Undiluted (1:1)	1:2	1:4	1:8	1:16	1:32	
R	W	N	N	N	N	Reactive, undiluted
R	R	W	N	N	N	Reactive, 1:2 dilution
R	R	R	W	N	N	Reactive, 1:4 dilution
W	W	R	R	W	N	Reactive, 1:8 dilution
N (rough)	W	R	R	R	N	Reactive, 1:16 dilution
W	N	N	N	N	N	Weakly reactive, undiluted

If reactive results are obtained through dilution 1:32, prepare further twofold serial dilutions in 0.9% saline (1:64, 1:128 and 1:256) and retest using the quantitative test procedure.

### Interpretation

1. The results of the serum USR Test must be confirmed by a treponemal test.
2. The diagnosis of syphilis depends on the results of the USR Test, treponemal confirmatory test, clinical signs and symptoms, and risk factors.
3. A Reactive USR Test may indicate past or present infection with a pathogenic treponeme. However, it may be a false-positive reaction. A false positive is determined if the confirmatory treponemal test is negative.
4. A Nonreactive USR Test with clinical evidence of syphilis may indicate early, primary syphilis, a prozone reaction in secondary syphilis, or late syphilis.
5. A Nonreactive USR Test with no clinical evidence of syphilis indicates no current infection or an effectively treated infection.
6. A quantitative USR Test detects changes in reagin titer. Therefore, a serum specimen showing a fourfold increase in titer on a repeat specimen may indicate an infection, a reinfection or a treatment failure. Likewise, a fourfold decrease during treatment indicates adequate syphilis therapy.

### Limitations of the Procedure

1. A prozone reaction may occur in which reactivity with true positive undiluted serum is inhibited. The prozone phenomenon often gives Weakly Reactive or “rough” Nonreactive results in the qualitative test. Specimens with such nonreactive results must be quantitatively tested.
2. Biological false-positive reactions can occur with nontreponemal tests in persons who abuse drugs, have diseases such as lupus erythematosus, mononucleosis, malaria, leprosy or viral pneumonia, or who have recently been immunized.<sup>1</sup>
3. Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.
4. If the temperature of the testing area, specimens, or reagents is less than 23°C, test reactivity is decreased. If the temperature is greater than 29°C, test reactivity is increased.<sup>1</sup>
5. Test results are unpredictable when testing hemolyzed, contaminated or extremely turbid serum specimens.
6. Test results may be erroneous if the speed and time of rotation are not correct.
7. Positive results obtained by using USR Antigen should not be considered as conclusive evidence that the patient is syphilitic. Conversely, a nonreactive USR Test, by itself, does not rule out the diagnosis of syphilis.

### References

1. **Larsen, S. A., E. F. Hunter, and S. J. Kraus.** 1990. A manual of tests for syphilis. American Public Health Association.
2. **Creighton, E. T.** 1990. Dark field microscopy for the detection and identification of *Treponema pallidum*, p. 49-61. In S. A. Larsen, E. F. Hunter, and S. J. Kraus (ed.), Manual of tests for syphilis, 8th ed. American Public Health Association, Washington, D.C.

3. **Janda, W. M. (ed.)**. 1994. Immunology, p. 9.7.1-9.7.20. In H. D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 2. American Society for Microbiology, Washington, D.C.
4. **Norris, S. J., and S. A. Larsen.** 1995. *Treponema* and other host-associated spirochetes, p. 636-651. In P. R. Murray, E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
5. **Perine, P. L., A. L. Wallace, J. H. Blount, and S. T. Brown.** 1981. Syphilis, p. 631-673. In A. Ballows and W. J. Hausler, Diagnostic procedures for bacterial, mycotic and parasitic infections. American Public Health Association, Washington, D.C.
6. **Matthews, H. M., T. K. Yang, and H. M. Jenkin.** 1979. Unique lipid composition of *Treponema pallidum* (Nichols virulent strain). Infect. Immun. **24**:713-719.
7. **Portnoy, J., and W. Garson.** 1960. New and improved antigen suspension for rapid reagin test for syphilis. Public Health Rep. **75**:985-988.
8. **Portnoy, J., H. W. Bossak, V. H. Falcone, and A. Harris.** 1961. A rapid reagin test with unheated serum and new improved antigen suspension. Public Health Rep. **76**:933-935.
9. **U. S. Department of Health and Human Services.** 1988. Biosafety in microbiological and biomedical laboratories, 2nd ed. U. S. Department of Health and Human Services publication no. 88-8395. U. S. Government Printing Office, Washington, D.C.
10. **Centers for Disease Control.** 1988. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. Morbidity and Mortality Weekly Reports **37**:377-382, 387-388.
11. **Occupational Safety and Health Administration, U.S. Department of Labor.** 1991. 29 CFR, part 1910. Occupational exposure to bloodborne pathogens, final rule. Federal Register **56**:64175-64182.
12. **Miller, J. M., and H. T. Holmes.** 1995. Specimen collection, transport and storage. In P. R. Murray, E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
13. **Pettit, D. E., S. A. Larsen, V. Pope, M. D. Perryman, and M. R. Adams.** 1982. Unheated serum reagin test as a quantitative test for syphilis. J. Clin Microbiol. **15**:238-242.

### Packaging

USR Antigen	6 x 3 ml	2405-46
USR Test Control Serum Set	1 set	3516-32
Contains:		
Nontreponemal Antigen		
Reactive Serum	3 ml	
USR Weakly Reactive Serum	3 ml	
Nontreponemal Antigen		
Nonreactive Serum	3 ml	
Aliquant Vials	3 vials	

# Bacto® VDRL Antigen

## Bacto VDRL Test Control Serum Set

### Intended Use

Bacto VDRL Antigen with Bacto VDRL Buffered Saline is used in the Venereal Disease Research Laboratory (VDRL)<sup>1</sup> Test for detecting reagin, an antibody-like substance, by the qualitative and quantitative slide procedures.

Bacto VDRL Test Control Serum Set is used for controlling the VDRL Test.

### Summary and Explanation

*Treponema pallidum* is the causative agent of syphilis. Syphilis is a chronic infection with many clinical manifestations which occur in distinct stages. Specific laboratory tests are recommended for the detection of each stage of the disease.<sup>2-4</sup>

The VDRL Antigen is a nontreponemal antigen composed of cardiolipin, cholesterol, and lecithin. The nontreponemal tests measure antilipid antibodies (IgG and IgM, which are formed by the host in response to lipoidal material released from damaged host cells early in infection with *T. pallidum*) and lipid like material from the treponemal cell

surface.<sup>5</sup> During infection with syphilis, an antibody-like substance called reagin can be detected in the patient's serum. In syphilis infection of the central nervous system, reagin can be detected in the cerebrospinal fluid (CSF).

Reactive nontreponemal tests confirm the diagnosis in the presence of early or late lesion syphilis. They offer a clue in latent subclinical syphilis, and are effective tools for detecting cases in epidemiological investigations. Nontreponemal tests are superior to the treponemal test for following the response to therapy.<sup>3</sup>

Nontreponemal antigen tests are not entirely specific for syphilis, nor do they have satisfactory sensitivity in all stages of syphilis. Whenever the results of a nontreponemal antigen test disagree with the clinical impression, a treponemal antigen test such as the FTA-ABS<sup>2,3</sup> should be performed. Nontreponemal tests such as the VDRL are used to screen patient serum, while treponemal tests such as the FTA-ABS are used for confirmation. The likelihood of obtaining a reactive VDRL test result in various stages of untreated syphilis has been reported as follows<sup>3</sup>:

STAGES OF UNTREATED SYPHILIS	% REACTIVE VDRL TEST
Primary	78
Secondary	100
Latent	96
Late	71

### Principles of the Procedure

In the VDRL Test procedure, the patient's serum is heat-inactivated and then mixed with a buffered saline suspension of VDRL Antigen containing cardiolipin, lecithin and cholesterol. The combination of reagin and VDRL Antigen forms microscopic clumping called flocculation. With certain modification, the serum test procedure can be used for testing CSF.

### Reagents

**VDRL Antigen** is 0.03% cardiolipin and 0.9% cholesterol dissolved in absolute alcohol with sufficient lecithin (approximately 0.18-0.2%) to produce standard reactivity.

It is prepared according to the modifications of Harris, Rosenberg and Riedel.<sup>6</sup> Cardiolipin and lecithin are prepared according to the directions of Pangborn.<sup>7,8,9</sup>

**VDRL Buffered Saline** is a 1% NaCl solution at pH 6.0 ± 0.1. It is packaged with VDRL Antigen and used to prepare the VDRL Antigen suspension.

**Nontreponemal Antigen Reactive Serum** is a lyophilized human serum standardized to provide a reactive reading when tested according to theUSR or VDRL test procedure.

**VDRL Weakly Reactive Serum** is a lyophilized human serum standardized to provide a weakly reactive reading when tested according to the VDRL test procedure.

**Nontreponemal Antigen Nonreactive Serum** is a lyophilized human serum standardized to provide a nonreactive reading when tested according to theUSR or VDRL test procedure.

### User Quality Control

#### Identity Specifications

##### VDRL Antigen

Appearance: Clear, colorless solution.

##### VDRL Buffered Saline

Appearance: Clear, colorless solution.

##### Nontreponemal Antigen Reactive Serum

Lyophilized Appearance: White to cream colored, button to powdered cake.

Rehydrated Appearance: Light gold to light amber, clear to slightly opalescent.

##### VDRL Weakly Reactive Serum

Lyophilized Appearance: White to cream colored, button to powdered cake.

Rehydrated Appearance: Light gold to light amber, clear to slightly opalescent.

##### Nontreponemal Antigen Nonreactive Serum

Lyophilized Appearance: White to cream colored, button to powdered cake.

Rehydrated Appearance: Light gold to light amber, clear to slightly opalescent.

#### Performance Response

Rehydrate the sera contained in the VDRL Test Control Serum Set per label directions. Perform the VDRL Slide Test according to the Test Procedure. Each serum in the VDRL Test Control Serum Set should yield appropriate reactions when tested with the VDRL Antigen.

Use an antigen suspension only if it produces the expected reactivity with the control sera.

## Precautions

1. For In Vitro Diagnostic Use.
2. **WARNING! POTENTIALLY BIOHAZARDOUS REAGENTS.** Each donor unit used in preparation of VDRL Antigen and VDRL Test Control Serum Set was tested by an FDA approved method for the presence of the antibody to human immunodeficiency virus (HIV) and for hepatitis B surface antigen and found negative (were not repeatedly reactive).  
Because no test method can offer complete assurance that HIV, hepatitis B virus or other infectious agents are absent, these reagents should be handled at the Biosafety Level 2 as recommended for any potentially infectious human serum or blood specimen.<sup>10</sup>
3. Observe universal blood and body fluid precautions in handling and disposing of specimens.<sup>11,12</sup>
4. **VDRL Antigen**  
**HIGHLY FLAMMABLE. IRRITATING TO EYES, RESPIRATORY SYSTEM AND SKIN.**<sup>us</sup> POSSIBLE RISK OF IRREVERSIBLE EFFECTS.<sup>us</sup> POSSIBLE RISK OF HARM TO THE UNBORN CHILD.<sup>us</sup> Avoid contact with skin and eyes. Do not breathe mist. Wear suitable protective clothing. Keep container tightly closed. Keep away from sources of ignition. No smoking. Target Organs: Blood, Intestines, Liver, Muscles, Nerves.  
**FIRST AID:** In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. If inhaled, remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Seek medical advice. If swallowed seek medical advice immediately and show this container or label.
5. **VDRL Test Control Serum Set**  
The Packaging of This Product Contains Dry Natural Rubber.
6. Follow proper established laboratory procedure in handling and disposing of infectious materials.

## Storage

Store VDRL Antigen at 15-30°C in the dark.

Store VDRL Buffered Saline at 15-30°C. After the bottle is opened, store at 2-8°C.

Store the lyophilized control sera in the VDRL Test Control Serum Set at 2-8°C. Store the rehydrated control sera at 2-8°C or divide into aliquots sufficient for one day of testing and store at -20°C. Do not thaw and refreeze.

## Expiration Date

The expiration date applies to the product in its intact container when stored as directed. Do not use a product if it fails to meet specifications for identity and performance.

## Procedure

### Materials Provided

VDRL Antigen with VDRL Buffered Saline  
VDRL Test Control Serum Set

### Materials Required But Not Provided

0.9% saline  
Nondisposable syringe, 1 cc

Nondisposable calibrated needles without bevel:

Serum test: 18 gauge

CSF test: 21 or 22 gauge

Bottles, 30 ml, round, narrow-mouthed, 35 mm in diameter with glass stoppers and a flat inner-bottom surface

Micropipettor, 50 µl

Pipettes, serological, graduated to tip:

1.0 ml, graduated in 1/100 ml

5.0 ml, graduated in 1/10 ml

10.0 ml, graduated in 1/10 ml

Slides:

Serum test: 2 x 3 inches with paraffin or ceramic rings approximately 14 mm in diameter and high enough to prevent spillage during rotation

CSF test: Kline concavity slides, 3 x 2 3 inches x 3 mm thick, 12 concavities measuring 16 mm in diameter and 1.75 mm in depth

Slide holder for 2 x 4 inch slides

Mechanical rotator adjustable to 180 ± 2 rpm, circumscribing a circle 19 mm in diameter on a horizontal plane

Waterbath, 56°C

Light microscope with 10X ocular and 10X objective

Sterile distilled or deionized water

Absolute alcohol

Acetone

Timer

### Reagent Preparation

VDRL Antigen and VDRL Buffered Saline are ready to use in preparing VDRL Antigen suspension.

Equilibrate all materials to room temperature before performing the tests. Ensure that all glassware and pipettes are clean and free of detergent residues.

**VDRL Test Control Serum Set:** To rehydrate the control sera, add 3 ml sterile distilled or deionized water each and rotate gently to completely dissolve the contents.

### Specimen Collection and Preparation

Collect a blood specimen by aseptic venipuncture. After the specimen has clotted, centrifuge to obtain serum. Store serum specimens at room temperature for no longer than 4 hours; for prolonged storage, keep at 2-8°C for up to 5 days or maintain below -20°C. Serum specimens must be clear, free of hemolysis and show no visible evidence of bacterial contamination such as turbidity, hemolysis or particulate matter. Consult appropriate references for more information on collection of specimens.<sup>1,3,13</sup>

Before testing, heat the test sera at 56°C for 30 minutes. Specimens that are not tested within four hours must be reheated for 10 minutes at 56°C.

### Preparation of Specific Glassware

Syringes with needles and emulsion bottles:

1. Pre-rinse with tap water.
2. Soak and hand wash thoroughly in a glassware detergent solution.
3. Rinse with tap water 6-8 times.
4. Rinse with unused distilled or deionized water.
5. Rinse with absolute alcohol.
6. Rinse with acetone.

7. Air dry until the acetone odor is completely eliminated.
8. Remove needles from syringes for storage.

Ceramic-ringed slides:

1. Pre-rinse with tap water.
2. Wash with a glassware detergent solution. Avoid prolonged soaking in detergent solution because the ceramic rings will become brittle and flake off.
3. Rinse with tap water 3-4 times.
4. Rinse with unused distilled or deionized water.
5. Wipe dry with a clean lint-free cloth. If, after cleaning, the slides do not allow serum to spread evenly within the inner surface of the circle, proceed as follows.
6. Scrub the slides with a nonscratching cleanser.
7. Rinse, dry and polish with a clean lint-free cloth.

### Prepare the Antigen Suspension

Check the pH of VDRL Buffered Saline before preparing the VDRL Antigen emulsion. VDRL Buffered Saline outside the range of pH  $6.0 \pm 0.1$  should be discarded.

Allow VDRL Antigen and VDR Buffered Saline to reach 23-29°C before preparing the VDRL Antigen emulsion.

Use only emulsion bottles with flat inner-bottom surfaces that allow the initial VDR Buffered Saline to evenly cover the inner-bottom surface of the bottle. If the VDRL Buffered Saline beads or does not spread evenly to cover the bottom of the bottle, rewash the bottle.

For reproducible results, the VDRL Antigen emulsion must be checked daily for proper reactivity by testing with VDRL Test Control Serum Set. Only those VDRL emulsions producing the established reactivity pattern of the control serum should be used.

1. Prepare a fresh VDRL Antigen suspension each testing day. The temperature of the VDRL Buffered Saline, VDRL Antigen and equipment should be between 23-29°C at the time the antigen suspension is prepared.
2. Pipette 0.4 ml of VDRL Buffered Saline to the bottom of a round, 30 ml glass-stoppered bottle with a flat inner-bottom surface. Gently tilt the bottle so that the VDR Buffered Saline will cover the entire inner-bottom surface of the bottle.
3. Add 0.5 ml of VDRL Antigen (from the lower half of a 1.0 ml pipette graduated to the tip) directly into the saline while continuously but gently rotating the bottle on a flat surface. Add antigen drop by drop at a rate allowing approximately 6 seconds for each 0.5 ml of antigen. Keep the pipette tip in the upper third of the bottle. Do not splash saline onto the pipette. The proper speed of rotation is obtained when the center of the bottle circumscribes a 2-inch diameter circle approximately three times per second.
4. Expel the last drop of antigen from the pipette without touching the pipette to the saline and continue rotating the bottle for 10 seconds.
5. Add 4.1 ml of buffered saline from a 5 ml pipette. Do not drop the saline directly onto the antigen; allow it to flow down the side of the bottle.
6. Cap the bottle and shake it from bottom to top and back approximately 30 times in 10 seconds. Let the VDRL Antigen emulsion stand without further disturbance for 10 minutes. The antigen suspension is ready for use and may be used during 1 day (8 hours).
7. Mix the VDRL Antigen suspension by gently swirling it each time it is used. Do not mix the suspension by forcing it back and forth

through the syringe and needle, since this may cause breakdown of particles and loss of reactivity.

### Testing the Accuracy of the Antigen Suspension Needle

1. The accuracy of the test depends on the amount of antigen suspension used. Check the calibration of the needle periodically to ensure delivery of the correct volume of VDRL Antigen suspension.
2. For the qualitative and quantitative tests on serum, dispense the antigen suspension from a syringe fitted with an 18-gauge needle without bevel that will deliver  $60 \pm 2$  drops of antigen suspension per ml when held vertically.
3. Place the needle on a 1 ml syringe. Fill the syringe with VDRL Antigen suspension. Holding the syringe in a vertical position, count the number of drops delivered in 0.5 ml. The needle is correctly calibrated if  $30 \pm 1$  drops are delivered in 0.5 ml.
4. Adjust or replace the needle if it does not meet this specification. Repeat calibration of the new needle.

### Testing and Storing the VDRL Antigen Suspension

1. Prepare a fresh antigen suspension each testing day. Once prepared, it should be used within 8 hours.
2. Store the antigen suspension at 23-29°C.
3. Test the reactivity of the antigen suspension with the Reactive, Weakly Reactive and Nonreactive control sera. Test the serum dilutions within 1 hour after inactivation.
4. Use the antigen suspension only if it produces the expected reactivity with the control sera (Reactive, Weakly Reactive and Nonreactive).
5. After each day of use, clean the dispensing needle, bottle and syringe by rinsing with water, alcohol and acetone, as described above. Remove the needle from the syringe after cleaning.

### VDRL Qualitative Slide Test on Serum

1. Slide flocculation tests for syphilis are affected by the temperature of the room. For reliable and reproducible test results, the VDRL Antigen suspension, controls and test specimens must be at 23-29°C when tests are performed.
2. Pipette 50 µl of serum into one ring of a paraffin or ceramic-ringed slide using a safety pipetting device. Do not use a glass slide with concavities, wells or glass rings. Spread the serum with a circular motion of the pipette tip so that the serum covers the entire inner surface of the paraffin or ceramic ring. Use only clean slides that allow the serum to evenly cover the entire surface within the ceramic or paraffin ring.
3. Gently Resuspend the VDRL Antigen suspension.
4. Holding the VDRL Antigen suspension dispensing needle and syringe in a vertical position, dispense several drops to clear the needle of air. Then add exactly 1 free-falling drop (17 µl) of antigen suspension to each circle containing serum. Do not allow the needle to touch the serum.
5. Place the slide on the mechanical rotator. Rotate the slide for 4 minutes at  $180 \pm 2$  rpm. When performing the test in a dry climate, cover the slides with a moist, humidifying cover during rotation to prevent excessive evaporation.
6. Immediately after rotating the slide, remove it from the rotator and read the test results microscopically using a 10X ocular and a 10X objective.

### Results - Qualitative Slide Test

- Read and record results as follows:  
**Reactive (R)** - Medium to large clumps  
**Weakly reactive (WR)** - Small clumps  
**Nonreactive (N)** - No clumping or very slight roughness
- Verify that the control sera results are as expected. If reactions are not as expected, the test is invalid and results cannot be reported.
- Perform a quantitative test on all serum specimens that produce Reactive, Weakly Reactive or "rough" Nonreactive results, since prozone reactions are encountered occasionally.

### VDRL Quantitative Test on Serum

- To quantitate serum samples to an endpoint titer, prepare serum dilutions on the slide at 1:1, 1:2, 1:4 and 1:8, as follows.
- Dispense 50 µl of 0.9% saline in circles 2-4. Do not spread the saline.
- Dispense 50 µl of serum in circles 1 and 2.
- Mix the saline and the serum in circle 2 by drawing the mixture up and down in the pipette 8 times. Avoid forming bubbles.
- Transfer 50 µl from circle 2 (1:2) to circle 3 (1:4) and mix.
- Transfer 50 µl from circle 3 (1:4) to circle 4 (1:8), mix, and then discard 50 µl from circle 4.
- Gently resuspend the antigen suspension.
- Holding the antigen suspension dispensing needle and syringe in a vertical position, dispense several drops to clear the needle of air. Then add exactly 1 free-falling drop (17 µl) of antigen suspension to each circle.
- Place the slide on the mechanical rotator. Rotate the slide for 4 minutes at 180 ±2 rpm. When performing the test in a dry climate, place the slides under a moist, humidifying cover during rotation to prevent excessive evaporation.
- Immediately after rotation, read the test results microscopically using a 10X ocular and a 10X objective.
- If the highest dilution tested (1:8) is reactive, prepare a 1:8 dilution of the test specimen by adding 0.1 ml of serum to 0.7 ml of 0.9% saline. Mix thoroughly. Retest as in steps 1-10, above.

### Results - Quantitative Test

Report the titer as the highest dilution that produces a Reactive (not Weakly Reactive) result.

**Table 1.** Sample quantitative VDRL Test results.

Undiluted (1:1)	SERUM DILUTIONS					REPORT
	1:2	1:4	1:8	1:16	1:32	
R	W	N	N	N	N	Reactive, undiluted
R	R	W	N	N	N	Reactive, 1:2 dilution
R	R	R	W	N	N	Reactive, 1:4 dilution
W	W	R	R	W	N	Reactive, 1:8 dilution
N (rough)	W	R	R	R	N	Reactive, 1:16 dilution
W	N	N	N	N	N	Weakly reactive, undiluted

If reactive results are obtained through dilution 1:32, prepare further twofold serial dilutions in 0.9% saline (1:64, 1:128 and 1:256) and retest using the quantitative test procedure.

### Interpretation

- The results of the serum VDRL Test must be confirmed by a treponemal test.
- The diagnosis of syphilis depends on the results of the VDRL test, the treponemal confirmatory test, clinical signs and symptoms, and risk factors.
- A reactive VDRL Test may indicate past or present infection with a pathogenic treponeme. However, it may be a false-positive reaction. A false positive is determined if the confirmatory treponemal test is negative.
- A nonreactive VDRL Test with clinical evidence of syphilis may indicate early, primary syphilis, a prozone reaction in secondary syphilis, or late syphilis.
- A nonreactive VDRL Test with no clinical evidence of syphilis indicates no current infection or an effectively treated infection.
- A quantitative VDRL Test detects changes in reagin titer. Therefore a serum specimen showing a fourfold increase in titer on a repeat specimen may indicate an infection, a reinfection or a treatment failure. Likewise, a fourfold decrease during treatment indicates adequate syphilis therapy.

### VDRL Test on Spinal Fluid

Consult an appropriate reference for the procedure to use when testing spinal fluids by the VDRL Test.<sup>1</sup>

### Limitations of the Procedure

- A prozone reaction may occur in which reactivity with undiluted serum is inhibited. The prozone phenomenon often gives Weakly Reactive or "rough" Nonreactive results in the qualitative test. Specimens with such results must be quantitatively tested.
- Biological false-positive reactions can occur with nontreponemal tests in persons who abuse drugs, have diseases such as lupus erythematosus, mononucleosis, malaria, leprosy or viral pneumonia, or who have recently been immunized.<sup>1</sup>
- During manufacturing, VDRL Antigen with VDRL Buffered Saline is tested only with serum. To modify the serum test products and procedures for testing CSF, consult the appropriate reference.<sup>1</sup> The user is responsible for modifying the products and procedures and for the required quality control standards according to this manual.
- Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.
- VDRL Buffered Saline showing turbidity or mold growth should be discarded.
- If the temperature of the testing area, specimens or reagents is less than 23°C, test reactivity is decreased. If the temperature is greater than 29°C, test reactivity is increased.<sup>1</sup>
- Test results are unpredictable when testing hemolyzed, contaminated or extremely turbid serum specimens.
- For correct test results, adhere strictly to the correct speed and length of time for rotating the specimens and antigen.

## References

1. **Larsen, S. A., E. F. Hunter, and S. J. Kraus.** 1990. A manual of tests for syphilis. American Public Health Association.
2. **Creighton, E. T.** 1990. Dark field microscopy for the detection and identification of *Treponema pallidum*, p. 49-61. In S. A. Larsen, E. F. Hunter, and S. J. Kraus (ed.), Manual of tests for syphilis, 8th ed. American Public Health Association, Washington, D.C.
3. **Janda, W. M. (ed.).** 1992. Immunology, p. 9.7.1-9.7.20. In H. D. Isenberg (ed.), Clinical microbiology procedures handbook, vol. 2. American Society for Microbiology, Washington, D.C.
4. **Norris, S. J., and S. A. Larsen.** 1995. *Treponema* and other host-associated spirochetes, p. 636-651. In P. R. Murray, E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
5. **Matthews, H. M., T. K. Yang, and H. M. Jenkin.** 1979. Unique lipid composition of *Treponema pallidum* (Nichols virulent strain). Infect. Immun. **24**:713-719.
6. **Harris, A., A. A. Rosenberg, and L. M. Riedel.** 1946. A microfloculation test for syphilis using cardioliipin antigen. J. Ven. Dis. Infor. **27**:169-174.
7. **Pangborn, M. C.** 1941. A new serologically active phospholipid from beef heart. Proc. Soc. Exp. Biol. and Med. **48**:484-486.
8. **Pangborn, M. C.** 1944. Acid cardioliipin and an improved method for the preparation of cardioliipin from beef heart. J. Biol. Chem. **153**:343-348.
9. **Pangborn, M. C.** 1945. A simplified preparation of cardioliipin, with a note on purification of lecithin for serologic use. J. Biol. Chem. **161**:71-82.
10. **U. S. Department of Health and Human Services.** 1988. Biosafety in microbiological and biomedical laboratories, 2nd ed. U. S. Department of Health and Human Services publication no. 88-8395. U. S. Government Printing Office, Washington, D.C.
11. **Centers for Disease Control.** 1988. Update: universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus, and other bloodborne pathogens in health-care settings. Morbidity and Mortality Weekly Reports **37**:377-382, 387-388.
12. **Occupational Safety and Health Administration, U.S. Department of Labor.** 1991. 29 CFR, part 1910. Occupational exposure to bloodborne pathogens, final rule. Federal Register **56**:64175-64182.
13. **Miller, J. M., and H. T. Holmes.** 1995. Specimen collection, transport and storage. In P. R. Murray, E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.

## Packaging

VDRL Antigen w/Buffered Saline	5 ml	0388-56
	10 x 0.5 ml	0388-49
VDRL Test Control Serum Set	1 set	3520-32
Contains:		
Nontreponemal Antigen Reactive Serum	3 ml	
VDRL Weakly Reactive Serum	3 ml	
Nontreponemal Antigen Nonreactive Serum	3 ml	
Aliquant Vials	3 vials	

## Bacto® Vibrio Cholerae Antisera

### Vibrio Cholerae Antiserum Inaba · Vibrio Cholerae Antiserum Ogawa · Vibrio Cholerae Antiserum Poly

#### User Quality Control

##### Identity Specifications

**Vibrio Cholerae Antiserum Inaba**  
**Vibrio Cholerae Antiserum Ogawa**  
**Vibrio Cholerae Antiserum Poly**

Lyophilized appearance: Light gold to amber button to powdered cake

Rehydrated appearance: Light gold to amber, clear liquid

##### Performance Response

Rehydrate Vibrio Cholerae Antisera per label directions. Perform the slide agglutination test using appropriate known cultures of *Vibrio cholerae* as positive and negative controls.

#### Intended Use

Bacto Vibrio Cholerae Antisera are used for serotyping *Vibrio cholerae* in slide agglutination tests.

#### Summary and Explanation

*V. cholerae* are facultative anaerobes that in microscopic morphology are gram-negative, curved or straight bacilli. The microorganisms either require sodium chloride or grow best in its presence. Various media designed to select and cultivate the growth of this microorganism are used in the clinical laboratory, the food industry and environmental testing.

*V. cholerae* is the causative agent of a secretory diarrhea known as cholera. Two biotypes, El Tor and Classical, are associated with human disease.<sup>1</sup> These two biotypes cannot be differentiated serologically. The spread of the disease is primarily through contaminated

water and by the fecal-oral route. Infections with *V. cholerae* may be asymptomatic, mild or severe. If not treated, patients with severe cholera may die within five hours from massive fluid and electrolyte loss.<sup>1</sup> Because *Vibrio* species are natural inhabitants of aquatic environments, food products such as uncooked or incorrectly handled seafood can spread infection.

Gastrointestinal symptoms of cholera include “rice water stools” caused by a potent enterotoxin. The primary human specimen is feces. Seafood products are frequently tested as vehicles of human infection.

In 1935, Gardner and Ventkatraman<sup>2</sup> published a classification scheme for *Vibrio cholerae*. *V. cholerae* isolated from cholera patients was classified as O1. All other strains were designated non-O1. Non-O1 *V. cholerae* causes both gastroenteritis and systemic infections. Some strains produce cholera enterotoxin. Non-O1 *V. cholerae* has been isolated from blood, wounds, ears, sputum, cerebro-spinal fluid and urine.<sup>3,4</sup>

*V. cholerae* of the O1 serogroup are divided into the serotypes Ogawa, Inaba and Hikojima.<sup>5,6</sup> The antigenic factors are:

SEROTYPE	O ANTIGEN FACTORS
Ogawa	AB
Inaba	AC
Hikojima	ABC

Cholera can be diagnosed retrospectively. A fourfold rise in titer between acute-phase serum and that collected 10-14 days later is considered diagnostic.<sup>7</sup>

## Principles of the Procedure

Identification of *V. cholerae* includes the isolation of the microorganism as well as biochemical identification and serological confirmation.

Serological confirmation requires that the microorganism (antigen) react with its corresponding antibody. This *in vitro* reaction produces macroscopic clumping called agglutination. The desired homologous reaction is rapid, does not dissociate (high avidity) and binds (high affinity).

Because a microorganism may agglutinate with an antibody produced in response to another species, heterologous reactions are possible. These are characterized as weak in strength or slow in formation. Such unexpected and unpredictable reactions may lead to some confusion in serological identification. A positive homologous agglutination reaction should support the morphological and biochemical identification of the microorganism.

Homologous reactions occur rapidly and are strong. Heterologous reactions form slowly and are weak.

## Reagents

**Vibrio Cholerae Antisera** are lyophilized, polyclonal rabbit *Vibrio cholerae* O1 antisera containing approximately 0.04% Thimerosal as a preservative. *Vibrio Cholerae* Antisera Inaba and Ogawa are monospecific absorbed antisera.

*Vibrio Cholerae* Antisera possess the following antibodies:

ANTISERUM	O ANTIBODIES
<i>Vibrio Cholerae</i> Antiserum Inaba	C
<i>Vibrio Cholerae</i> Antiserum Ogawa	B
<i>Vibrio Cholerae</i> Antiserum Poly	ABC

Each vial of *Vibrio Cholerae* Antiserum contains sufficient reagent for 20 slide tests.

## Precautions

1. For In Vitro Diagnostic Use.
2. **Vibrio Cholerae Antiserum Inaba**  
**Vibrio Cholerae Antiserum Ogawa**  
**Vibrio Cholerae Antiserum Poly**  
The Packaging of This Product Contains Dry Natural Rubber.
3. Follow proper established laboratory procedure in handling and disposing of infectious materials.

## Storage

Store lyophilized and rehydrated *Vibrio Cholerae* Antisera at 2-8°C.

## Expiration Date

The expiration date applies to product in its intact container when stored as directed. Do not use a product if it fails to meet specifications for identity and performance.

## Procedure

### Materials Provided

*Vibrio Cholerae* Antiserum Inaba  
*Vibrio Cholerae* Antiserum Ogawa  
*Vibrio Cholerae* Antiserum Poly

### Materials Required But Not Provided

Agglutination slides  
Applicator sticks  
Sterile 0.85% NaCl solution

### Reagent Preparation

Equilibrate all materials to room temperature before performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergent.

**Vibrio Cholerae Antisera:** To rehydrate, add 3 ml sterile 0.85% NaCl solution and rotate gently to completely dissolve the contents. The rehydrated antiserum is considered a 1:2 working dilution.

### Specimen Collection and Preparation

*V. cholerae* can be recovered from clinical specimens on selective media such as thiosulfate-citrate-bile salts-sucrose (TCBS) agar. For specific recommendations, consult appropriate references.<sup>1,8,9</sup>

*V. cholerae* can be recovered from various types of foods, particularly seafood. Samples are processed to prevent overgrowth of competing microorganisms and selective media are used to enhance the growth of the microorganism. Some isolation media have been specifically developed for the food industry. Consult appropriate references for recommended procedures for isolating *V. cholerae* from food.<sup>10,11</sup>

The isolate for serological testing should be subcultured from selective media to a nonselective agar such as Nutrient Agar. Consult standard protocols in appropriate references for the type of specimen and the appropriate medium.<sup>1,8-11</sup>

Having followed an established protocol, determine that a pure culture of the microorganism has been obtained and that biochemical test reactions are consistent with the identification of the organism as *V. cholerae*. After these criteria are met, serological identification can be performed.

**Testing the isolate for autoagglutination:**

1. From the test culture, transfer a loopful of growth to a drop (35 µl) of sterile 0.85% NaCl solution on a clean slide and emulsify the organism.
2. Rotate the slide for one minute, then observe for agglutination.
3. If autoagglutination occurs, the culture is rough and cannot be tested. Subculture to a nonselective medium, incubate and test the organism again as described in steps 1 and 2.  
If no agglutination occurs, proceed with testing the organism as described below.

**Test Procedure**

Use Vibrio Cholerae Antiserum Poly to screen possible *V. cholerae* isolates. Continue testing with Vibrio Cholerae Antisera Inaba and Ogawa. Include known positive and negative control cultures.

1. **Vibrio Cholerae Antiserum:** Dispense a drop of the antiserum to be tested on an agglutination slide.
2. **Test organism:** Transfer a loopful of growth to the drop of antiserum and mix thoroughly.
3. Rotate the slide for one minute and read for agglutination.

**Results**

1. Read and record results as follows:
  - 4+ 100% agglutination; background is clear to slightly hazy.
  - 3+ 75% agglutination; background is slightly cloudy.
  - 2+ 50% agglutination; background is moderately cloudy.
  - 1+ 25% agglutination; background is cloudy.
  - No agglutination.
1. **Positive control:** Should produce 3+ or greater agglutination.
2. **Negative control:** Should show no agglutination.
3. Test isolates: 3+ or greater agglutination within one minute is a positive result.

**Interpretation**

Agglutination of the monospecific antiserum, when used, provides preliminary presumptive identification of the serotype.

0.85% NaCl SOLUTION CONTROL	POLY	OGAWA	INABA	INTERPRETATION
–	–			Not <i>V. cholerae</i>
–	+			Presumptively <i>V. cholerae</i>
–	+	+	–	<i>V. cholerae</i> serotype Ogawa
–	+	–	+	<i>V. cholerae</i> serotype Inaba
–	+	+	+	<i>V. cholerae</i> serotype Hikojima
+	Any	Any	Any	Autoagglutination. Unsuitable test culture.

+ agglutination  
– no agglutination

1. Positive agglutination using Vibrio Cholerae Antiserum Poly with typical biochemical test results gives presumptive identification of *V. cholerae* O1.
2. Cultures with positive agglutination in Vibrio Cholerae Antiserum Poly may be serotyped using the Vibrio Cholerae Antiserum Ogawa and Vibrio Cholerae Antiserum Inaba. Positive agglutination in both antisera is rare and, when it occurs, is usually interpreted as

identifying *V. cholerae* serotype Hikojima.<sup>6</sup> *V. cholerae* serotype Hikojima is a rare serotype and should be sent to a reference laboratory for further study.

3. Positive agglutination will be immediate and strong. The strongest and most rapid reaction should be used to identify the serotype. *V. cholerae* O1 strains frequently cross-react slowly or weakly in monospecific antiserum for the other serotype.
4. Isolates that weakly or slowly agglutinate with Vibrio Cholerae Antiserum Poly and not with Vibrio Cholerae Antiserum Inaba or Vibrio Cholerae Antiserum Ogawa are usually considered *V. cholerae* non-O1. The isolate may be sent to a reference laboratory for further study.

**Limitations of the Procedure**

1. Correct interpretation of serological reactions depends on culture purity as well as morphological characteristics and biochemical reactions that are consistent with identification of the microorganism as *V. cholerae*.
2. Serological methods alone cannot identify the isolate as *V. cholerae*.
3. Excessive heat from external sources (hot bacteriological loop, burner flame, light source, etc.) may prevent a smooth suspension of the microorganism or cause evaporation or precipitation of the test mixture. False-positive reactions may occur.
4. Rough culture isolates do occur and will agglutinate spontaneously, causing agglutination of the negative control reaction (autoagglutination). Smooth colonies must be selected and tested in serological procedures.
5. Vibrio Cholerae Antisera have been tested using undiluted cultures taken from agar media. These antisera have not been tested using antigen suspensions in NaCl solution or other diluents. If the user applies variations to the recommended procedure, each lot of antiserum must be tested with known control cultures to verify expected reactions under the modified procedure.
6. Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.
7. Rehydrated Vibrio Cholerae Antiserum that is cloudy or has a precipitate at any time during its use should be discarded.

**References**

1. **McLaughlin, J. C.** 1995. *Vibrio*, p. 465-476. In P. R. Murray, E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover (ed.), Manual of clinical microbiology, 6th ed. American Society for Microbiology, Washington, D.C.
2. **Gardner, A. D., and K. V. Venkatraman.** 1935. The antigens of the cholera group of vibrios. *J. Hyg.* **35**:262-282.
3. **Florman, A. L., A. H. Cushing, T. Byers, and S. Popejoy.** 1990. *Vibrio cholerae* bacteremia in a 22-month-old New Mexican child. *Pediatr. Infect. Dis. J.* **9**:63-65. (Letter).
4. **Morris, J. G., Jr.** 1990. Non-O group 1 *Vibrio cholerae*: a look at the epidemiology of an occasional pathogen. *Epidemiol. Rev.* **12**:179-191.
5. **Nobecki, K.** 1923. Contributions to the knowledge of *V. cholerae*. 3. Immunological studies upon the types of *V. cholerae*. *Sci. Repr. Govt. Inst. Infect. Dis.* **2**:43-87.

6. **Kelly, M. T., F. W. Hickman-Brenner, and J. J. Farmer III.** 1991. *Vibrio*, p. 384-395. In A. Balows, W. J. Hausler, K. L. Herrmann, H. D. Isenberg, and H. J. Shadomy (ed.), Manual of clinical microbiology, 5th ed. American Society for Microbiology, Washington, D.C.
7. **Smith, H. L.** 1981. *Vibrio* Infections, p. 715-722. In A. Balows, and W. J. Hausler (ed.), Bacterial, mycotic and parasitic infections, 6th ed. American Public Health Association, Washington, D.C.
8. **Pezzlo, M.** 1992. Aerobic bacteriology, p. 1.0.1.-1.20.47. In H. D. Isenberg (ed.). Clinical microbiology procedures handbook, vol. 1. American Society for Microbiology, Washington, D.C.
9. **Baron, E. J., L. R. Peterson, and S. M. Finegold.** 1994. Bailey & Scott's diagnostic microbiology, 9th ed. Mosby-Year Book, Inc., St. Louis, MO.
10. **Kaysner, C. A., M. L. Tamplin, and R. M. Twedt.** 1992. *Vibrio*, p. 451-473. In C. Vanderzant and D. F. Splittstoesser, (ed.), Compendium of methods for the microbiological examination of foods, 3rd ed. American Public Health Association, Washington, D.C.
11. **Elliot, E. L., C. A. Kaysner, L. Jackson, and M. L. Tamplin.** 1995. *Vibrio cholerae*, *V. parahaemolyticus*, *V. vulnificus*, and other *Vibrio* sp., p. 9.01-9.27. In Food and Drug Administration, Bacteriological analytical manual, 8th ed. AOAC International, Gaithersburg, MD.

### Packaging

Vibrio Cholera Antiserum Inaba	3 ml	2430-47
Vibrio Cholera Antiserum Ogawa	3 ml	2431-47
Vibrio Cholera Antiserum Poly	3 ml	2432-47

## Bacto® Widal Antigen Set

**CONTAINS: Salmonella O Antigen · Salmonella H Antigens  
Salmonella Vi Antigen · Febrile Positive Control Polyvalent  
Febrile Negative Control**

**ALSO AVAILABLE: Salmonella O Antigens · Salmonella H Antigen**

### Intended Use

Bacto Widal Antigen Set, which contains four *Salmonella* antigens, is used in detecting *Salmonella* antibodies by slide and tube agglutination tests.

Bacto Salmonella H, O and Vi Antigens, available individually, are used in detecting *Salmonella* antibodies by slide and tube agglutination tests.

### Summary and Explanation

*Salmonellae* represent many species of pathogenic microorganisms that, upon invasion, produce a fever in their host. Consequently, they are often called "Febrile Antigens". *Salmonella* species cause a variety of human diseases called salmonellosis. The range of disease is from mild self-limiting gastroenteritis to a more severe form, possibly with bacteremia or typhoid fever, which can be severe and life-threatening. Severe disease and bacteremia are associated primarily with *S. Choleraesuis*, *S. Paratyphoid A* and *S. Typhi*, while most of the other 2,300 or more strains are associated with gastroenteritis. The severity of the diarrheal disease depends on the virulence of the strain and the condition of the human host.

*Salmonella* is found in nature and occurs in the intestinal tract of many animals, both wild and domestic. The microorganism can spread to man from contact with the environment or from eating contaminated meat or vegetable food products.

The genus *Salmonella* is in the family *Enterobacteriaceae*. Salmonellae are facultatively anaerobic, gram-negative bacilli that typically are

oxidase negative, non-lactose fermenting, H<sub>2</sub>S positive and produce gas. Serotypes of *Salmonella* are defined based on antigenic structure, both somatic or cell wall (O) antigens and flagellar (H) antigens. The antigenic formula lists the O antigen(s) followed by the H antigen(s).

In 1896, Widal introduced techniques for testing patient serum for antibodies in cases of typhoid fever.<sup>1</sup> The Widal test was used diagnostically in two ways. First, it was considered diagnostic when a single high antibody titer occurred during the first week of illness. Further, it was diagnostic if a greater than fourfold rise in titer existed in serum samples taken 1 to 2 weeks apart.<sup>2,3,4</sup> The Widal test was developed to include *Salmonella* Typhi and other species of *Salmonella* detected by a variety of O and H antigens. *S. Typhi* and *S. Paratyphi A* and *B* are the major pathogens in this group that can produce clinically distinct systemic illness. In areas such as developing countries where typhoid is highly endemic, the diagnostic value of the Widal test has been well documented.<sup>5,6,7</sup> The Widal test for antibodies to the O antigens of *Salmonella* serotypes most likely to cause typhoid fever (usually *S. Typhi* and *S. Paratyphi A* and *B*) can be useful in helping diagnose typhoid fever when other methods have failed.<sup>8</sup>

Diagnosis of the cause of febrile disease cannot be based solely on the analysis of serum samples for antibody response. Many factors may affect measurable antibody levels. For example, the patient's immune response can be affected by age, immune status, general state of health and previous immunizations. Patients with known typhoid fever have developed diagnostic titers of antibodies that are low.<sup>9,10</sup> Also, patients treated with antibiotics early in their disease may not develop a significant titer rise.<sup>2</sup>

The various species of *Salmonella* contain multiple antigens that are cross-reactive. This prevents using increased antibody levels, alone, to identify infecting organisms by species or serotype. Other non-typhoidal febrile illness or unrelated immunological disorders may produce significant elevations of antibody titers.<sup>2</sup>

Certain organisms may share cross-reacting antigens leading to the production of heterologous antibodies. These heterologous antibodies may react with one or more antigens in an antibody test procedure, resulting in low-level antibody titers that may not, singly, suggest disease.

The rapid slide procedure is a screening test designed to detect agglutinins. The macroscopic tube test<sup>11</sup> is a confirmatory procedure designed to quantitate agglutinin compositions. Any positive results obtained in the screening (slide) test of specimens must be confirmed by a tube test.

The rapid slide test is the most widely used procedure employing the Widal antigens because of the simplicity with which the results may be reported. Negative slide test reactions can usually be reported as such if all five serum dilutions have been used. Although the slide test is not quantitative, running the series of dilutions is necessary to detect agglutinins that might be overlooked with a “prozone phenomenon”. This often occurs in serum containing a high titer of typhoid agglutinins where higher concentrations of the serum may yield negative results but a dilution of the serum is positive.

## Principles of The Procedure

Agglutination tests involving the use of *Salmonella* antigens determine the presence of antibodies that react with the test antigen. The serological procedure involves serially diluting the patient serum and then adding

a standard volume of antigen. The end point of the test is the last dilution of the serum that shows a specific amount of agglutination. The end point, reported as a dilution of the serum, is called the patient’s antibody “titer.”

## Reagents

### Antigens

1. **Salmonella Antigens** are ready-to-use suspensions of the *Salmonella* organisms listed below. Salmonella O Antigens contain 20% glycerin.

#### Widal Antigen Set contents:

Salmonella O Antigen Group D - *Salmonella* Typhi O901, factors 9,12 (selected strain)

Salmonella H Antigen a - *Salmonella* Paratyphi A.

Salmonella H Antigen b - *Salmonella* Paratyphi B.

Salmonella H Antigen d - *Salmonella* Typhi H901.

**Salmonella O and H Antigens**, available individually, are prepared from selected strains containing the following group-specific antigens:

Salmonella H Antigen c - flagellar antigen c

Salmonella O Antigen Group A - factors 1, 2, 12

Salmonella O Antigen Group B - factors 1, 4, 5, 12

Salmonella O Antigen Group C - factors 6, 7, (8), 20

Salmonella Vi Antigen

When used as described, each vial contains sufficient reagent for 20 slide tests or 25 tube tests.

2. **Concentration of Antigen:** Salmonella O, H and Vi Antigens are used undiluted for the slide test and diluted 1:20 for the tube test.
3. **Antigen Density:** Salmonella O, H and Vi Antigens are adjusted to a density approximating 20 times a McFarland Barium Sulfate Standard No. 3 ( $1.8 \times 10^{10}$  organisms per ml).

Because antigen density may vary, it is adjusted for optimum performance when standardized with hyperimmune sera obtained from laboratory animals.

Variation in color intensity of the antigen is normal and will not affect the outcome of the test.

4. Salmonella Antigens contain the following preservatives:
  - Salmonella O Antigens:** 0.5% phenol, and approximately 0.002% crystal violet and 0.005% brilliant green.
  - Salmonella H Antigens:** 0.5% formaldehyde, and approximately 0.002% crystal violet and 0.005% brilliant green.
  - Salmonella Vi Antigen:** 0.5% phenol, and approximately 0.002% crystal violet and 0.005% brilliant green.

### Antisera:

1. **Febrile Positive Control Polyvalent** is lyophilized, polyclonal, polyvalent goat antisera containing approximately 0.04% Thimerosal as a preservative. This reagent contains antibodies at a titer of 1:80 or greater for the Salmonella O and H Antigens in the Widal Antigen Set.

Each vial of Febrile Positive Control Polyvalent contains sufficient reagent for 32 slide tests or 50 tube tests using the four antigens contained in the Widal Antigen Set. When using the Salmonella O, H and Vi Antigens separately, there is sufficient reagent for 20 slide or 25 tube tests.

## User Quality Control

### Identity Specifications

#### Salmonella O Antigens

#### Salmonella H Antigens

#### Salmonella Vi Antigen

Appearance: Turquoise-blue-violet suspension.

#### Febrile Positive Control Polyvalent

Lyophilized Appearance: Light gold to amber, button to powdered cake.

Rehydrated Appearance: Light gold to amber, clear liquid.

#### Febrile Negative Control

Lyophilized Appearance: Colorless to light gold, button to powdered cake.

Rehydrated Appearance: Colorless to light gold, clear liquid.

### Performance Response

Rehydrate Febrile Positive Control Polyvalent and Febrile Negative Control per label directions. Perform the slide or tube agglutination test using Salmonella O, H or Vi Antigens and positive and negative controls diluted in the same proportion as a patient serum.

A Salmonella Antigen is considered satisfactory if it does not agglutinate with the negative control and gives a 2+ or greater reaction at a 1:80 dilution with the positive control.

2. **Febrile Negative Control** is a lyophilized, standard protein solution containing approximately 0.02% Thimerosal as a preservative.

Each vial of Febrile Negative Control contains sufficient reagent for 32 slide tests when using the four antigens contained in the Widal Antigen Set. When using the Salmonella O, H and Vi Antigens separately, there is sufficient reagent for 20 slide tests.

## Precautions

1. For In Vitro Diagnostic Use.
2. Observe universal blood and body fluid precautions in the handling and disposing of specimens.<sup>12,13</sup>
3. **Salmonella H Antigen a**  
**Salmonella H Antigen b**  
**Salmonella H Antigen c**  
**Salmonella H Antigen d**  
POSSIBLE RISK OF IRREVERSIBLE EFFECTS. (US) Avoid contact with skin and eyes. Do not breathe mist. Wear suitable protective clothing. Keep container tightly closed. TARGET ORGAN(S): Eyes, Kidneys, Lungs, Skin.  
FIRST AID: In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. If inhaled, remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Seek medical advice. If swallowed seek medical advice immediately and show this container or label.
4. Follow proper established laboratory procedure in handling and disposing of infectious materials.
5. Salmonella O, H and Vi Antigens are not intended for use in the immunization of humans or animals.

## Storage

Store Salmonella O, H and Vi Antigens at 2-8°C.

Store lyophilized and rehydrated Febrile Positive Control Polyvalent at 2-8°C.

Store lyophilized and rehydrated Febrile Negative Control at 2-8°C.

## Expiration Date

The expiration date applies to the product in its intact container when stored as directed. Do not use a product if it fails to meet specifications for identity and performance.

## Procedure

### Materials Provided

#### Widal Antigen Set:

Salmonella O Antigen Group D  
Salmonella H Antigen Group a  
Salmonella H Antigen Group b  
Salmonella H Antigen Group d  
Febrile Positive Control Polyvalent  
Febrile Negative Control

#### Available separately:

Salmonella O Antigens  
Salmonella H Antigens  
Salmonella Vi Antigen  
Febrile Positive Control Polyvalent

Febrile Negative Control  
Droppers, 0.03 ml per drop (supplied)

## Materials Required but not Provided

### Slide test

Agglutination slides, 5 squares, 1" each  
Applicator sticks  
Sterile distilled or deionized water  
Serological pipettes, 0.2 ml

### Tube Test

Culture tubes 12 x 75 mm and rack  
Waterbath, 50 ± 2°C  
Serological pipettes, 1 ml and 5 ml  
Sterile 0.85% NaCl solution

## Reagent Preparation

**Salmonella O, H and Vi Antigens** are ready to use in the slide test. Salmonella O and H antigens must be diluted 1:20 for the tube tests. (See Test Procedure for preparation instructions).

Equilibrate all materials to room temperature before performing the tests. Ensure that all glassware and pipettes are clean and free of residues such as detergent.

**Febrile Positive Control Polyvalent:** To rehydrate, add 5 ml sterile distilled or deionized water and rotate gently to completely dissolve the contents.

**Febrile Negative Control:** To rehydrate, add 5 ml sterile distilled or deionized water and rotate gently to completely dissolve the contents.

## Specimen Collection and Preparation

Collect a blood specimen by aseptic venipuncture. Serum is required for the test. Store serum specimens at room temperature for no longer than 4 hours; for prolonged storage, keep at 2-8°C for up to 5 days or maintain at or below -20°C. Serum specimens must be clear, free of hemolysis and show no visible evidence of bacterial contamination (turbidity, hemolysis or particulate matter). Refer to appropriate references for more information on collection of specimens.<sup>14,15</sup> Serum specimens must not be heated. Heat may inactivate or destroy certain antibodies.

Because changes in titer over a period of time are the best indicators of active infection, and because the accuracy and precision of the tests can be affected not only by test conditions but also by the subjectivity of the person reading the endpoint, the following protocol is recommended.

A preliminary test using either the rapid slide test and/or the macroscopic tube test may be performed on the initial serum specimen and reported to the physician at that time. An aliquot of the serum should be transferred to a sterile test tube, sealed tightly, and kept in the freezer. When the second serum is obtained, it should be run in parallel with the original specimen. In this manner, the original serum will serve as a control. Any difference in titer will be more credible, since the bias associated with the performance of the test and determining the endpoint will be reduced.

## Test Procedure

### Rapid Slide Test

Use the slide test only as a screening test. Confirm positive results with the tube test.

- Test serum:** Using a 0.2 ml serological pipette, dispense 0.08, 0.04, 0.02, 0.01 and 0.005 ml of each test serum into a row of squares on an agglutination slide.
- Positive control:** Using a 0.2 ml serological pipette, dispense 0.08, 0.04, 0.02, 0.01 and 0.005 ml of Febrile Positive Control Polyvalent into a row of squares on the agglutination slide.
- Negative control:** Using a 0.2 ml serological pipette, dispense 0.08, 0.04, 0.02, 0.01 and 0.005 ml of Febrile Negative Control into a row of squares on the agglutination slide.
- Salmonella Antigen:** Shake the vial of antigen well to ensure a smooth, uniform suspension. Place one drop (35 µl) of antigen suspension in each drop of test serum, positive control and negative control.
- Mix each row of test sera and control sera, using a separate applicator stick for each row. Start with the most dilute mixture (0.005 ml) and work to the most concentrated (0.08 ml).
- Rotate the slide for 1 minute and read for agglutination.
- The final dilutions in squares 1-5 correspond with tube dilutions of 1:20, 1:40, 1:80, 1:160, 1:320, respectively.

### Results

- Read and record results as follows:
  - 4+ 100% agglutination; background is clear to slightly hazy.
  - 3+ 75% agglutination; background is slightly cloudy.
  - 2+ 50% agglutination; background is moderately cloudy.
  - 1+ 25% agglutination; background is cloudy.
  - No agglutination.
- Positive control:** Should show 2+ or greater agglutination at the 1:80 dilution.
- Negative control:** Should show no agglutination.
- If results for either the positive or negative control are not as described, the test is invalid and results cannot be reported.
- Test specimens:** The serum titer is that dilution that shows a 2+ or greater agglutination. See Table 1.
- The slide test is a screening test, only, and results must be confirmed with the tube test.

**Table 1.** Sample Rapid Slide Test reactions.

ml SERUM	CORRELATED DILUTION	REACTIONS		
		SPECIMEN 1	SPECIMEN 2	SPECIMEN 3
0.08	1:20	3+	4+	4+
0.04	1:40	2+	3+	3+
0.02	1:80	1+	3+	2+
0.01	1:160	–	2+	+
0.005	1:320	–	1+	–
<b>Serum titer</b>		<b>1:40</b>	<b>1:160</b>	<b>1:80</b>

### Macro Tube Test


Prepare a 1:20 dilution of the Salmonella O and H Antigens by adding 1 part of the antigen to 19 parts of sterile 0.85% NaCl solution.

- Prepare a row of eight culture tubes (12 x 75 ml) for each test serum, including a row for the Febrile Positive Control Polyvalent.


- Sterile 0.85% NaCl solution:** Dispense 0.9 ml in the first tube of each row and 0.5 ml in the remaining tubes.
- Test serum:** Using a 1 ml serological pipette, dispense 0.1 ml of test serum in the first tube in the row and mix thoroughly. Transfer 0.5 ml from tube 1 to tube 2 and mix thoroughly. Similarly, continue transferring 0.5 ml through tube 7, discarding 0.5 ml from tube 7 after mixing. Tube 8 is the antigen control tube and contains only sterile 0.85% NaCl solution.
- Positive control:** Using a 1 ml serological pipette, dispense 0.1 ml of Febrile Positive Control Polyvalent in the first tube in the row and mix thoroughly. Transfer 0.5 ml from tube 1 to tube 2 and mix thoroughly. Similarly, continue transferring 0.5 ml through tube 7, discarding 0.5 ml from tube 7 after mixing. Tube 8 is the antigen control tube and contains only sterile 0.85% NaCl solution.
- Salmonella Antigen:** Add 0.5 ml of the diluted antigen to all 8 tubes in each row and shake the rack to mix the suspensions.
- The final dilutions in tubes 1-7 are 1:20, 1:40, 1:80, 1:160, 1:320, 1:640 and 1:1280, respectively.
- H antigen tests:** Incubate in a waterbath at  $50 \pm 2^\circ\text{C}$  for 1 hour.  
**O antigen tests:** Incubate in a waterbath at  $50 \pm 2^\circ\text{C}$  for  $17 \pm 1$  hours.
- Remove from the waterbath. Avoid excessive shaking before reading the reactions, either when the tubes are in the waterbath or when removing them from the waterbath.
- Read and record the results.

### Results

- Tube agglutination reactions detect antibodies to either somatic (O) antigens or flagellar (H) antigens and these antibodies give two different reactions. An O antigen and the corresponding antibody give a coarse, compact agglutination which may be difficult to disperse. An H antigen and its corresponding antibody give a loose, flocculent agglutination. Do not vigorously shake tubes containing H antigens. Characteristic O and H agglutination is shown in the following diagrams.



**Somatic "O" Agglutination**



**Flagellar "H" Agglutination**

- Read and record results as follows:
  - 4+ 100% agglutination; background is clear to slightly hazy.
  - 3+ 75% agglutination; background is slightly cloudy.
  - 2+ 50% agglutination; background is moderately cloudy.
  - 1+ 25% agglutination; background is cloudy.
  - No agglutination.
- Positive control:** Should show 2+ or greater agglutination at the 1:80 dilution.

4. **Antigen control:** Tube #8 of each row should show no agglutination.
5. If results of the positive control and antigen control are not as described, the test is invalid and results cannot be reported.
6. **Test serum:** The serum titer is that dilution which shows 2+ or greater agglutination. See Table 2.

**Table 2.** Sample Macroscopic Tube Test reactions.

REACTIONS			
SERUM DILUTION	SPECIMEN 1	SPECIMEN 2	SPECIMEN 3
1:20	4+	3+	4+
1:40	4+	2+	4+
1:80	3+	1+	4+
1:160	2+	–	4+
1:320	1+	–	3+
1:640	–	–	2+
1:1280	–	–	1+
<b>Serum Titer</b>	<b>1:160</b>	<b>1:40</b>	<b>1:640</b>

### Interpretation

For a single serum specimen, a titer of 1:80 suggests infection.

A pair of serum specimens (acute and convalescent) showing a two-dilution increase in titer is significant and suggests infection. A one dilution difference is within the limits of laboratory error.

Table 3 presents data that will be helpful in interpreting serological tests with the *Salmonella* antigens. The values tabulated will vary in certain cases.

**Table 3.** Disease states and associated *Salmonella* Antigens.

SALMONELLA ANTIGEN	SUGGESTED PATHOLOGY	TIME TO MAXIMUM TITER	SIGNIFICANT TITER
<i>Salmonella</i> H Antigen d (Typhoid H)	Typhoid Fever	4-5 weeks	1:80
<i>Salmonella</i> O Antigen Group D (Typhoid O)	Typhoid Fever	3-5 weeks	1:80
<i>Salmonella</i> H Antigen a (Paratyphoid A)	Paratyphoid Fever	3-5 weeks	1:80
<i>Salmonella</i> H Antigen b (Paratyphoid B)	Paratyphoid Fever	3-5 weeks	1:80

### Limitations of the Procedure


1. The slide test is intended for screening only and should be confirmed by the tube test. Slide test dilutions are made to detect a prozone reaction and do not represent true quantitation of the antibody. A serum specimen with a prozone reaction shows no agglutination because of excessively high antibody concentrations. To avoid this occurrence, all five serum dilutions of the slide test should be run.
2. Detection of antibodies in serum specimens may complete the clinical picture of a patient having infection. However, isolation of the causative agent from patient specimens may be required. A definitive diagnosis must be made by a physician and must be based on patient history, physical examination and data from all laboratory tests.

3. In some cases of typhoid fever, sera may show a prozone reaction, the inability of an antigen to react in higher serum antibody concentrations. It is advised that all five serum dilutions be run in the rapid slide test, rather than just one dilution, to eliminate the possibility of missing positive reactions due to the prozone phenomenon.
4. Cross-reacting heterologous antibodies are responsible for many low titer reactions. Infections with other organisms, vaccinations and a history of disease may result in a low level of antibody titer. Antimicrobial therapy may suppress antibody production.
5. Previous immunization with typhoid vaccine or previous infection with *Salmonella* species sharing common antigens with *S. Typhi* can cause elevated antibody titers for prolonged periods. Other non-typhoid febrile illnesses may cause elevation of cross-reacting antibodies.
6. While a single serum specimen showing a titer of 1:80 suggests infection, it is not diagnostic.
7. Nonspecific agglutination has been noted with *Salmonella* O Group D antigen in the sera of patients with influenza.
8. Sera from narcotic addicts appear to contain broad nonspecific activity to the Widal antigen.<sup>16</sup>
9. Sera from patients with chronic active liver disease may show high agglutinin titers.<sup>17</sup>
10. To test for a significant rise in antibody titer, at least two specimens are necessary: an acute specimen (obtained at the time of initial symptoms) and a convalescent specimen (obtained 7 to 14 days later). A two-dilution difference in the titer is a significant increase in antibody level and suggests infection.
11. Antibody titers may be clinically useful in detecting chronic carriers of *S. Typhi* and in diagnosing typhoid fever in endemic areas. Detection of antibodies to heat-labile *Salmonella* envelope antigen (Vi) may be useful for the detection of a chronic carrier state for *S. Typhi*. The presence of capsular Vi antigen may mask somatic (O) antigen activity.
12. In *Salmonella* infections, the stage of the patient's disease is important and may affect test results. The optimum time for detecting peak titers may be missed if symptoms do not correlate with increased antibody levels. Antibodies to O antigens appear earlier and disappear earlier than antibodies to H antigens.
12. Prolonged exposure of reagents to temperatures other than those specified is detrimental to the products.
13. Exposure to temperatures below 2°C can cause autoagglutination. Antigens must be smooth uniform suspensions; before use, examine antigen vials for agglutination. Suspensions with agglutination are not usable and should be discarded.
14. Discard rehydrated Febrile Positive Control Polyvalent or Febrile Negative Control that is cloudy or has a precipitate anytime during its period of use.

### References

1. **Widal, F.** 1896. Serodiagnostic de la fièvre typhoïde. *Sem. Med.* **16**:259.
2. **Sack, B. R.** 1986. Serologic tests for the diagnosis of enterobacterial infections, p. 359-362. *In* N. R. Rose, H. Friedman, and J. L. Fahey (eds.), *Manual of clinical laboratory immunology*, 3rd ed. American Society for Microbiology, Washington, D.C.





# Reference Guide

## Agar Selection Guide

APPLICATIONS	BACTERIOLOGICAL USES OF AGARS		
	Bacto® Agar	Agar, Granulated	Agar, Technical
Auxotrophic studies	++	+	
Bacteriology, research	++	+	
Bacteriology, general purpose	++	+	
Bacteriophage studies	++		
Biotechnology production	+	++	+/-
General microbial production	+	++	+/-
Growth of fastidious organisms	++		
Identification of pathogenic organisms	++	+	+/-
Microaerophilic studies	++		
Molecular genetics	++	+	
Prepared plate manufacture	+	++	+
Quality control, production	++	+	+/-
Quality control, environmental	+	++	+/-
Transformation of bacteria	++		
Transformation of yeast	++		

APPLICATIONS	NON-BACTERIOLOGICAL USES OF AGARS		
	Bacto® Agar	Agar, Granulated	Agar, Technical
Immunodiffusion		++	
Electrophoresis		++	
Tissue Culture, mammalian		++	
Tissue Culture, plant	++	+	
Histology, tissue embedding	+	++	
Histology, bone marrow embedding	++	+	
Insect growth substrate			++

### Key

++ Recommended

+ Suitable

+/- Marginal

## Anaerobes - General

ANAEROBES - GENERAL
Anaerobe Broth MIC
Anaerobic Agar
Blood Culture Bottles Columbia Broth w/CO <sub>2</sub>
Blood Culture Bottles Columbia Broth w/SPS and CO <sub>2</sub>
Blood Culture Bottles Fluid Thioglycollate Medium w/SPS and CO <sub>2</sub>
Blood Culture Bottles Thioglycollate w/CO <sub>2</sub>
Blood Culture Bottles Thioglycollate w/SPS and CO <sub>2</sub>
Brain Heart Infusion Agar
Brewer Anaerobic Agar
Brewer Thioglycollate Medium
CHO Medium Base
Clostridium Difficile Antimicrobial Supplement CC
Cooked Meat Medium
ESP Anaerobic Broth
Differential Reinforced Clostridial Agar
Fluid Thioglycollate Medium
Liver Veal Agar
McClung Toabe Agar Base
NIH Thioglycollate Broth
Reinforced Clostridial Medium
Schaedler Agar
Schaedler Broth
SFP Agar Base
Sterility Bottles with Septum Fluid Thioglycollate Medium
Sterility Bottles with Screw Cap Fluid Thioglycollate Medium
SPS Agar
Sulfite Agar
Thioglycollate Medium USP Alternative
Thioglycollate Medium w/o Dextrose
Thioglycollate Medium w/o Dextrose or Indicator
Thioglycollate Medium w/o Indicator
Wilkins-Chalgren Agar

## Antimicrobial Selective Agents For Culture Media

AGENT	CONCENTRATION PER VIAL	PRODUCT NAME
Cefsulodin-novobiocin	4 mg/2.5 mg	Yersinia Antimicrobial Supplement CN
Ceftozidime	40 mg	Palcam Antimicrobial Supplement
Chloramphenicol	0.05 g	Rose Bengal Antimicrobial Supplement C
Chlortetracycline (Aureomycin®)	25 mg	Antimicrobial Vial A
Colistin Sulfate - Nystatin - Vancomycin	7,500 mcg/12,500 units/ 3,000 mcg	Antimicrobial Vial CNV
Colistin Sulfate - Nystatin - Vancomycin Trimethoprim	7,500 mcg/12,500 units/ 3,000 mcg/5,000 mcg (10 ml vial)	Antimicrobial Vial CNVT
Cyclohexamide - Colistin Sulfate - Acridavine - Cefotetan - Fosfomycin	400 mg/20 mg/5 mg/2 mg/ 10 mg	Oxford Antimicrobial Supplement
Cycloserine - Cefoxitin	125 mg/ 5 mg	Clostridium Difficile Antimicrobial Supplement CC
Kanamycin	25,000 mcg	Antimicrobial Vial K
Moxalactam	20 mg	Moxalactam Antimicrobial Supplement
Moxalactam - Colistin Sulfate	20 mg/10 mg	Modified Oxford Antimicrobial Supplement
Novobiocin	20 mg	Novobiocin Antimicrobial Supplement
Oxytetracycline	100 mg	Antimicrobial Vial Oxytetracycline
Polymyxin B	30,000 units	Antimicrobial Vial P
Potassium Tellurite Solution 1%	1%	Chapman Tellurite Solution
Potassium Tellurite Solution 3.5%	3.5%	Potassium Tellurite Solution 3.5%
Sodium 7- ethyl - 2 - methyl - 4 - undecyl sulfate		XLT4 Supplement

# Cosmetic Testing

PRODUCTS	APPLICATIONS					
	Environmental <small>(See also Environmental Sampling and Diagnostic Testing)</small>	Gram-Negative Screening	Pseudomonas Isolation	Staphylococcus	Sterility Testing	Yeast & Mold Isolation
AC Broth					✓	
AC Broth Medium w/o Dextrose <sup>†</sup>					✓	
Cetrimide Agar Base/PSEUDOSEL Agar <sup>†</sup>			✓			
RODAC™ Plates	✓					
D/E Neutralizing Agar <sup>†</sup>	✓					
D/E Neutralizing Broth <sup>†</sup>	✓					
EMB Agar/Eosin Methylene Blue Agar Modified <sup>†</sup>		✓				
Fluid Thioglycollate Medium <sup>†</sup>					✓	
HC Agar Base <sup>†</sup>						✓
HYcheck	✓					
Lethen Agar <sup>†</sup>	✓					
Lethen Broth/Lethen Broth AOAC <sup>†</sup>	✓					
m Staphylococcus Broth				✓		
MacConkey Agar <sup>†</sup>		✓				
Malt Agar <sup>†</sup>						✓
Malt Extract Agar						✓
Malt Extract Broth <sup>†</sup>						✓
Mannitol Salt Agar <sup>†</sup>				✓		
Microbial Content Test Agar/ TRYPTICASE™ Soy w/Lec. poly. <sup>†</sup>	✓					
Modified Lethen Agar					✓	
Modified Lethen Broth					✓	
Mycological Agar/MYCOPHIL™ Agar <sup>†</sup>						✓
Mycobiotic Agar/MYCOPHIL Agar <sup>†</sup>						✓
Neutralizing Buffer <sup>†</sup>	✓					
Phenylethanol Agar/Phenylethyl Alcohol Agar <sup>†</sup>				✓		
Potato Dextrose Agar <sup>†</sup>						✓
Pseudomonas Agar F/Flo Agar <sup>†</sup>			✓			
Pseudomonas Agar P/Tech Agar <sup>†</sup>			✓			
Pseudomonas Isolation Agar/Pseudomonas ISO <sup>†</sup>			✓			
Sabouraud Dextrose Agar <sup>†</sup>						✓
Staph Latex Test Kit/STAPHYLOSLIDE™ Test Kit <sup>†</sup>				✓		
Staphylococcus Medium 110/Staphylococcus Agar #10 <sup>†</sup>				✓		
Sterility Test Bottles, Prepared <sup>†</sup>					✓	
TAT Broth Base <sup>†</sup>					✓	
Tryptic Soy Broth/ TRYPTICASE™ Soy Broth <sup>†</sup>					✓	
VJ Agar/Vogel & Johnson Agar <sup>†</sup>				✓		

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Dairy Product Testing - Products and Applications

PRODUCTS	APPLICATIONS								
	Beta-Lactams in Milk	Clostridium	Coliform Analysis	Environmental	Lipolytic Microorganisms	Listeria Analysis	Standard Plate Count	Streptococcus Detection	Yeast & Mold Analysis
Antibiotic Medium 1 <sup>†</sup>	✓								
Antibiotic Medium 4 <sup>†</sup>	✓								
Azide Blood Agar Base <sup>†</sup>							✓		
Brilliant Green Bile 2%/Brilliant Green Bile Broth <sup>†</sup>			✓						
Bryant and Burkey Medium		✓							
Concentration Disks 1/2" Penase	✓								
Concentration Disks 1/2" Penicillin 0.05 Unit/PM Discs, 0.05, 1/4" Taxo™ <sup>†</sup>	✓								
Concentration Disks 1/2" Sterile Blanks, Antibiotic Detect Disc 1/2" Taxo <sup>†</sup>	✓								
Desoxycholate Lactose Agar <sup>†</sup>			✓						
D/E Neutralizing Agar <sup>†</sup>				✓					
D/E Neutralizing Broth <sup>†</sup>				✓					
EC Medium w/MUG/EC Broth w/ MUG <sup>†</sup>			✓						
m FC Agar/M-FC Agar <sup>†</sup>			✓						
Fraser Broth Base and/Fraser Broth Base Modified <sup>†</sup>						✓			
Fraser Broth Supplement/Fraser Broth Additive <sup>†</sup>						✓			
HYcheck™				✓					
LPM Agar Base <sup>†</sup>						✓			
Lactose Broth <sup>†</sup>			✓						
Lauryl Tryptose Broth w/MUG/ Lauryl Sulfate Broth w/ MUG <sup>†</sup>			✓						
Lethen Agar <sup>†</sup>				✓					
Lethen Broth/Lethen Broth AOAC <sup>†</sup>				✓					
Levine EMB Agar/Levine Eosin Methylene Blue Agar <sup>†</sup>			✓						
Listeria Antisera						✓			
M17 Broth <sup>†</sup>							✓		
Malt Extract Agar <sup>†</sup>								✓	
Malt Extract Broth <sup>†</sup>								✓	
McBride Listeria Agar <sup>†</sup>						✓			
Microbial Content Test Agar/ Trypticase™ Soy Agar w/Lec & Polysorbate <sup>†</sup>				✓					
Milk Agar IDF Formulation							✓		
Modified Listeria Enrichment Broth						✓			

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Dairy Product Testing - Products and Applications

PRODUCTS	APPLICATIONS								
	Beta-Lactams in Milk	Clostridium	Coliform Analysis	Environmental	Lipolytic Microorganisms	Listeria Analysis	Standard Plate Count	Streptococcus Detection	Yeast & Mold Analysis
Modified Oxford Antimicrobial Supplement						✓			
Motility Test Medium <sup>†</sup>						✓			
Moxalactam Antimicrobial Supplement						✓			
Neutralizing Buffer <sup>†</sup>				✓					
Nutrient Agar <sup>†</sup>				✓					
Nutrient Broth <sup>†</sup>		✓							
Oxford Antimicrobial Supplement						✓			
Oxford Medium Base/ Oxford Agar Base Modified <sup>†</sup>						✓			
PALCAM Medium Base and PALCAM Antimicrobial Supplement						✓			
PM Indicator Agar	✓								
PM Negative Control	✓								
PM Positive Control	✓								
Penase Concentrate/Penicillinase Concentrate <sup>†</sup>	✓								
Phenylethanol Agar/Phenylethanol Alcohol Agar <sup>†</sup>								✓	
Plate Count Agar/Standard Methods Agar <sup>†</sup>							✓		
Potato Dextrose Agar <sup>†</sup>									✓
Spirit Blue Agar w/Lipase Reagent					✓				
Strep Grouping Kit <sup>†</sup>								✓	
Subtilis Spore Suspension	✓								
Thermospore Suspension PM	✓								
Tryptic Soy Broth/ <b>TRYPTICASE</b> <sup>™</sup> Soy Broth <sup>†</sup>				✓					
Tryptone Glucose Extract Agar/ <b>TRYPTICASE</b> Glucose Extract Agar							✓		
UVM Modified Enrichment Broth/ UVM Mod. Listeria Enrichment Broth <sup>†</sup>						✓			
Universal Preenrichment Broth						✓			
Violet Red Bile Agar <sup>†</sup>			✓						
Violet Red Bile Agar w/MUG <sup>†</sup>			✓						
Yeast Extract Glucose Chloramphenicol Agar									✓

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

## Environmental Sampling and Disinfectant Testing

CULTURE MEDIA	HYGIENE CONTACT SLIDES
Aseptic Commissioning Medium Bushnell-Haas Broth D/E Neutralizing Agar <sup>†</sup> D/E Neutralizing Broth <sup>†</sup> Disinfectant Test Broth AOAC Egg Meat Medium Lethen Agar <sup>†</sup> Lethen Broth/Lethen Broth AOAC <sup>†</sup> Microbial Content Test Agar /Trypticase™ Soy Agar w/ Lec. & Polysorbate <sup>†</sup> Neutralizing Buffer <sup>†</sup> Synthetic Broth AOAC	HYcheck™ D/E Neutralizing Agar HYcheck for Disinfection Control HYcheck for Enterobacteriaceae HYcheck for Plate Count Agar with TTC HYcheck for Total Count HYcheck for Yeasts and Molds HYcheck for Yeasts and Molds with TTC
CONTACT PLATES	SETTLING PLATES (TRIPLE WRAPPED)
<b>DOUBLE WRAPPED (IRRADIATED)</b> D/E Neutralizing Agar Sabouraud Dextrose Agar w/Lecithin and Polysorbate 80 Standard Methods Agar Tryptic Soy Agar w/Lecithin and Polysorbate 80	<b>STERILE 100mm OR 150mm</b> Sabouraud Dextrose Agar Tryptic Soy Agar Tryptic Soy Agar w/Lecithin and Polysorbate 80

## Fermentation Products

FERMENTATION*		
Beef Extract, Desiccated <sup>†</sup> Brain Heart Infusion <sup>†</sup> Brucella Broth <sup>†</sup> Casamino Acids/Select Casamino Acids <sup>†</sup> Casamino Acids, Technical <sup>†</sup> Casitone Eugon Broth/ <b>EUGONBROTH™</b> Gelatone/Gelysate Peptone <sup>†</sup> Heart Infusion Broth <sup>†</sup> M Broth <sup>†</sup>	Malt Extract <sup>†</sup> Neopeptone Peptamin Peptone, Bacto® Peptone Bacteriological, Technical Proteose Peptone/Meat Peptone <sup>†</sup> Proteose Peptone No. 2 Proteose Peptone No. 3 Soytone/ <b>PHYTONE™</b> Peptone <sup>†</sup> TC Lactalbumin Hydrolysate TC Yeastolate/Yeastolate, TC <sup>†</sup>	Todd Hewitt Broth <sup>†</sup> Tryptic Soy Broth/ <b>TRYPTICASE™</b> Soy Broth Tryptic Soy Broth w/o Dextrose/ <b>TRYPTICASE</b> Soy Broth w/o Dextrose <sup>†</sup> Tryptone/Select <b>TRYPTICASE</b> Peptone <sup>†</sup> Tryptose/ <b>POLYPEPTONE™</b> Peptone <sup>†</sup> Tryptose Phosphate Broth <sup>†</sup> Yeast Extract/Select Yeast Extract <sup>†</sup> Yeast Extract, Technical <sup>†</sup>

\*Custom formulations and packaging are also available.

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Escherichia coli - Products and Applications

PRODUCTS	APPLICATIONS					
	Nonselective Media	Purity Plate	Selective-Differential Media	Serology	Specimen Collection and Shipment	Rapid Test
Brain Heart Infusion Agar <sup>†</sup>	✓					
BAGG Broth			✓			
Brilliant Green Bile Agar <sup>†</sup>			✓			
CULTURESWAB™ Transport System Amies Medium					✓	
CULTURESWAB Transport System Amies Medium w/o Charcoal/CULTURETTE™ Amies w/o Charcoal <sup>†</sup>					✓	
CULTURESWAB Transport System Cary-Blair Medium/Anaerobic CULTURETTE Cary-Blair Single <sup>†</sup>					✓	
CULTURESWAB Transport System Stuart's Medium Modified/CULTURETTE Modified Stuart's Medium <sup>†</sup>					✓	
EC Medium w/MUG <sup>†</sup>			✓			
EMB Agar/Eosin Methylene Blue Agar Modified <sup>†</sup>			✓			
E. coli H Antiserum H7 <sup>†</sup>				✓		
E. coli O Antiserum O157 <sup>†</sup>				✓		
m Endo Agar LES/M-Enda Agar LES <sup>†</sup>			✓			
m Endo Broth MF <sup>®</sup> /M-Endo Broth			✓			
EZ Coli™ Rapid Detection System for E. coli 0157						✓
Lauryl Tryptose Broth w/MUG/Lauryl Sulfate Broth w/MUG <sup>†</sup>			✓			
MacConkey Agar <sup>†</sup>			✓			
MacConkey Sorbitol Agar/MacConkey II Agar w/Sorbitol <sup>†</sup>			✓			
Nutrient Agar <sup>†</sup>	✓					
Nutrient Agar w/MUG <sup>†</sup>	✓					
Transport Medium Amies					✓	
Transport Medium Stuart <sup>†</sup>					✓	
Transport Medium Amies w/o Charcoal					✓	
Tryptic Soy Agar/TRYPICASE™ Soy Agar <sup>†</sup>	✓	✓				
TSA Blood Agar Base	✓	✓				
Veal Infusion Agar <sup>†</sup>	✓	✓				
Violet Red Bile Agar <sup>†</sup>			✓			
Violet Red Bile Agar w/MUG <sup>†</sup>			✓			

BIOCHEMICAL TESTS			
Purple Agar/Broth Base		Ammonium Citrate	-
Adonitol	-	Gelatin	-
Dulcitol	V	H <sub>2</sub> S	-
Glucose	+	Indole	+
Inositol	-	KCN	-
Lactose	V	Methyl Red	+
Mannitol	+	Phenylalanine	-
Salicin	V	Sodium Malonate	-
Sucrose	V	Tryptone Water	+
		Urease	-
		Voges-Proskauer	-
		<b>Key</b>	
		- Negative	
		+ Positive	
		d Delayed	
		V Variable +, d or -	

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Food and Beverage Testing - Products and Applications

PRODUCTS	APPLICATIONS													
	Bacillus	Beverage Analysis	Campylobacter Analysis	Clostridium	Coliform Analysis	Environmental	Lactobacillus	Listeria Analysis	Salmonella/Shigella spp.	Staphylococcus	Total Plate Count	Vibrio	Yeast & Mold Analysis	Yersinia spp.
A-1 Medium/A-1 Broth <sup>†</sup>				✓										
APT Agar <sup>†</sup>						✓								
APT Broth <sup>†</sup>						✓								
Baird-Parker Agar Base w/EY Tellurite Enrichment <sup>†</sup>									✓					
Bismuth Sulfite Agar <sup>†</sup>								✓						
Brilliant Green Agar <sup>†</sup>								✓						
Brilliant Green Agar Modified (Edel-Kampelmacher)								✓						
Brilliant Green Bile 2%/Brilliant Green Bile Broth 2% <sup>†</sup>				✓										
Brucella Agar <sup>†</sup>		✓												
Brucella Broth <sup>†</sup>		✓												
Bryant and Burkey Medium			✓											
Buffered Peptone Water <sup>†</sup>								✓						
Campylobacter Agar Kit Blaser/Campylobacter Agar w/5 Antimicrobics and 10% Sheep Blood <sup>†</sup>		✓												
Campylobacter Agar Kit Skirrow		✓												
Campylobacter Agar Base		✓												
Coagulase Plasma (Rabbit) <sup>†</sup>									✓					
Coagulase Plasma EDTA (Rabbit)/ Coagulase Plasma, Rabbit w/ EDTA <sup>†</sup>									✓					
Cooke Rose Bengal Agar												✓		
Cooked Meat Medium <sup>†</sup>			✓											
D/E Neutralizing Agar <sup>†</sup>					✓									
D/E Neutralizing Broth <sup>†</sup>					✓									
DNase Test Agar <sup>†</sup>									✓					
DNase Test Agar w/Methyl Green/ DNase Test Agar w/Toluidine Blue <sup>†</sup>									✓					
DRBC Agar												✓		
Demi-Fraser Broth Base							✓							
Desoxycholate Citrate Agar Hynes <sup>†</sup>								✓						
Differential Reinforced Clostridial Agar			✓											
EC Medium/EC Broth <sup>†</sup>				✓										
EC Medium with MUG/EC Broth w/MUG <sup>†</sup>				✓										
Elliker Broth <sup>†</sup>						✓								
m Endo Agar LES <sup>†</sup>				✓										
m Endo Broth MF <sup>®</sup> / mEndo Broth, ALPHA <sup>†</sup>				✓										
EZ Coli™ Rapid Detection System				✓										
Fluid Thioglycollate Medium <sup>†</sup>			✓											
Fraser Broth Base/Fraser Broth Base Modified <sup>†</sup>							✓							
Fraser Broth Supplement/Fraser Broth Base Supplement <sup>†</sup>							✓							
Hektoen Enteric Agar <sup>†</sup>								✓						
HYcheck™ D/E Neutralizing Agar					✓									
HYcheck for Disinfection Control					✓									
HYcheck for Enterobacteriaceae					✓									
HYcheck for Total Count					✓									
HYcheck for Yeasts and Molds					✓									
HYcheck for Yeasts and Molds w/TTC					✓									
LPM Agar Base <sup>†</sup>							✓							

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Food and Beverage Testing - Products and Applications

PRODUCTS	APPLICATIONS													
	Bacillus	Beverage Analysis	Campylobacter Analysis	Clostridium	Coliform Analysis	Environmental	Lactobacillus	Listeria Analysis	Salmonella/Shigella	Staphylococcus	Total Plate Count	Vibrio	Yeast & Mold Analysis	Yersinia sp.
Lactobacilli MRS Agar/MRS Agar <sup>†</sup>						✓								
Lactobacilli MRS Broth/MRS Broth <sup>†</sup>						✓								
Lactose Broth <sup>†</sup>							✓							
Lauryl Tryptose Broth/Lauryl Sulfate Broth <sup>†</sup>				✓										
Lauryl Tryptose Broth with MUG/Lauryl Sulfate Broth w/MUG <sup>†</sup>				✓										
Lethen Agar <sup>†</sup>					✓									
Lethen Broth/Lethen Broth AOAC <sup>†</sup>					✓									
Levine EMB Agar/Levine Eosin Methylene Blue Agar <sup>†</sup>				✓										
Listeria Enrichment Broth <sup>†</sup>							✓							
Listeria O Antisera							✓							
Liver Veal Agar <sup>†</sup>			✓											
Lysine Medium		✓											✓	
MYP Agar	✓													
M Broth <sup>†</sup>								✓						
MacConkey Agar <sup>†</sup>				✓				✓						
MacConkey Sorbitol Agar/MacConkey II Agar w/Sorbitol <sup>†</sup>				✓										
Malt Agar <sup>†</sup>													✓	
Malt Extract Agar <sup>†</sup>													✓	
Malt Extract Broth <sup>†</sup>													✓	
Mannitol Salt Agar <sup>†</sup>									✓					
McBride Listeria Agar <sup>†</sup>							✓							
McClung Toabe Agar Base <sup>†</sup>			✓											
Microbial Content Test Agar/ TRYPTICASE <sup>™</sup> Soy Agar w/Lec. & Polysorbate 80 <sup>†</sup>					✓									
Minerals Modified Glutamate Agar				✓										
Modified Listeria Enrichment Broth							✓							
Modified Oxford Antimicrobial Supplement							✓							
Modified EC Medium w/Indicator and Novobiocin Antimicrobial Supplement				✓										
Motility Test Medium <sup>†</sup>							✓							
Moxalactam Antimicrobial Supplement <sup>†</sup>							✓							
Muller Kauffmann Tetrathionate Broth Base								✓						
Mycological Agar/MYCOPHIL <sup>™</sup> Agar <sup>†</sup>													✓	
Neutralizing Buffer <sup>†</sup>					✓									
OGYE Agar Base w/Antimicrobial Vial Oxytetracycline													✓	
Orange Serum Agar <sup>†</sup>	✓													
Orange Serum Broth Concentrate 10X	✓													
Oxford Antimicrobial Supplement							✓							
Oxford Medium Base/Oxford Agar Base Modified <sup>†</sup>							✓							
PALCAM Medium Base with PALCAM Antimicrobial Supplement							✓							
Plate Count Agar/Standard Methods Agar <sup>†</sup>									✓					
Potato Dextrose Agar <sup>†</sup>													✓	
Potato Dextrose Broth <sup>†</sup>													✓	
Rappaport-Vassiliadis (MSRV) Medium Semisolid Modification								✓						

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Food and Beverage Testing - Products and Applications

PRODUCTS	APPLICATIONS													
	Bacillus	Beverage Analysis	Campylobacter Analysis	Clostridium	Coliform Analysis	Environmental	Lactobacillus	Listeria Analysis	Salmonella/Shigella sp.	Staphylococcus	Total Plate Count	Vibrio	Yeast & Mold Analysis	Yersinia sp.
Rappaport-Vassiliadis R10 Broth								✓						
Raka-Ray No. 3 Broth	✓						✓							
Raka-Ray No. 3 Medium	✓						✓							
Reinforced Clostridial Agar†			✓											
Rogosa SL Agar†							✓							
Rogosa SL Broth†							✓							
Rose Bengal Agar Base												✓		
Rose Bengal Antimicrobial Supplement C/ Chloramphenicol Selective†												✓		
SFP Agar Base/TSN Agar†			✓											
SS Agar/Salmonella Shigella Agar†								✓						
Sabouraud Dextrose Agar†												✓		
Sabouraud Dextrose Broth												✓		
Salmonella H Antiserum Poly a-z†								✓						
Salmonella O Antisera †								✓						
Selenite Broth/Selenite F Broth†								✓						
Selenite Cystine Broth†								✓						
Special Yeast and Mold Medium												✓		
m Staphylococcus Broth									✓					
Staphylococcus Medium 110/Staphylococcus Agar 110†									✓					
Staph Latex Test/STAPHYLOSLIDE™ Test†									✓					
Sulfite Agar			✓											
TCBS Agar†											✓			
TT Broth Base, Hajna/Tetrathionate H†								✓						
Tetrathionate Broth Base†								✓						
Tomato Juice Agar†							✓							
Tryptic Soy Broth/TRYPICASE™ Soy Broth†					✓									
Tryptone Water				✓										
UBA Medium/Universal Beer Agar†	✓													
UVM Modified Listeria Enrichment Broth†								✓						
Universal Preenrichment Broth								✓	✓					
Vibrio Cholerae Antisera											✓			
Violet Red Bile Agar†				✓										
Violet Red Bile Agar w/MUG†				✓										
Violet Red Bile Glucose Agar				✓										
WL Differential Medium/WL Differential Agar†	✓													
WL Nutrient Broth	✓													
WL Nutrient Medium†	✓													
Wort Agar†	✓											✓		
XLD Agar†								✓						
XLT4 Agar Base and XLT4 Supplement								✓						
YM Agar												✓		
YM Broth												✓		
Yeast Extract Glucose Chloramphenicol Agar												✓		
Yersinia Selective Agar Base/CIN Agar Base†													✓	
Yersinia Antimicrobial Supplement CN/CN Inhibitor†														✓

† Available from Difco &amp; Becton Dickinson Microbiology Systems.

## McFarland Standard Preparation

"MCFARLAND"	SULFURIC ACID, 1% AQUEOUS SOLUTION ML	BARIUM CHLORIDE, 1% AQUEOUS SOLUTION ML	CORRESPONDING DENSITY OF BACTERIA -10 <sup>6</sup>	INTERNATIONAL UNITS (IU) OF OPACITY
1	9.9	0.1	300	3
2	9.8	0.2	600	7
3	9.7	0.3	900	10
4	9.6	0.4	1200	12
5	9.5	0.5	1500	15
6	9.4	0.6	1800	–
7	9.3	0.7	2100	20
8	9.2	0.8	2400	–
9	9.1	0.9	2700	–
10	9.0	1.0	3000	30

1. Prepare the tubes by mixing 1% sulfuric acid with 1% barium chloride according to the table.
2. Make sure the tubes are uniform in size and made of chemically resistant glass.
3. Plug the tubes with rubber stoppers and carefully seal with paraffin. Store the tubes upright.
4. To estimate bacterial cell density, compare the bacterial suspension with the standards.
5. The above set is used to determine bacterial density in saline suspension. To estimate bacterial density in broth, make the set by dissolving the sulfuric acid and barium chloride in sterile broth.

REFERENCE IN: Gradwohl's Clinical Laboratory Methods and Diagnosis. In. A.C. Sonnenwirth and L. Jarett (ed.). C.V. Mosby Company, 1980 p. 1363.

## Molecular Genetics - Media and Ingredients

LB MEDIA	NZ MEDIA AND INGREDIENTS	GENERAL MEDIA AND INGREDIENTS
LB Agar, Lennox/LB Agar (Lennox LAgar) <sup>†</sup>	Casein Digest/Casein Digest Peptone <sup>†</sup>	M9CA Medium
LB Agar, Miller (Luria-Bertani)/ Luria Agar <sup>†</sup>	NZCYM Broth <sup>†</sup>	M9 Minimal Salts, 5x/ M9 Minimal Salt <sup>†</sup>
LB Broth, Lennox/LB Broth <sup>†</sup>	NZM Broth	Minimal Agar Davis
LB Broth, Miller/Luria Broth <sup>†</sup>	NZYM Broth/NZY Broth <sup>†</sup>	Minimal Broth Davis w/o Dextrose
Luria Agar Base, Miller		SOB Medium
Luria Broth Base, Miller		Terrific Broth
		2xYT
		YPD Agar/YEPD Agar <sup>†</sup>
		YPD Broth/YEPD Broth <sup>†</sup>
		Yeast Nitrogen Base <sup>‡</sup>
		Yeast Nitrogen Base w/o Amino Acids <sup>†</sup>
		Yeast Nitrogen Base w/o Amino Acids & Ammonium Sulfate

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Mycobacteria

MYCOBACTERIA		
ATS Medium	ESP Myco GS	Middlebrook 7H9 Broth/Middlebrook 7H9 Broth Base <sup>†</sup>
Dubos Albumin Broth	ESP Myco PVNA	Middlebrook 7H10 Agar/Middlebrook & Cohn 7H10 Agar Base <sup>†</sup>
Dubos Broth Base	Lowenstein Medium Base/ Lowenstein Jensen Medium Base <sup>†</sup>	Middlebrook ADC Enrichment <sup>†</sup>
Dubos Medium Albumin	Lowenstein Medium Gruft	Middlebrook OADC Enrichment
Dubos Oleic Agar Base	Lowenstein Medium, Jensen	Middlebrook OADC Enrichment w/WR 1339
Dubos Oleic Albumin Complex/ Oleic Albumin Complex <sup>†</sup>	Lowenstein Medium, Jensen Deeps	Mycobacteria 7H11 Agar/Seven H11 Agar Base <sup>†</sup>
ESP Myco	Lowenstein Medium w/5% NaCl	Petragnani Medium

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Mycology

CLINICAL REAGENTS FOR DETECTION OF <i>CANDIDA ALBICANS</i>
BiGGY Agar <sup>†</sup>
Candida Isolation Agar
Candida Albicans Antiserum

MEDIA FOR CLASSIFICATION OF YEASTS
Yeast Carbon Base
Yeast Morphology Agar
Yeast Nitrogen Base <sup>†</sup>
Yeast Nitrogen Base w/o Amino Acids <sup>†</sup>
Yeast Nitrogen Base w/o Amino Acids and Ammonium Sulfate

MEDIA	NON-INHIBITORY MEDIA FOR FUNGI		
	General Use	Vitamin Assay	Candida and/or Dermatophytes
Blood Agar Base <sup>†</sup>	✓		
Brain Heart Infusion <sup>†</sup>	✓		
Brain Heart Infusion Agar <sup>†</sup>	✓		
Brain Heart Infusion w/PAB <sup>†</sup>	✓		
Corn Meal Agar <sup>†</sup>			✓
Fluid Sabouraud Medium	✓		
Malt Agar <sup>†</sup>	✓		
Malt Extract Agar <sup>†</sup>	✓		
Malt Extract Broth <sup>†</sup>	✓		
Mycological Agar/MYCOPHIL™ Agar <sup>†</sup>	✓		
Mycological Agar w/Low pH/MYCOPHIL™ w/ Low pH <sup>†</sup>	✓		
Oatmeal Agar	✓		
Pagano Levin Base and TTC			✓
Potato Dextrose Agar <sup>†</sup>	✓		
Potato Dextrose Broth <sup>†</sup>	✓		
Rice Extract Agar <sup>†</sup>			✓
SABHI Agar Base/SABHI Agar <sup>†</sup>	✓		
Sabouraud Agar Modified/ Sabouraud Dextrose Agar Emmons <sup>†</sup>	✓		
Sabouraud Dextrose Agar <sup>†</sup>	✓		
Sabouraud Dextrose Broth	✓		
Sabouraud Maltose Agar <sup>†</sup>	✓		
Special Yeast and Mold Medium	✓		
YM Agar	✓		
YM Broth	✓		

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Mycology

## Reagents for Direct Microscopic Detection of Fungi

PRODUCT	USES
Gram Stain Sets and Reagents <sup>†</sup>	General screening
3-Step™ Gram Stain Sets and Reagents	General screening
SpotTest™ Calcofluor White Reagent/Calcofluor White Droppers <sup>†</sup> A nonspecific fluorochrome stain which allows the rapid examination of clinical specimens for the presence of fungi under an FA microscope	Screening of cultures for presence of fungi
SpotTest India Ink/India Ink Droppers <sup>†</sup> For use in the direct microscopic examination of clinical material for the presence of encapsulated yeast cells	Screening CSF for <i>Cryptococcus neoformans</i> Negative staining
SpotTest KOH 10%/Potassium Hydroxide 10% Droppers <sup>†</sup> A 10% solution of potassium hydroxide	Used in wet mounts of clinical specimens to examine for the presence of fungi
SpotTest Lactophenol Cotton Blue Stain/ Lactophenol Cotton Blue Droppers <sup>†</sup> For use in the direct mounting and staining of yeast and molds	Teased mount method Slide culture method

## Selective and/or Differential Media for Fungi

GENERAL USE	CANDIDA AND/OR DERMATOPHYTES
Brain Heart CC Agar <sup>†</sup>	BiGGY Agar <sup>†</sup>
Cooke Rose Bengal Agar and Antimicrobial Vial A	Candida BCG Agar Base
DRBC Agar	Candida Isolation Agar
HC Agar Base <sup>†</sup>	DTM Agar
Mycobiotic Agar/Mycosel™ Agar <sup>†</sup>	Littman Oxgall Agar
Mycological Agar w/Low pH/Mycophil™ Agar w/Low pH <sup>†</sup>	Trichophyton Agars 1-7
Rose Bengal Agar Base and Rose Bengal Antimicrobial Supplement C	
Yeast Extract Glucose Chloramphenicol Agar	

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

## Pasco® Panel Contents - Selection Guide

ANTIMICROBIAL AGENT (MCG/ML)	MIC GRAM-NEG.	MIC/ID GRAM-NEG.	MIC/ID GRAM- NEG. SPECIAL	BREAKPOINT/ID GRAM-NEG.	MIC ORAL	MIC SUPP.	MIC GRAM-POS.	MIC/ID GRAM-POS.
AK Amikacin	32,24,16-1	32-4	32-4	32,24,16-4	—	—	—	—
A/C Amoxicillin/ Clavulanic Acid	—	—	16-4/8-2	—	16-0.5/8-0.25	16-0.25/8-0.12	—	—
AM Ampicillin	16-2	16-2	16-8	16-8	16-0.12	—	8-0.12	8-0.12
A/S Ampicillin/ Sulbactam	16-2/8-1	16-8/8-4	16-8/8-4	16-8/8-4	—	—	16-2/8-1	16-2/8-1
AZT Aztreonam	16-4	16-8	16-8	16-8	—	—	—	—
CB Carbenicillin	—	—	—	—	256-128	256-128	—	—
CCL Cefaclor	—	—	—	—	16-1	16-1	—	—
CFZ Cefazolin	16-2	16-2	16-8	16-8	—	—	16-2	16-8
FIX Cefixime	2-1	2-1	—	2-1	2-0.25	—	2-0.5	2-1
CPZ Cefoperazone	32-4	32-16	32-16	32-16	—	—	—	—
CTX Cefotaxime	32-8	—	32-8	32-8	—	16-0.12	32-4	32-8
CTN Cefotetan	32-8	32-8	32-16	32-16	—	—	—	—
CX Cefoxitin	16-2	16-4	16-8	16-8	—	—	—	—
TAZ Ceftazidime	16-2	16-4	16-4	16-8	—	—	—	—
CTZ Ceftizoxime	32-4	32-8	—	32-8	—	—	—	—
FRX Ceftriaxone	32-4	32-8	32-8	32-8	—	8-0.25	32-4	32-8
CFX Cefuroxime	16-2	16-2	16-8	16-8	16-0.5	4-0.25	16-4	16-8
CF Cephalothin	—	—	16-8	16-8	16-1	16-8	—	—
C Chloramphenicol	16-2	—	16-8	16-8	16-0.5	8-1	16-4	16-4
CIP Ciprofloxacin	2-0.25	2-0.25	2-1	2-1	2-0.06	—	2-0.25	2-0.25
CLM Clarithromycin	—	—	—	—	4-0.25	4-0.12	4-0.5	4-0.5
CD Clindamycin	—	—	—	—	2-0.25	2-0.06	2-0.25	2-0.25
EN Enoxacin	—	—	—	4-2	—	—	—	—
E Erythromycin	—	—	—	—	4-0.25	4-0.25	4-0.5	4-0.5
GM Gentamicin	8,6,4-0.25	8,6,4-0.5	8,6,4-0.5	8,6,4-1	—	—	500,8,6,4-0.25	500,8,6,4-1
IMI Imipenem	—	—	8-4	—	—	8-0.06	—	—
LOM Lomefloxacin	4-2	4-2	—	4-2	4-2	—	4-2	4-2
MZ Mezlocillin	64-8	64-8	64-16	—	—	—	—	—
NET Netilmicin	16-8	—	—	—	—	—	—	—
FD Nitrofurantoin	64-32	64-32	64-32	64-32	64-32	64-32	—	64-32
NOR Norfloxacin	8-4	—	8-4	8-4	8-4	8-2	—	—
OFX Ofloxacin	4-0.5	4-0.5	4-2	4-2	4-0.12	8-1	4-0.25	4-1
OX Oxacillin	—	—	—	—	6,4-0.5	—	6,4-0.5	6,4-1
P Penicillin	—	—	—	—	8-0.03	2-0.03	8-0.03	8-0.03
PIP Piperacillin	64-8	—	64-16	64-8	—	—	—	—
P/T Piperacillin/ Tazobactam	64-8/4	64-8/4	64-16/4	64-16/4	—	—	16-2/4	16-4/4
RI Rifampin	—	—	—	—	2-1	4-0.5	2-1	2-1
STR Streptomycin	—	—	—	—	—	—	1000	1000
SFX Sulfisoxazole	—	—	—	—	256	—	—	—
TE Tetracycline	—	—	8-4	8-4	8-0.5	8-1	8-2	—
TC Ticarcillin	—	—	—	64-8	—	—	—	—
T/C Ticarcillin/ Clavulanic Acid	—	—	64-16/2	—	—	64-16/2	—	—
TO Tobramycin	8,6,4-0.25	8,6,4-0.5	8,6,4-0.5	8,6,4-1	—	—	8,6,4-0.25	8,6,4-1
T/S Trimethoprim/ Sulfamethoxazole	2-1/38-19	2-1/38-19	2-38	2/38	2-0.5/38-9.5	2-0.5/38-9.5	2-0.5/38-9.5	2/38
VA Vancomycin	—	—	—	—	16-4	4-0.5	16-1	16-2
<b>Biochemical Substrates</b>		30*	30*	30*				18**

\* Gram-Negative Panels

\*\* Gram-Positive Panels

# Peptones & Hydrolysates Selection Guide

## Typical Analyses

PRODUCT	AMINO ACIDS - %																	
	Alanine	Arginine	Aspartic Acid	Cystine	Glutamic Acid	Glycine	Histidine	Isoleucine	Leucine	Lysine	Methionine	Phenylalanine	Proline	Serine	Threonine	Tryptophan	Tyrosine	Valine
Beef Extract	2.54	1.39	1.67	0.18	6.01	4.14	4.94	0.53	1.00	1.45	0.30	<0.01	2.16	0.90	0.67	0.05	1.99	0.86
Beef Extract, Dessicated	8.96	5.66	4.30	0.17	12.55	16.25	2.50	1.45	3.63	3.27	1.08	2.00	9.58	2.10	1.42	0.32	1.03	2.62
Casamino Acids	3.26	2.20	4.76	0.16	15.30	1.31	1.66	3.34	5.47	5.71	1.28	2.11	6.17	2.19	2.41	<0.01	0.47	4.30
Casamino Acids, Technical	1.64	1.77	3.42	0.34	10.97	1.09	1.31	0.13	2.54	2.14	1.19	5.23	4.44	2.64	1.99	0.01	1.33	3.24
Casein Digest	2.92	3.02	6.75	0.17	23.10	1.93	2.52	4.66	8.29	7.70	2.66	4.27	11.04	5.55	4.33	1.16	2.54	6.51
Casitone	3.01	3.76	6.61	0.02	20.03	1.97	2.17	4.16	8.74	13.62	1.71	4.02	8.57	4.82	3.74	0.14	2.09	4.06
Neopeptone	4.03	4.14	6.19	0.26	13.22	7.02	<0.01	0.36	3.65	5.16	2.00	8.67	6.73	4.22	3.69	0.96	4.21	4.96
Bacto Peptone	8.67	6.76	5.60	0.20	10.21	15.59	0.58	1.45	3.01	3.42	1.19	1.81	8.80	2.87	1.81	0.36	0.64	2.35
Proteose Peptone	6.50	5.12	7.28	0.87	11.95	9.68	2.01	3.04	5.66	5.33	1.97	2.86	5.93	3.49	3.14	0.60	2.35	3.76
Proteose Peptone No. 2	6.08	5.47	7.45	0.40	10.57	10.84	<0.01	1.00	3.57	5.22	1.51	7.94	5.31	4.64	3.90	0.94	1.92	4.73
Proteose Peptone No. 3	5.99	5.49	6.92	1.12	12.38	9.26	1.74	2.65	5.70	5.02	1.86	2.72	4.94	3.65	3.32	0.59	1.96	3.62
Soytone	2.46	3.82	7.27	1.45	12.76	2.51	1.24	2.37	4.03	3.45	0.86	2.46	2.92	2.87	2.17	0.47	1.93	2.65
TC Lactalbumin Hydrolysate	3.70	2.67	7.13	0.55	16.30	2.02	1.83	1.18	5.43	4.11	2.11	10.81	6.62	4.62	4.76	1.75	1.22	5.82
TC Yeastolate	4.84	2.99	5.58	0.45	10.53	4.02	<0.01	0.56	3.23	3.82	0.96	5.59	2.59	2.81	2.80	0.79	1.21	3.80
Tryptone	2.86	3.03	6.11	0.42	17.05	1.75	2.02	4.40	7.11	6.70	2.57	3.71	7.45	4.29	3.58	0.71	1.42	5.00
Tryptose	4.45	4.65	6.34	0.44	13.92	2.84	<0.01	0.34	3.67	4.64	1.92	7.52	6.33	4.09	3.55	0.62	2.21	1.93
Yeast Extract	5.36	3.02	6.69	0.74	14.20	3.25	1.20	3.23	4.69	5.15	1.05	2.53	2.60	2.84	2.95	1.36	1.20	3.79

PRODUCT	INORGANICS - %															CARBOHYDRATE Total (%)	
	Calcium	Chloride	Cobalt	Copper	Iron	Lead	Magnesium	Manganese	Phosphate	Potassium	Sodium	Sulfate	Sulfur	Tin	Zinc		
Beef Extract	0.068	1.284	<0.001	<0.001	<0.001	<0.001	0.239	<0.001	5.458	5.477	2.315	0.629	0.707	<0.001	<0.001		0.2
Beef Extract, Dessicated	0.018	1.576	<0.001	0.001	0.011	<0.001	0.022	<0.001	0.345	1.994	2.774	0.829	0.661	<0.001	0.002		<0.1
Casamino Acids	<0.001	7.400	<0.001	<0.001	<0.001	<0.001	0.002	<0.001	3.325	0.410	8.710	0.045	0.420	<0.001	<0.001		
Casamino Acids, Technical	0.019	21.212	<0.001	<0.001	<0.001	<0.001	0.008	<0.001	1.358	0.273	13.721	0.167	0.424	<0.001	0.002		3.4
Casein Digest	0.019	0.059	<0.001	<0.001	0.002	<0.001	0.018	<0.001	2.670	0.078	2.937	0.167	0.621	<0.001	0.002		4.2
Casitone	0.010	0.110	<0.001	<0.001	0.003	<0.001	0.019	<0.001	2.604	0.162	3.073	0.339	0.676	<0.001	0.004		0.2
Neopeptone	0.012	0.344	<0.001	<0.001	<0.001	<0.001	0.006	<0.001	2.209	0.149	2.057	0.340	0.657	<0.001	<0.001		0.8
Bacto Peptone	0.008	1.086	<0.001	<0.001	0.004	<0.001	0.007	<0.001	0.445	0.203	1.759	0.244	0.410	<0.001	0.001		6.9
Proteose Peptone	0.021	4.510	<0.001	<0.001	0.002	<0.001	0.027	<0.001	0.872	0.685	3.677	0.162	0.812	<0.001	0.002		<0.1
Proteose Peptone No. 2	0.024	3.644	<0.001	<0.001	<0.001	<0.001	0.024	<0.001	1.674	0.815	3.956	0.232	0.698	<0.001	0.003		1.3
Proteose Peptone No. 3	0.023	3.581	<0.001	<0.001	0.002	<0.001	0.027	<0.001	1.447	0.982	3.815	0.232	0.975	<0.001	0.007		1.4
Soytone	0.055	0.165	<0.001	<0.001	0.008	<0.001	0.161	<0.001	0.820	2.220	3.404	2.334	1.660	<0.001	0.001		24.0
TC Lactalbumin Hydrolysate	0.095	<0.010	<0.001	<0.001	0.002	<0.001	0.026	<0.001	1.309	0.932	1.357	0.379	0.750	<0.001	0.002		9.2
TC Yeastolate	0.002	0.299	<0.001	<0.001	<0.001	<0.001	0.025	<0.001	2.633	5.085	0.819	0.488	0.528	<0.001	0.007		10.3
Tryptone	0.013	0.186	<0.001	<0.001	<0.001	<0.001	0.017	<0.001	2.669	0.229	2.631	0.241	0.740	<0.001	0.003		7.7
Tryptose	0.001	1.886	<0.001	<0.001	0.002	<0.001	0.022	<0.001	2.144	0.679	3.410	0.308	0.737	<0.001	0.005		7.1
Yeast Extract	0.013	0.380	<0.001	<0.001	<0.001	<0.001	0.075	<0.001	3.270	3.195	1.490	0.091	0.634	<0.001	0.011		17.5

# Peptones & Hydrolysates Selection Guide

## Typical Analyses

PRODUCT	PHYSICAL CHARACTERISTICS					BIOLOGICAL TESTING - CFU/g					NITROGEN CONTENT		
	Ash (%)	Clarity 1% Soln (NTU)	Filterability (g/cm <sup>2</sup> )	Loss on Drying (%)	pH, 1% Soln	Coliform	Salmonella	Spore Count	Standard Plate Count	Thermophile Count	Total Nitrogen (%)	Amino Nitrogen	AN/TN (%)
Beef Extract	24.1	116.8	0.1	77.2*	5.4	neg	neg	299	117	33	11.2	3.8	33.8
Beef Extract, Dessicated	10.2	1.7	0.6	2.5	6.9	neg	neg	585	690	28	14.0	2.2	15.7
Casamino Acids	24.4	0.5	2.9	4.5	6.4	neg	neg	390	950	25	10.5	8.8	83.8
Casamino Acids, Technical	38.3	0.3	2.6	4.5	6.7	neg	neg	2375	2250	<50	8.1	6.4	79.0
Casein Digest	6.4	0.4	2.6	4.7	7.2	neg	neg	235	250	178	13.4	7.2	53.9
Casitone	7.0	0.6	1.7	3.7	7.2	neg	neg	300	1850	100	13.3	4.7	35.3
Neopeptone	7.0	1.2	0.3	3.2	7.4	neg	neg	175	400	75	13.7	3.3	23.8
Bacto Peptone	4.4	0.5	0.5	3.0	7.0	neg	neg	90	273	13	15.5	3.1	20.0
Proteose Peptone	11.1	1.4	0.9	3.1	7.2	neg	neg	393	443	73	14.0	2.9	20.7
Proteose Peptone No. 2	12.7	1.5	0.6	3.5	7.2	neg	neg	75	1450	<50	12.6	5.0	39.7
Proteose Peptone No. 3	11.4	2.2	0.5	4.0	7.2	neg	neg	890	915	25	13.2	3.5	26.5
Soytone	12.0	1.0	1.2	4.6	7.2	neg	neg	10	38	<3	9.4	3.1	33.0
TC Lactalbumin Hydrolysate	7.2	0.4	7.3	4.6	7.1	neg	neg	<50	300	<50	13.0	6.3	48.3
TC Yeastolate	13.0	1.4	4.5	3.6	6.9	neg	neg	175	175	<50	10.8	6.5	59.8
Tryptone	6.8	0.5	1.3	3.7	7.2	neg	neg	73	870	8	13.0	5.2	40.0
Tryptose	9.7	0.8	2.3	3.2	7.4	neg	neg	875	825	100	13.4	4.4	32.5
Yeast Extract	11.2	1.5	2.7	3.1	6.7	neg	neg	9	60	<5	10.9	6.0	55.0

PRODUCT	VITAMINS - µg/g											
	Biotin	Choline	Cyanocobalamin	Folic Acid	Inositol	Nicotinic Acid	PABA	Pantothenic Acid	Pyridoxine	Riboflavin	Thiamine	Thymidine
Beef Extract	0.1	1171.5	0.5	3.3	4113.2	774.7	20.0	91.0	7.3	0.4	<0.1	1093.4
Beef Extract, Dessicated	0.1	1300.0	<0.1	0.6	2100.0	138.1	40.5	8.7	2.8	<0.1	<0.1	111.3
Casamino Acids	<0.1	160.0	<0.1	<0.1	<100.0	<20.0	<5.0	<0.1	<0.1	1.8	1.2	<30.0
Casamino Acids, Technical	<0.1	<50.0	<0.1	<0.1	<38.0	<0.1	9.8	<0.6	<0.1	0.2	0.2	<14.0
Casein Digest	0.1	<40.0	<0.1	1.0	490.0	14.1	6.1	6.7	0.4	<0.1	1.5	8297.0
Casitone	0.2	550.0	<0.1	0.8	980.0	20.3	15.9	7.7	1.3	0.4	<0.1	342.9
Neopeptone	0.2	3100.0	<0.1	0.4	3600.0	52.2	2.9	16.0	2.3	1.3	<0.1	<14.0
Bacto Peptone	0.2	2000.0	<0.1	0.3	2400.0	21.9	<0.5	5.9	1.7	3.9	<0.1	413.0
Proteose Peptone	0.1	2300.0	<0.1	0.4	5000.0	79.9	4.2	20.0	1.1	<0.1	1.2	99.7
Proteose Peptone No. 2	0.3	4500.0	<0.1	0.5	4700.0	157.1	1.2	47.0	4.0	6.4	1.6	1319.0
Proteose Peptone No. 3	0.4	3700.0	<0.1	0.3	8900.0	124.2	<0.5	20.0	1.3	6.8	0.1	659.6
Soytone	0.2	2200.0	<0.1	3.0	2100.0	19.1	9.0	13.0	11.0	<0.1	1.2	113.2
TC Lactalbumin Hydrolysate	<0.1	280.0	<0.1	0.2	360.0	<0.1	11.3	4.3	1.0	8.5	0.5	<14.0
TC Yeastolate	6.7	3400.0	<0.1	25.2	1900.0	945.0	96.6	300.0	77.5	21.8	54.3	2975.0
Tryptone	0.1	350.0	<0.1	0.3	1400.0	97.8	3.7	5.3	0.6	<0.1	0.4	93.4
Tryptose	0.2	2700.0	<0.1	0.4	5400.0	47.4	11.4	16.0	1.4	4.3	0.1	769.0
Yeast Extract	3.3	300.0	<0.1	1.5	1400.0	597.9	763.0	273.7	43.2	116.5	529.9	217.5

\*Represents Total Solids value rather than loss on drying.

## Pharmaceutical Testing - Products and Applications

PRODUCTS	APPLICATIONS		
	Sterility Testing	Antibiotic Assay	Vitamin/Amino Acid Assay
Agar Medium No. F	✓		
Antibiotic Medium 1 <sup>†</sup>		✓	
Antibiotic Medium 2/Base Agar, Penicillin Assay <sup>†</sup>		✓	
Antibiotic Medium 3 <sup>†</sup>		✓	
Antibiotic Medium 4/Yeast Beef Agar <sup>†</sup>		✓	
Antibiotic Medium 5/ Streptomycin Assay Agar w/YE <sup>†</sup>		✓	
Antibiotic Medium 8/Base Agr w/Low pH <sup>†</sup>		✓	
Antibiotic Medium 9/Polymixin Base Agar <sup>†</sup>		✓	
Antibiotic Medium 10/Polymixin Seed Agar <sup>†</sup>		✓	
Antibiotic Medium 11/Neomycin Assay Agar <sup>†</sup>		✓	
Antibiotic Medium 12		✓	
Antibiotic Medium 19/Nystatin Assay Agar <sup>†</sup>		✓	
APT Agar <sup>†</sup>			✓
APT Broth <sup>†</sup>			✓
Aseptic Commissioning Medium	✓		
B12 Assay Medium USP			✓
B12 Culture Agar USP			✓
B12 Inoculum Broth USP			✓
Biotin Assay Medium			✓
Cystine Assay Medium			✓
Fluid Thioglycollate Medium <sup>†</sup>	✓		
Fluid Thioglycollate Medium w/Beef Extract <sup>†</sup>	✓		
Fluid Thioglycollate Medium w/K Agar	✓		
Folic Acid Assay Medium			✓
Folic Acid Casei Medium			✓
Folic AOAC Medium			✓
Folic Buffer A, Dried			✓
Inositol Assay Medium			✓
Lactobacilli Agar AOAC			✓
Lactobacilli Broth AOAC			✓
Lysine Assay Medium			✓
Methionine Assay Medium			✓

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Pharmaceutical Testing - Products and Applications

PRODUCTS	APPLICATIONS		
	Sterility Testing	Antibiotic Assay	Vitamin/Amino Acid Assay
Micro Assay Culture Agar			✓
Micro Inoculum Broth			✓
Neurospora Culture Agar			✓
Niacin Assay Medium			✓
NIH Thioglycollate Medium <sup>†</sup>	✓		
Panthenol Assay Medium			✓
Pantothenate Assay Medium			✓
Pantothenate Medium AOAC USP			✓
Pyridoxine Y Medium			✓
Riboflavin Assay Medium			✓
Sterility Bottles w/ Screw Cap Fluid Thioglycollate Medium/ Fluid Thioglycollate Medium <sup>†</sup>	✓		
Sterility Bottles w/ Screw Cap Tryptic Soy Broth/Trypticase™ Soy Broth <sup>†</sup>	✓		
Sterility Bottles w/ Septum Fluid A/ Peptone Water (0.1%) <sup>†</sup>	✓		
Sterility Bottles w/ Septum Fluid D/ Peptone Water (0.1%) w/Polysorbate <sup>†</sup>	✓		
Sterility Bottles w/ Septum Fluid Thioglycollate Medium <sup>†</sup>	✓		
Sterility Bottles w/ Septum Tryptic Soy Broth/ Septum Fluid Thioglycollate Medium <sup>†</sup>	✓		
Thiamine Assay Medium LV			✓
Tryptic Soy Broth/Trypticase™ Soy Broth <sup>†</sup>	✓		
Vitamin B12 Assay Medium			✓

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

## Salmonella - Products and Applications

PRODUCTS	APPLICATIONS						
	Biochemical Tests	Differential Agar	Enrichment	Highly Selective Media	Primary Differential Media	Serological Identification	Streak lactose negative colonies on
Adonitol/Adonitol CP <sup>†</sup>	-						
Alginate	-						
Arginine	+ or (+)						
Bismuth Sulfite Agar <sup>†</sup>				✓			
BG Sulfa Agar				✓			
Brilliant Green Agar <sup>†</sup>				✓			
Brilliant Green Agar Modified (Edel-Kampelmacher)				✓			
Desoxycholate Agar <sup>†</sup>		✓					
Desoxycholate Citrate Agar <sup>†</sup>				✓			
Dulcitol	(+) d						
EMB Agar <sup>†</sup>		✓					✓
Erythritol	-						
Esculin	-						
Gelatin	-						
Glucose	+						
GN Broth Hajna <sup>†</sup>			✓				
Hektoen Enteric Agar <sup>†</sup>		✓					
H <sub>2</sub> S	+(-)						
Individual O and H Antisera						✓	
Indole	-						
Inositol	d						
KCN	-						
Lactose	-						
Lysine	+						
Lysine Iron Agar <sup>†</sup>					✓		
MacConkey Agar <sup>†</sup>		✓					✓
MacConkey Agar CS		✓					
Methyl Red	+						

### Key

- Negative
- + Positive
- d Delayed
- (+) Variable
- (-) Variable

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Salmonella - Products and Applications

PRODUCTS	APPLICATIONS						
	Biochemical Tests	Differential Agar	Enrichment	Highly Selective Media	Primary Differential Media	Serological Identification	Streak lactose negative colonies on
MIO Medium					✓		
Muller Kauffmann Tetrathionate Broth Base				✓			
Nitrate Reduction	+						
Ornithine	+ or (+)						
Oxidase	-						
Phenylalanine	-						
Polyvalent Antisera						✓	
Raffinose	-						
Salicin	-						
SBG Enrichment			✓				
SBG Sulfa Enrichment			✓				
Selenite Broth			✓				
Selenite Cystine Broth <sup>†</sup>			✓				
SIM Medium <sup>†</sup>					✓		
Simmons Citrate <sup>†</sup>	(+) d						
Sodium Malonate	-						
SS Agar <sup>†</sup>				✓			
Sucrose	-						
m Tetrathionate Broth Base			✓				
Tetrathionate Broth Base <sup>†</sup>			✓				
Triple Sugar Iron Agar <sup>†</sup>					✓		
Universal Preenrichment Broth			✓				
Urea Agar Base <sup>†</sup>					✓		
Urease	-						
VP	-						
XLD Agar <sup>†</sup>		✓					
XLT4 Agar Base and XLT4 Supplement		✓					

## Key

- Negative
- + Positive
- d Delayed
- (+) Variable
- (-) Variable

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

## Shigella and Alkalescens-Dispar Group - Products and Applications

PRODUCTS	APPLICATIONS						
	Shigella	A-D Group	Enrichment	Biochemical Tests	Primary Differential Media	Primary Plating Media	Streak lactose negative colonies on
Acetate Differential Agar				✓			
Adonitol	-	-					
Arabinose	d	+					
Arginine Dihydrolase	d	d					
Christensen Citrate	-	d					
Desoxycholate Agar					✓		
Desoxycholate Citrate Agar					✓		
Dulcitol	d	d					
EMB Agar <sup>†</sup>					✓	✓	
Gas from Glucose	-(1)	-					
Gelatin (22°C)	-	-					
GN Broth Hajna <sup>†</sup>			✓				
Hektoen Enteric Agar <sup>†</sup>					✓		
H <sub>2</sub> S	-	-					
Indole	- or+	+					
Inositol	-	-					
KCN	-	-					
Lactose	-(1)	d					
Lysine Decarboxylase	-	d					
Lysine Iron Agar <sup>†</sup>				✓			
MacConkey Agar <sup>†</sup>					✓	✓	
MacConkey Agar CS					✓		
Malonate	-	-					
Maltose	d	+					
Mannitol	+ or-	+					
Methyl Red	+	+					
Motility	-	-					

### Key

- Negative
- + Positive
- d Different reactions

(1) Certain biotypes of *S. flexneri* produce gas; cultures of *S. sonnei* ferment lactose and sucrose slowly and decarboxylate ornithine.

(-) Variable

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Shigella and Alkalescens-Dispar Group - Products and Applications

PRODUCTS	APPLICATIONS					
	Shigella	A-D Group	Enrichment	Primary Differential Media	Primary Plating Media	Streak lactose negative colonies on
Mucate	-	d				
Ornithine Decarboxylase	d	d				
Phenylalanine	-	-				
Raffinose	d	d				
Rhamnose	d	d				
Salicin	-	-				
Simmons Citrate Agar <sup>†</sup>	-	-				
Sodium Acetate	-	+(+)				
Sucrose	-(1)	d				
Triple Sugar Iron Agar <sup>†</sup>				✓		
Urease	-	-				
Voges-Proskauer	-	-				
XLD Agar <sup>†</sup>					✓	
Xylose	d	+				

## Key

- Negative

+ Positive

d Different reactions

(1) Certain biotypes of *S. flexneri* produce gas; cultures of *S. sonnei* ferment lactose and sucrose slowly and decarboxylate ornithine.

(-) Variable

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

## Sterility Testing

### STERILITY TESTING, MANUAL

AC Broth	Sterility Bottles w/Screw Cap Tryptic Soy Broth
AC Medium	Sterility Bottles w/Septum Fluid A
Agar Medium No. F	Sterility Bottles w/Septum Fluid D
Brewer Thioglycollate Medium	Sterility Bottles w/Septum Fluid Thioglycollate Medium
Fluid Thioglycollate Medium <sup>†</sup>	Sterility Bottles w/Septum Tryptic Soy Broth
Fluid Thioglycollate Medium w/Beef Extract <sup>†</sup>	Thioglycollate Medium w/o Dextrose
Fluid Thioglycollate Medium w/K Agar	Thioglycollate Medium w/o Indicator
NIH Thioglycollate Medium	Thioglycollate Medium w/o Dextrose or Indicator
Sterility Bottles w/Screw Cap Fluid Thioglycollate Medium	Tryptic Soy Broth

See also: Environmental Sampling Section of the Reference Guide

### STERILITY TESTING, AUTOMATED

See ESP, Industrial Applications

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Veterinary Testing - Products and Applications

PRODUCTS	APPLICATIONS								
	Collection/Transport Bacterial	Collection/Transport, Viral	General Purpose Bacterial Culture Media	Brucella	Clostridium	Francisella	Gram Negative Enteric Bacteria (General)	Leptospira	Pseudomonas
Brain Heart Infusion†			✓						
Brain Heart Infusion Agar†			✓						
Brain Heart Infusion Agar and Clostridium Difficile Antimicrobial Supplement CC†					✓				
Brucella Agar†				✓					
Brucella Broth†				✓					
Brucella Antisera				✓					
Cellmatics™ Viral Transport Pack/ CULTURETTE™ Viral Single Swab†		✓							
Cetrimide Agar Base/PSEUDOSEL™ Agar†									✓
Cooked Meat Medium†					✓				
CULTURESWAB™ Amies Medium	✓								
CULTURESWAB Amies Medium w/o Charcoal/ CULTURETTE Amies w/o Charcoal†	✓								
CULTURESWAB Cary-Blair Medium/Anaero CULTURETTE Cary-Blair Single†	✓								
CULTURESWAB Perinasal Swab w/Amies Medium	✓								
CULTURESWAB Urethral Swab w/Amies Medium	✓								
CULTURESWAB Stuart's Medium Modified/ CULTURETTE Modified Stuart's Medium†	✓								
Cystine Heart Agar					✓				
EMB Agar/Eosin Methylene Blue Agar Modified†						✓			
Eugon Agar/EUGONAGAR™†				✓					
Eugon Broth/EUGONBROTH™†				✓					
Fletcher Medium Base								✓	
Fluid Thioglycollate Medium†		✓							
Francisella Tularensis Antigen					✓				
Francisella Tularensis Antiserum†					✓				
Leptospira Enrichment EMJH								✓	
Leptospira Medium Base EMJH								✓	
Levine EMB Agar/Levine Eosin Methylene Blue Agar†						✓			
MacConkey Agar						✓			
McClung Toabe Agar Base†					✓				
Pseudomonas Agar F/Flo Agar†									✓
Pseudomonas Agar P/Tech Agar†									✓
Pseudomonas Isolation Agar/Pseudomonas ISO†									✓
Reinforced Clostridium Medium/ Reinforced Clostridium Agar†					✓				

† Available from Difco & Becton Dickinson Microbiology Systems.

## Veterinary Testing - Products and Applications

PRODUCTS	APPLICATIONS								
	Collection/Transport Bacterial	Collection/Transport, Viral	General Purpose Bacterial Culture Media	Brucella	Clostridium	Francisella	Gram Negative Enteric Bacteria (General)	Leptospira	Pseudomonas
SFP Agar Base/TSN Agar <sup>†</sup>					✓				
Simmons Citrate Agar <sup>†</sup>							✓		
SPS Agar <sup>†</sup>					✓				
Sterile Swab/Aerobic Collection & Trans. System w/o Agar <sup>†</sup>	✓								
Triple Sugar Iron Agar <sup>†</sup>							✓		
Tryptic Soy Agar/ <b>TRYPTICASE</b> <sup>™</sup> Soy Agar <sup>†</sup>			✓						
Tryptic Soy Broth/ <b>TRYPTICASE</b> <sup>™</sup> Soy Broth <sup>†</sup>			✓						
Tryptose Phosphate Broth <sup>†</sup>			✓						

PRODUCTS	APPLICATIONS							
	Staphylococcus	Oxidation-Fermentation (O-F) Test	Mycological Media	Blood Culture	Salmonella/Shigella	Streptococcus	Vibrio	Susceptibility Testing
Azide Blood Agar Base <sup>†</sup>						✓		
Baird-Parker Agar Base w/ EY Tellurite Enrichment <sup>†</sup>	✓							
Bismuth Sulfite Agar <sup>†</sup>					✓			
Blood Culture Bottles Brain Heart Infusion w/ PAB, SPS + CO <sub>2</sub>					✓			
Blood Culture Bottles Columbia Broth w/CO <sub>2</sub>					✓			
Blood Culture Bottles Columbia Broth w/SPS + CO <sub>2</sub>					✓			
Blood Culture Bottles Fluid Thioglycollate Medium w/SPS + CO <sub>2</sub>					✓			
Blood Culture BottlesThiol Broth w/CO <sub>2</sub>					✓			
Blood Culture BottlesThiol Broth w/SPS + CO <sub>2</sub>					✓			
Blood Culture Bottles Tryptic Soy Broth w/CO <sub>2</sub>					✓			
Blood Culture Bottles Tryptic Soy Broth w/SPS + CO <sub>2</sub> / <b>SEPTIC-CHEK</b> <sup>™</sup> TSB w/ SPS + CO <sub>2</sub> <sup>†</sup>					✓			
Brain Heart Infusion Agar <sup>†</sup>			✓					
Brilliant Green Agar <sup>†</sup>						✓		
Coagulase Plasma (Rabbit)/ Coagulase Plasma, Rabbit <sup>†</sup>	✓							

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Veterinary Testing - Products and Applications

PRODUCTS	APPLICATIONS							
	Staphylococcus	Oxidation-Fermentation (O-F) Test	Mycological Media	Blood Culture	Salmonella/Shigella	Streptococcus	Vibrio	Susceptibility Testing
Coagulase Plasma EDTA (Rabbit)/ Coagulase Plasma, Rabbit w/EDTA†	✓							
Columbia CNA Agar†					✓			
Corn Meal Agar†			✓					
Decarboxylase Medium Base		✓						
Dispens-O-Disc Susceptibility Disks: Apramycin 15 mcg								✓
Ceftiofur 30 mcg								✓
Enrofloxacin 5 mcg								✓
Tilmicosin 15 mcg								✓
DNase Test Agar†	✓							
DNase Test Agar w/ Methyl Green/DNase Test Agar w/Toluidine Blue†	✓							
DTM Agar			✓					
Lysine Iron Agar†					✓			
Mannitol Salt Agar†	✓							
MIO Medium					✓			
MR-VP Medium/MR-VP Broth†					✓			
Mueller Hinton Medium/Mueller Hinton II Agar†								✓
Mycological Agar/MYCOPHII™ Agar†			✓					
OF Basal Medium†		✓						
Phenol Red Broth Base†		✓						
Phenylethanol Agar/Phenylethyl Alcohol Agar†	✓					✓		
Purple Agar Base		✓						
Purple Broth Base†		✓						
Sabouraud Dextrose Agar†			✓					
Sabouraud Dextrose Broth			✓					
Selenite Broth/Selenite F Broth					✓			
Selenite Cystine Broth†					✓			
SS Agar/Salmonella Shigella Agar†					✓			
Staph Latex Test/STAPHYLOSLIDE™ Test Kit†	✓							
Strep Grouping Kit†						✓		
Tetrathionate Broth Base†					✓			
TCBS Agar†							✓	
Urea Agar Base†			✓					
Vibrio Cholerae Antisera							✓	
VJ Agar (Vogel-Johnson Medium)†	✓							
XLD Agar†					✓			
XLT4 Agar Base and XLT4 Supplement					✓			

† Available from Difco & Becton Dickinson Microbiology Systems.

# Water/Wastewater Testing - Products and Applications

PRODUCTS	APPLICATIONS						
	Standard Plate Count	Total Coliforms	Fecal Coliforms / E. coli	Fecal Streptococcus	Salmonella	Stressed Organisms	Staphylococcus
A-1 Medium/A-1 Broth <sup>†</sup>			✓				
Azide Dextrose Broth <sup>†</sup>				✓			
BAGG Broth			✓				
Baird-Parker Agar Base w/ EY Tellurite Enrichment/Egg Yolk Tellurite Solution <sup>†</sup>							✓
Bile Esculin Agar <sup>†</sup>				✓			
Bile Esculin Azide Agar/ENTEROCOCCOSEL™ Agar <sup>†</sup>				✓			
Bismuth Sulfite Agar <sup>†</sup>					✓		
Brilliant Green Agar <sup>†</sup>					✓		
Brilliant Green Agar Modified (Edel-Kempelmacher)					✓		
Brilliant Green Bile 2%/Brilliant Green Bile Broth 2% <sup>†</sup>	✓						
m Brilliant Green Broth					✓		
Desoxycholate Lactose Agar <sup>†</sup>	✓						
m E Agar/M-E Agar Base <sup>†</sup>				✓			
EC Medium/EC Broth <sup>†</sup>			✓				
EC Medium with MUG/EC Broth w/ MUG <sup>†</sup>			✓				
m Endo Agar LES/m Endo Agar LES <sup>†</sup>	✓						
m Endo Broth MF®/m Endo Broth <sup>†</sup>	✓						
m Enterococcus Agar/m Enterococcus Agar <sup>†</sup>				✓			
Esculin Iron Agar				✓			
EVA Broth/Ethyl Violet Azide Broth <sup>†</sup>				✓			
m FC Agar/ m FC Agar <sup>†</sup>			✓				
m FC Basal Medium			✓				
m FC Broth Base/m FC Broth <sup>†</sup>			✓				
m HPC Agar/m HPC Agar Base <sup>†</sup>	✓						
Lactose Peptone Broth		✓					
Lauryl Tryptose Broth/Lauryl Sulfate Broth <sup>†</sup>		✓					
Lauryl Tryptose Broth with MUG/ Lauryl Sulfate Broth w/ MUG <sup>†</sup>		✓					
Levine EMB Agar/Levine Eosin Methylene Blue Agar <sup>†</sup>		✓					
Minerals Modified Glutamate Broth		✓					
Muller Kauffmann Tetrathionate Broth Base					✓		

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.

# Water/Wastewater Testing - Products and Applications

PRODUCTS	APPLICATIONS						
	Standard Plate Count	Total Coliforms	Fecal Coliforms / E. coli	Fecal Streptococcus	Salmonella	Stressed Organisms	Staphylococcus
Nutrient Agar w/ MUG <sup>†</sup>			✓				
Plate Count Agar/Standard Methods Agar <sup>†</sup>	✓						
Presence-Absence Broth <sup>†</sup>		✓					
R2A Agar <sup>†</sup>						✓	
Rosolic Acid			✓				
Selenite Broth/Selenite F Broth <sup>†</sup>					✓		
Selenite Cystine Broth <sup>†</sup>					✓		
m Staphylococcus Broth							✓
m T7 Agar						✓	
m TEC Agar						✓	
m Tetrathionate Broth Base					✓		
Tetrathionate Broth Base <sup>†</sup>					✓		
Tryptone Water			✓				
XLD Agar <sup>†</sup>					✓		

<sup>†</sup> Available from Difco & Becton Dickinson Microbiology Systems.